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TABLE OF CONTENTS

SECTION	PAGE
1. DEVICE DESCRIPTION	18
2. INDICATIONS FOR USE	18
3. CONTRAINDICATIONS	19
4. WARNINGS AND PRECAUTIONS	19
5. HOW THE PRODUCT IS SUPPLIED	19
6. DEVICE SELECTION	20
7. CASE PRE-PLANNING/ INDIVIDUALIZATION OF TREATMENT	21
8. EQUIPMENT REQUIREMENTS	21
9. IMPLANTATION PROCEDURE	21
10. MAGNETIC RESONANCE SAFETY INFORMATION	23
11. ADVERSE EVENTS	24
12. AVAILABLE DEVICE CONFIGURATIONS	25
13. SYMBOLS AND DEFINITIONS	25

1. DEVICE DESCRIPTION

The RelayPro Thoracic Stent-Graft System (NBS Configuration) is an endovascular device intended to treat thoracic aortic pathologies. Once placed in the aorta, the RelayPro Stent-Graft (NBS Configuration) provides an alternative conduit for blood flow while excluding the lesion. The system consists of an implantable stent-graft and delivery system.

1.1. Stent-Graft

The RelayPro Stent-Graft (NBS Configuration) is comprised of self-expanding nitinol stents sutured to polyester graft fabric. The stent scaffold is a series of serpentine stents stacked in a tubular configuration and spaced along the length of the graft fabric. Longitudinal support for the stent-graft is provided by a curved nitinol wire called the spiral support strut.

The RelayPro Stent-Graft (NBS Configuration) has one proximal end configuration. The proximal end configuration consists of a sinusoidal nitinol stent plus a crown-shaped nitinol stent that are both covered with fabric. These covered wires, when in position, expand to the vessel wall anchoring the device in place and aid in creating a seal zone oriented to the vasculature.

The RelayPro Stent-Graft (NBS Configuration) offers a single distal configuration. The straight end configuration consists of nitinol stents placed evenly over the fabric around the circumference of the stent-graft.

All stent-grafts have platinum/iridium radiopaque marker bands that indicate the fabric edge and also serve as guides for positioning the spiral support strut.

The RelayPro Stent-Graft (NBS Configuration) does not contain any natural rubber latex; however, it may incidentally come into contact with latex during the manufacturing process. Table 1 lists the RelayPro Stent-Graft materials.

1.2. Delivery System

The RelayPro Delivery System consists of a series of coaxially-arranged sheaths and catheters, along with a tubular handle control system. The tapered tip and introducer sheath have a lubricious hydrophilic coating. The delivery system is provided in outer diameters ranging from 19 to 23 French, depending on the corresponding stent-graft diameter, with a working length of 90 cm. The delivery system is intended for single use and is disposable.

The delivery procedure consists of two stages. The first stage consists of a hydrophilically coated introducer (outer sheath), which is used to advance and track over a guidewire through the access vessels to the distal intended landing zone. Within the first sheath is the second sheath. During the second stage this flexible sheath maintains the compressed stent-graft. The flexibility of the second sheath facilitates tracking through tortuous and curved portions of the thoracic aorta. The larger diameters of the non-bare stent configurations have a delivery system that has two heat shaped nitinol wires, called support wires, attached to the delivery system catheter. The distal end of the support wires present atraumatic tips that are tethered to the inferior proximal portion of the graft. The support wires are designed to control the expansion of the inferior portion of the graft and ensure proper apposition against the anatomical inner curvature. These support wires are only present in devices with stent-grafts proximal ends that are 32 mm in diameter and greater.

2. INDICATIONS FOR USE

The RelayPro Thoracic Stent-Graft System (NBS Configuration) is intended for the treatment of thoracic aortic pathologies such as aneurysms, pseudoaneurysms, dissections, penetrating ulcers and intramural hematoma, in adult patients (as defined by local statutes). Anatomical/sizing features specified in Tables 1, 2 and 3 should be observed.

Table 1: The RelayPro Stent-Graft (NBS Configuration) Materials

Component	Material
Stents	Nitinol wire
Spiral Support Strut	Nitinol wire
Graft fabric	Polyester
Suture	PTFE-impregnated polyester
Radiopaque markers	Platinum-iridium wire

3. CONTRAINDICATIONS

The RelayPro Thoracic Stent-Graft System (NBS Configuration) is contraindicated when patients present with any of the following characteristics/conditions:

- Pregnancy / lactation
- Aneurysm / lesion location not accessible to the delivery system and stent placement
- Arterial access size insufficient for delivery system entrance
- Treatment of lesion that would require a delivery system with usable length greater than 90 cm
- Excessive arterial disease precluding delivery system entrance or passages
- Untreatable allergy or history of allergic reaction to radiographic contrast media (X ray dye)
- Untreatable allergy or history of allergic reaction to anticoagulants
- Systemic infection
- Arterial tortuosity not allowing passage of the delivery system
- Arterial or aneurysm / lesion size incompatible with stent-graft
- Congenital connective tissue disease
- Mycotic aneurysm / lesions
- Aortic inner diameter that cannot accommodate the expanded inner sheath outer diameter of approximately 10 mm
- Hypersensitivity to polyester or nitinol
- Massive thrombus
- Bleeding diathesis
- Aortic angulations (radius) less than 15 mm at intended proximal landing zone

4. WARNINGS AND PRECAUTIONS

- Placement of the RelayPro Stent-Graft (NBS Configuration) in the thoracic aorta often requires proximity to the great vessels perfusing the brain, increasing the possibility of thrombus or embolization proximally. Care should be taken to ensure air has been purged from the system before its introduction into the artery.
- Proximal and distal required landing zones vary with stent-graft size.
- Excessive aortic tortuosity may result in inability to properly position the stent-graft, or stent-graft kinking with thrombus formation.
- If balloon modelling is desired, use a compliant balloon equal in size to the largest target vessel's diameter. Balloon inflation should not exceed 1 atm.

- Do not use power / pressure injections through the delivery system.
- Do not use the system if defects are noted.
- Care should be taken with respect to occlusion of intercostals / spinal cord arteries.
- Care should be taken when treating morbidly obese patients, as visualization of imaging may be comprised.
- Careful consideration should be given to treating patients where tracking through a previously placed endovascular or surgical prosthesis is required.
- Given the contraindication for pregnant women, care should be taken when treating women of childbearing potential.
- Consideration should be given to treating patients with significant comorbidities / high risk for open surgical repair.
- The minimum recommended amount of overlap between devices is three overlapping covered stents (approximately 50 mm). Less than this amount of overlap may result in endoleak (with or without component separation). Device lengths should be selected accordingly.

5. HOW THE PRODUCT IS SUPPLIED

The RelayPro Thoracic Stent-Graft System (NBS Configuration) is provided STERILE.

DO NOT RESTERILIZE. FOR SINGLE USE ONLY.

Re-sterilization of the device for re-use will result in loss of component integrity (e.g., reduction of stent-graft radial force, component cracking or discoloration, etc.).

Do not reuse, reprocess, clean/disinfect, or re-sterilize. Reuse, reprocessing, cleaning/disinfection, or re-sterilization may compromise the structural integrity of the device and/or lead to device failure which, in turn, may result in deterioration of health or death of patients or users. Reuse, reprocessing or re-sterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the end-user. Also, each single-use device carries specific labeling instructions relative to storage, use and handling to minimize exposure to conditions, which could compromise the product, the patient or the user. These conditions cannot be assured once the packaging is opened and discarded.

5. HOW THE PRODUCT IS SUPPLIED (cont.)

The RelayPro Stent-Grafts (NBS Configuration) are available in the following approximate lengths:

90 mm	190 mm
100 mm	200 mm
110 mm	210 mm
130 mm	230 mm
150 mm	250 mm
170 mm	

All stent-grafts are available in 2 mm incremental diameters ranging from 22 mm to 46 mm. Product sizes can be found in Section 12: AVAILABLE DEVICE CONFIGURATIONS

The product is supplied with the following model designation identified on the label. As an example, reference number 28N4-46-259-46S can be decoded as follows:

Product ID	Device Type	Design Modification Number	Stent Diameter Proximal (mm)	Stent Covered Length (mm)	Stent Diameter Distal (mm)	Device Designation
28	N: NBS	4	46	259	46	S: CE-Marked Product

6. DEVICE SELECTION

Table 2 addresses the recommended healthy landing zone length depending on the stent-graft diameter selected. Tables 3 addresses the selection of the appropriate stent-graft diameters for the RelayPro Stent-Graft (NBS Configuration) based on the vessel size.

Table 2: Target Landing Zone

Stent-Graft Diameter (mm)	Proximal Length (mm)	Stent-Graft Diameter (mm)	Distal Length (mm)
22 – 28	22	22 – 38	25
30 – 38	25	40 – 46	30
40 – 46	30		

Table 3: Standard Stent-Graft Length and Oversizing Recommendations

Stent-Graft Size (mm)	Thoracic Vessel Size (mm)	% Graft Oversizing	Delivery System French Size (O.D.)
22	19	16%	19
24	20-21	14-20%	19
26	22-23	13-18%	19
28	24-25	12-17%	20
30	26-27	11-15%	20
32	28-29	10-14%	21
34	30-31	9-13%	21
36	32-33	9-12%	22
38	34	11%	22
40	35-36	11-14%	22
42	37-38	10-13%	23
44	39-40	10-13%	23
46	41-42	9-12%	23

Note: For tapered devices, the delivery system sheath size is driven by the largest diameter of the stent-graft. Oversizing recommendations are the same for proximal and distal diameters. Product sizes can be found in Section 12: Available Device Configurations

7. CASE PRE-PLANNING/ INDIVIDUALIZATION OF TREATMENT

Practitioners using the RelayPro Thoracic Stent-Graft System (NBS Configuration) should have a thorough understanding of endovascular procedures and techniques. In particular, the RelayPro System should only be used by physicians and teams with experience and training in vascular interventional techniques, including, but not limited to, training on the use of RelayPro. This will include practitioners with formal education/training in vascular surgery, interventional radiology, cardiothoracic surgery and interventional cardiology. Selecting the proper stent-graft with the appropriate length and diameter is paramount to the successful exclusion of the indicated thoracic aortic pathologies. Therefore, measure all parameters needed for proper sizing of the stent-graft carefully. Bolton Medical recommends evaluation of all imaging studies available (e.g., angiograms, CT scans, MRI scans, MRA scans, and plain radiographs). Each imaging modality offers additional information to the sizing process. The physical characteristics of the vessel should be evaluated in addition to its size. Factors such as stenosis, atherosclerotic disease, ectasia and tortuosity may affect stent-graft selection and placement strategy. Potential vessel's morphology changes between pre-op CT and time of the procedure should be considered. The final stent-graft selection will be the responsibility of the physician.

Practitioner must ensure that the access vessel diameter is compatible with the selected delivery system's Outer Sheath French size. The physician and patient (and/or family) should review the risks and benefits when discussing this endovascular device and the need for compliance with follow-ups. Any pertinent actions to be avoided or precautions to be taken should also be discussed.

8. EQUIPMENT REQUIREMENTS

Fluoroscopic equipment, including a high-resolution image intensifier on a freely angled C-arm (which can be ceiling, pedestal-mounted or portable), will be needed for the procedure. Preferably, the image intensifier will have a complete range of motion to achieve AP projections to lateral projections.

Its capabilities should include:

- Digital subtraction angiography
- High-resolution angiography
- Roadmapping

Supportive/supplementary equipment:

- 0.035" (0.89 mm) by 300 cm Meier guidewire
- 0.035" (0.89 mm) by 260 cm or 300 cm Lunderquist guidewire

- Guidewire torque device
- Inflation device with pressure gauge
- Compliant stent-graft modeling balloons of the appropriate size
- Arterial puncture needles 18 G or 19 G
- Nitinol retrieval device/snare (10-15 mm diameter)
- Assortment of vascular stents
- Assortment of angiographic and graduated catheters

9. IMPLANTATION PROCEDURE

PREPARATION

(Procedural Steps 1-13)

Position the patient on the surgical table where standard aseptic preparation of the surgical site is conducted. Drape the patient with sterile surgical drapes leaving exposed the bilateral groin access sites.

Anticoagulation and antiplatelet therapies are performed at the discretion of the physician. Similarly, blood pressure adjustment and spinal cord protection measures are also at the discretion of the physician.

1. Verify devices are correct for the patient.
2. Access may be obtained either via an arteriotomy of the access artery or through percutaneous approach at the access site that will be used to introduce the device.
3. Determine a suitable secondary access site for diagnostic and imaging purposes.
4. Under fluoroscopy, advance a 0.035 in (0.89 mm) guidewire from the access site and a graduated angiographic catheter via the secondary access site.
5. Place the C-Arm DSA system into a left anterior oblique position in preparation for the initial angiogram. Perform an angiogram to confirm the preoperative case planning and mark the target area.
6. Inspect the system packaging for visible tears, breaks, or openings.
7. Take the Delivery System from the sterile packaging and bring it to the surgical table. Examine the delivery system for structural integrity. DO NOT USE the system if defects are noted.
8. Ensure that the Controller is in the "1" position. If it is not change it to the "1" position to prevent premature deployment of the stent-graft.
9. Verify that the delivery system Tip is properly seated in the outer sheath and that the tip side hole is not covered. If not, correct by rotating the deployment grip until the delivery system tip is properly seated.
10. Flush the delivery system with heparinized saline through the flush port to purge air

9. IMPLANTATION PROCEDURE (cont.)

from the coaxially situated sheaths. Ensure that a continuous stream of saline exits the tip side hole. It may be necessary to elevate the distal end of the system to different positions to bring air to the highest point for purging. The flush port valve must be closed under pressure to prevent air from re-entering the system. Visually inspect it for remaining air and repeat if necessary.

11. Verify that the Apex Holder Knob [ITEM 13] is securely engaged in the V-shaped notch of the Guidewire Luer [ITEM 14] (see Detail A from Figure 1). DO NOT USE the system if the Apex Holder Knob is not engaged in the V-shaped notch provided by the Guidewire Luer. Also, do not attempt to re-engage the Apex Holder Knob with the Guidewire Luer.
12. Flush the delivery system with heparinized saline through the gray guidewire luer. Attention should be paid to not rotate the apex holder knob during this step (Figure 2).
13. Activate the hydrophilic coating by wetting the tip and the outer sheath with saline.

INTRODUCTION/ADVANCEMENT OF THE OUTER SHEATH (Procedural Steps 14-15)

14. Advance the outer sheath into the artery over the guidewire. The guidewire should always remain in the delivery system while inside the patient.
15. Under fluoroscopic control, advance the outer sheath until the delivery system tip is just below the intended distal landing zone. If the descending aorta presents tight tortuosity, the tip should be advanced past the tight curvature(s) to facilitate navigation of the inner sheath.

NOTE: Do not advance the outer sheath into the thoracic arch.

NOTE: If the outer sheath cannot be advanced beyond the region of tight curvatures the delivery system should be removed from the patient and an alternate procedure should be considered.

ADVANCEMENT OF THE INNER SHEATH (Procedural Steps 16-21)

Once the inner sheath is advanced, the user will be committed to implant the stent-graft.

CAUTION: The controller must be in the “1” position.

NOTE: The handle body can be rotated to position the top of the device facing the user. The handle body should not be rotated more than half a turn.

16. While holding the gray grip so that the Main Body remains stationary, rotate the deployment grip clockwise to advance the inner sheath. Ensure that while the deployment grip is being turned, it is being pushed toward the gray stationary grip (Figure 3).

If preferred, the mechanical advantage can be bypassed by pressing on the “disengagement button” while advancing the deployment grip (Figure 4).

17. Advance the inner sheath until the stent-graft proximal markers reach the proximal landing zone. Upon reaching the intended proximal landing zone, visually confirm that the stent-graft’s distal marker bands can be seen approximately 2 cm outside of the outer sheath. If the stent-graft’s distal marker bands do not appear to have exited the outer sheath, while in position 1, press the disengagement button and hold the deployment grip in place while pulling back on the gray stationary grip until the stent-graft’s distal marker bands have exited the outer sheath by approximately 2 cm. The arrow mark on the main handle body should be used to approximate when the inner sheath and therefore the stent-graft have completely exited the outer sheath.
18. As the inner sheath is advanced out of the outer sheath, note the alignment of the spiral support strut by locating the spiral support strut markers under fluoroscopy.
19. If the device is to be implanted in a curved section of the aorta, verify that the D-shaped marker on the inner sheath and the spiral support strut marker(s) face the greatest curvature.

If radial adjustment is needed, retract the deployment grip (disengage mechanical advantage by pressing on the disengagement button) to bring the stent-graft to a straight portion of the vessel. When retracting the deployment grip, ensure that the distal end of the stent-graft is not pulled into the outer sheath (the arrow marker can be used as a reference). It may be necessary to retract the whole device a few centimeters to bring the stent-graft to a straight position. After the stent-graft is in the straight position, while holding the front nose cap, rotate the whole handle body to manually align the spiral support strut markers toward the greatest curvature of the aorta (Figure 5). The D-shaped marker can be used to aid in this placement. If the round portion of the D-shaped marker is facing the greater curvature, the handle body should be turned clockwise. If the round portion is facing the lesser curvature the turn should be counterclockwise. One to three handle revolutions may be required before the

9. IMPLANTATION PROCEDURE (cont.)

stent-graft begins rotating. Once alignment is confirmed, re-advance the stent graft into the desired position.

20. Perform an angiogram of the area of interest to confirm proper position of the device in preparation for deployment.
21. Finalize the longitudinal placement of the stent-graft in relation to the proximal landing zone by adjusting the deployment grip as necessary. Confirm the position of the proximal and distal marker bands as well as the spiral support strut markers.

DEPLOYMENT OF STENT-GRAFT (Procedural Steps 22-26)

22. With the stent-graft in the desired deployment position, turn the controller to the "2" position (Figure 6).
23. While holding the gray grip fixed, rotate the deployment grip counterclockwise to pull down the inner sheath and expose the first covered stent. Stop as soon as you see the D-shaped marker is aligned with the spiral support strut marker.

NOTE: The inner sheath has the D-Shaped radiopaque marker located near the tip that can be used to visualize its movement under fluoroscopy.

CAUTION: DO NOT pass the overlapping of the two markers with the rotation of the deployment grip since it could lead to distal migration of the device or to retroflex effect.

24. Make any final linear position adjustments (proximally or distally).
25. Deploy the stent-graft by holding the gray grip fixed and, while pressing the disengagement button down, retract the deployment grip with one continuous motion without stopping until the stent-graft is fully deployed and the inner sheath completely retracted.

CAUTION: Failure to retract the deployment grip with one continuous motion without stopping may cause blood pressure to increase and may result in distal migration of the device during deployment.

CAUTION: Attachment of the proximal end of the stent graft is maintained by the Apex Holder. As such, the handle body should not be moved until the apex holder is released.

26. Under fluoroscopic control, release the proximal stent clasp by rotating the apex holder knob and sliding it toward the guidewire luer. The stent graft is now fully released (Figure 7).

SYSTEM REMOVAL

(Procedural Steps 27-34)

27. Place the controller in the "4" position (Figure 8).
28. Under fluoroscopic control, retract the stainless steel rod allowing the tip to rejoin the outer sheath (Figure 9). Monitor the travel of the delivery system tip through the deployed stent-graft to ensure that the stent-graft's position is not affected. If the tip does not rejoin easily, apply slightly greater force until the tip rejoins with the outer sheath.
29. Withdraw the entire system from the patient.
30. Perform a final angiogram to assess for any endoleaks and/or migration. Confirm successful aneurysm/lesion exclusion.
31. If an endoleak is detected, consider balloon modeling to correct the leak.

CAUTION: Do not exceed 1 atm. balloon pressure. Balloon only within the covered portion of the stent-graft. Ballooning outside the covered portion could cause aortic rupture, atherosclerotic plaque embolization, or other complications. Always recheck position of the stent-graft following ballooning.
32. Straighten the angiographic pigtail catheter and remove the catheter and sheath from access site.
33. Perform standard closure of the access site.
34. Assess blood flow to the distal extremities.

10. MAGNETIC RESONANCE SAFETY INFORMATION

The RelayPro Stent-Graft (NBS Configuration) was determined to be MR-conditional. Specifically, when present in a patient undergoing a magnetic resonance imaging (MRI) procedure at 3-Tesla or less will not create an additional hazard or risk to the patient under the conditions used for the testing.

Magnetic resonance imaging (MRI) procedures must be performed according to the following guidelines:

- A patient with the RelayPro Stent-Graft (NBS Configuration) may safely undergo an MRI procedure using an MR system with a static magnetic field of 3-Tesla or less. The RelayPro Stent-Graft does not exhibit magnetic field interactions with respect to translational attraction (tested at a maximum spatial gradient, 3.3-Tesla/meter) and showed no torque during exposure to a 3-Tesla MR system. Therefore, there is no additional risk to a patient with the RelayPro Stent-Graft (NBS configuration)

10. MAGNETIC RESONANCE SAFETY INFORMATION (cont.)

with respect to movement or dislodgement using an MR system with a static magnetic field of 3-Tesla or less. In addition, because of the lack of magnetic field interactions at 3-Tesla, an MRI procedure may be performed on a patient immediately after implantation of the RelayPro Stent-Graft.

- MRI procedures must not exceed exposures to radiofrequency (RF) fields greater than a whole-body averaged specific absorption

rate (SAR) of 2.0-W/kg for 15 minutes at 3-Tesla in a patient with the RelayPro Stent-Graft (NBS Configuration).

- Safety information for magnetic (MRI) procedures pertain to use of MR systems with static magnetic fields of 3-Tesla or less (maximum spatial gradient 3.3-Tesla/m) and a whole body averaged specific absorption rate (SAR) of 2.0-W/kg for 15 minutes of MR imaging. The effects of performing MRI procedures using MR systems with static magnetic fields greater than 3-Tesla and other conditions have not been determined.

11. ADVERSE EVENTS

Adverse events that may occur in conjunction with endovascular procedures include, but are not limited to those listed below.

Fever	Transient Ischemic Attack
Post implantation syndrome	Cerebral vascular accident (stroke)
Hematoma	Congestive heart failure
Blood Loss	Paralysis/Paresthesia/Paraparesis
Hemorrhage	Death
Infection	Incision site complications
Cardiac events	Limb ischemia
Anaphylaxis	Delivery system failure
Vessel dissection	Access failure
Vessel occlusion/Thrombosis	Stent graft misplacement
Vessel damage	Endoleak
Emboli	Stent graft migration
Hepatic failure	Stent graft failure
Ischemia (spinal cord, perfusion pathways)	Wire form fractures
Renal failure or complications	Suture fracture
Aneurysm/Lesion rupture	Deployment failure
Arteriovenous fistula/Aorto-esophageal fistula	Perforation
Radiation overexposure or reaction	Device dehiscence
Pseudoaneurysm	Sten-graft tearing/wear

In the event of surgical removal or postmortem examination, please contact Bolton Medical for guidance on removal and disposal of implant.


12. AVAILABLE DEVICE CONFIGURATIONS

All available configurations are listed on the table beginning on the next page. All device configurations are listed with the following format:

Internal Bolton #	Internal Bolton #	Proximal Diameter (mm)	Family Length (mm)*	Distal Diameter (mm)	Internal Bolton	French Size
28	N4	PD	LEN	DD	S	(XX Fr)

*Family lengths listed. Final lengths to be listed on product labeling (Tolerance of ± 10 mm)

13. SYMBOLS AND DEFINITIONS

	Consult Instructions For Use	REF	Catalogue Number
	Do Not Reuse	LOT	Batch code
	Do Not Resterilize		Do Not Use If Package Is Damaged
STERILE R	Sterilized using Irradiation		Use By
	Manufacturer		MR-Conditional
	Date of Manufacture		Keep dry
	Temperature limit		Caution
NBS	Non-Bare Stent Configuration	Bare Stent	Bare Stent Configuration
EC REP	European Community Authorized Representative		

Bolton Medical, Inc.
799 International Parkway
Sunrise, FL 33325 U.S.A.
Tel: +1 954 838-9699

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