RG34M121-01

2013 - 04

 **TERUMO[®]**

GUIDE WIRE

GT

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TERUMO CORPORATION

44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO
151-0072, JAPAN MADE IN JAPAN

TERUMO EUROPE N.V.

INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM
www.terumo-europe.com

Terumo Australia Pty Ltd

Level 4 Building B 11 Talavera Rd
Macquarie Park NSW 2113 Australia
T: +61 2 9878 5122 E: cs_australia@terumo.co.jp

©TERUMO CORPORATION 2013-04

TERUMO®
GUIDE WIRE
GT
RADIFOCUS® GUIDE WIRE GT with Gold Coil

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TERUMO CORPORATION

44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO
151-0072, JAPAN MADE IN JAPAN

TERUMO EUROPE N.V.

INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM
www.terumo-europe.com

Terumo Australia Pty Ltd

Level 4 Building B 11 Talavera Rd
Macquarie Park NSW 2113 Australia
T: +61 2 9878 5122 E: cs_australia@terumo.co.jp

ENGLISH

RADIFOCUS® GUIDE WIRE GT with Gold Coil **Guide wire with Gold Coil**

INSTRUCTION FOR USE

Read the following warnings, precautions and directions for use carefully.

INDICATION FOR USE

The Radifocus GUIDE WIRE GT is designed to direct a catheter to the desired anatomical location in the vascular system during diagnostic or interventional procedures.

CONTRAINDICATIONS

No-contraindications are known when using guide wires in vascular procedures. Carefully read the list of Contraindications in the instructions for use accompanying the other interventional devices to be used with the guide wire.

COMPLICATIONS

Complications that may be related to the use of guide wire can include but may not be limited to the following:

- Tissue trauma
- Infection
- Vessel perforation
- Hematoma
- Embolism
- Allergic reaction
- Thrombus formation
- Hemorrhage

Carefully read the list of Complications in the instructions for use accompanying the other interventional devices to be used with the guide wire.

WARNINGS

Failure to abide by the following warnings might result in damage to the vessel, shearing of the GUIDE WIRE GT, and release of plastic fragments from the GUIDE WIRE GT. Such pieces or fragments from the wire may have to be removed from the vessel.

- Do not manipulate or withdraw the GUIDE WIRE GT through a metal entry needle, a metal dilator, or a metal guide wire inserter. Manipulation and/or withdrawal through a metal entry needle, a metal dilator, or a metal guide wire inserter may result in destruction and/or separation of the outer polyurethane coating requiring retrieval. A plastic entry needle is recommended when using this wire for initial placement.
- Do not use the GUIDE WIRE GT with devices which contain metal parts such as atherectomy catheters, laser catheter, or metal introduction devices as they may cause the GUIDE WIRE GT plastic coating to shear and/or sever the wire.
- Do not reshape the GUIDE WIRE GT by any means. Attempting to reshape the wire may cause damage, resulting in the release of wire fragments into the vessel.
- When exchanging or withdrawing a catheter over the GUIDE WIRE GT, secure and maintain the guide wire in place under fluoroscopy to avoid unexpected guide wire advancement, otherwise damage to the vessel wall by the wire's tip may occur.
- A retrieving device, such as a gripper or basket forceps, can only be used after the GUIDE WIRE GT has been removed from the patient's vessel. Using a retrieving device while the GUIDE WIRE GT is in the vessel may cause the GUIDE WIRE GT to break.
- Manipulate the GUIDE WIRE GT slowly and carefully in the vessel while confirming the behaviour and location of the wire's tip under fluoroscopy. Improper manipulation of the GUIDE WIRE GT without fluoroscopic confirmation may result in vessel perforation.
- Do not apply repetitive bending force to one specific point of the device as this may cause damage to the GUIDE WIRE GT.
- If any resistance is felt or if the tip's behaviour and/or location seems improper, stop manipulating the GUID WIRE GT and/or the catheter and determine the cause by fluoroscopy. Continuing to manipulate or rotate the GUIDE WIRE GT or failure to exercise proper caution may result in bending, kinking, separation of the guide wire's tip, damage to the catheter, or damage to the vessel.
- Do not attempt to use the GUIDE WIRE GT if it has been bent, kinked or damaged. Use of a damaged wire may result in damage to the vessel or the release of wire fragments into the vessel.
- Consider the use of systemic heparinization to prevent or reduce the possibility of thrombus formation on the surface of the GUIDE WIRE GT.

PRECAUTIONS

- The GUIDE WIRE GT should be used by a physician, who is well trained in manipulation and observation of guide wires under fluoroscopy.
- Sterile in an unopened and undamaged unit package. Do not use if the unit package or the guide wire is broken or soiled. The GUIDE WIRE GT should be used immediately after opening the package and be disposed of safely and properly after use, following local regulations for medical waste management.
- When using a drug or a device concurrently with the GUIDE WIRE GT, the operator should have a full understanding of the properties/characteristics of the drug or device so as to avoid damage to the GUIDE WIRE GT. For example when using the GUIDE WIRE GT with any device that emits energy (laser, pressure, ultrasound, etc.) confirm that the GUIDE WIRE GT is retracted into a position where it will not be impacted by the energy.
- Consider the use of systemic heparinization.
- The surface of the GUIDE WIRE GT is not lubricous unless it is wet. Before taking it out of its holder and inserting it through a catheter, fill the holder and the catheter with heparinised physiological saline solution.
- When reinserting the GUIDE WIRE GT back into the holder, take care not to damage the wire's hydrophilic polymer coating with the edge of the holder.
- Use the enclosed RADIFOCUS TORQUE DEVICE to handle the GUIDE WIRE GT. Do not use a metal torque device with the GUIDE WIRE GT. Use of a metal torque device may result in damage to the wire.
- Do not slip a tightened up torque device or Y-connector over the wire, as this may result in damage to the wire.
- Due to variations of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the catheter, it is advisable to stop using such catheters.
- Do not manipulate the GUIDE WIRE GT through a tightened up rotating haemostasis valve, as this may result in damage to the wire.
- After removal from the patient's vessel, and prior to reinserting it into the same patient during the same catheterisation, the GUIDE WIRE GT should be rinsed in a bowl full of heparinised physiological saline solution. Any blood residues still adhering to the wire can be removed by wiping once with a gauze moistened with heparinised physiological saline solution. Use of alcohol, antiseptic solutions or other solvents must be avoided, because they may adversely affect the surface of the GUIDE WIRE GT.
- The GUIDE WIRE GT contains a metallic core, do not use with any inappropriate equipment (e.g. MRI).
- The entire operation should be carried out aseptically.
- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Avoid exposure to water, direct sunlight, extreme temperatures and high humidity during storage.

DIRECTIONS FOR USE

1. Remove the GUIDE WIRE GT and the holder together from the tray.
2. Fill the holder with heparinised physiological saline solution through the hub of the holder using a syringe.
3. Remove the GUIDE WIRE GT from the holder and inspect the GUIDE WIRE GT prior to use, to verify if it is lubricated. If the GUIDE WIRE GT can not be easily removed from the holder, inject more heparinised physiological saline solution into the holder and try again.
4. Prior to use prime the catheter with heparinised physiological saline solution to ensure smooth movement of the GUIDE WIRE GT within it.
5. The GUIDE WIRE GT may slide entirely into the catheter or slip out of the catheter because of its low sliding friction.
6. Keep at least 5 cm of the wire extended out of the hub of the catheter during introduction.
7. A RADIFOCUS TORQUE DEVICE, which is included in the package, is recommended for easier handling/manipulation of the wire.

FRANÇAIS / FRENCH

RADIFOCUS® GUIDE WIRE GT with Gold Coil

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