
# IMUGARD® III S-PL

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

© TERUMO CORPORATION November 2009 09K04

## ENGLISH

### IMUGARD<sup>®</sup> III S-PL For Laboratory use, with platelet storage bag

#### DESCRIPTION/DESIGNATION

- This device is intended to obtain/store an equivalent of 10 leucodepleted platelet concentrates obtained from buffy coats, platelet rich plasma or single donor platelets (apheresis).
- The integrated platelet storage bag is suitable for storage of platelets at 22±2°C up to 5 days on a platelet agitator, if prepared in a closed system.
- Sterilised by ethylene oxide. Sterile in an unopened and undamaged unit package.
- Non-pyrogenic fluid path.

#### PRECAUTIONS

- This device contains DEHP. Based on animal studies, significant exposure to DEHP may interfere with the normal development of the male reproductive tract. For children, pregnant and nursing women, alternative devices may be appropriate.
- Dispose safely following the local procedures for the disposal of medical waste.
- Do not use if the unit package or the product has been damaged or soiled.
- Do not store at extreme temperatures and humidity. Avoid direct sunlight.
- Read the following instructions carefully before use.
- For single use only. Do not reuse. Do not resterilize. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.

#### INSTRUCTIONS FOR USE

**WARNING** Open system processing does not allow preservation of platelets for more than 24 hours. USE CLOSED SYSTEM PROCESSING.

- CAUTION**
  - A strict aseptic technique should be respected at all times. Do not vent.
  - If an automated component extractor is used, refer to the manufacturer's instructions.

During and after filtration, do not squeeze or apply pressure on the filter (to prevent this, Terumo is offering the filter holder ref. ME-05TA69) until the filter is separated from the platelet storage bag.

#### 1. Attaching the unit to be filtered

- Open the packaging and close all clamps on the set.
- Use a sterile connecting device to connect the filter with the unit to be filtered.

**CAUTION** Follow the instructions for use given by the manufacturer of the sterile connecting device.

#### 2. Filtration

**CAUTION** In case of the pooled buffy-coats, proceed the paragraph << Instructions for pooled buffy-coats >>. In case of the pooled platelet concentrates or the single donor platelets, proceed the paragraph << Instructions for pooled platelet concentrates / single donor platelets >>.

<< Instructions for pooled buffy-coats >>

- Hang the unit to be filtered on an I.V. stand.
- Position all the clamps in the upper level of the cup and avoid direct contact with the filled bag.
- In case the sampling bag is attached to the platelet storage bag, envelop the sampling bag and the soft platelet filter with the platelet storage bag.

Apply soft centrifugation conditions according to the local protocol.

Carefully remove the centrifuged bag set from the centrifugation cup and place the filled blood bag into an extractor to apply pressure (Fig. 1).

If an automated component extractor is used, refer to the manufacturer's instructions.

Invert the filter (outlet upward positioning) and open the YELLOW and the WHITE clamps for priming the filter.

When the air is expelled from the filter (platelet concentrate reaches the outlet), return it to its upright position.

When the layer of the remaining buffy-coat is detected in the tube, the transfer of the platelet concentrate is stopped by closing the YELLOW clamp.

Allow the platelets contained in the filter housing to drain to the platelet storage bag. When the filter housing collapses, close the WHITE clamp.

<< Instructions for pooled platelet concentrates / single donor platelets >>

- Hang the unit to be filtered on an I.V. stand.
- Place the platelet storage bag on the table.
- Invert the filter (outlet upward positioning) and open the YELLOW and the WHITE clamps for priming the filter.
- When the air is expelled from the filter (platelet concentrate reaches the outlet), return it to its upright position.
- When the unit has emptied, close the YELLOW clamp.
- Allow the platelets contained in the filter housing to drain to the platelet storage bag.
- When the filter housing collapses, close the WHITE clamp.
- Seal off the tube below the filter and dispose the used filter and bag.

**3. Air removal procedure from the platelet storage bag and sampling collection for quality control**

Agitate the platelet storage bag to ensure adequate mixing.

Slightly incline the platelet storage bag to draw the air in the corner closest to the sampling bag.

Open the BLUE clamp and gently squeeze the platelet storage bag to remove the residual air and the desired volume of platelet concentrate sample into the sampling bag (Fig. 2). Close the BLUE clamp.

If a sample of platelet concentrate is not required, close the BLUE clamp when the residual air is completely expelled from the platelet storage bag.

Seal off the tube below platelet storage bag and sampling bag.

Obtain platelet samples for bacterial analysis and/or other QC sampling by using the sampling port. By twisting the sampling port, the upper part is detached, releasing a female luer that can be connected to a standard male luer (re. VENOJECT<sup>®</sup> Luer Adapter) (Fig. 3).

**CAUTION** For bacterial testing sample under laminar flow conditions.

#### CAUTION

When transfusing the platelets to a patient, make use of an administration set with a filter.

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