REF T401231, T401231LS, T401231U vô trùng

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Instructions for Use

Titan™ Stabilizer Attachments

Attachments for Hercules™ 3/360 Universal Stabilizer Arms
Catalogue Numbers T401231, T401231LS, T401231U

Sterile

CAUTION: *Before use, please read carefully the instructions for use of Hercules™ 3/360 Universal Stabilizer Arm in addition to these instructions for use.*

PRODUCT DESCRIPTION

The Titan™ Stabilizer Attachment, Titan™ 360 Stabilizer Attachment, and Titan™ Stabilizer Attachment, Long Shaft are sterile single use accessory devices for use with the Hercules™ 3 or 360 Universal Stabilizer Arms (Product Reference T401161 or T401160). The attachment is provided in a double barrier sterile package and is intended for single use only. The device provides stabilization and coronary artery isolation. The malleable metal design with textured silicone coating allows the surgeon to bend the attachment to conform to anatomic shape in beating heart CABG (coronary artery bypass graft) procedures. The Titan Stabilizer Attachment, Long Shaft has a long attachment shaft with two malleable zones to meet additional positional requirements. The Titan Stabilizer Attachments include a vacuum option to maintain the position of the device.

Included in the sterile package is:

Titan consisting of the following parts (see Figure 1):

- The Titan attachment assembly with an approximately 22" length of ¼" outer diameter tubing with female luer connector on one end to connect to stopcock.
- Hi-flow yellow 4 way stopcock

Accessory vacuum tubing consisting of the following parts (see Figure 2):

- A 12' length of ¼" internal diameter tubing with a male luer connector on one end and vacuum source connector on the other end.
- A 3' segment of ¼" internal diameter tubing with female connectors on both sides to connect the canister to the vacuum regulator if required.

Figure 1 – Titan™ Stabilizer Attachment with tubing and stopcock

Figure 2 – Titan™ Stabilizer Attachment accessory vacuum tubing

Refer to Table 1 for available versions of the Titan Stabilizer Attachments.

Table 1

Product Reference Number	Type of Quick Connect	Interfacing Stabilization Platforms
T401231	Shaft Quick Connect	Hercules™ 3 Universal Stabilizer Arm (Product Reference: T401161)
T401231LS		
T401231U	Ball Quick Connect	Hercules™ 360 Universal Stabilizer Arm (Product Reference: T401160)

INTENDED USE

The Titan Stabilizer Attachments are intended for use in beating heart CABG (coronary artery bypass graft) procedures to stabilize the anastomotic site (coronary artery isolation).

CONTRAINDICATIONS

The Titan Stabilizer Attachments are contraindicated for use in patients with local or systemic infection.

Do not position the Titan Stabilizer Attachments on newly infarcted tissue, aneurysmal tissue, fragile tissue, or directly over a coronary artery.

WARNINGS

- This device is sterilized using gamma irradiation and is intended **FOR SINGLE USE ONLY. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT REPROCESS. REPROCESSING MAY COMPROMISE THE STERILITY, BIOCOMPATIBILITY AND FUNCTIONAL INTEGRITY OF THE DEVICE.**
- If the packaging is damaged, do not use, as sterility of the device may have been compromised.

- Inspect the product thoroughly for damage prior to use. Do not use damaged product.
- Proper surgical procedure and technique is the responsibility of the attending medical professional (physician). Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- Repeated bending the malleable zone of the shaft will eventually cause the shaft to break. Slowly bending the malleable zone will increase the allowable range of deformation of the shaft. The integrity of the shaft should be visually verified after bending.
- Excessive flexing of the malleable zone may cause breakage.
- **Take care not to occlude the vacuum lumen prior to obtaining proper positioning.**
- **Use of excessive negative vacuum pressure may cause tissue damage.**
- **Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue.**

CAUTION

United States Federal Law restricts this device to sale by or on the order of a physician or other licensed practitioner.

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the stabilizer attachment and the accessory vacuum tubing from its packaging and pass to the surgical field for attachment to the compatible Hercules Universal Stabilizer Arm (see Table 1).
 - a. For connection to the Hercules 3 Universal Stabilizer Arm:

Slide the sleeve of the Hercules 3 Universal Stabilizer Arm forward and insert the shaft of the attachment into the quick connect with the detent on the shaft (of the attachment) facing the flat area with logo on the sleeve of the Hercules 3 Universal Stabilizer Arm (see Figure 3). Return the sleeve to its original position to lock the accessory device into the quick connect. Check for secure connection.


Figure 3 – Attachment of the Titan™ Stabilizer Attachment to the Hercules™ 3 Universal Stabilizer Arm

- b. For connection to the Hercules 360 Universal Stabilizer Arm: Grasp the sleeve of the Hercules Universal Stabilizer Arm and press the ball of the attachment into the collet of the Hercules Universal Stabilizer Arm (see Figure 4). Check for secure connection.

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


Figure 4 – Attachment of the Titan™ 360 stabilizer attachment to the Hercules™ 360 Universal Stabilizer Arm

- The Titan Stabilizer Attachments can be bent to conform to anatomic requirements as shown in the pictures below indicating the maximum deflection ranges and the proper technique to be used during a procedure.

- The malleable zone of the shaft can be bent to a maximum of 45° in the downward and reversed upward positions as shown below. One bend cycle for the shaft is 0° to 45°.

WARNING

- Flexing and bending of the shaft and foot should not exceed 5 cycles during use. Visually verify the integrity of the shaft and foot after bending.
 - Repeated bending of the malleable zone of the shaft will eventually cause the shaft to break. Slowly bending the malleable zone will increase the allowable range of deformation of the shaft.
 - Excessive flexing of the malleable zone may cause breakage.
 - When bending devices with vacuum option, take care not to occlude the vacuum lumen.
3. Position the stabilizer attachment as required. Turn the main handle of the Hercules 3/360 Universal Stabilizer Arm clockwise until the gap between the compression washers closes. For further instructions refer to the instructions for use of the Hercules 3/360 Universal Stabilizer Arm.

Use of vacuum option with Titan Stabilizer Attachments:

1. Attach male luer connection of accessory vacuum tubing to the female luer port of the vacuum control stopcock, making sure the vacuum source is switched OFF. Attach the vacuum source connector of the accessory vacuum tubing to the vacuum source; turn the vacuum control stopcock OFF to the vacuum source.
2. Turn the vacuum regulator of the vacuum source on to sufficient negative pressure to affix the Titan Stabilizer Attachment to tissue. Due to variations in OR equipment and utilities, use the least amount of vacuum necessary for desired product performance. The regulator should not exceed -400 mmHg.

NOTE: The limit of -400 mmHg is clearly indicated on the yellow labels that are found on the accessory vacuum tubing.

WARNING

Use of excessive negative vacuum pressure may cause tissue damage.

3. Turn the vacuum control stopcock OPEN to the vacuum source and attachment once proper position has been obtained. Verify attachment of the stabilizer to the tissue.
4. To release attachment from tissue, turn the vacuum control stopcock OPEN to the atmosphere prior to removal of the stabilizer attachment.

WARNING

Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue.

STORAGE

The Titan Stabilizer Attachments storage conditions:

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DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

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English

Instructions for Use

Atlas™ Positioner Attachments Attachments for Hercules™ 3/360 Universal Stabilizer Arms Catalogue Numbers T401261 & T401261U Sterile

CAUTION: *Before use, please read carefully the instructions for use of the Hercules™ 3/360 Universal Stabilizer Arms in addition to these instructions for use.*

PRODUCT DESCRIPTION

The Atlas™ Positioner Attachment and Atlas™ 360 Positioner Attachment are sterile single use accessory devices for use with the Hercules™ 3 or 360 Universal Stabilizer Arms (Product Reference T401161 or T401160). The attachment is provided in a double barrier sterile package and is intended for single use only. This attachment is used to position the heart to provide exposure and access during beating heart CABG (coronary artery bypass graft) procedures. The positioner attachment is affixed to the tissue by means of negative pressure.

Included in the sterile package is:

Atlas consisting of the following parts (see Figure 1):

- The Cup/Spring/Elbow Assembly with an approximately 22" length of ¼" outer diameter tubing with female luer connector on one end to connect to stopcock.
- Hi-flow blue 4 way stopcock

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- A 12' length of ¼" internal diameter tubing with a male luer connector on one end and vacuum source connector on the other end
- A 3' segment of ¼" internal diameter tubing with female connectors on both sides to connect the canister to the vacuum regulator if required.


Figure 1 – Atlas™ Positioner Attachment with tubing and stopcock

Figure 2 – Atlas™ Positioner Attachment accessory vacuum tubing

Refer to Table 1 for available versions of the Atlas Positioner Attachments.

Table 1

Product Reference Number	Type of Quick Connect	Interfacing Stabilization Platforms
T401261	Shaft Quick Connect	Hercules™ 3 Universal Stabilizer Arm (Product Reference: T401161)
T401261U	Ball Quick Connect	Hercules™ 360 Universal Stabilizer Arm (Product Reference: T401160)

INTENDED USE

The Atlas Positioner Attachments are intended for use to position the heart to provide exposure and access during beating heart CABG (coronary artery bypass graft) procedures.

CONTRAINDICATIONS

The Atlas Positioner Attachments are contraindicated for use in patients with local or systemic infection.

Do not position the attachment on newly infarcted tissue, aneurysmal tissue, fragile tissue, or directly over a coronary artery.

WARNINGS

- This device is sterilized using gamma irradiation and is intended **FOR SINGLE USE ONLY. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT REPROCESS. REPROCESSING MAY COMPROMISE THE STERILITY, BIOCOMPATIBILITY AND FUNCTIONAL INTEGRITY OF THE DEVICE.**
- If the packaging is damaged, do not use, as sterility of the device may have been compromised.
- Inspect the product thoroughly for damage prior to use. Do not use damaged product.
- Proper surgical procedure and technique is the responsibility of the attending medical professional (physician). Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- Repeated bending of the malleable zone of the shaft will eventually cause the shaft to break. Slowly bending the malleable zone will increase the allowable range of deformation of the shaft. The integrity of the shaft should be visually verified after bending.
- Excessive flexing of the malleable zone may cause breakage.
- **Take care not to occlude the vacuum lumen prior to obtaining proper positioning.**

- Use of excessive negative vacuum pressure may cause tissue damage.
- Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue.

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INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the positioner attachment and the accessory vacuum tubing from its packaging and pass to the surgical field for attachment to the compatible Hercules Universal Stabilizer Arm (see Table 1).
 - a. For connection to the Hercules 3 Universal Stabilizer Arm:

Slide the sleeve of the Hercules Universal Stabilizer Arm forward and insert the shaft of the attachment into the quick connect with the detent on the shaft (of the attachment) facing the flat area with logo on the sleeve of the Hercules Universal Stabilizer Arm (see Figure 3). Return the sleeve to its original position to lock the accessory device into the quick connect. Check for secure connection.


Figure 3 – Attachment of the Atlas™ Positioner Attachment to the Hercules™ 3 Universal Stabilizer Arm

- b. For connection to the Hercules 360 Universal Stabilizer Arm:

Grasp the sleeve of the Hercules Universal Stabilizer Arm and press the ball of the attachment into the collet of the Hercules Universal Stabilizer Arm (see Figure 4). Check for secure connection.

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
Figure 4 – Attachment of the Atlas™ 360 Positioner Attachment to the Hercules™ 360 Universal Stabilizer Arm

2. Attach male luer connection of accessory vacuum tubing to the female luer port of the vacuum control stopcock, making sure the vacuum source is switched OFF. Attach the vacuum source connector of the accessory vacuum tubing to the vacuum source and turn the vacuum control stopcock OFF to the vacuum source.
3. Turn the vacuum regulator of the vacuum source on to sufficient negative pressure to affix the Atlas Positioner Attachment to tissue. Due to variations in OR equipment and utilities, use the least amount of vacuum necessary for desired product performance. The regulator should not exceed -250 mmHg.

NOTE: The limit of -250 mmHg is clearly indicated on the blue labels that are found on the accessory vacuum tubing.

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4. Turn the vacuum control stopcock OPEN to the vacuum source and positioner attachment. Verify attachment of the positioner to the tissue.

WARNING

Take care not to occlude the vacuum lumen prior to obtaining positive positioning.

5. Position the positioner attachment as required. Turn the main handle of the Hercules 3/360 Universal Stabilizer Arm clockwise until the gap between the compression washers closes. For further instructions refer to the instructions for use of the Hercules 3/360 Universal Stabilizer Arms.
6. To release the tissue from the positioner attachment, turn the vacuum control stopcock OPEN to the atmosphere prior to removal of the positioner attachment.

WARNING

Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue

STORAGE

The Atlas Positioner Attachments storage conditions:

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Hercules™ 3/360 Universal Stabilizer Arm

Catalogue Number T401160 & T401161

CAUTION: Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these directions. Failure to do so may result in patient complications.

PRODUCT DESCRIPTION

The Hercules™ 360 and Hercules™ 3 Universal Stabilizer Arms are reusable and are supplied non-sterile. They consist of a distal connection mechanism (either standard Terumo® Ball Quick Connect or Terumo® Shaft Quick Connect), an articulating arm with links, a handle for providing cable tension, and a mechanism to mount the device to a sternal retractor. Refer to Table 1 for available versions of the Hercules Universal Stabilizer Arms. See Figure 1 for a diagram of the Hercules 3 Universal Stabilizer Arm.

WARNING Device must be cleaned, lubricated, and sterilized prior to first use and before subsequent uses.


Figure 1 – Hercules™ 3 Universal Stabilizer Arm

Table 1

Model	Type of distal connection mechanism	Length	# of Links
Hercules™ 360 Universal Stabilizer Arm – Product Reference T401160	Terumo Ball Quick Connect	6.50 in (16.51 cm)	19
Hercules™ 3 Universal Stabilizer Arm – Product Reference T401161	Terumo Shaft Quick Connect	6.35 in (16.13 cm)	18

Note: The packaging configuration consists of one Hercules™ Universal Stabilizer Arm unit (1) per box.

A wide variety of single-use and reusable attachments are available from Terumo Cardiovascular Group intended for use with the Hercules Universal Stabilizer Arm.

See Table 2 for potential attachments for use with Hercules™ 360 Universal Stabilizer Arm (Product Reference: T401160).

Table 2

Product Reference	Description	
T401261U	Atlas™ 360 Positioner Attachment	Single use
T401231U	Titan™ 360 Stabilizer Attachment	Single use

See Table 3 for potential attachments for use with Hercules™ 3 Universal Stabilizer Arm (Product Reference: T401161).

Table 3

Product Reference	Description	
T401261	Atlas™ Positioner Attachment	Single use
T401231	Titan™ Stabilizer Attachment	Single use
T401231LS	Titan™ Stabilizer Attachment, Long Shaft	Single use
T401244	StableSoft™ II Stabilizer	Single use

T401245	StableSoft™ II Long Shaft	Single use
T401209S	Assistant with StableSoft™ 3L	Single use
T401210S	Assistant with StableSoft™ 3S	Single use
T401211S	Assistant with StableSoft™ 2L	Single use
T401212S	Assistant with StableSoft™ 2S	Single use

INTENDED USE

The Hercules 3 Universal Stabilizer Arm and Hercules 360 Universal Stabilizer Arm are intended for use in cardiothoracic surgery to support stabilization and positioning of various anatomical structures during various procedures.

CONTRAINDICATIONS

The Hercules 360 and Hercules 3 Universal Stabilizer Arms are contraindicated for use in patients with local or systemic infection.

WARNINGS AND PRECAUTIONS

WARNINGS

- The Hercules Universal Stabilizer Arm should only be operated by experienced physicians properly trained in beating heart surgery and the assistance of an experienced anesthetist is required.
- Proper surgical procedure and technique is the responsibility of the attending medical professional (physician). Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- Excessive torque might cause cable to fray, snap or break. If the Hercules Universal Stabilizer Arm does not hold position, it should not be used.
- DO NOT reposition the arm when the Hercules Universal Stabilizer Arm is tensioned. It will cause the cable in the flexible Universal Stabilizer Arm to fray and possibly break.
- Heavy or improper use of devices may cause damage requiring replacement of unit.
- Visually inspect cable for fray. Do not use fingers to check for cable fray.

- This unit is not to be reprocessed using hydrogen peroxide gas plasma sterilization such as STERRAD® sterilizer which can cause damage or breakage to the Nylon bearing component. Please contact the manufacturer of the automated washer/disinfector prior to use to determine if the washer/disinfector processing parameters are compatible with the Hercules Universal Stabilizer Arm.
- Strictly follow the manufacturer's instructions for concentration of the detergent solution and proper use to avoid high acidity or alkaline pH balances, which may cause corrosion and result in breakage.

PRECAUTION

- Cable fatigue and wear can be easily detected, and failure prevented, by examining the cable before each use. Prior to beginning each procedure, while the arm is fully loosened, look at the arm to determine if the cable shows signs of fraying (see Figure 3 for fraying example). Do not use if signs of fraying are present as a cable may suddenly fail during a procedure.
- Do not use chlorine based solutions (such as bleach) on the Hercules Universal Stabilizer Arm as this can cause corrosion.

CAUTION

- United States Federal Law restricts this device to sale by or on the order of a physician or other licensed practitioner.

INSTRUCTIONS FOR USE

I. PRE-USE INSPECTION

1. Prior to use, the unit should be fully loosened and visually inspected. The entire cable should be inspected by examining the spaces between links. No signs of wear or fraying of the cable should be present. If any wear is found, the device should not be used. See Figure 2 and Figure 3 for examples of normal and fraying cables.

Figure 2 - Appearance of a normal cable

Figure 3 - Appearance of a fraying cable – DO NOT USE

2. Verify function of the device by tightening and loosening the cable by turning the main handle.

II. PROCEDURE

WARNING Device must be cleaned and sterilized prior to first use and before subsequent uses (Refer to chapter on **REPROCESSING INSTRUCTIONS**).

1. Attach an accessory device to the distal connection mechanism. There are two different types of attachment interfaces that can be used.
 - a. For connection of a standard Terumo SHAFT QUICK CONNECT attachment, slide the sleeve of the Hercules Universal Stabilizer Arm forward and insert the shaft of the attachment into the quick connect with the detent on the shaft (of the attachment) facing the flat area with logo on the sleeve of the Hercules Universal Stabilizer Arm (see Figure 4). Return the sleeve to its original position to lock the accessory device into the quick connect. Check for secure connection.


Figure 4 – Attachment of accessory device to the Hercules™ 3 Universal Stabilizer Arm (Product Reference T401161)

- b. For connection of a Terumo BALL QUICK CONNECT attachment, grasp the sleeve of the Hercules Universal Stabilizer Arm and press the ball of the attachment into the collet of the Hercules Universal Stabilizer Arm (see Figure 5). Check for secure connection.

NOTE: Ensure that the sleeve is grasped during attachment insertion. Failure to grasp the sleeve may result in an inability to insert the attachment.


Figure 5 – Attachment of accessory device to the Hercules™ 360 Universal Stabilizer Arm (Product Reference T401160)

2. Attach the retractor clamp of the Hercules 360 Universal Stabilizer Arm to the retractor by sliding the rear jaw backward and placing the retractor clamp onto the retractor arm. Allow the rear jaw to spring closed against the retractor arm before rotating the retractor clamp handle clockwise to clamp the Hercules Universal Stabilizer Arm.

3. Position accessory device as required. Turn the main handle clockwise until the gap between the compression washers closes as depicted in Figure 6. The closing of the gap between the washers indicates that the recommended cable tension is reached as shown in the following picture.

CAUTION Do not reposition the arm without loosening the main handle.

Compression Washers


Figure 6 – Example of closed compression washers

III. POST-USE INSPECTION

After every use, the arm should be fully loosened. A visual inspection should be conducted between the last link and the clamp assembly for cable fray. Cable fatigue and fraying typically occurs in the same location, at the back end of the arm, nearest the retractor mounting clamp, making it very easy to detect early. There should be no signs of wear or fraying of the cable. See Figure 7 with a cable demonstrating fatigue and fraying. In addition, check for broken or missing nylon (red) washer depicted in Figure 8 and Figure 9.

If wear is observed, do not continue use, and contact your local Terumo Account Representative or call Customer Service for replacement. Again, the unit should never be used in the conditions shown in Figure 7, Figure 8, and Figure 9.

Figure 7 - Cable demonstrating fatigue and fraying


Figure 8 - Broken nylon washer


Figure 9 - Missing nylon washer

REPROCESSING INSTRUCTIONS

WARNING

Should cleaning and sterilization methods other than the recommended be performed, validation of the effectiveness of these methods must be performed, as well as verifying the effects on the Hercules Universal Stabilizer Arm materials and functions.

CAUTION It is recommended that the device is reprocessed as soon as is reasonably practical following use (recommended within one hour).

I. POINT OF USE

Remove excess soil with disposable cloth/paper wipe.

II. PREPARATION FOR CLEANING

Ensure the main handle is fully loosened prior to cleaning, and that any attachments are disconnected from the distal end of the Hercules Universal Stabilizer Arm. Perform cleaning and decontamination using either the manual or automated methods described below, followed by inspection, lubrication and sterilization.

III. CLEANING AND DECONTAMINATION

MANUAL METHOD

1. Rinse thoroughly under reverse osmosis/deionized (RO/DI) water for a minimum of 1 minute to remove gross soil and debris, taking into account the gaps between the links.
2. Immerse the device in a 0.8% enzymatic detergent instrument cleaner for a minimum of 5 minutes. Terumo Cardiovascular Group recommends using ENZOL/Cidezyme, or a suitable equivalent. ENZOL/Cidezyme is a brand of protein resolving cleaner to remove all traces of blood and debris.
3. Actuate and tap the device to insure penetration of the cleaner and release of trapped air bubbles.

4. Scrub each link and cable thoroughly using a soft brush or lint-free cloth under 0.8% enzymatic cleaner (ENZOL/Cidezyme). Remove all traces of blood and debris. Make sure all moving parts are cleaned thoroughly to prevent debris from interfering with movement.
5. Rinse the device thoroughly for a minimum of 1 minute in RO/DI water to remove dislodged surgical debris and the detergent solution, taking into account the gaps between the links.
6. Flush any hard to reach areas under water with a spray gun. Thoroughly flush the cable between each link for 5 seconds. Flush the distal connection mechanism for 15 seconds.
7. Sonicate device in 0.8% ENZOL/Cidezyme cleaning solution at 104°F/40°C for 10 minutes. Prior to sonication, actuate device in the cleaning solution to insure penetration of the cleaner and release of trapped air bubbles.
8. Rinse the device thoroughly for a minimum of 1 minute in RO/DI water taking into account the gaps between the links.
9. Flush any hard to reach areas under water with a spray gun. Thoroughly flush the cable between each link for 5 seconds. Flush the distal connection mechanism for 15 seconds.
10. Dry device completely using compressed air.

CAUTION Do not use chlorine based solutions (such as bleach) on the Hercules Universal Stabilizer Arm as this can cause corrosion.

CAUTION Incomplete rinsing after cleaning can cause enzymatic cleaner to form a residue on the Hercules Universal Stabilizer Arm that may result in binding of the links.

AUTOMATED CLEANING METHOD

MANUAL PRE-CLEAN

1. Rinse thoroughly under reverse osmosis/deionized (RO/DI) water for a minimum of 1 minute to remove gross soil and debris taking into account the gaps between the links.
2. Immerse the device in a 0.8% enzymatic detergent instrument cleaner for a minimum of 5 minutes. Terumo Cardiovascular Group recommends using ENZOL/Cidezyme, or a suitable equivalent. ENZOL/Cidezyme is a brand of protein resolving cleaner to remove all traces of blood and debris.
3. Actuate and tap the device to insure penetration of the cleaner and release of trapped air bubbles.

4. Scrub each link and cable thoroughly using a soft brush or lint-free cloth under 0.8% enzymatic cleaner (ENZOL/Cidezyme). Remove all traces of blood and debris. Make sure all moving parts are cleaned thoroughly to prevent debris from interfering with movement.
5. Flush any hard to reach areas under water with a spray gun. Thoroughly flush the cable between each link for 5 seconds. Flush the distal connection mechanism for 15 seconds.
6. Sonicate device in 0.8% ENZOL/Cidezyme cleaning solution at 104°F/40°C for 10 minutes. Prior to sonication, actuate device in the cleaning solution to insure penetration of the cleaner and release of trapped air bubbles.
7. Flush any hard to reach areas under water with a spray gun. Thoroughly flush the cable between each link for 5 seconds. Flush the distal connection mechanism for 15 seconds.

AUTOMATED PROCESS

1. Two minute pre-clean with cold tap water. Drain.
2. Five minute clean with deionized water at 55°C using 0.5 % of Neodisher Mediclean forte (Dr. Weigert, Hamburg). Drain.
3. Three minute rinsing and neutralization with deionized water. Drain.
4. Three minute thermal disinfection at 93°C with deionized water. Drain.
5. Ten minute drying.

CAUTION Do not use chlorine based solutions (such as bleach) on the Hercules Universal Stabilizer Arm as this can cause corrosion.

IV. INSPECTION, MAINTENANCE AND TESTING

Visually inspect the entire Hercules Universal Stabilizer Arm for cleanliness to assure residue has been removed. If any tissue, blood, pus or soil is still present, repeat cleaning process.

Ensure the arm is fully loosened, and inspect the entire cable for damage. There should be no signs of wear or fraying of the cable. If wear is observed, return to Terumo Cardiovascular Group for replacement. Cable fatigue and wear is a gradual event that occurs over many arm rigidifying cycles. As the cable begins to fatigue, the outer strands will begin to fray and break. The user can recognize onset of fatigue as continued revolutions of the handle do not rigidify the arm and an audible popping is heard. This sound indicates the individual cable wires breaking. Once cable fraying begins, the cable will eventually sever after numerous continued revolutions of the handle.

In the rare and unfortunate event where a cable completely severs and the articulating links disperse into the field, it is imperative to completely account for the distal assembly and all links per Table 1 of this IFU.

Check for smooth movement of retractor clamp handle, rear jaw, and main handle and articulating arm. Inspect the mated stainless steel and Nylon thrust bearing washers near the main handle to verify all are intact. The instrument should be returned for replacement if any defects are found.

V. LUBRICATION

For optimal device performance and to maximize useful life, lubrication of the device using Steris Hinge-Free instrument lubricant (or equivalent) per manufacturer's instructions is recommended prior to each device sterilization cycle. Ensure the handle is fully loosened to expose the screw prior to lubrication.

VI. STERILIZATION

The device should be sterilized by autoclave sterilization method using the following processing parameters shown in Table 4. Ensure the arm is not tightened during sterilization.

Table 4

Description	Temperature	Exposure Time	Min. Drying Time
Gravity Displacement (wrapped)	132°C to 138°C	15-20 minutes	30 minutes
Prevacuum (wrapped)	132°C to 138°C	4-20 minutes	20 minutes
Prevacuum (wrapped)	134°C to 138°C	3-20 minutes	20 minutes

STORAGE

The Hercules Universal Stabilizer Arm storage conditions:

- A clean, dust-free location out of direct sunlight and free from risk of high temperature, excessive humidity, and water leakage.
- Do not store the device near chemicals.

LIMITATIONS ON REPROCESSING

Repeated reprocessing has no significant influence on the length of life of the device.

END OF LIFE

Inspect the device prior to each use for evidence of unusual wear and/or damage. If such wear is observed, do not use the device. This device has a shelf life of 48 months from the date of manufacture.

DISPOSAL

At end of life, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

Manufacturer

TERUMO CARDIOVASCULAR SYSTEMS CORPORATION

125 Blue Ball Road
Elkton, Maryland 21921, USA
Tel: (410) 398-8500
Fax: (410) 392-7171
www.terumo-cvs.com

TERUMO EUROPE N.V.

Interleuvenlaan 40, 3001 Leuven, Belgium
www.terumo-europe.com

AU Sponser **TERUMO AUSTRALIA PTY LTD**
Macquarie Park NSW 2113 Australia

TERUMO CORPORATION

Tokyo 151-0072, Japan

Terumo® is a registered trademark of TERUMO CORPORATION
Hercules™, Titan™, StableSoft™, Assistant with StableSoft™, Atlas™ are
trademarks of Terumo Cardiovascular Systems Corporation
U.S. Patent No. 6,860,668

ENZOL® is a registered trademark of Advanced Sterilization Products,
Division of Ethicon US, LLC.

Hinge-Free® is a registered trademark of STERIS Corporation

STERRAD® is a registered trademark of Advanced Sterilization Products.
Advanced Sterilization Products (ASP), a Division of Ethicon, Inc. a
Johnson & Johnson company.

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