TERUFUSION®

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**TERUMO VIETNAM MEDICAL
EQUIPMENT Co., Ltd.**
14th Floor, Geleximco Building,
36 Hoang Cau Street, O Cho Dua Ward,
Dong Da District, Hanoi, Vietnam
Tel.: +84 4 3936 1643
Fax: +84 4 3936 1641

**TERUMO VIETNAM MEDICAL
EQUIPMENT Co., Ltd.
HO CHI MINH CITY BRANCH**
Room 505, 5th floor,
HOANG ANH SAFOMEC Office Building,
7/1 Thanh Thai Str., Ward14, District 10,
Ho Chi Minh City, Vietnam
Tel.: +84 8 3866 9263
Fax: +84 8 3866 9261

P.T. TERUMO INDONESIA
Wisma KEMAI 5th Floor JL.Jend.Sudirman
KAV.3
Jakarta 10220, Indonesia
Tel.: +62 21 572 4071
Fax: +62 21 572 4072

TERUMO MALAYSIA Sdn. Bhd.
Suite C405, 4th Floor, Centre Tower
Wisma Consplant 1
No. 2, Jalan SS 16/4 47500 Subang Jaya
Selangor Darul Ehsan, Malaysia
Tel.: +60 3 5880 8888
Fax: +60 3 5880 8891

**TERUMO MARKETING PHILIPPINES
INC.**
Unit 3203 A & B, West Tower, Philippine
Stock Exchange Centre, Exchange Road,
Ortigas Centre, Pasig City, Philippines
Tel.: +63 02 632 1674
Fax: +63 02 632 7966

TERUMO TAIWAN MEDICAL Co., Ltd.
7C, No. 170 Tun-Hwa North Road
Taipei, Taiwan, R.O.C.
Tel.: +886 2 2545 1250
Fax: +886 2 2545 1251

TERUMO CHINA (HONG KONG) Ltd.
Room 607-608, 6/F, Harcourt House,
39 Gloucester Road, Wanchai,
Hong Kong
Tel.: +852 2866 0811
Fax: +852 2529 0451

**TERUMO MEDICAL (SHANGHAI) Co.,
Ltd.**
Rm.901, Zhongshan Expo Plaza,
666 Huihai Road West,
Changning, Shanghai 200052, China
Tel.: +86 21 6237 1155
Fax: +86 21 6237 1150

TERUMO KOREA CORPORATION
6th Fl. Shinwon bldg., 21, Teheran-ro,
8-gil, Gangnam-gu, Seoul, 135-033
Korea
Tel.: +82-2-565-9225
Fax: +82-2-565-9224

TERUMO AUSTRALIA Pty LIMITED
Macquarie Park NSW 2113, Australia
Tel.: +61 2 9878 5122

**TERUMO CORPORATION CHENNAI
BRANCH**
Alexander Square, 2nd Floor,
No.34 & 35 Sardar Patel Road, Guindy,
Chennai 600 032, India
Tel.: +91 44 2230 0634
Fax: +91 44 2230 0622

Terumo Middle East FZE
Building BWB - Office#317
Dubai Airport Free Zone (DAFZ)
P.O. Box 54814 Dubai
United Arab Emirates
Tel.: +971 4 292 0200
Fax: +971 4 292 0222

As of August, 2017

Original instructions

TE_RS700N_E_600_002

TERUMO

TERUMO CORPORATION
44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN
TERUMO EUROPE N.V.
INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM

MADE IN JAPAN

TERUMO is a trademark of TERUMO CORPORATION.

® Registered trademark ©TERUMO CORPORATION 2017-08-01 17108

2017-09

TERUFUSION®

Standard Rack System TE-RS700

Instruction Manual

Store this instruction manual in a convenient location for future reference whenever necessary. Read and carefully follow the instructions described in this instruction manual before using this product. For safe and long-term use, perform periodic maintenance and inspection.

* For details on handling the pumps, see the instruction manual for each device.

Table of Contents

Purpose of the Product	1	After use	11
Features of the Product	1	Cleaning	11
Parts Description	2	Storage	12
Exterior Diagram	2	Maintenance and Inspections	12
Standard Accessory	2	Users Inspection	12
Precautions	3	Waste and Recycle	13
Use	3	Troubleshooting	13
Operation Procedure	5	Specifications	14
Attaching the TERUFUSION Standard Rack System	5	For Medical Staff	16
Connection Procedure	6	EMC Technical Information	15
Connecting the TERUFUSION Standard Rack System	8	Symbols	19
Attaching/Detaching a Pump	9	FOR INFORMATION ABOUT TERUMO PRODUCTS	20

Purpose of the Product

The TERUFUSION Standard Rack System is a modular system designed to meet the request of a multi-infusion system which combines the advantages of a stacking system (compact and extendable) with those of a racking system (the pumps can be removed or added independently). This modular system allows the use of multiple TERUFUSION Syringe Pump and TERUFUSION Infusion Pump for the same patient with only one power cable.

• Compatible pumps

Product name	Model	Catalogue number
TERUFUSION Infusion Pump Type LF/ TERUFUSION Infusion Pump Type LF3	TE-LF600/ TE-LF630	TE-LF6xxxxx TE*LF6xxxxx
TERUFUSION Infusion Pump Type LM/ TERUFUSION Infusion Pump Type LM3	TE-LM700/ TE-LM730	TE-LM7xxxxx TE*LM7xxxxx
	TE-LM800/ TE-LM830	TE-LM8xxxxx TE*LM8xxxxx
TERUFUSION Syringe Pump Type SS/ TERUFUSION Syringe Pump Type SS3	TE-SS700/ TE-SS730	TE-SS7xxxxx TE*SS7xxxxx
	TE-SS800/ TE-SS830	TE-SS8xxxxx TE*SS8xxxxx
TERUFUSION Syringe Pump Type SS3TCI	TE-SS830T	TE-SS8xxxxx TE*SS8xxxxx


Alphanumeric characters are entered for x in the Catalogue number.

Features of the Product

The product enables up to two additional TERUFUSION Standard Rack System to be combined, and up to nine pumps to be attached and powered.

Parts Description

Exterior Diagram




Note

• Up to two TERUFUSION Standard Rack System can be added, enabling up to nine pumps to be attached on the racks and supplied with power.

e.g.: When using one rack

e.g.: When using two combined racks

e.g.: When using three combined racks

Standard Accessory

Standard accessory supplied with this product.

Instruction Manual

Precautions

To ensure safe and correct use of this system, please be sure to observe all precautions.
Non-compliance with precautions and incorrect system use may result in damage or injury.
The following are signs used in this manual and their meanings:

Warning

This label preceding a precaution indicates that there will be a possible risk of death or personal injury if the precaution given is not complied with.

Cautions

This label preceding a precaution indicates that there will be a possible risk of personal injury or property damage if the precaution given is not complied with.

Use

Warnings

- Since this product does not have an airtight structure, it should not be used or stored in an active gas environment (including sterilizer gas), nebulizer-sprayed environment, high-humidity environment, etc. It should not be submerged into water. [The electronic components inside the device are affected, and any subsequent damage and time degradation will cause failure of this product.]
- Periodically check that no malfunction is observed in the operation of this product and the devices being used. If any malfunction is observed, take appropriate actions, such as immediately stopping the operation. [This product is not equipped with the function that indicates that an alarm has been issued at the pumps.]
- The product should not be used or stored in a flammable environment. [This may lead to ignition or explosion of the product.]
- Attach this product to the IV pole with the pole being same diameter. [Otherwise this product may be damaged or dropped.]

Cautions

- When attaching this product to the IV pole, check the withstand load and the stability. Do not release hands until securely fixed and confirmed. Also, when attaching a pump to this product, be sure to fix firmly. [Otherwise this product or the pump may drop or the IV pole may be damaged.]
- When this product is fixed to IV pole, ensure that the IV pole is on a stable surface to prevent it from falling down.
- Do not use this product placed or laid on the floor or a desk. [The product may fall. Also, the drip-proof function cannot be guaranteed.]
- Do not connect more than three racks. [Otherwise the electrical safety of this product cannot be guaranteed.]
- When connecting the racks, start attaching from the bottom. A connection cannot be made below a rack attached to the IV pole.
- When connecting the racks, check that AC inlet and AC outlet engage each other. [Power will not be supplied to this product if AC inlet and AC outlet are not properly engaged.]
- When removing connected rack systems, start removing from the top. [Loosening the pole clamp of the bottom rack first may damage the AC inlet and AC outlet, resulting in drop of the rack.]
- Communication Rack System (TE-RS800) and Communication Rack System (Extension) (TE-RS811) cannot be connected to this product.
- Disconnect the AC power cable from this product to remove it from the AC power source. Do not place this product in an area with obstacles that block the AC power cable from being disconnected.
- Do not use any pump other than that specified for this product. Read the instruction manual of the specified pump. [Using a pump other than that specified may not allow this product to work properly.]
- This product is not equipped with communication functions, and does not allow for communication even if a pump equipped with the infrared function is attached.
- Do not use, store, carry or place on the floor this product. When moving this product, attach this product to the IV pole.
- Use the AC power cable supplied with the pump. Connect to an AC outlet that is earthed at all times. [Using an AC power cable other than one supplied with the pump may result in failure of this product. In addition, if used without an earth connection, the electrical safety of this product cannot be guaranteed.]
- Attach the provided pole clamps to the pump before attaching it to the product. [Otherwise the pump cannot be properly fixed to the product and may drop.]

Cautions

- When attaching the pump to this product, check that fingers, the tube attached to the pump and the cord of the drip sensor are not caught between this product and the pump.
- After attaching the pump to this product, check that the tube of the pump is not caught or tangled. [Otherwise the tube may be deformed or occluded, thus causing a flow rate error or Occlusion alarm.]
- When attaching/detaching a pump to/from the product, a sudden drop or oscillation may occur on the pump. Be careful about this. [A sudden drop or oscillation affects the accuracy of the flow rate of the pump and the alarm function.]
- When attaching/detaching a pump to/from the product, check that the pump does not contact another pump attached to the product. [Otherwise the pump may be damaged.]
- When attaching/detaching a syringe pump to/from the product, be very careful not to touch the syringe. [Touching the syringe may cause an unexpected operation of the syringe pump (alarm issuance, flow rate error, etc.).]
- When attaching/detaching a pump to/from the product, do not touch the LCD and the switches. [Otherwise the LCD of the pump may be damaged or an unexpected operation (power on/off, stop, start, rapid infusion, etc.) may occur.]
- After attaching a pump on this product, check that the AC icon is lit on LCD of the pump. [The AC icon that is not lit indicates that the pump is operating on the battery, possibly causing the pump to become unavailable in case of emergency.]
- When detaching the pump with the infusion set from this product, check that the tube is not caught by the tube holder or other parts of the product. [If caught, the tube is pulled, possibly causing the IV pole to fall down or the tube to be damaged.]
- This product should not be brought into a control area for radiation devices/MRI or inside a hyperbaric oxygen therapy room. If this product is accidentally brought into such an environment, immediately stop using it and report to TERUMO trained service technicians.
- This product should not be used if it has received any impact (drop to floor, falling the IV pole, violent shock). Even though no fault is observed in the product appearance, the original functionality or performance (pump holding or power supply) of this product may not be achieved due to internal damage. In such cases, inspection and checking are required.
- Before use, check the instruction manual of the medical supplies, medical equipment and combined equipment. (The IV pole etc.)
- This product should only be operated by skilled personnel.
- Prior-to-use inspections must be conducted for this product. If any fault is observed, do not use this product and request an inspection and repair.
- In order to use the product safely, regularly conduct maintenance and inspection. [If any fault is found during any of the checkups, immediately stop using the product.]
- Attention is required not to apply any strong static electricity. [It may result in failure or malfunction.]
- Hold the handle of the rack to transport it. When the racks are combined, hold the back side or the bottom. [Otherwise the rack(s) may be damaged or dropped.]
- Even within the use conditions, this product should not be used under conditions that cause a sharp temperature change. [Condensation inside the device results in damage and time degradation, and thus the original functionality or performance of this product may not be achieved.]
- AC power cable used for this product should not be pinched with forceps or equipment, e.g. caster, or punctured with a needle. In addition, when the cables are laid on the floor, they should not be stepped on by a caster, etc. [Damaging the cables may lead to electric shock or fire. In addition, the original functionality or performance of this product will not be achieved.]
- This product should not be used in a place where vibration, dust, mist or corrosive gas occurs or in a place where the product is sprayed with liquid. If this product is sprayed with liquid, thoroughly wipe with a dry soft cloth.
- Since dripping drug solution may result in a short circuit by wetting the AC outlet, make sure that the connecting section is not wet when connecting. If moisture is present, make sure that the AC power cable is plugged out, and then thoroughly wipe it off with a dry cloth immediately. [Since this product does not have a waterproof structure, the components inside may be affected, resulting in device failure.]
- When moving the product fixed to the IV pole, be careful with slopes and uneven surfaces on the floor, and make sure not to run over any cables. [The product may drop or the IV pole may fall down.]
- Do not disassemble, make alterations to, or repair this product (including actions that interfere with the functionality or performance such as taping a movable part). [It may result in failure, damage or device performance degradation for this product.]
- This product should be used only after the power supply is confirmed to be sufficient. [If the power supply is not sufficient, the pumps operate on internal battery, resulting in making it unavailable in emergency situations. In the case that the supply capability is exceeded by connecting this product, it may affect other equipment.]
- When using devices (mobile phones, radio devices, radio knives, defibrillators, etc.) that emit an electromagnetic wave within the area of this product, they should be used as far away as possible and you should verify the normal operation of this product in the configuration it is used.


Cautions

- When used within the vicinity area of a radio knife: Medical radio knives are surgical equipment for incision and coagulation by high energy radiofrequency current. When this product is used in combination within the vicinity area of a radio knife, check the following before use.
 - (1) Radio knives have a different level of high-frequency noise emission depending on the types, and using old models (vacuum tube gap type) in combination should be particularly avoided as the noise levels from those are higher.
 - (2) The distance from the radio knife cord (knife holder, knife cord and return electrode cord) and radio knife body to this product should be kept 25 cm or more.
 - (3) The radio knife and this product should be operated by a power supply from a different system, and both should be securely earthed.

Operation Procedure

- Below, this product (Standard Rack System) will be referred as rack.

Attaching the TERUFUSION Standard Rack System

- **Fix the rack firmly to the IV pole by using all the pole clamps on the back side of the product.**

Warning

- Attach this product to the IV pole with the pole being same diameter. [Otherwise this product may be damaged or dropped.]


Caution

- Check the withstand load and the stability of the IV pole, and do not release hands until securely fixed and checked.


Connection Procedure

Caution


- When connecting the racks, start attaching from the bottom. A connection cannot be made below a rack attached to the IV pole.

- 1) **Attach a rack to the IV pole.**


- 2) **Remove the cover fixing screw of this product and lift the connector cover to remove it.**


- 3) **Connect the second rack on top of the rack in step 2).**

Check

- AC inlet and AC outlet engage each other and the rear sides of the racks on the first and second rows are aligned.

4) Tighten the pole clamp.

Caution

- Check the withstand load and the stability of the IV pole and do not release hands until securely fixed and checked.

Note

- Incorrect connection causes the pole clamp to be misaligned from the pole, not allowing you to close it.


Note

- To connect the third rack, repeat the steps 2) to 4).


Cautions

- When connecting the racks, check that AC inlet and AC outlet engage each other. Power will not be supplied to this product if AC inlet and AC outlet are not properly engaged.
- Do not connect more than three racks. [Otherwise the electrical safety of this product cannot be guaranteed.]
- When removing connected rack systems, start removing from the top. [Loosening the pole clamp of the bottom rack first may damage the AC inlet and AC outlet, resulting in drop of the rack.]
- Communication Rack System (TE-RS800) and Communication Rack System (Extension) (TE-RS811) cannot be connected to this product.
- This product is not equipped with communication functions, and does not allow for communication even if a pump equipped with the infrared function is attached.

Connecting the TERUFUSION Standard Rack System

- 1) Firmly connect the AC power cable supplied with the pump into the AC Inlet on the back of this product.


- 2) Insert the plug into the AC outlet with an earth connection.

Caution

- Use the AC power cable supplied with the pump. Connect to an AC outlet that is earthed at all times. [Using an AC power cable other than one supplied with the pump may result in failure of this product. In addition, if used without an earth connection, the electrical safety of this product cannot be guaranteed.]

Attaching/Detaching a Pump

Attaching a Pump


- 1) Remove the cover from the AC outlet.

Caution

- Keep the supplied cover attached to the AC outlet when not in use.


- 2) Hold the pump with both hands and align the pole clamp of the pump with the pump attaching guide.

- 3) Push the pump into the holder.

Checks


- The pump is pressed until it makes a click and is firmly fixed.
- The AC icon of the pump is lit.
- The Attaching check window is complete green.


Note

- A "red" Attaching check window indicates that the pump is not firmly fixed. Repeat the procedure from step 2) to firmly attach the pump.

Detaching a Pump

- 1) Pulling the release lever makes the pump come out forward.

- 2) Hold the pump with both hands and detach it.

Cautions

- Use both hands to attach/detach a pump to/from the product.
- After attaching the pump on this product, check that the AC icon is lit on LCD of the pump. [The AC icon that is not lit indicates that the pump is operating on the battery, possibly causing the pump to become unavailable in the case of an emergency.]
- When attaching the pump to this product, check that fingers, the tube attached to the pump or the cord of the drip sensor are not caught between this product and the pump. [Otherwise the pump may drop from the product. Also, the cord, etc. may be damaged or the tube may be damaged or deformed, and therefore the accuracy of the flow rate of the pump and the alarm functions cannot be guaranteed.]
- Attach the provided pole clamps to the pump before attaching it to this product. [Otherwise the pump cannot be properly fixed to the product and may drop.]
- After attaching the pump to this product, check that the tube of the pump is not caught or tangled. [Otherwise the tube may be deformed or occluded, thus causing a flow rate error or Occlusion alarm.]
- When attaching/detaching a pump to/from the product, a sudden drop or oscillation may occur on the pump. Be careful about this. [A sudden drop or oscillation affects the accuracy of the flow rate of the pump and the alarm function.]
- When attaching/detaching a syringe pump to/from the product, be very careful not to touch the syringe. [Touching the syringe may cause an unexpected operation of the syringe pump (alarm issuance, flow rate error, etc.).]
- When attaching/detaching a pump to/from the product, do not touch the LCD and the switches. [Otherwise the LCD of the pump may be damaged or an unexpected operation (power on/off, stop, start, rapid infusion, etc.) may occur.]
- When attaching/detaching a pump to/from the product, check that the pump does not contact another pump attached to the product. [Otherwise the pump may be damaged.]
- When detaching the pump with the infusion set from this product, check that the tube is not caught by the tube holder or other parts of the product. [If caught, the tube is pulled, possibly causing the IV pole to fall down or the tube to be damaged.]
- Since dripping drug solution may result in a short circuit by wetting the AC outlet, make sure that the connecting section is not wet when connecting. If moisture is present, make sure that the AC power cable is plugged out, and then thoroughly wipe it off with a dry cloth immediately. [Since this product does not have a waterproof structure, the components inside may be affected, resulting in device failure.]

After use


Cleaning

- Clean this product by regularly. When disinfecting, use a gauze (damp with antiseptic solution) to wipe off the product, and wipe off the antiseptic solution with a gauze (damp with cold/luke warm water), and then thoroughly wipe off any moisture with a dry soft cloth. Follow the instruction manual for the dilution ratio. Examples of disinfectants (ingredient names) which can be used are listed below.

Examples of Cleaning Disinfectants (Ingredient Names)

Ingredient Name	Dilution (e.g.)
Chlorhexidine gluconate	5%
Benzalkonium chloride	10%
Ethanol	76.9 - 81.4vol%

When using disinfectant, follow the instruction manual of each disinfectant (regarding degree of dilution, etc.).

Cautions

- When disinfecting, do not use a sterilizer. Use a gauze (damp with antiseptic solution) to wipe off the product, and wipe off the antiseptic solution with a gauze (damp with cold/luke warm water), and then thoroughly wipe off any moisture with a dry soft cloth. For the dilution ratio of the antiseptic solution, refer to the instructions for use of the antiseptic solution. Examples of the available antiseptic solutions (component name) are as follows:
Chlorhexidine gluconate/Benzalkonium chloride/Ethanol
- Before cleaning, always disconnect the AC power cable and pump. [Failure to do so may lead to failure of this product or electric shock etc.]
- If any adhered drug solution is present, the movable parts, such as the release lever, may not properly operate. If any drug solution is present, immediately conduct cleaning such as wiping off the dirt with a cotton swab. (See figure above.)
- Do not wipe with undesignated alcohol, other organic solvents such as thinner, or povidone-iodine. Using organic solvent or any antiseptic solution other than those permitted for use may result in damage to or failure of this product.
- Do not use any replacement parts other than those specified. [The original functionality or performance of this product may not be achieved.]
- Do not use an electric dryer to dry this product. [This product may be damaged.]
- Do not wash this product with running water or immerse in water. [Since this product does not have a waterproof structure, it may result in damage or failure.]
- Before wiping off dirt, remove the pumps and the AC power cable from the product. [If drug solution or other liquid adheres to the AC outlet, a short circuit may occur.]

Storage

Warning

- The product should not be used or stored in a flammable environment. [This may lead to ignition or explosion of the product.]

Cautions

- Do not store this product in a place with a high level of vibration, dust, mist or corrosive gas.
- Do not expose this product to the sunlight or ultraviolet irradiation for a long time. [The exterior may experience colour change, deformation or deterioration.]
- Do not store this product in a place where an environment of atmospheric pressure, temperature, humidity, ventilation, salt content or sulfur content may cause an adverse effect.
- Do not store this product in a storage area for chemicals or in a place that generates gases.
- When transporting this product, avoid any shock, vibration, dust and/or high temperatures and high humidity.
Transport conditions: Temperature: -20 to 60°C, Humidity: 10 to 95%RH (non-condensing), Atmospheric pressure: 50 to 106kPa (500 to 1060hPa)

Maintenance and Inspections

In order to ensure safe use of the product for an extended time, regularly conduct maintenance and inspection. If any fault is found during any of the checkups, immediately stop using the product and contact TERUMO trained service technicians.

Cautions

- This product should not be used if it has received any impact (drop to floor, falling the IV pole, violent shock). Even though no fault is observed in the product appearance, the original functionality or performance (pump holding or power supply) of this product may not be achieved due to internal damage. In such cases, inspection and checking are required.
- In order to use the product safely, regularly conduct maintenance and inspection. If any fault is found during any of the checkups, immediately stop using the product.
- Do not disassemble, make alterations to, or repair this product (including actions that interfere with the functionality or performance such as taping the operation panel or a movable part). [It may result in failure, damage or device performance degradation for this product.]

Users Inspection

Pre-use Inspection

- 1) The appearance of the body is not damaged.
- 2) The movable parts including the AC outlet move smoothly.
(Check when no current is being supplied to the product.)
- 3) When the release lever is not operated, the colour of the Attaching check window changes to "green"; when operated, to "red".
- 4) This product can be firmly fixed to the IV pole with no rattling.
- 5) The pumps can be firmly fixed to this product with no rattling.

Waste and Recycle

Electrical and electronic equipment (EEE) contain materials, components and substances which can be dangerous to the environment and harmful to human health if waste electrical and electronic equipment (WEEE) are not disposed of correctly.

Waste electrical and electronic equipment must not be disposed of with the remainder of unseparated waste, but should instead be collected separately. In this way, the environmental impact associated with disposal of WEEE is reduced and there will be more opportunity for reusing, recycling and recovering WEEE.

At the end of life, please dispose of this equipment according to your local regulations. Contact your local distributor or municipality for details of the available collection schemes.

Troubleshooting

Keyword	Occurrence	Cause	Action
Power supply	Power is not supplied to the pump.	The AC power cable is not properly connected to the AC inlet.	Check the connection of the AC power cable.
		The AC power cable is damaged.	Replace the AC power cable.
		The pump is not securely attached.	Attach the pump again by following the procedure (see page 9).
		The AC inlet of the pump is not firmly connected to the AC outlet on the rack.	Attach the pump again by following the procedure (see page 9).
Attaching of the pump	The pump cannot be attached.	The pole clamp of the pump is not correctly attached.	Re-attach the pole clamp correctly.
		A drug solution is adhered to the movable part.	Clean by following page 11 "Cleaning".

Specifications

Product name		TERUFUSION Standard Rack System
Model		TE-RS700
Catalogue number		TE•RS700N
Functions		Power output: Supplies AC power source to the attached pumps.
Operating Conditions		Temperature: 5 to 40°C, Humidity: 20 to 90%RH (no condensation), Atmospheric pressure: 70 to 106kPa (700 to 1060hPa)
Storage conditions		Temperature: -20 to 45°C, Humidity: 10 to 95%RH (no condensation), Atmospheric pressure: 50 to 106kPa (500 to 1060hPa)
Transport conditions		Temperature: -20 to 60°C, Humidity: 10 to 95%RH (no condensation), Atmospheric pressure: 50 to 106kPa (500 to 1060hPa)
Power source		AC 100-240 V, 50-60 Hz
Rated voltage/Rated current/ Rated frequency		AC 100-240 V, 1.8-1.3A, 50-60 Hz (Up to nine pumps can be supplied with power*)
Power consumption		Standard Rack System x1: Max.102VA (When three pumps are attached*) Standard Rack System x2: Max.204VA (When six pumps are attached*) Standard Rack System x3: Max.306VA (When nine pumps are attached*)
Classification	Electric shock protection	Class I equipment
	Mode of operation	Continuous operation
	Protection against harmful ingress of water or particulate matter	IP22
Dimensions		Standard Rack System x1: 220 mm (W) x 458 mm(H) x 132 mm(D) (Approx.) Standard Rack System x2: 220 mm (W) x 845 mm(H) x 132 mm(D) (Approx.) Standard Rack System x3: 220 mm (W) x 1232 mm(H) x 132 mm(D) (Approx.)
Weight		Standard Rack System x1: Approx. 3.3 kg Standard Rack System x2: Approx. 6.5 kg Standard Rack System x3: Approx. 9.7 kg
Standard accessory		Instruction manual

* Up to three racks can be combined and up to nine pumps can be attached at a time.

• This product is compatible with EMC (electromagnetic compatibility) standard IEC 60601-1-2:2007 (EN 60601-1-2:2007) (CISPR group classification and class classification are Group 1 and Class B). It is also compatible with the EMC level required by IEC 60601-2-24:2012 (EN 60601-2-24:2015).

• Conformity standard and directive
 IEC 60601-1:2005+A1:2012 (EN 60601-1:2006+A1:2013)
 IEC 60601-1-2:2007 (EN 60601-1-2:2007)
 IEC 60601-1-6:2010+A1:2013 (EN 60601-1-6:2010+A1:2015)
 IEC 60601-1-8:2006+A1:2012 (EN 60601-1-8:2007+A1:2013)
 IEC 60601-2-24:2012 (EN 60601-2-24:2015)
 MDD (Medical Device Directive) 93/42/EEC (Class IIb)

For Medical Staff

EMC Technical Information

Medical electrical equipment requires particular care in regards to EMC, and it is necessary to install and use the equipment in accordance with the following EMC information.

Caution

- Make sure to use the AC power cable supplied with the pump. If they are not, this product's emissions may increase and its immunity may decrease.

Table 1 – Guidance and manufacturer's declaration – electromagnetic emissions – for all EQUIPMENT and SYSTEMS (see 5.2.2.1 c))

Guidance and manufacturer's declaration – electromagnetic emissions		
The TERUFUSION Standard Rack System is intended for use in the electromagnetic environment specified below. The customer or the user of the TERUFUSION Standard Rack System should assure that it is used in such an environment.		
Emissions test	Compliance	Electromagnetic environment-guidance
RF emissions CISPR 11	Group 1	The TERUFUSION Standard Rack System uses RF energy only for its internal function. Therefore, its RF emissions are very low and are not likely to cause any interference in nearby electronic equipment.
RF emissions CISPR 11	Class [B]	The TERUFUSION Standard Rack System is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class [A]	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Complies	

Table 2 – Guidance and manufacturer's declaration – electromagnetic immunity – for all EQUIPMENT and SYSTEMS (see 5.2.2.1 f))

Guidance and manufacturer's declaration – electromagnetic immunity			
The TERUFUSION Standard Rack System is intended for use in the electromagnetic environment specified below. The customer or the user of the TERUFUSION Standard Rack System should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Electrostatic discharge (ESD) IEC 61000-4-2	± 6kV contact ± 8kV air	±8 kV contact ±15 kV air (*)	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	±2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV line(s) to line(s) ±2 kV line(s) to earth	±1 kV line(s) to line(s) ±2 kV line(s) to earth	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply input lines IEC 61000-4-11	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s	<5% U_T (>95% dip in U_T) for 0.5 cycle 40% U_T (60% dip in U_T) for 5 cycles 70% U_T (30% dip in U_T) for 25 cycles <5% U_T (>95% dip in U_T) for 5 s	Mains power quality should be that of a typical commercial or hospital environment. If the user of the TERUFUSION Standard Rack System requires continued operation during power mains interruptions, it is recommended that the TERUFUSION Standard Rack System be powered from an uninterruptible power supply or a battery.
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Power frequency (50-60 Hz) magnetic field IEC 61000-4-8	3 A/m	3 A/m (*)	Power frequency magnetic fields should be at levels characteristic of a typical location in a typical commercial or hospital environment.
<p>Note</p> <ul style="list-style-type: none"> • U_T is the a.c. mains voltage prior to application of the test levels. 			

* The TERUFUSION Standard Rack System complies with the more stringent levels of IEC 60601-2-24:2012 (EN 60601-2-24:2015).

Table 3 – Guidance and manufacturer's declaration – electromagnetic immunity – for LIFE-SUPPORTING EQUIPMENT and SYSTEMS (see 5.2.2.2)

Guidance and manufacturer's declaration – electromagnetic Immunity			
The TERUFUSION Standard Rack System is intended for use in the electromagnetic environment specified below. The customer or the user of the TERUFUSION Standard Rack System should assure that it is used in such an environment.			
Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment-guidance
Conducted RF IEC 61000-4-6	3 Vrms 150 kHz to 80 MHz outside ISM bands ^a	3 Vrms	Portable and mobile RF communications equipment should be used no closer to any part of the TERUFUSION Standard Rack System, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance $d=1.2 \sqrt{P}$ $d=1.2 \sqrt{P}$
	10 Vrms 150 kHz to 80 MHz in ISM bands ^a	10 Vrms	
Radiated RF IEC 61000-4-3	10 V/m 80 MHz to 2.5 GHz	10 V/m	$d=1.2 \sqrt{P}$ 80 MHz to 800 MHz $d=2.3 \sqrt{P}$ 800 MHz to 2.5 GHz where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in metres (m). ^b Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey ^c , should be less than the compliance level in each frequency range ^d . Interference may occur in the vicinity of equipment marked with the following symbol: 
Notes			
1. At 80 MHz and 800 MHz, the higher frequency range applies.			
2. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.			
a. The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.			
b. The compliance levels in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz are intended to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas. For this reason, an additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in these frequency ranges.			
c. Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the TERUFUSION Standard Rack System is used exceeds the applicable RF compliance level above, the TERUFUSION Standard Rack System should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as reorienting or relocating the TERUFUSION Standard Rack System.			
d. Over the frequency range 150 kHz to 80 MHz, field strengths should be less than 3 V/m.			

Table 5 –Recommended separation distances between portable and mobile RF communications equipment and the EQUIPMENT or SYSTEM—for LIFESUPPORTING EQUIPMENT and SYSTEMS (see 5.2.2.2)

Recommended separation distances between portable and mobile RF communications equipment and the TERUFUSION Standard Rack System				
The TERUFUSION Standard Rack System is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the TERUFUSION Standard Rack System can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the TERUFUSION Standard Rack System as recommended below, according to the maximum output power of the communications equipment.				
Separation distance according to frequency of transmitter				
m				
Rated maximum output power of transmitter W	150 kHz to 80 MHz outside ISM bands $d=1.2 \sqrt{P}$	150 kHz to 80 MHz in ISM bands $d=1.2 \sqrt{P}$	80 MHz to 800 MHz $d=1.2 \sqrt{P}$	800 MHz to 2.5 GHz $d=2.3 \sqrt{P}$
0.01	0.12	0.12	0.12	0.23
0.1	0.38	0.38	0.38	0.73
1	1.2	1.2	1.2	2.3
10	3.8	3.8	3.8	7.3
100	12	12	12	23
For transmitters rated at a maximum output power not listed above, the recommended separation distance d in metres (m) can be determined using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.				
Notes				
1. At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.				
2. The ISM (industrial, scientific and medical) bands between 150 kHz and 80 MHz are 6.765 MHz to 6.795 MHz; 13.553 MHz to 13.567 MHz; 26.957 MHz to 27.283 MHz; and 40.66 MHz to 40.70 MHz.				
3. An additional factor of 10/3 is used in calculating the recommended separation distance for transmitters in the ISM frequency bands between 150 kHz and 80 MHz and in the frequency range 80 MHz to 2.5 GHz to decrease the likelihood that mobile/portable communications equipment could cause interference if it is inadvertently brought into patient areas.				
4. These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.				

Symbols

Symbol	Description of Symbols	Symbol	Description of Symbols
	CE means compliance with the Medical Device Directive, 93/42/EEC.		Warning: Ensure that the IV pole is on a stable surface to prevent it from falling.
	Authorized Representative in the European Community		Warning: Do not exceed the withstand load of the IV pole to prevent it from falling.
	Serial number		Warning: Fix the pole clamp securely to prevent the rack from falling.
	Catalogue number		Up to three rack systems can be connected
	Date of Manufacture		This way up
	Manufacturer		Fragile, handle with care
IP22	IEC 60529 Degrees of protection provided by enclosures (IP code)		Keep dry
	Consult instructions for use		Keep away from sunlight
	Follow instructions for use		Temperature limit
	Warning		Humidity limitation
	Alternating current		Atmospheric pressure limitation
	DC power supply		Stacking limit by number
	Separate collection of electrical and electronic equipment, European Community		Contents

FOR INFORMATION ABOUT TERUMO PRODUCTS

• If this product should fail to perform as intended, immediately stop using the product and contact the nearest branch or sales office of TERUMO.

<p>TERUMO EUROPE N.V. Interleuvenlaan 40 3001 Leuven, Belgium Tel.: +32 16 38 12 11 Fax: +32 16 40 02 49</p>	<p>TERUMO MEDICAL CORPORATION 2101 Coltonail Lane Somerset, New Jersey 08873, U.S.A. Tel.: +1 732 302 4600 Fax: +1 732 302 3083 Toll free: +1 800 283 7866</p>	<p>TERUMO VIETNAM MEDICAL EQUIPMENT Co., Ltd. 14th Floor, Geleximco Building, 36 Hoang Cau Street, O Cho Dua Ward, Dong Da District, Hanoi, Vietnam Tel.: +84 4 3936 1643 Fax: +84 4 3936 1641</p>	<p>TERUMO MEDICAL (SHANGHAI) Co., Ltd. Rm.901, Zhongshan Expo Plaza, 666 Huailin Road West, Changning, Shanghai 200052, China Tel.: +86 21 6237 1155 Fax: +86 21 6237 1150</p>
<p>TERUMO EUROPE N.V. BENELUX SALES DIVISION Interleuvenlaan 40 3001 Leuven, Belgium (Belgium) Tel.: +32 16 39 25 80 Fax: +32 16 39 25 99 (Netherlands) Tel.: +0800 0220396 Fax: +0800 0220414</p>	<p>TERUMO LATIN AMERICA CORPORATION Doral Corporate Center I 8750 NW 36th street, Suite 600 Miami, FL 33178, U.S.A. Tel.: +1 305 477 4822 Fax: +1 305 477 4872</p>	<p>TERUMO VIETNAM MEDICAL EQUIPMENT Co., Ltd. HO CHI MINH CITY BRANCH Room 505, 5th floor, HOANG ANH SAFOMECH Office Building, 7/1 Thanh Thai Str., Ward14, District 10, Ho Chi Minh City, Vietnam Tel.: +84 8 3866 9263 Fax: +84 8 3866 9261</p>	<p>TERUMO KOREA CORPORATION 6th Fl. Shinwon bldg., 21, Teheran-ro, 8-gil, Gangnam-gu, Seoul, 135-933 Korea Tel.: +82-2-565-9225 Fax: +82-2-565-9224</p>
<p>TERUMO FRANCE S.A.S. Bâtiment Renaissance 3 rond-point des Saules 78280 GUYANCOURT France Tel.: +33 1 30 96 13 00 Fax: +33 1 30 43 60 85</p>	<p>TERUMO MEXICO, S.A. de C.V. Av. Insurgentes Sur # 1647 piso 11 Colonia San Jose Insurgentes C.P. 03900 Delegacion Benito Juarez Mexico D.F. Tel.: +52 55 10 85 07 70 Fax: +52 55 10 85 07 71</p>	<p>P.T. TERUMO INDONESIA Wisma KEIAI 5th Floor JL.Jend.Sudirman KAV.3 Jakarta 10220, Indonesia Tel.: +62 21 572 4071 Fax: +62 21 572 4072</p>	<p>TERUMO AUSTRALIA Pty LIMITED Macquarie Park NSW 2113, Australia Tel.: +61 2 9878 5122</p>
<p>TERUMO DEUTSCHLAND GmbH Ludwig-Erhard-Strasse 6 65760 Eschborn Germany Tel.: +49 6196 80 230 Fax: +49 6196 80 23 200</p>	<p>TERUMO MEDICAL do BRASIL Ltda. Praça General Gentil Falcão, 108 cj 91 e 92 - São Paulo - SP - Brasil - CEP: 04571 - 150 Tel.: +55 11 3594 3800 Fax: +55 11 3594 3801</p>	<p>TERUMO MALAYSIA Sdn. Bhd. Suite C405, 4th Floor, Centre Tower Wisma Consplant 1 No. 2, Jalan SS 16/4 47500 Subang Jaya Selangor Darul Ehsan, Malaysia Tel.: +60 3 5880 8898 Fax: +60 3 5880 8891</p>	<p>TERUMO CORPORATION CHENNAI BRANCH Alexander Square, 2nd Floor, No.34 & 35 Sardar Patel Road, Guindy, Chennai 600 032, India Tel.: +91 44 2230 0634 Fax: +91 44 2230 0622</p>
<p>TERUMO EUROPE ESPANA SL Avda. Juan Carlos I, N° 13 - 7ª Planta, Edificio Torre La Garena 28806 Alcalá de Henares (Madrid), Spain Tel.: +34 902 10 12 98 Fax: +34 902 10 13 58</p>	<p>TERUMO CHILE Ltda. Carmencita 25, Oficina 22 piso 2, Edificio Central Park, Las Condes, Santiago 755-0000, Chile Tel.: +56 2 2480 9600 Fax: +56 2 2480 9608</p>	<p>TERUMO PHILIPPINES INC. Unit 3203 A & B, West Tower, Philippine Stock Exchange Centre, Exchange Road, Ortigas Centre, Pasig City, Philippines Tel.: +63 02 632 1674 Fax: +63 02 632 7966</p>	<p>TERUMO MIDDLE EAST FZE Building #WB - Office#317 Dubai Airport Free Zone (DAFZ) P.O. Box 54614 Dubai United Arab Emirates Tel.: +971 4 292 0200 Fax: +971 4 292 0222</p>
<p>TERUMO ITALIA SRL Via Paolo di Dono 73 00142 Roma Italy Tel.: +39 06 94 80 28 00 Fax: +39 06 60 51 32 33</p>	<p>TERUMO SINGAPORE Pte. Ltd. 300 Beach Road, #33-06 The Concourse Singapore 199555 Tel.: +65 6 291 3603 Fax: +65 6 291 2696</p>	<p>TERUMO TAIWAN MEDICAL Co., Ltd. 7C, No. 170 Tun-Hwa North Road Taipei, Taiwan, R.O.C. Tel.: +886 2 2545 1250 Fax: +886 2 2545 1251</p>	<p>TERUMO THAILAND Co., Ltd. #1206, 12th Floor 54 B.B. Building Asoket Road (Sukhumvit 21), Kwaeng Klongloeyua, Khel Wattana, Bangkok 10110, Thailand Tel.: +66 2 260 7020 Fax: +66 2 260 7018</p>
<p>TERUMO SWEDEN AB Nya Varvet, Byggnad 80 Sven Kållfells Gata 18 426 71 Västra Frölunda, Sweden Tel.: +46 3174 85 880 Fax: +46 3174 85 890</p>	<p>TERUMO UK Ltd. Otium House 2 Freemantle Road, Bagshot Surrey, GU19 5LL United Kingdom Tel.: +44 1276 480440 Fax: +44 1276 480465</p>	<p>TERUMO CHINA (HONG KONG) Ltd. Room 607-608, 8/F, Harcourt House, 39 Gloucester Road, Wanchai, Hong Kong Tel.: +852 2866 0811 Fax: +852 2529 0451</p>	

As of August, 2017

Original Instructions

TE_RS700N_E_500_002

TERUMO

TERUMO CORPORATION
 44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN
TERUMO EUROPE N.V.
 INTERLEUVENLAAN 40, 3001 LELUVEN, BELGIUM
 MADE IN JAPAN

TERUMO is a trademark of TERUMO CORPORATION.

© Registered trademark © TERUMO CORPORATION 2017-09-01 17108