REF T401261 & Y401261U Tiệt trùng

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This page is intentionally left blank.

Atlas™ Positioner Attachments

Table of Contents

English	5
اللغة العربية	10
Česky	14
Dansk	19
Nederlands	24
Suomi	29
Français (France)	34
Deutsch	39
Ελληνικά	44
Magyar	49
Italiano	54
日本語	59
Norsk	64
Polski	69
Português (Brasil)	74
Русский язык	79
简体中文	84
Español (España)	89
Svenska	94
Türkçe	99
Graphical Symbols & Definitions	104

This page is intentionally left blank.

English

Instructions for Use

Atlas™ Positioner Attachments Attachments for Hercules™ 3/360 Universal Stabilizer Arms Catalogue Numbers T401261 & T401261U Sterile

CAUTION: *Before use, please read carefully the instructions for use of the Hercules™ 3/360 Universal Stabilizer Arms in addition to these instructions for use.*

PRODUCT DESCRIPTION

The Atlas™ Positioner Attachment and Atlas™ 360 Positioner Attachment are sterile single use accessory devices for use with the Hercules™ 3 or 360 Universal Stabilizer Arms (Product Reference T401161 or T401160). The attachment is provided in a double barrier sterile package and is intended for single use only. This attachment is used to position the heart to provide exposure and access during beating heart CABG (coronary artery bypass graft) procedures. The positioner attachment is affixed to the tissue by means of negative pressure.

Included in the sterile package is:

Atlas consisting of the following parts (see Figure 1):

- The Cup/Spring/Elbow Assembly with an approximately 22" length of ¼" outer diameter tubing with female luer connector on one end to connect to stopcock.
- Hi-flow blue 4 way stopcock

Accessory vacuum tubing consisting of the following parts (see Figure 2):

- A 12' length of ¼" internal diameter tubing with a male luer connector on one end and vacuum source connector on the other end
- A 3' segment of ¼" internal diameter tubing with female connectors on both sides to connect the canister to the vacuum regulator if required.


Figure 1 – Atlas™ Positioner Attachment with tubing and stopcock

Figure 2 – Atlas™ Positioner Attachment accessory vacuum tubing

Refer to Table 1 for available versions of the Atlas Positioner Attachments.

Table 1

Product Reference Number	Type of Quick Connect	Interfacing Stabilization Platforms
T401261	Shaft Quick Connect	Hercules™ 3 Universal Stabilizer Arm (Product Reference: T401161)
T401261U	Ball Quick Connect	Hercules™ 360 Universal Stabilizer Arm (Product Reference: T401160)

INTENDED USE

The Atlas Positioner Attachments are intended for use to position the heart to provide exposure and access during beating heart CABG (coronary artery bypass graft) procedures.

CONTRAINDICATIONS

The Atlas Positioner Attachments are contraindicated for use in patients with local or systemic infection.

Do not position the attachment on newly infarcted tissue, aneurysmal tissue, fragile tissue, or directly over a coronary artery.

WARNINGS

- This device is sterilized using gamma irradiation and is intended **FOR SINGLE USE ONLY. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT REPROCESS. REPROCESSING MAY COMPROMISE THE STERILITY, BIOCOMPATIBILITY AND FUNCTIONAL INTEGRITY OF THE DEVICE.**
- If the packaging is damaged, do not use, as sterility of the device may have been compromised.
- Inspect the product thoroughly for damage prior to use. Do not use damaged product.
- Proper surgical procedure and technique is the responsibility of the attending medical professional (physician). Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- Repeated bending of the malleable zone of the shaft will eventually cause the shaft to break. Slowly bending the malleable zone will increase the allowable range of deformation of the shaft. The integrity of the shaft should be visually verified after bending.
- Excessive flexing of the malleable zone may cause breakage.
- **Take care not to occlude the vacuum lumen prior to obtaining proper positioning.**

- Use of excessive negative vacuum pressure may cause tissue damage.
- Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue.

CAUTION

United States Federal Law restricts this device to sale by or on the order of a physician or other licensed practitioner.

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the positioner attachment and the accessory vacuum tubing from its packaging and pass to the surgical field for attachment to the compatible Hercules Universal Stabilizer Arm (see Table 1).
 - a. For connection to the Hercules 3 Universal Stabilizer Arm:

Slide the sleeve of the Hercules Universal Stabilizer Arm forward and insert the shaft of the attachment into the quick connect with the detent on the shaft (of the attachment) facing the flat area with logo on the sleeve of the Hercules Universal Stabilizer Arm (see Figure 3). Return the sleeve to its original position to lock the accessory device into the quick connect. Check for secure connection.


Figure 3 – Attachment of the Atlas™ Positioner Attachment to the Hercules™ 3 Universal Stabilizer Arm

- b. For connection to the Hercules 360 Universal Stabilizer Arm:

Grasp the sleeve of the Hercules Universal Stabilizer Arm and press the ball of the attachment into the collet of the Hercules Universal Stabilizer Arm (see Figure 4). Check for secure connection.

NOTE: Ensure that the sleeve is grasped during attachment insertion. Failure to grasp the sleeve may result in an inability to insert the attachment.


Figure 4 – Attachment of the Atlas™ 360 Positioner Attachment to the Hercules™ 360 Universal Stabilizer Arm

2. Attach male luer connection of accessory vacuum tubing to the female luer port of the vacuum control stopcock, making sure the vacuum source is switched OFF. Attach the vacuum source connector of the accessory vacuum tubing to the vacuum source and turn the vacuum control stopcock OFF to the vacuum source.
3. Turn the vacuum regulator of the vacuum source on to sufficient negative pressure to affix the Atlas Positioner Attachment to tissue. Due to variations in OR equipment and utilities, use the least amount of vacuum necessary for desired product performance. The regulator should not exceed -250 mmHg.

NOTE: The limit of -250 mmHg is clearly indicated on the blue labels that are found on the accessory vacuum tubing.

WARNING

Use of excessive negative vacuum pressure may cause tissue damage.

4. Turn the vacuum control stopcock OPEN to the vacuum source and positioner attachment. Verify attachment of the positioner to the tissue.

WARNING

Take care not to occlude the vacuum lumen prior to obtaining positive positioning.

5. Position the positioner attachment as required. Turn the main handle of the Hercules 3/360 Universal Stabilizer Arm clockwise until the gap between the compression washers closes. For further instructions refer to the instructions for use of the Hercules 3/360 Universal Stabilizer Arms.
6. To release the tissue from the positioner attachment, turn the vacuum control stopcock OPEN to the atmosphere prior to removal of the positioner attachment.

WARNING

Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue

STORAGE

The Atlas Positioner Attachments storage conditions:

- A clean, dust-free location out of direct sunlight and free from risk of high temperature, excessive humidity, and water leakage.
- Do not store the device near chemicals.

DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

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Terumo® is a registered trademark of TERUMO CORPORATION
Atlas™ and Hercules™ are trademarks of Terumo
Cardiovascular Systems Corporation