REF T401231, T401231LS, T401231U vô trùng

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Titan™ Stabilizör Eklentileri

This page is intentionally left blank.

Titan™ Stabilizer Attachments

Table of Contents

English	5
اللغة العربية	11
Česky	16
Dansk	22
Nederlands.....	28
Suomi	34
Français (France).....	40
Deutsch	46
Ελληνικά.....	52
Magyar	58
Italiano	64
日本語.....	70
Norsk.....	76
Polski	82
Português (Brasil)	88
Русский язык.....	94
简体中文.....	100
Español (España)	105
Svenska	111
Türkçe	117
Graphical Symbols & Definitions.....	123

This page is intentionally left blank.

English

Instructions for Use

Titan™ Stabilizer Attachments

Attachments for Hercules™ 3/360 Universal Stabilizer Arms
Catalogue Numbers T401231, T401231LS, T401231U

Sterile

CAUTION: *Before use, please read carefully the instructions for use of Hercules™ 3/360 Universal Stabilizer Arm in addition to these instructions for use.*

PRODUCT DESCRIPTION

The Titan™ Stabilizer Attachment, Titan™ 360 Stabilizer Attachment, and Titan™ Stabilizer Attachment, Long Shaft are sterile single use accessory devices for use with the Hercules™ 3 or 360 Universal Stabilizer Arms (Product Reference T401161 or T401160). The attachment is provided in a double barrier sterile package and is intended for single use only. The device provides stabilization and coronary artery isolation. The malleable metal design with textured silicone coating allows the surgeon to bend the attachment to conform to anatomic shape in beating heart CABG (coronary artery bypass graft) procedures. The Titan Stabilizer Attachment, Long Shaft has a long attachment shaft with two malleable zones to meet additional positional requirements. The Titan Stabilizer Attachments include a vacuum option to maintain the position of the device.

Included in the sterile package is:

Titan consisting of the following parts (see Figure 1):

- The Titan attachment assembly with an approximately 22" length of ¼" outer diameter tubing with female luer connector on one end to connect to stopcock.
- Hi-flow yellow 4 way stopcock

Accessory vacuum tubing consisting of the following parts (see Figure 2):

- A 12' length of ¼" internal diameter tubing with a male luer connector on one end and vacuum source connector on the other end.
- A 3' segment of ¼" internal diameter tubing with female connectors on both sides to connect the canister to the vacuum regulator if required.


Figure 1 – Titan™ Stabilizer Attachment with tubing and stopcock

Figure 2 – Titan™ Stabilizer Attachment accessory vacuum tubing

Refer to Table 1 for available versions of the Titan Stabilizer Attachments.

Table 1

Product Reference Number	Type of Quick Connect	Interfacing Stabilization Platforms
T401231	Shaft Quick Connect	Hercules™ 3 Universal Stabilizer Arm (Product Reference: T401161)
T401231LS		
T401231U	Ball Quick Connect	Hercules™ 360 Universal Stabilizer Arm (Product Reference: T401160)

INTENDED USE

The Titan Stabilizer Attachments are intended for use in beating heart CABG (coronary artery bypass graft) procedures to stabilize the anastomotic site (coronary artery isolation).

CONTRAINDICATIONS

The Titan Stabilizer Attachments are contraindicated for use in patients with local or systemic infection.

Do not position the Titan Stabilizer Attachments on newly infarcted tissue, aneurysmal tissue, fragile tissue, or directly over a coronary artery.

WARNINGS

- This device is sterilized using gamma irradiation and is intended **FOR SINGLE USE ONLY. DO NOT REUSE. DO NOT RESTERILIZE. DO NOT REPROCESS. REPROCESSING MAY COMPROMISE THE STERILITY, BIOCOMPATIBILITY AND FUNCTIONAL INTEGRITY OF THE DEVICE.**
- If the packaging is damaged, do not use, as sterility of the device may have been compromised.

- Inspect the product thoroughly for damage prior to use. Do not use damaged product.
- Proper surgical procedure and technique is the responsibility of the attending medical professional (physician). Each surgeon must evaluate the appropriate use of this device, case by case, based upon medical training, experience, and the type of surgical procedure employed.
- Repeated bending the malleable zone of the shaft will eventually cause the shaft to break. Slowly bending the malleable zone will increase the allowable range of deformation of the shaft. The integrity of the shaft should be visually verified after bending.
- Excessive flexing of the malleable zone may cause breakage.
- **Take care not to occlude the vacuum lumen prior to obtaining proper positioning.**
- **Use of excessive negative vacuum pressure may cause tissue damage.**
- **Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue.**

CAUTION

United States Federal Law restricts this device to sale by or on the order of a physician or other licensed practitioner.

INSTRUCTIONS FOR USE

1. Using aseptic technique, remove the stabilizer attachment and the accessory vacuum tubing from its packaging and pass to the surgical field for attachment to the compatible Hercules Universal Stabilizer Arm (see Table 1).
 - a. For connection to the Hercules 3 Universal Stabilizer Arm:

Slide the sleeve of the Hercules 3 Universal Stabilizer Arm forward and insert the shaft of the attachment into the quick connect with the detent on the shaft (of the attachment) facing the flat area with logo on the sleeve of the Hercules 3 Universal Stabilizer Arm (see Figure 3). Return the sleeve to its original position to lock the accessory device into the quick connect. Check for secure connection.


Figure 3 – Attachment of the Titan™ Stabilizer Attachment to the Hercules™ 3 Universal Stabilizer Arm

- b. For connection to the Hercules 360 Universal Stabilizer Arm: Grasp the sleeve of the Hercules Universal Stabilizer Arm and press the ball of the attachment into the collet of the Hercules Universal Stabilizer Arm (see Figure 4). Check for secure connection.

NOTE: Ensure that the sleeve is grasped during attachment insertion. Failure to grasp the sleeve may result in an inability to insert the attachment.




Figure 4 – Attachment of the Titan™ 360 stabilizer attachment to the Hercules™ 360 Universal Stabilizer Arm

2. The Titan Stabilizer Attachments can be bent to conform to anatomic requirements as shown in the pictures below indicating the maximum deflection ranges and the proper technique to be used during a procedure.

- The malleable zone of the shaft can be bent to a maximum of 45° in the downward and reversed upward positions as shown below. One bend cycle for the shaft is 0° to 45°.

WARNING

- Flexing and bending of the shaft and foot should not exceed 5 cycles during use. Visually verify the integrity of the shaft and foot after bending.
 - Repeated bending of the malleable zone of the shaft will eventually cause the shaft to break. Slowly bending the malleable zone will increase the allowable range of deformation of the shaft.
 - Excessive flexing of the malleable zone may cause breakage.
 - When bending devices with vacuum option, take care not to occlude the vacuum lumen.
3. Position the stabilizer attachment as required. Turn the main handle of the Hercules 3/360 Universal Stabilizer Arm clockwise until the gap between the compression washers closes. For further instructions refer to the instructions for use of the Hercules 3/360 Universal Stabilizer Arm.

Use of vacuum option with Titan Stabilizer Attachments:

1. Attach male luer connection of accessory vacuum tubing to the female luer port of the vacuum control stopcock, making sure the vacuum source is switched OFF. Attach the vacuum source connector of the accessory vacuum tubing to the vacuum source; turn the vacuum control stopcock OFF to the vacuum source.
2. Turn the vacuum regulator of the vacuum source on to sufficient negative pressure to affix the Titan Stabilizer Attachment to tissue. Due to variations in OR equipment and utilities, use the least amount of vacuum necessary for desired product performance. The regulator should not exceed -400 mmHg.

NOTE: The limit of -400 mmHg is clearly indicated on the yellow labels that are found on the accessory vacuum tubing.

WARNING

Use of excessive negative vacuum pressure may cause tissue damage.

3. Turn the vacuum control stopcock OPEN to the vacuum source and attachment once proper position has been obtained. Verify attachment of the stabilizer to the tissue.
4. To release attachment from tissue, turn the vacuum control stopcock OPEN to the atmosphere prior to removal of the stabilizer attachment.

WARNING

Place vacuum control stopcock OPEN to atmosphere and OFF to the vacuum source prior to removal of attachment from tissue.

STORAGE

The Titan Stabilizer Attachments storage conditions:

- A clean, dust-free location out of direct sunlight and free from risk of high temperature, excessive humidity, and water leakage.
- Do not store the device near chemicals.

DISPOSAL

After use, dispose of product and packaging in accordance with hospital, administrative and/or local, state, federal and international laws and regulations.

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Terumo® is a registered trademark of TERUMO CORPORATION
Titan™ and Hercules™ are trademarks of Terumo Cardiovascular Systems Corporation