
2019-XX

CAPIOX[®] EBS[™] Circuit with X coating[™]

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The logo consists of a stylized, curved arrow-like shape pointing to the right, followed by the word "TERUMO" in a bold, sans-serif font.

 **TERUMO CORPORATION**
44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN
MADE IN JAPAN


® : Registered Trademark  **TERUMO**, EBS and X coating are trademarks of TERUMO CORPORATION. ©TERUMO CORPORATION 2019-XX

2019-XX

CAPIOX[®] EBS[™] Circuit with X coating[™]

Instructions for Use

Please read this Instructions for Use carefully before use, and use this product according to the instructions described here.

 **TERUMO**

Symbols on the packages of this product

Symbols	Instructions
	Catalogue number
	Batch code
	Use-by date
	Manufacturer
	Consult instructions for use
	Do not re-use
	Do not re-sterilize
	Sterilized using ethylene oxide
	Keep away from sunlight
	Keep dry
	Fragile, handle with care
	Do not use if package is damaged
	Temperature limit
	This way up
	Non-pyrogenic

Contents/Description of Product

Contents

Contents/Description of Product

Contents 3

Description of Product 3

Structure 4

Warnings/Precautions

Precautions for operation 5 - 6

Precaution for storage 6

Instructions for Use

Set up 7 - 9

Priming procedure 10 - 12

Connection 13

Initiation of BYPASS 14

During perfusion 14 - 16

Termination of BYPASS 16

Circuit replacement 17

Data 18

Specifications 19

Limited warranty 20

Description of Product


This Product is a kit consisting of membrane oxygenator with microporous polymethylpentene hollow fibers, a centrifugal pump to send blood, and blood tubing. It is intended for use during extracorporeal circulation for cardiovascular surgery and ECLS (Extracorporeal Life Support). Blood and gas flow on the outside and inside of fiber respectively in the membrane oxygenator. Straight path design employed in the centrifugal pump reduces disorder in blood flow and minimizes blood component damage.

The blood contacting surfaces are coated with Xcoating.

This product is used with CAPIOX SP Pump Console SP-101 or CAPIOX Centrifugal Pump Controller SP-200 as a driving instrument of the centrifugal pump.

Structure

Structure

WARNINGS/PRECAUTIONS

Precautions for operation

WARNINGS

• **Adequate heparinization of the blood is required to prevent blood coagulation in the system.**

- Never clamp the venous line (blue) during circulation. Otherwise, pressure inside of the pump will turn into negative, which is a potential hazard of producing gas bubbles.
- Always close the branch line of the venous line (blue) during circulation. Otherwise, pressure in this part will turn into negative, which is a potential hazard of intermixing bubbles.
- **DO NOT OBSTRUCT GAS OUTLET PORT.**
Avoid build-up of excess pressure in the gas phase to prevent gaseous emboli entering the blood phase.
- The oxygenator must be positioned lower than the patient. Pressure in the blood phase should always be higher than that in the gas phase to prevent gaseous emboli entering the blood phase.
- The gas flow should not exceed 20L/min. Excessive gas flow will bring about pressure increase in the gas phase, allowing gaseous emboli to enter the blood phase.
- Do not exceed a 15 °C (27 °F) temperature difference between blood and water in heat exchanger to prevent gas dissolved in blood from forming bubbles.
- This product should be used under monitoring. If abnormalities such as product-damage and deteriorated performance are observed, replace the product with a new one properly after ascertaining the state of patient.
There are possibilities that increasing blood cell damage or blood coagulation is formed, or the circulation is obstructed.

PRECAUTIONS

- This product should be used by a physician or under the supervision and/or direction of a physician.
- This product is sterilized by ethylene oxide gas.
- For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.

- Do not use if the unit package or the product has been damaged or soiled. Use it immediately after opening the unit package and dispose of it safely after use.
- Spare pack of CAPIOX EBS Circuit should be available for emergency exchange.
- Do not allow solvents such as alcohol, ether, acetone, etc. to come in contact with the device as it may cause damage to the device.
- Use CAPIOX SP Pump Console SP-101 and CAPIOX Centrifugal Pump Controller SP-200 according to the instructions described in the Instructions for Use.
- Set this product correctly with a holder or a drive motor.
- The centrifugal pump of this product contains a magnet. Therefore, keep away from metals and magnetic tapes.
- Use this product at flow rate lower than the maximum 7 L/min.
- Do not exceed the maximum pump speed of 3000rpm. Otherwise, it would increase blood damage and cause breakage of this product and leakage.
- Do not install or remove the pump until the motor has completely stopped.
- All procedures must be aseptically performed.
- Perform priming procedure using crystalloid solution which contains no blood, plasma and/or blood derived products.
- For monitoring temperature, Y.S.I. 400 Series manufactured by Measurement Specialties Inc. can be used. Use it with its exclusive cable.
- In case of using CAPIOX EBS Cannula Kit, follow instructions described in the Instructions for Use.
- Water pressure at the heat exchanger inlet should not exceed 196kPa (2kgf/cm²). Pressure greater than 196kPa (2kgf/cm²) may cause leaks or damage to the device.

Precaution for storage

- PRECAUTION** • Pay attention not to wet the device. Avoid direct sunlight, extreme temperature and humidity for storage.

Instructions for Use

Thoroughly read “Warnings/Precautions” (Page 4 - 5) before use.

Instructions for use

Set Up

1. Set a holder and a drive motor.
 - Fix the drive motor to the holder tightly with screws for the drive motor (Figure 1).
2. Fix the holder to the pole with a pole clamp of the holder (Figure 2).
3. Take the product out of the unit package and confirm there is nothing wrong.

- CAUTION**
- Do not use if the unit package or the product has been damaged or soiled.
 - When using CX*XSA0171, CX*XSA0173 or CX*XSA0251, do not hold the product by the tube between oxygenator and centrifugal pump. This may cause heavy load at the connection parts of tubing and lead to leakage.

4. Confirm that all the screws of 3-way stopcocks, caps and connectors are fixed tightly, and that the clamps of the arterial line (red), the venous line (blue) and the priming line are open.

5. Set the centrifugal pump to the drive motor fixed to the holder.
(For CX*XSA: Remove the cover of the centrifugal pump.)

- ① Remove the magnet cover from the drive motor (Figure 3).
- ② Insert the rim of the centrifugal pump into the fixed hook of the drive motor (Figure 4).
- ③ Insert the rim of back side of the oxygenator into the fixed hook of the holder (XX*EB04), and push forward and set the oxygenator (Figure 5).


Figure 1


Figure 2


Figure 3


Figure 4


Figure 5

- CAUTIONS**
- Make sure that the underside of the centrifugal pump is in close contact with the pump receptacle and that the slide hook is in position (Figure 6).
 - Make sure that the rim of back side of the oxygenator and the hook of the exclusive holder are connected to the slide hook (Figure 7).
 - Only when the display of the motor rpm indication shows "0", the centrifugal pump can be attached or detached .


Figure 6


Figure 7

6. Connect the water lines (inner diameter of the tube; 12.7mm (1/2") or couplers) to the water ports of CAPIOX EBS circuit (Figure 8-a, 8-b).
7. Start water circulation through heat exchanger and circulate for at least 5 minutes. Check for leaks.

CAUTION Do not use an oxygenator that leaks. Replace it with another EBS circuit.

8. Connect the thermistor probe and the temperature monitor cable (Figure 8-c).

CAUTION For a monitor, Y.S.I. 400 Series manufactured by Measurement Specialties Inc. can be used. Use it with its exclusive cable.

9. Connect the gas line (inner diameter of the tube; 6.4mm (1/4")) to the gas port (GAS IN) (Figure 8-d).

WARNING Do not block the gas port (GAS OUT).


Figure 8

10. After the circuit was assembled, confirm all parts are connected properly, and the tubes are not folded or kinked.

11. When using with CAPIOX SP Console SP-101:

Attach the CAPIOX SP Pump Flow Sensor 3/8" to the flow connector on the arterial line following the Instructions for Use of CAPIOX SP Pump Console SP-101.

CAUTIONS • Ultrasonic sensor gel must be fully applied to the flow connector.

- When a flow connector is connected with a flow sensor, attach recommended ultrasonic gel (AQUASONIC 100 made by Parker) adequately between the flow sensor and flow connector. Do not use gel which solvents are included such as alcohol, thinner and amine, etc. These solvents may cause damage the sensor and connector.

12. When using with CAPIOX Centrifugal Pump Controller SP-200:

Attach the ABD/Flow Sensor to 3/8" tubing, following the Instruction for Use of CAPIOX Centrifugal Pump Controller SP-200.

CAUTION • When attaching the ABD/Flow Sensor, apply enough vaseline to evenly coat the detection window of the clamp-on ABD/Flow sensor. If the amount of vaseline applied is insufficient, the Flow Signal Unstable Alarm may be issued, and air bubbles may not be detected.

13. After the circuit was assembled, confirm all parts are connected properly, and the tubes are not folded or kinked.

Priming procedure

- CAUTIONS**
- This product should be primed with crystalloid solution which contains no blood, plasma and/or blood derived products.
 - In case of adding blood or plasma, etc. in the priming solution, after priming, follow the procedure below.
 - Infuse blood or plasma, etc. from the priming line.
 - At the same time, remove the clear priming solution from the sampling line to replace in the circuit.

1. Set the circuit at the lower position than the oxygenator.
2. Remove the protection cap from the plastic spike after closing the clamp on the priming line. Then connect the plastic spike to a bag or a bottle filled with priming solution.
3. After confirming that the 3-way stopcock on the priming line is open, release the clamp.
4. Fill the tubing circuit except the centrifugal pump and the oxygenator with priming solution.

CAUTION It may take more time for priming if a large quantity of bubbles remains in the circuit.

5. Fill the oxygenator and the centrifugal pump with priming solution.
6. After the circuit is fully primed or almost fully primed, pull out the handle of the exclusive holder (XX*EB04) (Figure 9).
7. Unlock the rotation of the exclusive holder (XX*EB04) (Figure 10) and rotate the whole holder 90 degrees (Figure 11).

- CAUTIONS**
- After the rotation, lock the rotation, and make sure that the holder is fixed.
 - When rotating 90 degrees, the tube must be put on the holder for rotation not to fold the circuit (Figure 12).
 - Pay attention to prevent the tubing from tucking in during rotation.


Figure 9


Figure 10


Figure 11


Figure 12

8. Press the Auto-Priming button or the icon of the applicable centrifugal pump controller.

The Auto-Priming function works on the centrifugal pump to turn intermittently.

The procedure to expel bubbles in the circuit through microporous of the oxygenator is automatically performed.

CAUTION Do not operate the centrifugal pump without filling it with priming solution. Otherwise, the centrifugal pump may get damaged.

9. Confirm that bubbles in the circuit are completely removed while the centrifugal pump turns intermittently. At this time, remove the cap (with air filter) of the sampling port on the sampling line of the oxygenator and open the 3-way stopcock, and then remove bubbles from the sampling line using a syringe.

CAUTION When bubbles remain at the junction parts of connectors etc., remove them by tapping the parts.

10. Hold the handle for rotation of the exclusive holder (XX*EB04) and unlock the rotation. Then rotate the holder 90 degrees and replace it at the previous position (Figure 13). Put away the handle for rotation depending on the situation.

CAUTION After the rotation, lock the rotation, and make sure that the holder is fixed.


Figure 13

11. Press the Auto-Priming button or the icon again to complete auto priming.
12. Confirm that the motor rpm adjustment knob of the centrifugal pump console is at the position of "0".
13. With the pump outlet clamped, set the pump speed to 3000rpm (maximum speed) for approximately 30 seconds. Since this condition makes internal pressure of the centrifugal pump maximum, check it for something wrong such as leakage from EBS circuit.

CAUTIONS

- Do not operate the centrifugal pump for a long time with its outlet clamped, which will raise the temperature of the priming solution.
- Do not use when something wrong is noted and replace the EBS circuit with a new one.

14. Turn the motor rpm adjustment knob to “0” and stop the rotation of the pump.
15. Pulling up the slide hook of the drive motor, remove the centrifugal pump from the drive motor (Figure 14).
After confirming that the centrifugal pump has been thoroughly primed and there is nothing wrong with it, set the pump to the drive motor again.

CAUTIONS

- Only when the display of the motor rpm indication shows “0”, the centrifugal pump can be attached or detached.
- In case that the bubbles still remain in the centrifugal pump, set the centrifugal pump to the drive motor again, and remove them circulating priming solution.
- After the bubbles were completely removed from the circuit, until the extracorporeal circulation starts, stop the centrifugal pump with its pump speed "0". Do not continue circulation excessively.

16. Close the clamp and the 3-way stopcock on the priming line.

CAUTION

If circulation starts without closing the clamp and the 3-way stopcock, air will enter the circuit. Therefore, be sure to close them.


Figure 14

Connection

1. Close clamps of the arterial line (red) and venous line (blue) from over the clean cover (Figure 15).


Figure 15

2. Remove the clean cover pulling from the corner (Figure 16).


Figure 16

3. Loosen each screw wing for the arterial line and the venous line to remove the bypass connector (Figure 17).


Figure 17

4. Connect the inserted arterial cannula (red) to the arterial line (red) and the inserted venous cannula (blue) to the venous line (blue) respectively (Figure 18) (Figure 19).

CAUTION When the arterial and venous cannulas are connected to corresponding lines, prevent them from kinking. Otherwise, flow rate will not be secured, and it may lead to damage of cannulas.


Figure 18

5. Take off the cap (with air filter) from the 3-way stopcock for removing bubbles, and remove bubbles remaining in the circuit using a syringe.

CAUTION The 3-way stopcock must be closed after completing removal of bubbles.


Figure 19

6. Fix the circuit to the patient or the bed.

CAUTION Pay attention not to fold or kink the circuit.

Initiation of BYPASS

- WARNINGS**
- Prior to starting circulation, confirm that priming line is closed as well as 3-way stopcocks for removing bubbles and those of the sampling port. If the 3-way stopcocks are not closed, bubbles may enter the circuit as the pressure inside the venous line (blue) is negative during circulation.
 - Prior to delivering gas, confirm that the gas port (GAS OUT) is not blocked up. Otherwise, the pressure of the gas side in the oxygenator will increase, and bubbles will flow into the blood side.
 - Make sure that the underside of the centrifugal pump is in close contact with the pump receptacle and that the slide hook is in position.
1. Open the clamp of the venous line (blue) and forceps on the tube of the venous cannula (blue).

WARNING Do not operate the pump with the venous line (blue) clamped. Otherwise, negative pressure would be generated inside the pump and bubbles would be formed in the blood.

2. Start the rotation of the centrifugal pump and set the pump speed at around 1000rpm.
3. Gradually release the clamp on the arterial line (red) to start circulation (Figure 20).

CAUTION After-load pressure higher than the pump's outlet pressure causes back flow of the blood. To prevent this, do not unclamp the outlet of the pump at the start of circulation, until an adequate outlet pressure has been reached. (Refer to the data on Pump performance (19 page).)

4. Start gas supply.

CAUTIONS

- Start gas supply only after blood circulation is initiated.
- Start gas supply with $V/Q = 1$ and $FiO_2 = 100\%$.

5. Adjust the pump speed of the centrifugal pump console for appropriate flow rate.

During perfusion

WARNING Never clamp the venous line (blue) during circulation. Otherwise, negative pressure would be generated inside the pump and bubbles would be formed in the blood.


Figure 20

CAUTIONS

- Carefully adjust the pump speed to control the blood flow, because the blood flow varies depending on the load applied to the pump (the patient's arterial pressure, resistance of the circuit).
- Control blood flow by adjusting the pump speed. Partially clamping the outlet of the pump to control the blood flow may increase the blood damage.
- Entry of massive air into the pump will cause the pump to deprime and the blood flow to stop. When resuming circulation, stop the pump once and remove the air.
- After-load pressure higher than the pump's outlet pressure causes back flow. When decreasing the pump speed, adjust the pump speed taking care not to cause back flow.

1. Collect blood from the sampling line of the oxygenator (Figure 21).

WARNING

Never collect blood from branch lines of the arterial line (red) and the venous line (blue). Otherwise, bubbles may enter the circuit.


Figure 21

2. Measure blood gases and make necessary adjustments as follows:
 - a) Control PaO₂ by changing concentration of oxygen in ventilating gas using gas blender.
 - To decrease PaO₂, decrease FiO₂.
 - To increase PaO₂, increase FiO₂.
 - b) Control PaCO₂ by changing the total gas flow.
 - To decrease PaCO₂, increase total gas flow.
 - To increase PaCO₂, decrease total gas flow.

CAUTION

- A minimum of 0.5L/min oxygen gas flow is needed when blood is circulated. Less than 0.5L/min oxygen gas flow may result in inadequate gas exchange.

- During circulation, always carefully monitor the patient blood condition and the oxygenator gas exchange performance. If gas exchange performance deteriorates, temporarily raise the gas flow rate to flush the inside of the fiber and try to recover performance.

The gas flow rate for flushing is 20 L / min, the time is 10 seconds. If the performance does not recover immediately, it is considered to be caused by plasma leak. Therefore, replace this product with a new one without repeating the flush. A phenomenon called wet lung may occur when water condensation occurs inside fibers of microporous membrane oxygenators with blood flowing exterior to the fibers. This may occur when oxygenators are used for a longer period of time.

Termination of BYPASS

1. At the termination of circulation, clamp the arterial line (red) by forceps etc. taking care that no back flow occurs. After that, stop the pump by adjusting pump speed to 0.
2. Immediately after stopping the pump, stop gas supply.

- | | |
|-----------------|---|
| CAUTIONS | <ul style="list-style-type: none">▪ Do not supply gas while the circulation stops. When gas is supplied in case of recirculation, check blood gas pressure.
Excessive gas supply may cause low PaCO₂, alkalosis and blood damage.▪ When stopping circulation, be sure to clamp the outlet of the pump and then stop the motor. Do not clamp the venous line (blue) first.▪ Do not operate the pump for a long time with its outlet clamped. The temperature within the pump may rise and increase blood damage.▪ After-load pressure higher than the pump's outlet pressure causes back flow. While circulation is terminated, clamp both the arterial line (red) and the venous line (blue). |
|-----------------|---|

Circuit replacement

CAUTIONS

- Always prepare the spare circuit for emergency replacement.
- Exchange all of the circuit when a problem occurs in the oxygenator or the centrifugal pump.
- All procedures must be aseptically performed.


1. Perform priming procedure with the spare applicable centrifugal pump controller or the backup controller. Then clamp the arterial line (red) and the venous line (blue).
2. Disconnect the lines from the cannulas while taking care not to spill the blood.
Method in detail: Clamp the clamping tube of the cannulas and close the clamps on the lines. Loosen the screw rings of the arterial line and the venous line and then disconnect the lines.
3. Connect the cannulas with a newly primed EBS circuit.
4. Open the clamp on the venous line (blue) and the forceps on the venous cannula (blue).
5. Start the rotation of the centrifugal pump and set the pump speed at around 1000rpm.
6. Gradually release the forceps and the clamp on the arterial cannula (red) and the arterial line (red) to start circulation.

CAUTION

After-load pressure higher than the pump's outlet pressure causes back flow. To prevent this, do not unclamp the outlet of the pump at the start of circulation, until an adequate outlet pressure has been reached.


Data

O₂ transfer


<Conditions>	
Blood	:Bovine Blood
Hb	:12±1g/dL
pH	:7.4
SvO ₂	:65±5%
PvCO ₂	:45±5mmHg
B.E.	:0±5mEq/L
V/Q	:1.0
Blood Temperature	:37±1°C

CO₂ transfer

<Conditions>	
Blood	:Bovine Blood
Hb	:12±1g/dL
pH	:7.4
SvO ₂	:65±5%
PvCO ₂	:45±5mmHg
B.E.	:0±5mEq/L
V/Q	:1.0
Blood Temperature	:37±1°C


Heat exchanger performance factor


<Conditions>	
Blood	:Bovine Blood
Hb	:12±1g/dL
T _{bi}	:30±1°C
T _{wi}	:40±1°C
Water flow rate	:15 L/min

Pump performance

SP45

SL45

—	1000 rpm
—	1500 rpm
—	2000 rpm
—	2500 rpm
—	3000 rpm

Specifications

Specifications

Code number	CX*XSA	CX*XSA0171	CX*XSA0173	CX*XSA0251
Sterilization	Ethylene oxide gas			
Priming volume	550 mL	620mL	620mL	740mL
Centrifugal pump	SP45	SL45	SL45	SL45
Maximum blood flow rate	7 L/min			
Minimum blood flow rate	0.5 L/min			
Maximum pressure	1000 mmHg (133 kPa)			
Connector for venous cannula	Lock connector			
Connector for arterial cannula	Lock connector			

Items	Specifications
Oxygenator	
Fiber surface area/material	Approx. 2.5m ² / Poly methyl pentene
Heat exchanger surface area/material	Approx. 0.22m ² / Stainless steel
Purge port	Luer port
Gas port (Gas in)	1/4" (6.4 mm)
Gas port (Gas out)	1/4" (6.4 mm)
Water port	1/2" (12.7mm) Couples connect fitting
Thermistor probe	Blood outlet port
Centrifugal pump	
Maximum pump speed	3000rpm
Maximum outlet pressure	800mmHg (107kPa)
Applicable console	CAPIOX SP Pump Console SP-101 (Code number: ME*SP101K, ME*SP101K2, ME*SP101K3, ME*SP101S, ME*SP101S2, ME*SP101S3)
	CAPIOX Centrifugal Pump Controller SP-200 (Code number: ME*SP200C, ME*SP200C2, ME*SP200C3)
Applicable drive motor	CAPIOX SP Pump Drive Motor (Code number) ME*SP101M, ME*SP200M
Applicable flow sensor	CAPIOX SP Pump Flow Sensor 3/8" (9.5 mm) (Code number: ME*FS38) ABD/Flow Sensor (Code number: ME*SPFAS01, ME*SPFAS02)

Accessories (Sold separately)

Applicable holder	CAPIOX EBS Circuit Holder (Code number: XX*EB04)
Applicable cannula	CAPIOX EBS Cannula KIT For the code number, refer to the Instructions for Use of the Cannula Kit.

Applicable temperature monitor	YSI temperature YSI400 Series (Measurement Specialties Inc.)
--------------------------------	--

Disclosable Information

- How to sterilize
- List of raw materials used for blood path
- Data on plasma leakage of semipermeable membrane
- Pressure drop of the blood path at the range of blood flow rates specified by the manufacturer for clinical use
- Pressure drop of the gas path at the maximum blood and gas flow rate specified by the manufacturer for clinical use
- Information on blood cell damage
- Information on particulate emissions from the membrane oxygenator with microporous Poly methyl pentene hollow fiber
- Tolerance range of data provided

Limited Warranty

TERUMO Corporation warrants that the product has been manufactured and packaged with reasonable care, and will be free from defects under normal and proper use. TERUMO Corporation will not be liable for any direct or indirect, incidental, special, general or consequential loss, damage, or expense from the use of this product. TERUMO Corporation's sole obligation under this warranty shall be to repair or replace, at its option, any product that it believes was defective at the time of shipment if notice thereof is received before the expiration date indicated on the package.

THIS WARRANTY IS THE SOLE AND ONLY WARRANTY OF TERUMO CORPORATION AND IS EXPRESSLY IN LIEU OF AND EXCLUDES ANY OTHER EXPRESSED OR IMPLIED WARRANTY OF ANY KIND, INCLUDING THE WARRANTY OF MERCHANTABILITY OR FITNESS FOR ANY PARTICULAR PURPOSE. IN ADDITION, THE REMEDIES SET FORTH HEREIN ARE THE SOLE REMEDIES AVAILABLE TO ANY PERSON OR ENTITY.

No agent, employee or representative of TERUMO Corporation has any Authority to change, Amend or alter the foregoing or assume or bind TERUMO Corporation to any additional liability or responsibility in connection with this product.

The logo features the word "TERUMO" in a bold, sans-serif font. A curved line above the letters "E" and "R" suggests a stylized shape or motion.

 **TERUMO CORPORATION**
44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN
MADE IN JAPAN

® : Registered Trademark  **TERUMO**, EBS and X coating are trademarks of TERUMO CORPORATION. ©TERUMO CORPORATION 2019-XX