TERUMOBCT

41910

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
2019-11
Part No. 777021-001B

TERUMOBCT

This product may be covered by one or more patents or pending patent applications.
See TERUMOBCT.COM/patents for details.

CE 2797

Terumo BCT, Inc.
10811 W. Collins Avenue
Lakewood, Colorado 80215
USA
USA Phone: +1.877.339.4228
Phone: +1.303.231.4357
USA Fax: +1.866.715.6768
Fax: +1.303.542.5215

Terumo BCT Europe N.V.
Ikaraoslaan 41
1930 Zaventem
Belgium
Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

TERUMOBCT.COM

©2019 Terumo BCT, Inc.

41910

Platelet Pooling Set
Instructions for Use
Trombocytpoolingsæt
Brugsanvisning til
Thrombozyten-Pooling-Set
Gebrauchsanweisung
Σετ συγκέντρωσης αιμοπεταλίων
Οδηγίες χρήσης
Equipo de mezcla de plaquetas
Instrucciones de uso
Kit de mélange plaquettaire
Mode d'emploi
Circuito per pooling delle piastrine
Istruzioni per l'uso
Poolingsett for blodplater
Bruksanvisning
Zestaw do zlewania płytek krwi
Instrukcja używania
Conjunto de Mistura de Plaquetas
Instruções de Uso
Conjunto de Mistura de Plaquetas
Instruções de Utilização
Комплект для пулирования тромбоцитов
Инструкция по эксплуатации
Trombocytpoolingsset
Bruksanvisning
Trombosit Havuzlama Seti
Kullanım Yönergeleri


2019-11

Part No. 777021-001B

Intended Use

The Platelet Pooling Set is intended for use with whole-blood-derived platelets, including pooling and filtration, resulting in a platelet product for transfusion. The pooling kit has only been evaluated on the Atreus system and the Reveos system. Therefore, if it will be used with other systems an evaluation will need to be performed according to the institution's SOP.

Set Diagram

- 1 End cap
- 2 Platelet medium tubing clamp
- 3 Blue-striped platelet medium tubing
- 4 Pooling tubes
- 5 Pooling tubing clamp
- 6 Filter
- 7 Red-striped restrictive tubing
- 8 Pooled platelet storage bag
- 9 Sample tubing clamp
- 10 Sample bag
- 11 Needle tubing clamp
- 12 Needle and needle cap

Storage Conditions

- Long-term storage temperature: 1 °C to 35 °C
- Permitted temperature excursions:
 - -29 °C to 1 °C for up to 72 hours
 - Up to 50 °C for up to 6 weeks
- Storage humidity: 0% to 75% relative humidity (RH), non-condensing
- Permitted humidity excursions:
 - Up to 85% (\pm 5%) RH for up to 72 hours
- Protect platelet pooling sets from freezing and excessive heat during storage.

Warnings

1. **Do Not Reuse/Not for Reuse:** Terumo BCT, Inc. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

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Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
 - Fluid leaks
 - Parts that are warped or deformed
 - Plastics that are brittle and discolored
 - Filters that have reduced filtration capabilities
- Exposure to excessive ethylene oxide (EO) residuals
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

2. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.
3. Inspect the packaging and the platelet pooling set prior to use. Do not use the set if any of the following conditions are present:
 - There are tears or holes in the outer packaging.
 - The tubing has severe kinks.
 - The platelet pooling set is incorrectly assembled.
 - The platelet pooling set is defective or damaged.
 - Any clamps are closed.
 - The end cap is not in place.
 - The needle cap is not in place.
4. Seal the tubing between the sample bag and the pooled platelet storage bag before you open the clamp on the needle tubing. Failure to do so may result in contamination of the pooled Platelet product.

Using the platelet pooling set under these conditions may result in inadequate leukoreduction or contamination of the pooled Platelet product.

Cautions

1. When filtering platelets, ensure that the tubing is free of kinks or other obstructions, because these may result in inadequate leukoreduction of the pooled Platelet product.
2. Do not squeeze the IPU bags to increase the filtration rate. Squeezing the IPU bags may result in platelet activation and inadequate leukoreduction.
3. The wall thickness of the red-striped restrictive tubing is greater than the wall thickness of the other tubing on the platelet pooling set. Do not seal the red-striped restrictive tubing without first ensuring that the tubing meets the specifications of the tube sealing devices used in your facility.

Notes

1. The blood and fluid pathways of the platelet pooling set are sterilized with ethylene oxide and are non-pyrogenic.
2. The time ranges used to define “overnight,” “same day as collection,” and “day after collection” are approximate values. You will need to validate specific time ranges according to your facility’s standard operating procedures (SOPs).
3. When making a sterile connection to the platelet pooling set, always follow the instructions for use provided by the manufacturer of the sterile connection device.

Procedures for Use

This document provides platelet pooling and filtration instructions specific to the platelet pooling set. For warnings, cautions, and instructions on processing whole blood with either the Atrius system or the Reveos system, see the system operator’s manual.

The process for producing a pooled Platelet product from whole-blood-derived interim platelet units (IPUs) consists of the following parts:

1. Resting and agitating the IPUs
2. Pooling and filtering four to six IPUs using the platelet pooling set, according to one of the following methods:
 - Pooling and filtering four to six IPUs with a platelet medium
Either platelet additive solution (PAS) or an additional unit of plasma may be used as a platelet medium.
 - Pooling and filtering four to six plasma-suspended IPUs that contain sufficient plasma for storage when pooled and do not require the addition of a platelet medium

Note: Adjusting the volume of plasma in the IPUs requires the assistance of a Terumo BCT implementation specialist.

3. Storing the pooled Platelet product

The platelet pooling set includes a sample bag that can be used to collect samples of the pooled Platelet product after pooling and filtration are complete. Follow your institution’s SOP for collecting samples from the pooled Platelet product.

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Recommended Conditions for Resting and Agitating the IPU's

Before you begin the pooling and filtration process, rest and agitate the individual IPU's. The recommended conditions for resting and agitating the IPU's depend on when the unit of whole blood was processed.

- If the unit of whole blood was processed the same day as it was collected (between 2 hours and 14 hours after collection), follow the recommendations in Table 1.
- If the unit of whole blood was processed the day after it was collected (between 14 hours and 24 hours after collection), follow the recommendations in Table 2.

Table 1: Conditions for resting and agitating IPU's from whole blood units processed on the same day as collection

Action	Temperature Range	Time
Rest the IPU	18 °C to 28 °C	1 hour
Agitate the IPU	20 °C to 24 °C	Overnight
Pool and filter the IPU	18 °C to 28 °C	Within 36 hours of processing the whole blood unit on the device

Table 2: Conditions for resting and agitating IPU's from whole blood units processed the day after collection

Action	Temperature Range	Time
Rest the IPU	18 °C to 28 °C	1 hour
Agitate the IPU	20 °C to 24 °C	At least 2 hours
Pool and filter the IPU	18 °C to 28 °C	Within 36 hours of processing the whole blood unit on the device

Pooling and Filtering the Interim Platelet Units

In the following instructions, steps highlighted with a gray background apply only when you are pooling and filtering IPU's with an additional platelet medium (PAS or plasma). If you are pooling plasma-suspended IPU's, the gray-highlighted steps do not apply.

1. Prepare a leukoreduction rack. Ensure that the head height of the leukoreduction rack is sufficient for the pooled platelet storage bag to hang vertically.
2. Close all the clamps on the platelet pooling set.
3. Use a sterile connection device to connect the platelet medium bag to the blue-striped platelet medium tubing. Ensure that the remaining blue-striped tubing on the set is long enough to allow the platelet medium bag to hang so that the fluid level of the platelet medium is higher than the fluid level of the IPU's.
4. Use a sterile connection device to connect four to six IPU's to the pooling tubes.
5. Hang all of the IPU's from the same filtration hook.

6. Hang the platelet medium bag from a separate filtration hook. Ensure that the fluid level of the platelet medium is higher than the fluid level of the IPU's.
7. Ensure that all tubing is free of kinks or other obstructions and that the leukoreduction filter hangs vertically.
8. After all air bubbles flow into the IPU bags, open the platelet medium tubing clamp and gently squeeze the platelet medium bag.
This helps remove any additional air from all the tubing and forces the platelet medium to flow into the interim platelet unit bags.
9. Allow the platelet medium to flow into the IPU bags for at least 30 seconds.
10. Ensure that there are no air bubbles larger than 1 cm in any of the tubing.
11. Open the pooling tubing clamp to allow the IPU's (and platelet medium, if used) to flow through the leukoreduction filter and into the pooled platelet storage bag.
12. After filtration is complete, close the platelet medium tubing clamp and the pooling tubing clamp.
13. Seal the tubing between the red-striped restrictive tubing and the pooled platelet storage bag.
14. Disconnect the pooled platelet storage bag from the filtration tubing. Do not remove the sample bag that is attached to the pooled platelet storage bag.
15. After use, dispose of the empty IPU bags, the empty platelet medium bag (if applicable), and the leukoreduction filter assembly safely according to your institution's SOP and/or local regulations.
16. Open the sample tubing clamp and express residual air from the pooled platelet storage bag into the sample bag, if necessary.
17. Close the sample tubing clamp and agitate the pooled platelet storage bag to ensure adequate mixing.
18. Store the pooled Platelet product according to your institution's SOP. Recommended storage conditions are provided in Table 3.

Table 3: Storage conditions for the pooled Platelet product in the pooled platelet storage bag

Storage Parameter	Recommended Ranges for Storage in PAS	Recommended Ranges for Storage in Plasma
Volume per bag	100 mL to 400 mL	100 mL to 400 mL
Platelet concentration	0.5×10^6 platelets/ μ L to 2.1×10^6 platelets/ μ L	0.5×10^6 platelets/ μ L to 2.1×10^6 platelets/ μ L
Maximum total platelets per bag ([volume \times concentration] $\times 10^3$)	4.5×10^{11} platelets	5.1×10^{11} platelets
Storage temperature	22 °C \pm 2 °C	22 °C \pm 2 °C

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Sampling the Pooled Platelet Product

- Open the sample tubing clamp, transfer the desired volume of pooled Platelet product into the sample bag, and mix the product in the sample bag according to your institution's SOP.
- Seal the tubing that connects the sample bag to the pooled platelet storage bag. Follow your institution's SOP for sampling the pooled Platelet product using the sample bag.
- Disconnect the sample bag from the pooled platelet storage bag.
- Dispose of the sample bag safely according to your institution's SOP and/or local regulations.

Symbols and Certification

Symbol	Definition	Symbol	Definition
	Indicates the product quantity when the quantity is placed in the square.		Indicates the date after which the medical device is not to be used.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.		Indicates the manufacturer's catalog number so that the medical device can be identified.
	Indicates a non-pyrogenic fluid pathway.		Indicates the presence of a sterile fluid path. The method of sterilization is ethylene oxide (EO).
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.		Indicates a medical device that is not to be re-sterilized.
	Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended.		Indicates the Authorized representative in the European Community.
	Indicates the medical device manufacturer.		Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.
	Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).		Indicates a medical device that needs protection from light sources.
	Indicates the range of temperatures to which the device can be safely exposed.		Indicates that the product should not be used if the package is damaged.
	Indicates a medical device that needs to be protected from moisture.		Indicates that the contents of the transport package are fragile; therefore, it shall be handled with care.
	Indicates the correct upright position of the transport package.		Indicates that the transport package is not allowed to be vertically stacked beyond the specified number.

Symbol	Definition	Symbol	Definition
	Indicates that the user should consult the instructions for use.		

Return Of Used Product

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department. **IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT. Please contact your local representative for information regarding returned goods and product complaints.**

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