TERUMOBCT

4NG456S0

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Terumo BCT Ltd.
Old Belfast Road
Millbrook
Larne
BT40 2SH
United Kingdom
Phone: +44.28.2827.3631
Fax: +44.28.2827.3719

TERUMOBCT.COM

©2018 Terumo BCT, Inc.

2019-11
Part No. GBB-IURVNGS01911

TERUMOBCT


GBB-IURVNGS01911

This product may be covered by one or more patents or pending patent applications.
See TERUMOBCT.COM/patents for details.

Made in Vietnam.

Terumo BCT Ltd.
Old Belfast Road
Millbrook
Larne
BT40 2SH
United Kingdom
Phone: +44.28.2827.3631
Fax: +44.28.2827.3719

Terumo BCT Europe N.V.
Ikarooslaan 41
1930 Zaventem
Belgium
Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

TERUMOBCT.COM

©2019 Terumo BCT, Inc.

4NG456S0

Reveos® NLR Set
CPD/SAG-M Blood Bag Set
Instructions for Use
Σετ Reveos® NLR
Σετ ασκών αίματος CPD/SAG-M
Οδηγίες χρήσης
Equipo Reveos® NLR
Conjunto de bolsas de sangre con CPD/SAG-M
Instrucciones de uso
Kit NLR Reveos®
Kit de poches de sang CPD/SAG-M
Mode d'emploi
Zestaw Reveos® NLR
Zestaw pojemników na krew z CPD/SAG-M
Instrukcja używania
Conjunto Reveos NLR
CPD/SAG-M - Conjunto de Bolsa de Sangue
Instruções de Uso
Комплект Reveos NLR (без функции лейкоредукции)
Набор мешков для крови, содержащий консервант ЦФД/САГМ
Инструкция по эксплуатации

2019-11


Part No. GBB-IURVNGS01911

Intended Use and Description

- The Reveos® NLR Set is intended to collect a unit of whole blood and to process the whole blood unit on the Reveos Automated Blood Processing System, producing blood components.
- The Reveos NLR Set has been evaluated for use only with the Reveos system.
- The blood and fluid pathways of the blood bag set are steam-sterilized and are non-pyrogenic.

CPD Anticoagulant Description	SAG-M Additive Solution Description
Each 100 mL of CPD (citrate phosphate dextrose) anticoagulant contains: <ul style="list-style-type: none"> • Acidum citricum (anhydr.) 0.299 g • Natrii citras (dihydr.) 2.63 g • Mononatrii phosphas (monohydr.) 0.222 g • Dextrosum (monohydr.) 2.55 g • Aqua ad iniectionem q.s. 100 mL 	Each 100 mL of SAG-M (saline, adenine, glucose, mannitol) additive solution contains: <ul style="list-style-type: none"> • Natrii chloridum 0.877 g • Adeninum 0.0169 g • Dextrosum (anhydr.) 0.818 g • Mannitolum 0.525 g • Aqua ad iniectionem q.s. 100 mL

Set Diagram

- | | | |
|----------------------------|-----------------------------------|--------------------------------------|
| 1 Whole blood bag | 6 Sample/diversion bag | A, C: Y connector |
| 2 Platelet bag | 7 Needle | B, G: Blue clamp |
| 3 Plasma bag | 8 Needle injury protector (NIP) | D, F, H: CLIKTIP frangible connector |
| 4 Red blood cell (RBC) bag | 9 Sample tube holder/luer adapter | E: White clamp |
| 5 Residual leukocyte bag | 10 Cross-connector | |

This page intentionally left blank.

This page intentionally left blank.

Reveos NLR Storage Conditions

- Long-term storage temperature: 1 °C to 30 °C
- Permitted temperature excursions:
 - -20 °C to 1 °C for up to 2 weeks
 - Up to 50 °C for up to 1 week
- Storage humidity: 35% to 75% relative humidity (RH), non-condensing
- Permitted humidity excursions:
 - 15% to 35% RH for up to 6 months
 - Up to 85% (\pm 5%) RH for up to 72 hours

Use blood bag sets within 28 days after you open the outer aluminum foil pouch. To store unused blood bag sets, return them to the outer aluminum foil pouch and reclose the pouch with tape or a clip. Once you open the transparent packing wrap, you must use the blood bag set within 7 days, not exceeding 28 days from when you opened the outer aluminum foil pouch. Each outer aluminum foil pouch contains a sachet that absorbs oxygen. Dispose of the sachet and the outer aluminum foil pouch with normal waste.

Warnings

1. **Do Not Reuse/Not for Reuse:** Terumo BCT Ltd. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
 - Fluid leaks
 - Parts that are warped or deformed
 - Plastics that are brittle and discolored
 - Filters that have reduced filtration capabilities
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

2. This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.
3. Inspect the packaging and the blood bag set prior to use. Do not use the set if any of the following conditions are present:
 - There are tears or holes in the outer aluminum foil packaging or in the individual transparent packaging wrap.
 - The tubing has severe kinks.

This page intentionally left blank.

- The blood bag set is incorrectly assembled.
- The blood bag set is defective or damaged, or there are any leaks from the fluid-filled components of the set.
- Any clamps are closed.
- The needle cap is not in place.
- The solutions are cloudy or discolored or contain particulates.

Using the blood bag set under these conditions may result in product contamination or poor performance during collection and/or processing.

Note: It is normal to have some condensation in the outer aluminum foil pouch and individual transparent packaging wrap due to sterilization.

4. The blood bag set is no longer sterile if any of the following conditions occur:
 - You disconnect the sample/diversion bag before you seal the sample/diversion line.
 - You remove blood samples before you seal the sample/diversion line.
 - The integrity of the set is compromised for any reason.

Manage the blood bag set according to your institution's standard operating procedure (SOP).

Cautions

1. Do not process whole blood less than 2 hours after collection. Processing blood too soon after collection may result in a reduced platelet yield.
2. Residual leukocytes are a by-product of the processing procedure and contain mostly white blood cells (WBC), with some plasma, platelets, and red blood cells (RBC). Residual leukocytes are not intended for transfusion.
3. Use aseptic technique during blood collection to ensure donor safety and product quality.
4. Do not vent the blood bag set.

Notes


1. Process whole blood collected in the Reveos NLR set within 8 hours of collection.

Procedures for Use

This document provides blood collection and post-processing instructions specific to the Reveos NLR set. For warnings, cautions, and instructions on processing whole blood with the Reveos system, see the Reveos system operator's manual.




Blood Collection

1. If it is part of your institution's standard operating procedure (SOP), make a loose knot in the collection tubing between white clamp [E] and the whole blood bag.
2. Separate the sample/diversion bag from the other bags.
3. Load the bags onto a scale and/or blood mixing device according to your institution's SOP. Ensure that the whole blood bag and the sample bag are lower than the donor's arm.
4. Close blue clamp [B].
5. Prepare the venipuncture site and perform the venipuncture according to your institution's SOP.

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6. Open blue clamp .
7. Position the sample/diversion bag with blue clamp  at the top and the sample tube holder/luer adapter at the bottom (see Figure 1). Allow the desired volume of whole blood to flow into the sample/diversion bag. The nominal volume of the sample bag is 60 mL. When the volume reaches the marking, approximately 40 mL has been collected.
8. Close blue clamp .
9. Break CLIKTIP  to allow the whole blood to flow into the whole blood bag. Bend the CLIKTIP in both directions to ensure that you break it completely.
10. Seal the sample/diversion tubing as near as possible to Y connector . Seal the tubing according to your institution's SOP.
11. Invert the sample/diversion bag and transfer donor blood samples from the sample/diversion bag, using evacuated blood collection tubes. Transfer samples as soon as possible after venipuncture to avoid possible clot formation in the sample/diversion bag.
 - a. Hold the sample/diversion bag and sample tube holder/luer adapter in one hand with blue clamp  positioned at the bottom (see Figure 2).
 - b. Using the other hand, insert a blood collection tube firmly into the sample tube holder/luer adapter. After filling the blood collection tube, remove the tube from the sample tube holder. Repeat the process to take additional samples.
12. Mix the whole blood and the anticoagulant during collection according to your institution's SOP.
13. Collect the target volume of blood $\pm 10\%$, as specified on the whole blood bag label.
14. Close white clamp  and perform one of the following steps:
 - If you made a knot in the collection tubing in step 1, tighten the knot firmly.
 - If you did not make a knot in step 1, seal the collection tubing near CLIKTIP .
15. Hold the needle hub and then slide the NIP partially over the needle hub (see Figure 3).
16. Gently pull on the tubing to withdraw the needle from the donor's arm and into the NIP until the needle locks securely (see Figure 4).
17. Seal the collection tubing near CLIKTIP  if you did not do so in step 14. Disconnect the sealed tubing that includes the needle and the sampling assembly. Dispose of the needle, the NIP, and the sample tube holder safely according to your institution's SOP and/or local regulations.
18. Immediately after collection is complete, invert the whole blood bag several times to thoroughly mix the whole blood and the anticoagulant.
19. If required, strip the blood in the collection tubing into the whole blood bag according to your institution's SOP.
20. Seal and disconnect the collection tubing between 0.5 cm and 3.0 cm away from the whole blood bag. Dispose of the collection tubing safely according to your institution's SOP and/or local regulations.
21. Place the whole blood unit into a temperature-controlled environment according to your institution's SOP.
22. Pack and transport the whole blood unit to the processing laboratory according to your institution's SOP.


Figure 1: Sample/diversion bag


Figure 2: Collecting a whole blood sample


Figure 3: Sliding the NIP over the needle hub


Figure 4: Withdrawing the needle into the NIP

Post-Processing

1. After processing, remove the plasma bag, the residual leukocyte bag, and the platelet bag from the organizer and handle them according to your institution's SOP.
2. Remove the RBC bag from the organizer. The RBC bag is still attached to the whole blood bag.
3. Hang the RBC bag.
4. Break CLIKTIP [H] and open blue clamp [G].
5. Drain all of the additive solution into the whole blood bag.
6. Close blue clamp [G].
7. Gently mix the additive solution with the RBC unit in the whole blood bag.
8. Hang the whole blood bag.
9. Ensure that the tubing is free of kinks or other obstructions.
10. Open blue clamp [G] and allow the mixture of red blood cells and additive solution to flow into the RBC storage bag.
11. Seal the RBC tubing below CLIKTIP [H], and disconnect the RBC bag.
12. Dispose of the whole blood bag safely according to your institution's SOP and/or local regulations.
13. Seal the line for sampling segments according to your institution's SOP.
14. Store the RBC product according to your institution's SOP.

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Symbols and Certification

Symbol	Definition	Symbol	Definition
	Indicates the product quantity.		Indicates the date after which the medical device is not to be used.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.		Indicates the manufacturer's catalog number so that the medical device can be identified.
	Indicates a non-pyrogenic fluid pathway.		Indicates the presence of a sterile fluid path. The method of sterilization is steam.
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.		Indicates a medical device that is not to be re-sterilized.
	Indicates the medical device manufacturer.		Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile).
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.		Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).
	Indicates a medical device that needs protection from light sources.		Indicates a medical device that needs to be protected from moisture.
	Indicates the range of temperatures to which the device can be safely exposed.		Indicates that the contents of the transport package are fragile; therefore, it shall be handled with care.
	Indicates the correct upright position of the transport package.		Indicates that the transport package is not allowed to be vertically stacked beyond the specified number.
	Indicates that the user should consult the instructions for use.		Indicates that the product should not be used if the product, its sterile barrier system, or its packaging is damaged or shows any signs of deterioration.
	Indicates that the product should not be used if the package is damaged.		Indicates that the product should not be vented. Do not remove the end caps from the set.
	Indicates a donor needle.		Indicates the length of time within which the product must be used after the packaging is opened.

Symbol	Definition	Symbol	Definition
	Indicates the nominal volume of the blood container.		Indicates a whole blood bag.
	Indicates the target whole blood collection volume.		Indicates a red blood cell bag.
	Indicates a plasma bag.		Indicates a platelet bag.
	Indicates anticoagulant solution.		Indicates additive solution.
	Indicates a sampling site.		Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended.
	Indicates the Authorized representative in the European Community.		

Return of Used Product

If for any reason this product must be returned to Terumo BCT Ltd., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department. **IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT. Please contact your local representative for information regarding returned goods and product complaints.**

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