TERUMOBCT

4FG456S0

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TERUMOBCT


GBB-IURVFGS01912

This product may be covered by one or more patents or pending patent applications. See TERUMOBCT.COM/patents for details.

Made in Vietnam.

Terumo BCT Ltd.
Old Belfast Road
Millbrook
Larne
BT40 2SH
United Kingdom
Phone: +44.28.2827.3631
Fax: +44.28.2827.3719

Terumo BCT Europe N.V.
Ikarostraat 41
1930 Zaventem
Belgium
Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

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This page intentionally left blank.


This page intentionally left blank.

Intended Use and Description

- The Reveos® LR Set is intended to collect a unit of whole blood and to process the whole blood unit on the Reveos Automated Blood Processing System, producing blood components.
- The Reveos LR Set has been evaluated for use only with the Reveos system.
- The blood and fluid pathways of the blood bag set are steam-sterilized and are non-pyrogenic.

CPD Anticoagulant Description	SAG-M Additive Solution Description
<p>Each 100 mL of CPD (citrate phosphate dextrose) anticoagulant contains:</p> <ul style="list-style-type: none"> • Acidum citricum (anhydr.) 0.299 g • Natrii citras (dihydr.) 2.63 g • Mononatrii phosphas (monohydr.) 0.222 g • Dextrosum (monohydr.) 2.55 g • Aqua ad iniectionem q.s. 100 mL 	<p>Each 100 mL of SAG-M (saline, adenine, glucose, mannitol) additive solution contains:</p> <ul style="list-style-type: none"> • Natrii chloridum 0.877 g • Adeninum 0.0169 g • Dextrosum (anhydr.) 0.818 g • Mannitolum 0.525 g • Aqua ad iniectionem q.s. 100 mL

Set Diagram

- | | | |
|-----------------------------------|-----------------------------------|---|
| 1 Whole blood bag | 6 Sample/diversion bag | 11 RBC leukoreduction filter |
| 2 Interim platelet unit (IPU) bag | 7 Needle | A, C: Y connector |
| 3 Plasma bag | 8 Needle injury protector (NIP) | B, G: Blue clamp |
| 4 Red blood cell (RBC) bag | 9 Sample tube holder/luer adapter | D, F, H, J: CLIKTIP frangible connector |
| 5 Residual leukocyte bag | 10 Cross-connector | E: White clamp |

Reveos LR Storage Conditions

- Long-term storage temperature: 1 °C to 30 °C
- Permitted temperature excursions:
 - -20 °C to 1 °C for up to 2 weeks
 - 30 °C to 50 °C for up to 1 week
- Storage humidity: 0% to 75% relative humidity (RH), non-condensing

Use blood bag sets within 28 days after you open the outer aluminum foil pouch. To store unused blood bag sets, return them to the outer aluminum foil pouch and re-close the pouch with tape or a clip. Once you open the transparent packing wrap, you must use the blood bag set within 7 days, not exceeding 28 days from when you opened the outer aluminum foil pouch.

Warnings



- Do Not Reuse/Not for Reuse:** Terumo BCT Ltd. products bearing the "Do Not Reuse" symbol are intended for single use only and are not intended to be reused or re-sterilized in any manner. Terumo BCT cannot ensure the functionality or sterility of the product if it is reused or re-sterilized.

Reuse of a single-use product could result in:

- Product performance issues due to a loss of product integrity, including but not limited to the following:
 - Fluid leaks
 - Parts that are warped or deformed
 - Plastics that are brittle and discolored
 - Filters that have reduced filtration capabilities
- Viral infections such as hepatitis or human immunodeficiency virus (HIV)
- Bacterial infections
- Cross-contamination

Any of these risks could result in serious injury or death. These risks are shared by product users, donors, patients, and recipients of end products of the device.

- This product contains phthalates (DEHP). Whole blood donors are not exposed to DEHP; the potential health risk to apheresis donors is low, because the time-averaged DEHP dose exposure is very low. Patient groups that include pregnant or nursing women and children are considered to be the most at risk to potential harmful effects of exposure to DEHP. However, regulatory bodies have noted that the benefit of doing a needed procedure is far greater than the risk associated with exposure to DEHP. It is the responsibility of the treating physician to balance this risk for the patient.
- Inspect the packaging and the blood bag set prior to use. Do not use the set if any of the following conditions are present:
 - There are tears or holes in the outer aluminum foil packaging or in the individual transparent packaging wrap.
 - The tubing has severe kinks.
 - The blood bag set is incorrectly assembled.
 - The blood bag set is defective or damaged, or there are any leaks from the fluid-filled components of the set.
 - Any clamps are closed.
 - The needle cap is not in place.
 - The solutions are cloudy or discolored or contain particulates.

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Using the blood bag set under these conditions may result in product contamination or poor performance during collection and/or processing.

Note: It is normal to have some condensation in the outer aluminum foil pouch and individual transparent packaging wrap due to sterilization.

4. The blood bag set is no longer sterile if any of the following conditions occur:
 - You disconnect the sample/diversion bag before you seal the sample/diversion line.
 - You remove blood samples before you seal the sample/diversion line.
 - The integrity of the set is compromised for any reason.

Manage the blood bag set according to your institution's standard operating procedure (SOP).

Cautions

1. Do not process whole blood less than 2 hours after collection. Processing blood too soon after collection may result in incomplete RBC leukoreduction and/or a reduced platelet yield.
2. Residual leukocytes are a by-product of the processing procedure and contain mostly white blood cells (WBC), with some plasma, platelets, and red blood cells (RBC). Residual leukocytes are not intended for transfusion.
3. Use aseptic technique during blood collection to ensure donor safety and product quality.
4. Do not vent the blood bag set.
5. Precautions for the oxygen absorber sachet:
 - Each aluminum foil pack contains a sachet that absorbs oxygen and generates heat on contact with air.
 - Handle the oxygen absorber sachet with caution.
 - Do not open the oxygen absorber sachet.
 - Do not dispose of the oxygen absorber sachet together with waste containing volatile or flammable material.

Procedures for Use

This document provides blood collection and post-processing instructions specific to the Reveos LR set. For warnings, cautions, and instructions on processing whole blood with the Reveos system, see the Reveos system operator's manual.

Blood Collection

1. If it is part of your institution's standard operating procedure (SOP), make a loose knot in the collection tubing between white clamp [E] and the whole blood bag.
2. Separate the sample/diversion bag from the other bags.
3. Load the bags onto a scale and/or blood mixing device according to your institution's SOP. Ensure that the whole blood bag and the sample bag are lower than the donor's arm.
4. Close blue clamp [B].
5. Prepare the venipuncture site and perform the venipuncture according to your institution's SOP.
6. Open blue clamp [B].

7. Position the sample/diversion bag with blue clamp [B] at the top and the sample tube holder/luer adapter at the bottom (see Figure 1). Allow the desired volume of whole blood to flow into the sample/diversion bag. The nominal volume of the sample bag is 60 mL. When the volume reaches the marking, approximately 40 mL has been collected.
8. Close blue clamp [B].
9. Break CLIKTIP [D] to allow the whole blood to flow into the whole blood bag. Bend the CLIKTIP in both directions to ensure that you break it completely.
10. Seal the sample/diversion tubing as near as possible to Y connector [A]. Seal the tubing according to your institution's SOP.
11. Invert the sample/diversion bag and transfer donor blood samples from the sample/diversion bag, using evacuated blood collection tubes. Transfer samples as soon as possible after venipuncture to avoid possible clot formation in the sample/diversion bag.
 - a. Hold the sample/diversion bag and sample tube holder/luer adapter in one hand with blue clamp [B] positioned at the bottom (see Figure 2).
 - b. Using the other hand, insert a blood collection tube firmly into the sample tube holder/luer adapter. After filling the blood collection tube, remove the tube from the sample tube holder. Repeat the process to take additional samples.
12. Mix the whole blood and the anticoagulant during collection according to your institution's SOP.
13. Collect the target volume of blood \pm 10%, as specified on the whole blood bag label.
14. Close white clamp [E] and perform one of the following steps:
 - If you made a knot in the collection tubing in step 1, tighten the knot firmly.
 - If you did not make a knot in step 1, seal the collection tubing near CLIKTIP [D].
15. Slide the NIP as close to the needle hub as possible.
16. With one hand, grip the needle hub and the front of the NIP. Pull gently to withdraw the needle from the donor's arm (see Figure 3).
17. Immediately after you withdraw the needle from the donor's arm, gently pull on the tubing behind the NIP to cover the needle with the NIP. Ensure that the needle locks securely in place (see Figure 4).
18. Seal the collection tubing near CLIKTIP [D] if you did not do so in step 14. Disconnect the sealed tubing that includes the needle and the sampling assembly. Dispose of the needle, the NIP, and the sample tube holder safely according to your institution's SOP and/or local regulations.
19. Immediately after collection is complete, invert the whole blood bag several times to thoroughly mix the whole blood and the anticoagulant.
20. If required, strip the blood in the collection tubing into the whole blood bag according to your institution's SOP.
21. Seal and disconnect the collection tubing between 0.5 cm and 3.0 cm away from the whole blood bag. Dispose of the collection tubing safely according to your institution's SOP and/or local regulations.
22. Place the whole blood unit into a temperature-controlled environment according to your institution's SOP.
23. Pack and transport the whole blood unit to the processing laboratory according to your institution's SOP.

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
Figure 1: Sample/diversion bag


Figure 2: Collecting a whole blood sample


Figure 3: Withdrawing the needle


Figure 4: Covering the needle with the NIP

Post-Processing and Red Blood Cell Leukoreduction

1. After processing, remove the plasma bag, the residual leukocyte bag, and the IPU bag from the organizer and handle them according to your institution's SOP.
 -  **Note:** You must rest and agitate the IPU prior to platelet pooling. For conditions for resting and agitating the IPU, see the platelet pooling set instructions for use.
2. Prepare a leukoreduction rack with a head height of at least 110 cm.
3. Remove the RBC bag from the organizer.
4. Hang the RBC bag from the leukoreduction rack so that the filter hangs vertically.
5. Open blue clamp [G].
6. Break CLIKTIP [H].
7. Break CLIKTIP [J].
8. Drain all of the additive solution through the filter and into the whole blood bag.
9. Close blue clamp [G].
10. Gently mix the additive solution with the RBC unit in the whole blood bag.
11. Hang the whole blood bag on the leukoreduction rack.
12. Ensure that the tubing is free of kinks or other obstructions.
13. Open blue clamp [G] and allow the mixture of red blood cells and additive solution to flow through the leukoreduction filter.

Do not squeeze the whole blood bag to increase the draining rate. Filtration is complete when the inlet side of the leukoreduction filter collapses and the red blood cells drain from the inlet side of the filter.
14. Seal the RBC tubing below the leukoreduction filter, and disconnect the RBC bag. Ensure that the filter remains in a vertical position while you seal the tubing, as this prevents the loss of red blood cells.

15. Dispose of the whole blood bag and filter assembly safely according to your institution's SOP and/or local regulations.
16. Seal the line for sampling segments according to your institution's SOP.
17. Store the RBC product according to your institution's SOP.

Symbols and Certification

Symbol	Definition	Symbol	Definition
	Indicates the product quantity.		Indicates the date after which the medical device is not to be used.
	Indicates the manufacturer's batch code so that the batch or lot can be identified.		Indicates the manufacturer's catalog number so that the medical device can be identified.
	Indicates a non-pyrogenic fluid pathway.		Indicates the presence of a sterile fluid path. The method of sterilization is steam.
	Indicates a medical device that is intended for one use, or for use on a single patient during a single procedure.		Indicates a medical device that is not to be re-sterilized.
	Indicates the medical device manufacturer.		Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile).
	Indicates that the product packaging complies with European Directive 94/62/EC for packaging and packaging waste.		Indicates that the product contains phthalates, specifically Di(2-ethylhexyl)phthalate (DEHP).
	Indicates a medical device that needs protection from light sources.		Indicates a medical device that needs to be protected from moisture.
	Indicates the range of temperatures to which the device can be safely exposed.		Indicates that the contents of the transport package are fragile; therefore, it shall be handled with care.
	Indicates the correct upright position of the transport package.		Indicates that the transport package is not allowed to be vertically stacked beyond the specified number.
	Indicates that the user should consult the instructions for use.		Indicates that the product should not be used if the product, its sterile barrier system, or its packaging is damaged or shows any signs of deterioration.
	Indicates that the product should not be used if the package is damaged.		Indicates that the product should not be vented. Do not remove the end caps from the set.

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Symbol	Definition	Symbol	Definition
	Indicates a donor needle.		Indicates the length of time within which the product must be used after the packaging is opened.
	Indicates the nominal volume of the blood container.		Indicates a whole blood bag.
	Indicates the target whole blood collection volume.		Indicates a red blood cell bag.
	Indicates a plasma bag.		Indicates a processing container.
	Indicates anticoagulant solution.		Indicates additive solution.
	Indicates a sampling site.		Indicates a leukoreduction filter.
	Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended.		Indicates the Authorized representative in the European Community.

Return of Used Product

If for any reason this product must be returned to Terumo BCT Ltd., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping. Instructions for cleaning and materials, including appropriate shipping containers, proper labeling, and an RGA number, may be obtained from the Terumo BCT Quality Assurance Department. **IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT. Please contact your local representative for information regarding returned goods and product complaints.**

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©2017 Terumo BCT, Inc.