2019-08

Progreat[®]

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CE 2797

 **TERUMO**

REF

Số catalog

LOT

Số lô

Hạn sử dụng

STERILE EO

Tiệt trùng bằng khí Ethylene oxide

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Nội dung

Nhà sản xuất

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
I.D.

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Điểm đánh dấu

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TERUMO CORPORATION 44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN
MADE IN JAPAN

TERUMO EUROPE N.V. INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM www.terumo-europe.com

AU Sponsor Terumo Australia Pty Ltd Macquarie Park NSW 2113 Australia
T: +61 2 9878 5122 E: cs_australia@terumo.co.jp

- ®: Registered Trademark
- ®: Marque Enregistrée
- ®: Eingetragenes Warenzeichen
- ®: Marca Registrada
- ®: Marchio Registrato
- ®: Geregistreerd Handelsmerk
- ®: Registrerat Varumärke
- ®: Tescilli Ticari Marka

 is a trademark of TERUMO CORPORATION.

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
DESCRIPTION OF COMPONENTS

The catheter assembled with the guide wire is for angiographic or intravascular therapy. The catheter has a hydrophilic polymer coating on the surface over its entire length except its proximal end. The coating gives it lubricity when it is wet. Furthermore, the guide wire consists of super elastic alloy core, polyurethane jacket, hydrophilic coating on its surface, and a gold coil distal tip marker can help to advance the catheter to target vessels. There are two types of guide wire available, pre-shaped type and shapeable type. Shapeable type guide wire tip can be reshaped. When infusing a contrast media through the catheter, a power injector can be used.

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
ENGLISH

ENGLISH

Please read all instructions prior to use. To avoid complications, observe all warnings and precautions throughout these instructions.

INDICATION

This MICRO CATHETER SYSTEM is intended for the infusion of contrast media into all blood vessels. The catheter is also intended for drug infusion in intra-arterial therapy and for the infusion of embolic materials for haemostasis. Do not use this product in coronary and cerebral blood vessels.

CAUTIONS FOR USE

1. CONTRAINDICATIONS

- Generally, angiography or intravascular therapy is contraindicated for, but not limited to, the patients listed below.
- Patients in the acute phase of myocardial infarction
- Patients with serious arrhythmia
- Patients with serious serum electrolyte imbalance
- Patients who in prior procedures have developed an adverse reaction to the injection of contrast media
- Patients with renal dysfunction
- Patients with coagulopathy or those whose blood has suffered a serious change in coagulation capability for some reasons
- Patients who cannot lie on their back on the operating table because of congestive heart failure or some respiratory disorder
- Patients with mental disease or those who are not expected to lie quietly during angiography
- Patients who are pregnant or are supposed to be pregnant

2. COMPLICATIONS

Angiography or intravascular therapy may be accompanied by, but not limited to, the following:

- Headache
- Fever and chill
- Blood pressure drifting
- Myocardial infarction
- Infection and pain at the puncture site
- Haemorrhage, haematoma, arterio-venous fistula and false aneurysm at the puncture site
- Spasm, artery perforation, dissecting aneurysm and false aneurysm with the use of a guide wire or catheter
- Inflammation with embolic material
- Bradycardia
- Behavior disorder
- Nausea and vomiting
- Abnormality in blood sampling tests
- Shock
- Renal failure
- Cerebral edema
- Cerebral infarction from peripheral artery occlusion
- Death

3. WARNINGS

- **Flush the lumen of the guiding catheter and the micro catheter system continuously with heparinized saline solution. Residual contrast media or blood clots on the micro catheter system surface reduce its lubricity, preventing smooth catheter movement. If flushing fails to restore surface lubricity, discontinue the use of the micro catheter system and remove it slowly and carefully together with the guiding catheter. Excessive force used in pulling the catheter may cause breakage/rupture/separation, which may necessitate retrieval of fragments.**
- **Do not pressurize the catheter or advance the guide wire through the catheter when the catheter is kinked or blocked. This may result in breakage of the catheter and damage to the vessels.**
- **Monitor the manipulation of the micro catheter system in the vessel by confirming the position of the catheter tip / guide wire under high resolution fluoroscopy and a digital subtraction angiography monitor (DSA). If any resistance is felt in the vessel, do not advance or withdraw the micro catheter system until the cause of resistance is determined under high resolution fluoroscopy and DSA. Manipulating the catheter and/or the guide wire against resistance may result in damaging the vessel, the catheter or the guide wire. If the situation is not solved, withdraw the entire system of the catheter or the guide wire with the guiding catheter.**
- **Do not advance the micro catheter system by forcing in extremely tortuous vessels. This may result in kink of the catheter or damage to the vessel.**
- **Do not pass the catheter through a stent strut.**


4. PRECAUTIONS

- The micro catheter system should be used by a physician who is familiar to the intended procedures.
- Manipulation of the micro catheter system should be monitored under high resolution fluoroscopy and DSA.
- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Do not use if the unit packaging or the product have been damaged or soiled.
- The product should be used immediately after opening the package and be disposed of safely and properly after use.
- The entire procedure should be carried out aseptically.
- The surface of the micro catheter system must be completely wet with heparinized saline solution to maintain a lubricious surface.

- Before use, consult the instructions for use of the drugs and devices to be used along with this catheter to determine compatibility and to prevent catheter damage.

DIRECTIONS FOR USE

1. Carefully remove the micro catheter system in its holder from the package.
2. Fill the holder with heparinized saline solution through the hub of the holder using a syringe, to thoroughly wet the surface of the catheter (Fig. 1).


3. Remove the micro catheter system slowly from its holder. If resistance is felt, do not try to remove it against the resistance, but inject heparinized saline solution into its holder again, and try once more.

CAUTIONS


- Do not use if the micro catheter system has been damaged or if any other abnormality is observed.
- When wet, the shaft of the micro catheter system is very lubricious. Hold the catheter by its hub during handling.

4. Make sure that the lock adapter is not loose. Inject heparinized saline solution into the guide wire hub using the attached 2.5 mL luer lock syringe (Fig. 2). In order to prime the catheter sufficiently, slowly inject at least 1 mL of the solution into the catheter until more than 10 drops of the solution appear out of its tip. To maintain surface lubricity, immerse the catheter and the guide wire assembly in a heparinized saline solution bath or put it into its holder filled with heparinized saline solution.

CAUTIONS


- Use immediately after opening the package. Once the catheter pouch is opened, the enclosed syringe package is no longer sterile.
- Prime the catheter and guide wire sufficiently. Manipulation of an insufficiently primed catheter may cause wrinkling, separation of the catheter, and/or abrasion of the hydrophilic coating on the guide wire.

5. To shape the tip of the re-shapable type guide wire, maintain its surface lubricity and coil it carefully around your fingertip or the attached inserter (Fig. 3).

CAUTIONS

- Do not shape the pre-shaped type guide wire. This may cause damage to the guide wire.
- Do not shape the guide wire using method other than that described above. Do not handle the wire when it is dry, do not heat, shape with forceps or fingernails. Do not bend tightly or back and forth nor bend the same part of the guide wire repeatedly. This may result in guide wire breakage/separation. Wiping without surface lubricity may result in abrasion of hydrophilic coating.


6. When shaping this catheter by steam, insert the enclosed shaping mandrel into the distal tip of the catheter and gently shape the tip to the desired angle. Then expose the tip to the steam for approximately 10 seconds (Fig. 4). Check the resulting shape after removing of the shaping mandrel.

CAUTIONS

- Do not rub or bend the catheter tip with too small radius, pinch by forceps or tweezers, which may result in the damage of the surface coating, collapse of the catheter shaft and/or deformation of catheter.
- Positioning the catheter tip closer than 2 cm from the steam source may result in the damage of the surface coating or the tip of the catheter.
- Excessively re-shaping the catheter may damage the surface coating or the tip of the catheter.
- When shaping with steam, take care not to burn yourself.
- Do not insert the enclosed shaping mandrel into the patient's body.
- Do not stretch the catheter tip tightly or bend excessively when shaping it not with the enclosed shaping mandrel but with your fingers. It may result in collapse of the catheter shaft and/or deformation of the catheter.

7. Insert the guiding catheter into the patient's vessel. Attach a rotating haemostatic valve (Tuhoy-Borst type) to the guiding catheter and continuously irrigate the catheter with heparinized saline solution. Insert the micro catheter system with guide wire assembly through the valve into the guiding catheter and advance to the distal end of the guiding catheter. In case of difficulty of the micro catheter system insertion, loosen the lock adapter, slowly pull the guide wire back approximately 20mm from the catheter tip in order to straighten the catheter tip, and then carefully insert the catheter, avoiding catheter kinking (Fig. 5-2). After the catheter tip has successfully negotiated the guiding catheter hub, slowly advance the guide wire and re-tighten the lock adapter.

WARNINGS

- Do not manipulate and/or withdraw the micro catheter system through a metal entry needle or a metal dilator. Manipulation and/or withdrawal through a metal entry needle or a metal dilator may result in abrasion of the hydrophilic coating, destruction and/or separation of the micro catheter system.
- If the guiding catheter is fitted with a stopcock, do not close the stopcock with the micro catheter system inside the guiding catheter. The micro catheter system may be broken.
- Make sure that the guiding catheter does not slip out of the vessel. If the guiding catheter should leave the vessel when the catheter and/or the guide wire is moved, this may result in the damage of the micro catheter system.

CAUTIONS


- Do not tighten the rotating haemostatic valve excessively on the micro catheter system, and/or manipulate the micro catheter system through a tightened valve. Damage to the micro catheter system may occur.
 - If resistance is felt, do not force the micro catheter system into the guiding catheter as this may result in the damage of the micro catheter system.
8. Carefully advance the micro catheter system through the guiding catheter until it reaches the desired site. At bifurcation, rotate the guide wire hub so that the angled tip of the guide wire points in the desired direction (Fig. 5-1). If complicated vessels selectivity is observed, slowly pull the guide wire back approximately 20 mm from the catheter's tip to make its tip straight (Fig. 5-2).

WARNINGS

- Monitor the manipulation of the micro catheter system in the vessel, by confirming the position of the catheter tip / guide wire under high-resolution fluoroscopy and DSA. If any resistance is felt in the vessel, do not advance or withdraw the catheter until the cause of resistance is determined. Manipulating the catheter and/or the guide wire against resistance may result in damaging the vessel, the catheter or the guide wire. If the situation is not solved, withdraw the entire system of the catheter or the guide wire with the guiding catheter.
 - When advancing the micro catheter system into the peripheral vessel, draw it back slightly under fluoroscopy each time it has been advanced, to make sure that the micro catheter system has not been advanced so far that it can not be drawn back.
 - If the micro catheter system is manipulated in the vessel without the guide wire, it may result in damaging the vessel. When re-inserting the guide wire into the catheter, carefully advance the guide wire while making sure the position of the guide wire under high resolution fluoroscopy and DSA. The quick and unreasonable movement may cause damage to the vessel.
9. When advancing the micro catheter system to the desired site, the wire protruding from the catheter tip can be adjusted and affixed with the attached wire stopper. Set the slit part of the wire stopper to the guide wire between the catheter hub and the guide wire hub. Remove the guide wire from the catheter hub and adjust the protruding tip length. Slide the wire stopper into the catheter hub to fix the guide wire (Fig. 6).

CAUTION

Do not move the guide wire while the wire stopper is fixed on the catheter hub. If the guide wire is moved in such a situation, the hydrophilic coating and the coating resin on the guide wire may be abraded.

10. When the desired site is reached, remove the guide wire from the catheter.

CAUTIONS

- If any resistance is felt while removing the guide wire, do not remove the guide wire by force. The resistance may indicate insufficient lubricity of the guide wire. Flush the catheter again with heparinized saline solution. If any resistance is felt while removing the guide wire, do not remove the guide wire by force. Drawing back the guide wire against resistance may cause the catheter kink. Carefully remove the guide wire together with the catheter.
 - Rinse residual blood from the removed guide wire in a heparinized saline solution bath. If the residual stains do not come off, wipe the guide wire once with gauze moistened with heparinized saline solution. Blood remaining on the guide wire can affect its surface lubricity and could cause resistance when inserted into the catheter.
 - To avoid damage to the surface coating, do not wipe the guide wire with gauze with disinfectant.
11. Before introducing an embolic material or other agent, slowly inject a small volume of contrast media into the catheter using a syringe and verify under high resolution fluoroscopy and DSA that the media come out of the catheter tip. When using embolic materials and drugs, use them following their instruction for use to check their compatibility with the catheter. When using many embolic materials, it is recommended to change the catheter each time.

WARNINGS

- If any increase of resistance is felt during infusion, replace the catheter with a new one. Injection against increased resistance may cause the catheter to break, resulting in damage to the vessel.
- If no contrast media come out, it indicates possible kinking of the catheter. If drawing back the catheter fails to correct the kink, replace the catheter with a new one. Do not try to correct the kink by inserting a guide wire or by pressurized infusion. Starting the introduction of embolic material or the agent without correcting the kink or attempts to correct the kink by inserting guide wire or by infusion may cause the catheter to break/rupture/separate and this may result in damage to the vessel.

CAUTIONS

- **Friction between the catheter wall and the embolic material may work to advance the catheter, resulting in perforation of the vessel wall. To prevent this, take up the slack of the catheter by drawing it back slightly and hold.**
- Increased resistance to infusion suggests that the catheter is blocked with the drug or contrast media being infused or with blood clots. Discontinue infusion immediately and replace the catheter with a new one.
- Do not use the attached 2.5 mL syringe to inject contrast medium, embolic materials or other medicines. Using the 2.5 mL syringe for high pressure injection of such agents may cause damage of the syringe or dislocation of the syringe from the catheter hub.
- When introducing an embolic material, do not use the attached guide wire as a supportive device. It may run onto the embolic material inside the lumen.
- When a power injector is to be used, follow the instructions given below under "Instruction For Using a Power Injector with the catheter".
- Using organic solvents may cause damage to the catheter.
- When introducing an embolic material, use material and loading devices not exceeding 0.018" (0.46 mm) in diameter.
- In case of use of embolic material and supportive device suitable for catheter with 0.016" (0.41 mm) or smaller I.D., the supportive device may run onto the embolic material inside the lumen because of clearance between their O.D. and the catheter I.D.

12. When inserting the catheter further more to the other vessel, sufficiently inject the heparinized saline solution into the catheter. If any resistance is felt during the insertion of the guide wire, discontinue to advance the guide wire and replace with a new one.

In case the guide wire is difficult to be inserted into the catheter hub, insert the guide wire tip by turning the guide wire or the catheter hub right and left.

WARNING

When re-inserting the guide wire into the catheter, carefully advance the guide wire while making sure of the position of the guide wire tip under high resolution fluoroscopy and DSA. The quick and unreasonable movement may cause the catheter to break/rupture/separate and may result in damage to the vessel.

13. When the procedure is completed, remove the micro catheter system together carefully.

WARNING

If any resistance is felt, do not remove the micro catheter system by force. Withdraw the catheter carefully together with the guiding catheter. Removing the catheter by force may result in the catheter breakage/separation, which may necessitate retrieval.

Instruction For Using a Power Injector with the catheter

A power injector can be used to infuse a contrast media through the catheter. Observe the warnings and cautions given below. The flow rate depends upon such factors as the viscosity of the contrast media, which varies with the type and temperature of the media, the model and setting of the power injector and how the injector is connected to the catheter. The observed flow rate values indicated below are for reference only.

WARNINGS

- **Do not use a power injector to infuse agents other than contrast media, as the catheter may become blocked.**
- **Setting of injection pressure must not exceed below listed maximum injection pressure that correspond to outer diameter of each catheter tip. Exceeding of injection pressure beyond the maximum injection pressure may cause catheter rupture.**

Catheter O.D.	Maximum injection pressure
2.4Fr., 2.7Fr.	5171 kPa (750psi)
2.8Fr.	6205 kPa (900psi)

- **Under high resolution fluoroscopy and DSA monitor, inject small amount of contrast media with syringe and confirm the flow of contrast media out of the catheter tip before using injector.**
- **If expansion of the catheter O.D. is observed during the injection, it may be an excess over the maximum pressure limit. In such case, stop injection immediately.**
- **When securing the catheter in position, secure it by the hub so that catheter shaft is not damaged. In securing, do not hold the catheter shaft with forceps, or this may result in catheter separation.**

CAUTIONS

- If the catheter has been kinked or bent sharply, replace it with a new one.
- Connect the power injector to the catheter using a pressure-resistant extension tube.
- When re-inserting the guide wire after completion of angiography, flush out the catheter lumen with heparinized saline solution.

REFERENCE DATA

- 1. Injector use MARK V, MARK 7 Arterion (MEDRAD)
- 2. Conditions and injector setting Contrast media temperature 37 °C
Injection pressure monitor/limit
.....4137 kPa (600 psi), 5171 kPa (750 psi), 6205 kPa (900 psi)
Flow scale mL/sec
Linear rise seconds..... 0.3 sec

3. Methods and results
 The value of injection volume was selected to be three times the set flow rate value. The set injection volume is not attained because of varying conditions, including the viscosity of the contrast media.

Catheter O.D.	Usable length of Catheter (cm)	Contrast media	Iodine content (mg/mL)	Viscosity (cP)	Set Condition		Actual Flow Rate (mL/sec)			Dead space volume (mL)
					Flow Rate (mL/sec)	Volume (mL)	4137 kPa (600 psi)	5171 kPa (750 psi)	6205 kPa (900 psi)	
2.4/2.9Fr. (0.80/ 0.97mm)	100	Iopamidol	300	4.4	3.0	10	2.2	2.5	–	0.36
			370	9.1	3.0	10	1.2	1.4	–	
	110	Iopamidol	300	4.4	3.0	10	2.1	2.4	–	0.38
			370	9.1	3.0	10	1.1	1.3	–	
	130	Iopamidol	300	4.4	3.0	10	1.8	2.2	–	0.43
			370	9.1	3.0	10	0.9	1.2	–	
	150	Iopamidol	300	4.4	3.0	10	1.5	1.9	–	0.47
			370	9.1	3.0	10	0.8	1.1	–	
2.7/2.9Fr. (0.90/ 0.97mm)	100	Iopamidol	300	4.4	6.0	20	2.8	3.5	–	0.44
			370	9.1	3.0	10	1.9	2.2	–	
	110	Iopamidol	300	4.4	6.0	20	2.6	3.4	–	0.46
			370	9.1	3.0	10	1.7	2.0	–	
	130	Iopamidol	300	4.4	6.0	20	2.3	3.0	–	0.53
			370	9.1	3.0	10	1.5	1.8	–	
	150	Iopamidol	300	4.4	6.0	20	2.1	2.7	–	0.59
			370	9.1	3.0	10	1.4	1.6	–	
2.8/3.0Fr. (0.93/ 1.00mm)	100	Iopamidol	300	4.4	6.0	20	3.5	3.9	4.4	0.48
			370	9.1	6.0	20	2.1	2.5	3.0	
	110	Iopamidol	300	4.4	6.0	20	3.1	3.6	4.0	0.53
			370	9.1	6.0	20	1.9	2.4	2.7	
	130	Iopamidol	300	4.4	6.0	20	2.9	3.3	3.8	0.59
			370	9.1	6.0	20	1.6	1.9	2.3	
	150	Iopamidol	300	4.4	6.0	20	2.7	3.1	3.6	0.66
			370	9.1	6.0	20	1.5	1.7	2.1	

PRECAUTION FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperature or humidity during storage.

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© TERUMO CORPORATION 2019-08
TERUMO is a trademark of TERUMO CORPORATION.