TERUMO PENPOL®

D601

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MỤC LỤC

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TERUMOBCT
Unlocking the Potential of Blood

Terumo BCT, Inc.
10811 West Collins Ave.
Lakewood, Colorado 80215-4440
USA

USA Phone: 1.877.339.4228
Phone: +1.303.231.4357
Fax: +1.303.542.5215

Terumo BCT Europe N.V.
Europe, Middle East and Africa
Ikaroslaan 41
1930 Zaventem
Belgium

Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

Terumo BCT Latin America S.A.
La Pampa 1517 – 12th Floor
C1428DZE
Buenos Aires
Argentina

Phone: +54.11.5530.5200
Fax: +54.11.5530.5201

TERUMOBCT.COM

Terumo BCT (Asia Pacific) Ltd.
Room 3903-3903A, 39/F
ACE Tower, Windsor House
311 Gloucester Road
Causeway Bay, Hong Kong

Phone +852.2283.0700
Fax: +852.2576.1311

Terumo BCT Japan, Inc.
Takanawa Park Tower 13F
20-14, 3-chome
Higashi Gotanda, Shinagawa-ku,
Tokyo 141-0022 Japan

Phone +81.3.6743.7890
Fax: +81.3.6743.9800

TERUMO PENPOL®

D601

BLOOD COLLECTION MONITOR

INSTRUCTIONS FOR USE

APRIL 2020

CONTENTS

1. INTRODUCTION	3
1.1 Mixing of Anticoagulant.....	3
1.2 Blood Collection Monitor D601.....	3
1.3 Features of D601.....	3
2. SYMBOLS AND DEFINITIONS	4
3. IMPORTANT SAFETY INSTRUCTIONS	5
4. SPECIFICATIONS	6
5. TECHNICAL DESCRIPTION	7
5.1 Internal Block Diagram	7
5.2 Functional Description	8
6. INSTALLATION DETAILS	13
6.1 Packing Details.....	13
6.2 Site Planning Guidelines.....	13
6.3 Installation Procedure.....	13
7. OPERATING PROCEDURE.....	14
7.1. Setting of Volume.....	14
7.2 Collection of Blood.....	14
7.3 Pause Collection.....	14
7.4 Abort Collection.....	14
7.5 Time Setting.....	15
7.6 Previous Collection History.....	15
8. BATTERY BACKUP.....	15
8.1 Battery Charging Details.....	15
8.2 Battery Replacement.....	15
8.3 Equipment Shipment.....	16
9. OPTIONAL CARRY CASE/CARRY BAG.....	16
9.1 Carry case.....	16
9.2 Placement in carry bag.....	18
10. PRECAUTIONS.....	18
11. CLEANING INSTRUCTIONS.....	18
12. TROUBLE SHOOTING	19
13. EMC INFORMATION.....	20
14. DISPOSAL (END OF LIFE).....	22
15. CONTACT ADDRESS.....	23
16. WARRANTY.....	25

1. INTRODUCTION

Transfusion of Blood and Blood components is an established standard way of treating patients who are deficient in one or more blood components. If proper process of collection, separation and storage are not employed an appreciable fraction of red cells and a major fraction of platelets and factor VIII are lost or have no useful function by the time it is transfused to the patient. Much of these losses are due to the “lesion of collection”.

1.1 Mixing of Anticoagulant:

To prevent the occurrence of clotting in the transfusion apparatus anticoagulant solution is mixed with whole blood when it is taken from the donor. The ratio of volume of blood to anticoagulant is a factor, which determines the component yield. Occasionally blood that is collected is unsuitable for components preparation because of over or under collection with reference to the available volume of anticoagulant solution in the Blood bag. The Blood while it is drawn from the patient should be mixed with the anticoagulant gently and periodically.

1.2 Blood Collection Monitor D601


The Blood Collection Monitor D601 monitors the process of blood collection. The yield of blood components from whole blood is determined by factors like the rate of collection, proper mixing of blood with anticoagulant solution, accuracy of collected blood volume etc. The Blood Collection Monitor helps in the process of blood collection by monitoring of critical parameters. The basic functions of a Blood Collection Monitor include weighing of blood collected, mixing it with the anticoagulant and stopping the blood inflow at the pre-programmed volume. The equipment can be programmed to collect any volume between 50ml and 500ml in steps of 5ml. Also provisions to pause the started process and to abort the current process are available. During blood collection, the programmed volume and collected volume are displayed in the display for user verification. The features also include flow rate calculation and indications for low and high rates of flow. Flow rate alarm indication is set in such a way that it activates when the flow goes below 20ml/min or above 180ml/min. This alarm can be acknowledged. At the end of the collection, the total time taken for the collection is indicated. The Blood Collection Monitor D601 is having a keypad by which the user can operate different functions of the equipment. Optional Carry bag /case (to be ordered separately) helps to carry the device for camps and other collection points.

1.3 Features of D601

- Comply with EN 60601-1.
- Donor safety by way of programmed blood volume collection.
- Power backup by Rechargeable Lithium ion battery.
- 12 hours of battery back-up time.

- In built battery charger to enhance the battery life.
- Battery low voltage indication and automatic cut-off.
- Continuous sensing of blood, so no variation in the collection of blood.
- Flow rate display, Low and High flow rate indication and alarms. Alarms are acknowledgeable.
- Automatic clamping when low flow rate sustain for more than two minute.
- Previous data in memory, so user can retrieve the data of previous bleeding.
- Motorized clamping, so more reliable and less power consuming
- Power failure alarm.
- Removable tray for easy cleaning.
- Real time clock.
- Time taken for collection to determine suitability for component separation.
- Displays collection time in every 3 seconds during collection.

2. SYMBOLS AND DEFINITIONS

	European Conformities		This way up
	Serial Number		Fragile; Handle with care
	Manufacturer		Keep dry
	Indicates that the user must read the instructions for use before application.		Temperature Limitation
	EC Representative		Stacking Limit by number
	General Warning Sign		Product complies with Waste electrical and electronic equipment directive 2002/96/EC (WEEE).
	Caution		Alternating Current
	Class II equipment		Type B equipment
	Direct Current		

3. IMPORTANT SAFETY INSTRUCTIONS

WARNING: Identifies conditions that could result in personal injury or loss of life

CAUTION: Identifies conditions that could result in damage to the equipment

If the equipment is not used following the intended use described in this manual, the manufacturer will not guarantee the results nor the safety of the device and shall not be liable.

Carefully read the entire content of this manual before using the D601.

WARNING:

- The equipment uses motorized clamp to pause or stop the donation process. So the tube should always be routed through the clamp. The equipment should not be used for collection without routing in any circumstance, which can lead to a safety hazard to the donor.
- After the battery low indication the equipment can be shut down without any further warning. Kindly charge the battery before proceeding for next collections.
- Whenever a flow rate warning is heard please examine the donor's condition and check the needle before proceeding
- **Remove the adaptor from the power socket while not charging and equipment not in use. There is no ON/OFF switch for charging. The equipment switch functions as a standby switch.**

CAUTION:

- The equipment is not suitable for use in the presence of a flammable anesthetics mixture with air or oxygen or nitrous oxide.
- The Blood Collection Monitor contains a sensitive weight sensor. Do not place weight above 1000gms on the tray.
- Avoid exposure to water, extreme temperature and humidity.
- Do not use aggressive or corrosive products, abrasive powder, steel wool, abrasive sponges or chemical solvents for cleaning.
- No user serviceable part inside.
- While carrying the device from the installed location to other places for collection, ensure to place them in dedicated carry case/carry bag.
- Use only accessories supplied with the equipment.
- Turn the power switch OFF when not in use.
- Do not place any objects that blocks the movement of the tray.
- If battery is drained please charge the equipment for 3.5 hours minimum for fully charging the battery. This will be indicated by the green LED on the equipment.

- Do not plug in any other adaptor to the power socket, other than the one supplied along with the equipment
- Use only with the rated input voltage
- Special care to be given while transporting the blood collection monitor from one place to another
- Do not attempt collection in "Weight Check mode". This is only meant for checking the weight in grams.

4. SPECIFICATIONS


Product code	D601
Power supply	18V DC 3.3A for D601 100V -240V AC, 50-60Hz for power adaptor
Compliance	IEC 60601-1
Electrical Safety	Protection against electrical shock-Class II
Operation	Continuous operation
Power	70 VA (Max)
Battery.....	Type: Li-ion Nominal Voltage/Capacity: 14.8V, 5.2Ah Nominal power: 76.96Wh
Working environment	Temperature: 10 ^o to 40 ^o C Relative humidity: 30 to 75% Atm Pressure:700 to1060hPa
Storing and transportation	Temperature: -20 ^o to 40 ^o C Relative Humidity: 20 to 90% Atm Pressure:600 to1060hPa
Tray movement details	12± 2 oscillations per minute
Clamping mode	Motor activated clamping
Volume setting ranges	50ml to 500ml insteps of 5ml 350ml/450ml quick selection Automatic recall of set volume
Display	20x4 alphanumeric display
Time measurement	Collection time display
Accuracy.....	± 2% of the programmed ml
Dimension in mm* (WDXH)	190mm x 500mm x 185mm
Equipment Weight(Approx.)	D601: 4.5 kg, Power Adaptor:550gm
LED indication and audible alarm	1.Low/High flow indication 2.End of collection 3.Power failure 4.Battery Low
Automatic Clamping	1.When programmed volume reaches 2.When flow rate is less than 20 ml/min for more than two minute

*Tolerance applicable

- Power Cord According to Country specification
- Accessories AC Powercord, IFU English, Warranty Card

5. TECHNICAL DESCRIPTION

5.1 INTERNAL BLOCK DIAGRAM

CPU Card

The CPU card is based on the Micro-controller and controls the operation of the equipment.

LCD Display

LCD display shows the status of collection, time, programmed volume, collected volume.

Power Card

This control card provides the supply to the CPU card and other circuits. It converts the input voltage to +12V, +5V, -5V levels which is required for powering the CPU card, load cell and peripheral circuits

It also consists of a Lithium Ion battery charger which regulates the current to charge the lithium ion battery

Li-Ion Battery

The equipment is provided with a power backup with rechargeable lithium ion battery

Clamp Motor and Magnetic Reed Switch

Clamp motor clamps the blood tube to stop the flow of blood at the event of collection completion and collection abortion.

Membrane Keypad

This board consists of tactile keys for user interface. All the keys are the inputs of Micro controller, which monitor each key whether it is pressed or not. According to the operation of keys D601 perform its specified task.

LED

LEDs in the membrane keypad indicates the status of D601. Red LED Indicates warning for emergency attention and Yellow LED indicates need for attention.

Rocking Motor

Rocking motor rotates the tray assembly when the blood is collected and ensures proper mixing of blood with anticoagulant during collection.

Load Cell

Load cell senses the weight on the tray.

Buzzer


The alarm conditions are alerted using a buzzer. The service person can adjust the buzzer volume to the required level between factory set minimum and maximum levels.

Power Adaptor(External)

The equipment is powered using an external power adaptor 18V, 3.3A.

5.2 FUNCTIONAL DESCRIPTION

PARTS OF D601

1-Tray

The blood bag is placed on the tray before the collection starts. On pressing the “start” the bag weight is tarred and the tray starts to oscillate. As blood flows into the bag the weight is monitored by the equipment

2-Clamp

The tube is routed through the clamp before the collection procedure. In any case of collection need to be stopped or in case of collection completed, the clamp operates down to squeeze the tube, thereby stopping the donation.

3-Display

Display is a 20X 4 alphanumeric LCD where the user can see the current status of donation. An option to view the last collection history is provided by pressing the “pause” key


4- Keypad and Indications

The membrane keypad provides the controls and LEDs for operating the equipment. Detailed operating procedure is provided below.

5-ON/OFF Switch

The illuminated ON/OFF switch is used for turning the equipment to operating/standby mode. **The charging function will not be terminated by turning OFF this switch.**

FRONT PANEL

FRONT PANEL CONTROLS

START: Used to start collection of blood and gentle mixing. Also used to acknowledge the Low/ High flow alarm and indication.

PAUSE: This key is used to pause the collection for some time. During this time the clamp is activated and gives audible alarm. Press the key again to release the clamp and continue collection

CLAMP DOWN: Used to close the clamp. Collection will be aborted when this key is pressed in the collection mode. Also used for decreasing the minute settings when in Time Setting Mode.

CLAMP UP: Used to up the clamp manually. Also used for increasing the minute settings when in Time Setting Mode.

VOLUME DOWN: This key is used to decrease the set volume. Pressing this key decreases the set volume by 5ml. The lowest programmable volume is 50ml.

VOLUME UP: This key is used to set the programmed volume. Pressing the key increases the set volume by 5ml. The highest programmable volume is 500ml.

350/450 SELECT: This key selects between 350ml or 450ml volume programming.

On long press of this key the equipment goes into Weight check modewhere the weight can be measured in grams. Again long pressing it will change to the donation mode.

FRONT PANEL INDICATIONS


CHARGING INDICATION: This LED turns ON when the power adaptor is connected. The LED glows in Amber colour when the battery is charging and turns Green colour when the battery is fully charged.

POWER FAILURE INDICATION: When power fails this LED turns ON, and when battery is in low condition it blinks.

FLOW RATE WARNING INDICATION: This LED indicates that the flow rate is below 20ml/min or above 180ml/min. A buzzer and a 'LOW' / 'HIGH' display is accompanied with the collected volume display to indicate this.

DISPLAY MESSAGES AND INDICATIONS

HOME SCREEN

When the equipment is powered ON the default screen appears and shows the set volume which is previous set and the collected volume as zero. The real time is also displayed. The third line is a rolling display showing the instructions to follow before collection

- 1.SET VOLUME
- 2.PLACE BAG
- 3.ROUTE TUBE
- 4.DO PHLEBOTOMY
- 5.PRESS "START"

COLLECTION SCREEN

```

SET      C.VOL  FLOW
500ml   000ml  RATE
09:39   00:06
          D601
  
```

On pressing the “start” key the donation starts and the collection screen appears on the display. The second line shows the set volume and the collected volume and the flow rate. The flow rate is updated every 30 seconds and whenever a high rate (above 180ml/min) or low rate (below 20ml/min) occurs the display shows “HIGH” or “LOW” respectively along with the buzzer beeps. This is an indication for the doctor/ technician to check the flow of blood to decide whether any activity to be initiated.

The third line shows the real time and the collection duration

COLLECTION COMPLETE SCREEN

```

SET      C.VOL  FLOW
          END OF BLEEDING
09:44   00:42
          D601
  
```

When the donation is completed the screen displays “END OF BLEEDING” and the buzzer beeps

```

350 =350
LIFT THE BAG
          D601
  
```

```

                                00m42s
LIFT THE BAG
          D601
  
```

The display toggles between the collection volume and the collection duration and instructs the user to lift the bag.

On lifting the bag the clamp opens and the screen returns back to the default home screen.


PAUSE SCREEN

```

SET      C.VOL  FLOW
500ml   000ml  RATE
          BLEEDING PAUSED
          D601
  
```

In any case when the bleeding is paused by pressing the “PAUSE” key, the display shows bleeding paused and can be restarted again by pressing the “PAUSE” key again

ABORT SCREEN


```

SET      C.VOL  FLOW
500ml   000ml  LOW
        BLEEDING ABORTED
        D601
  
```

By pressing the “CLAMP DOWN” in between a collection process, the bleeding is aborted and the by pressing “CLAMP UP” again the equipment returns backs to default state

PREVIOUS DATA SCREEN


```

PREVIOUS DATA
START TIME: 10:44:10
END TIME: 10:46:41
VOLUME: 349
  
```

From the equipment idle condition if the “PAUSE” key is pressed the display shows the history of the last collection which includes the start time, the end time and the collected volume.

TIME SET MODE SCREEN

```

TIME SET MODE
500ml 000ml 09:48
PRESS "CLAMP UP"
        D601
  
```

From the equipment idle condition if the “CLAMP DOWN” key is pressed the display shows to press either press the “PAUSE” key to enter the “TIME SET MODE” or just press “CLAMP UP” to go back to default state.

In time set mode the “CLAMP UP” key is used to increment minutes and “CLAMP DOWN” key is used to decrement minutes. The time can be saved by using the “START” key. Exit without saving by pressing the “PAUSE” key followed by the “CLAMP UP” key to return to the default screen and function.

WEIGHT CHECK SCREEN

```

WEIGHT CHECK MODE
        000gm
  
```

By long press in the “350/450” select key the equipment enters into the weight check mode. In the weight check mode the display shows the weight in grams of the item in the tray. Again a long press on the “350/450” select key the equipment goes back to the default state

CAUTION:

The weight displayed will blink on this screen as alert. The weight check mode is only for calibration of weights and to be performed by a service personnel only.

6 INSTALLATION DETAILS

- Inspect the carton for shipping damages. If you find any, notify Terumo Penpol Private Limited / Branch Office.
- The items received may be cross verified against the P.O and packing slip. Any missing or damaged item must be reported immediately to the carrier and Terumo Penpol Private Limited / Branch Office

6.1. PACKING DETAILS

The package contains the following items

- Blood Collection Monitor D601
- Tray
- Power cord
- Adaptor
- User manual
- Test Certificate
- Do's & Don'ts
- Carry bag/Carry case (Optional – to be ordered separately)

6.2. SITE PLANNING GUIDELINES

When selecting a place to set up the Blood Collection Monitor, please keep the following in mind.

- Place the D601 on a flat stable surface.
- Leave adequate space around the D601 for easy operation and proper ventilation.
- Avoid locations that are subject to direct sunlight, excessive heat, moisture and dust.
- Avoid electrical outlets on the same line with large motors or other heavy power consumption equipment that may cause fluctuations in line voltage. This may cause variation in the volume of blood being sensed since D601 contains a sensitive weight sensor.

6.3. INSTALLATION PROCEDURE

- Ensure proper earth of power point.
- Ensure that the mains voltage is within the range as in the specification.
- Ensure that the adaptor is supplied along with the D601 package
- Ensure that the battery is reconnected by authorized service personnel.
- Remove the tray packing carefully and fix the tray on the magnetic cup.
- Do not apply pressure on the tray while cleaning or fixing it on the magnetic cup to avoid damage of the sensitive load cell.
- The tray should not touch any part of the body other than its mounting plate.
- Connect the power adaptor.
- Switch ON the equipment.
- Verify whether the equipment is working according to the user manual.

7 OPERATING PROCEDURE

- Connect the power adaptor on the back side of the equipment and turn ON.
- Switch ON the equipment power switch.
- When the equipment is turned ON the LCD display indicates the volume set for previous collection.

7.1. SETTING OF VOLUME

- Press the volume up or volume down key to change the volume of blood to be collected. This key need not be used if the volume of blood to be collected is the same as that of the previous collection.
- For Preset Volume: Press and hold 350/450 switch to toggle the selection.
- Continue action following the message displayed in the second line of the LCD, i.e. Place the Blood bag on the tray and route the tube through the clamp and perform phlebotomy.

7.2. COLLECTION OF BLOOD

- Press start key, only when the blood has started flowing into the Blood Bag. The machine automatically tares the bag weight and the weight of the anticoagulant solution so that the display shows only the volume of blood being collected.
- Press start key to acknowledge the flow rate alarm and indications if any occurred during donation. The equipment gives gentle agitation to the bag.
- Press pause key, if user want to change the programmed volume or adjust the needle or tubing.
- When the collected volume reaches “programmed volume - 10ml” the rocking motor stops its agitation and
- when collected volume reaches “programmed volume - 3 ml” the clamp will be activated.
- At the end of collection Collected volume and duration of the bleeding are displayed in the LCD in the toggling mode.
- To release the clamp at the end of collection, pick up the bag from the tray.

7.3. PAUSE COLLECTION

- Press “Pause” key.
- In this mode, the clamp activates and the alarm will sound. Now the user can change the programmed volume by using volume up and down keys. During this time the agitation continues. To resume collection press pause key again

7.4. ABORT COLLECTION

- To abort collection before the machine achieves the programmed volume, press clamp down key.
- Release the clamp by pressing clamp up key

7.5. TIME SETTING

- Time setting can be done only in HOME SCREEN.
- Press clamp down key
- Press “Pause” key. This will allow you to edit the minutes of the time
- Press clamp down key to decrease the minutes and press clamp up for increasing the minutes.
- Press “Start” key to enter new time settings.
- Press clamp up key to go to HOME SCREEN.

7.6. PREVIOUS COLLECTION HISTORY

- When the equipment is in HOME SCREEN press the “PAUSE” key.
- The display will show last collection details for few seconds and automatically go back to the home screen.

8 BATTERY BACKUP

- On a full battery recharge, the system is provided with a backup of 12 hours.
- Li-Io Rechargeable Battery is used for backup
- The battery full charged is indicated by the charging LED which turns GREEN
- Once the battery level comes below 13.3V the low battery indication is provided by blinking LED.
- After low battery indication the equipment can turn OFF without any further alerts

8.1. BATTERY CHARGING DETAILS

- The full charged battery is indicated by a GREEN colour of charging LED
- An empty battery takes 3.5 hours to reach full charge if equipment is non-operational
- If the system is not in use please turn OFF the equipment while charging
- There is no ON/OFF switch for battery charging. Remove the adaptor from the power socket to terminate power

8.2. BATTERY REPLACEMENT

- Batteries should be replaced by an authorized technician only
- Do not short circuit the battery terminals
- Old batteries should be disposed properly
- Replace with specified type of battery only

8.3. EQUIPMENT SHIPMENT / BATTERY SHIPMENT

The D601 equipment consists of a Li-Ion battery with a nominal voltage 14.8V and 5.2Ah capacity and is subject to the requirements of the IATA Dangerous Goods regulations (DGR). If contained in the D601 equipment, it is classified in Class 9 – Miscellaneous Dangerous Goods as UN3481 (Li-Ion batteries contained in equipment). If battery only needs to be shipped, it is classified in Class 9- Miscellaneous Dangerous Goods as UN3480 (Li-Ion batteries).

CAUTION:

- If for some reason, the D601 equipment needs to be shipped with battery, it is the sole responsibility of the shipper to comply with all requirements of the DGR. The most recent information on DGR and packing instructions of Li batteries can be found on:

<http://www.iata.org/whatwedo/cargo/dgr/Pages/lithium-batteries.aspx>

- In case of any further clarification required on this topic please contact the factory
- **Please also note that the shipping of damaged or defective batteries is strictly forbidden. Thus, if a battery is suspected to be defective, it should not be shipped under any circumstances.**

9. OPTIONAL CARRY CASE /CARRY BAG

9.1. CARRY CASE

As an optional accessory, a D601 carry case/carry bag can be ordered. D601 with accessories can be placed inside the case/bag while carrying it for mobile collections.

9.1.1. PLACEMENT IN CARRY CASE

There is space for placing the power cord and adapter as shown below

Place the equipment above the case and close with the lid so that the slot in the lid align with the position of clamp in equipment.

For easier transportation, the cases can be stacked on top of each other.

WARNING:

Do not stack more than 2 cases on top of each other. Since it is a rigid case the D601 can be installed on top of the Carry case if required – details mentioned in section 9.1.2.

9.1.2. INSTALLATION ON TOP OF CARRY CASE

The case is designed to make installation for mobile collections as easy and as fast as possible. Proceed as follows to install a device:

- Put the case next to the place where D601 should be positioned during donation.
- Open the case and put the lid upside down where D601 should be placed.

- Take out the powercord and adaptor and put them at the place where they should be while donating.
- Put the bottom side of the case on top of the lid (which is still upside down). The lid will serve as a table for the device.

9.2. PLACEMENT IN CARRY BAG

Place the equipment inside the Carry bag as per below picture and secure it with velcro tapes. The tray and power cord can be secured within the velcro. Then close the top cover of the carry bag, as shown below

WARNING:

Do not stack the carry bags on top of each other. Also, not recommended to install D601 on top of carry bag like Carry case.

10 PRECAUTIONS

CAUTION:

The following precautions must be taken while operating the Blood Collection Monitor:

- Do not put a weight of more than 1000gm on the tray.
- Turn OFF the system when not in use.
- Do not press the tray of the Blood Collection Monitor while cleaning the machine.
- Do not lean on the equipment.
- The frequency of maintenance should be in every 6 months that shall be serviced by authorized service personnel or trained Biomedical Engineer of the Hospital.
- Do not open the equipment under any circumstances. Call the service engineer. Warranty / AMC get invalid if equipment is opened by unauthorized or untrained personnel.

11 CLEANING INSTRUCTIONS

- Before Cleaning, be sure to switch OFF the D601 and to disconnect the AC power cord
- Do not use aggressive or corrosive products, abrasive powder, steel wool, abrasive sponges or chemical solvents.
- Clean the BCM using clinical detergents. Do not immerse the unit in any liquid.
- Take out the tray and clean it separately.
- Do not press the tray holding mechanism of the Blood Collection Monitor while cleaning.
- Do not rub the stickers with Hospital spirit or Iso- Propyl Alcohol.
- Be sure to dry the unit, especially the AC power supply connector before using it.

Disinfection Procedure: In any case of fluid spill to the trays please disinfect the surface before using again. Use hospital spirits or Iso- Propyl Alcohol dipped into a soft cloth or cotton to slowly rub on the affected surface.

12 TROUBLESHOOTING

Problems Observed	Possible Cause	Possible Remedy
Equipment is not turning ON when adaptor is connected or in battery operated condition	No power input to the unit	Check the adaptor connected properly
	Battery disconnected	Check the battery connection
No oscillations and no change in display with increase in volume collected	Weight not being sensed	Switch OFF the machine and fix the tray properly. If the problem persists, call the Service Engineer
Clamp Activated before the programmed volume is achieved.	There could be some obstruction to the free movement of the tray	Remove the tray and ensure that there are no obstructions, which could hamper free movement. If the problem persist contact service engineer
Fluctuations in collected volume during Blood collection	<ol style="list-style-type: none"> 1. The blood bag and tubing may be touching some parts of the tray 2. Equipment directly under a fan 	<ol style="list-style-type: none"> 1. Take care to see that the blood bag and tubing are contained within the tray. 2. Change the location of the BCM
No flow rate display and alarms	Problem with internal memory	Reset the real time clock and enter new time. Then check again for the flow rate. If the problem persist contact service engineer

CAUTION:

- Other problems should be corrected by Authorized service personnel only.
- Do not make an attempt to open the equipment under any circumstances.
- Warranty will be void if the equipment is opened.

13 EMC INFORMATION

BCM D601 complies to the requirements for electromagnetic compatibility(EMC) for medical devices as defines in IEC 60601-1-2. The customer or the user should ensure that it is used in such an environment.

Electromagnetic emission		
BCM D601 is intended for use in the electromagnetic environment specified below. The customer or the user of BCM D601 should assure that it is used in such an environment		
Emission test	Compliance	Electromagnetic environment guidance
RF Conducted Emissions CISPR 11	Group A	BCM D601 uses no RF energyfor its internal function. Therefore, its RF emissions arevery low and are not likely tocause any interference in nearby electronic equipments
RF Radiated Emissions CISPR 11	Class A	
Harmonic Emissions IEC 61000-3-2	Class A	
Voltage fluctuations/ flicker emissions IEC 61000-3-3	Compliance	

Electromagnetic immunity			
BCM D601 is intended for use in the electromagnetic environment specified below. The customer or the user of BCM D601 should assure that it is used in such an environment			
Immunity test	IEC 60601 test level	Compliance Level	Electromagnetic environment guidance
Electrostatic discharge (ESD) IEC 61000-4-2	up to $\pm 8\text{kV}$ contact up to $\pm 15\text{kV}$ air	Criteria B	Floor should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30%
Electrical fast transient IEC 61000-4-4	$\pm 2.0\text{ kV}$, 5/50ns	Criteria A	Mains power quality should be that of a typical commercial or hospital environment
High energy Surge IEC 61000-4-5	$\pm 2.0\text{ kV}$, 1.2/50 μs for common mode $\pm 1\text{ kV}$, 1.2/50 μs for differential mode	Criteria A	Mains power quality should be that of a typical commercial or hospital environment

Voltage Dips and Interruption Immunity IEC 61000-4-11	0% short interruption for 5sec 0% Dips for 0.5 cycles 0% Dips for 1 cycles 75% Dips for 25 cycles	Criteria A	Mains power quality should be that of a typical commercial or hospital environment
Power frequency magnetic field IEC 61000-4-8	30 A/m for continuous 1 min.	Criteria A	Power frequency magnetic field should be at levels characteristics of a typical commercial or hospital environment

Radiated Susceptibility 61000-4-3	80 MHz- 2.7 GHz , 3 V/m, AM, 1 kHz, 80%AM	Criteria A	<p>Portable and mobile RF communication equipment should be used no closer to any part of BCM D601, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter. Recommended separation distance</p> $d = 3.5 \sqrt{P}$ $d = 1.2 \sqrt{P}, 80\text{MHz to } 800 \text{ MHz}$ $d = 2.3 \sqrt{P}, 800\text{MHz to } 2.5\text{GHz}$ <p>Where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and d is the recommended separation distance in meters (m). Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey should be less than the compliance level in each frequency range. Interference may occur in the vicinity of the equipment marked with the following symbol</p> 
Conducted RF IEC 61000-4-6	3 Vrms Range 150 KHz - 80MHz	Criteria A	

Note 1: At 80MHz, the higher frequency range applies

Note 2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and peoples

Performance Criteria

A	Normal performance within limits specified by the manufacturer, requester or purchaser.
B	Temporary loss of function or degradation of performance which ceases after the disturbance ceases, and from which the equipment under test recovers its normal performance, without operator intervention.
C	Temporary loss of function or degradation of performance, the correction of which requires operator intervention.
D	Loss of function or degradation of performance, which is not recoverable, owing to damage to hardware or software, or loss of data.

Recommended separation distance between portable and mobile RF communications equipment and BCM D601

The BCM D601 is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of the BCM D601 can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and the BCM D601 as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter(W)	Separation distance according to frequency of transmitter(m)		
	150 KHz to 800 MHz $d=1.17\sqrt{P}$	80MHz to 800 MHz $d=1.17\sqrt{P}$	800MHz to 2.5GHz $d=2.33\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.37	0.37	0.74
1	1.17	1.17	2.33
10	3.69	3.69	7.38
100	11.67	11.67	23.3

The transmitters rated data maximum output power not listed above, the recommended separation distance d in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where P is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

Note1: At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

Note2: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

14 DISPOSAL (END OF LIFE)

Electrical and electronic equipment (EEE) contains materials, components and substances which can be dangerous to the environment and harmful to human health if waste electrical and electronic equipment (WEEE) is not disposed of separately. At end of life, please dispose of this device according to your local regulations or contact Terumo Penpol Private Limited/Terumo BCT/Branch office

15 CONTACT ADDRESS

For other complaints, suggestions please contact your nearest authorized service personnel or please contact:

TERUMO PENPOL Private Limited

Medical Systems Group
T.C.27/373, Andoor Buildings
General Hospital Road,
Thiruvananthapuram 695035, India.
Tel: +91 471 3015800
Fax: +91 471 3015805
Email: dist.penpolservicesection@terumobct.com

Corporate Office:

House No. I - 2, Jawahar Nagar
Thiruvananthapuram -695003, Kerala,
India.
Tel: 91- 471-3015500 / 3015501.
Fax: 91- 471-2721519

TERUMO BCT EUROPE N.V.

Ikaroslaan 41
1930 Zaventem
Belgium
Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

In the event you require some technical clarifications / help, do contact our helpline at the Medical Systems Group. Our technical experts will be able to provide you with the required support. Do let us have your feedback regarding the Equipment. Continuing with our Vision of “Better Ways to Better Health Care” your feedback will help us in serving our Customers better.

Zonal Offices in India

Chennai : Tele : +91-44 -42054528/42054538
Fax : +91-44 -42054528

Ahmedabad : Tele : +91-79 - 26580369/26580370
Fax : +91-79 - 26580370

Kolkata : Tele : +91-33 - 25344609
Fax : +91-33 - 25344623

Mumbai : Tele : +91-22 -28938683/28949944
Fax : +91-22 - 28949944

Bengaluru : Tele : +91 - 80 - 26618858
Fax : +91-80-26522265

Delhi : Tele : +91-11- 25741637/42437170
Fax : +91-11- 25748298

Visit www.terumo.co.jp for TERUMO worldwide offices.

* Due to continuous development, specifications are likely to change without prior notice.

This product is manufactured by Terumo Penpol Private Limited ® Registered Trademark

© Terumo Penpol Private Limited. April 2020. Printed in India.

16 WARRANTY

WARRANTY

Dear Customer,

Please accept our hearty congratulations on choosing TERUMO PENPOL equipment. We are sure that this marks the beginning of a long and mutually satisfying relationship. Details of the warranty period and terms of validity are outlined in this document. We would request you to go through the same.

Apart from covering your installation through the warranty with which your equipment comes, we have post warranty service options available. Your nearest Terumo Branch/Distributor will be happy to give you the details. It is advisable to get in touch before warranty expires. This will enable us to work out an option that suits you best. Once again we thank you for placing your trust in TERUMO PENPOL Private Limited.

Best Regards,

Authorized Signatory

Medical Systems Group

TERUMO PENPOL Private Limited

TERMS OF WARRANTY

The warranty covers repairing/replacement of any defective hardware items from the configuration listed in the Warranty Card, contingent to these defect having arisen only out of faulty design or workmanship if brought to the notice of TERUMO PENPOL Private Limited, in writing, on or before the date at expiry of warranty mentioned on the Warranty Card. Warranty is valid for a period of two years from the date of installation or 27 months from the date of delivery which ever is earlier specifically agreed in writing by TERUMO Branch officials.

TERUMO PENPOL®

WARRANTY

Warranty Number:

Date:

Customer Code:

Customer Name:

Installation Address:

D601

Date of dispatch:

Date of Installation:

DC No.	Description	Eqpt. SL No.

TERUMO PENPOL Private Limited, undertake to repair/replace any defective items from the products listed above subject to the conditions mentioned overleaf.

Installation Remarks:

TERUMO PENPOL Private Limited, Medical Systems Group
TC 27/373, Andoor Buildings, General Hospital Road, Thiruvananthapuram, India
Registered Office: I-2, Jawahar Nagar, Thiruvananthapuram 695003

If this product fails to perform as intended, immediately stop use and contact Terumo Penpol Private Limited/ TERUMO BCT/Branch Office.

FOR INFORMATION ABOUT TERUMO BCT/TERUMO PENPOL PRODUCTS

TERUMOBCT
Unlocking the Potential of Blood

Terumo BCT, Inc.
10811 West Collins Ave.
Lakewood, Colorado 80215-4440
USA

USA Phone: 1.877.339.4228
Phone: +1.303.231.4357
Fax: +1.303.542.5215

Terumo BCT Europe N.V.
Europe, Middle East and Africa
Ikaroslaan 41
1930 Zaventem
Belgium

Phone: +32.2.715.05.90
Fax: +32.2.721.07.70

Terumo BCT Latin America S.A.
La Pampa 1517 – 12th Floor
C1428DZE
Buenos Aires
Argentina

Phone: +54.11.5530.5200
Fax: +54.11.5530.5201

TERUMOBCT.COM

Terumo BCT (Asia Pacific) Ltd.
Room 3903-3903A, 39/F
ACE Tower, Windsor House
311 Gloucester Road
Causeway Bay, Hong Kong

Phone +852.2283.0700
Fax: +852.2576.1311

Terumo BCT Japan, Inc.
Takanawa Park Tower 13F
20-14, 3-chome
Higashi Gotanda, Shinagawa-ku,
Tokyo 141-0022 Japan

Phone +81.3.6743.7890
Fax: +81.3.6743.9800