RADIFOCUS®


RADIFOCUS® HAEMOSTASIS VALVE II

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CE 0197

RS34M017-01

Mã hàng

Số lô

Hạn sử dụng

Nhà sản xuất

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Contents

Nội dung

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Tiệt trùng bằng khí ethylene
oxide

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TERUMO CORPORATION

44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO
151-0072, JAPAN

MADE IN JAPAN

TERUMO EUROPE N.V.

INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM
www.terumo-europe.com

TERUMO MEDICAL CORPORATION

2101 COTTONTAIL LANE, SOMERSET,
NEW JERSEY 08873, U.S.A.

Terumo Australia Pty Ltd

Level 4 Building B 11 Talavera Rd
Macquarie Park NSW 2113 Australia
T: +61 2 9878 5122 E: cs_australia@terumo.co.jp

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
TERUMO CORPORATION

44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO
151-0072, JAPAN MADE IN JAPAN

TERUMO EUROPE N.V.

INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM
www.terumo-europe.com

TERUMO MEDICAL CORPORATION

2101 COTTONTAIL LANE, SOMERSET,
NEW JERSEY 08873, U.S.A.

Terumo Australia Pty Ltd

Level 4 Building B 11 Talavera Rd
Macquarie Park NSW 2113 Australia
T: +61 2 9878 5122 E: cs_australia@terumo.co.jp

RADIFOCUS® HAEMOSTASIS VALVE II Haemostasis valve

DESCRIPTION/DESIGNATION

The RADIFOCUS HAEMOSTASIS VALVE II prevents blood reflux when a catheter is inserted in a vessel for angiographic or other purposes.

FEATURES

- The valve locks securely onto the catheter hub when attached correctly.
- The valve effectively prevents blood reflux when the catheter is inserted in a vessel, whilst allowing manipulation of a catheter or guide wire inserted through the valve.

WARNING / PRECAUTIONS

WARNINGS

- **Use of alcohol, antiseptic solutions or other solvents must be avoided, as they may adversely affect the device.**
- **Apply appropriate anticoagulant therapy to the patient.**

PRECAUTIONS

- This valve shall only be used by a trained physician.
- Sterile and non-pyrogenic in an unopened and undamaged unit package. Do not use if the unit package or the product have been damaged or soiled.
- The product should be used immediately after opening the package and be disposed of safely and properly after use.
- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- The maximum device size the HAEMOSTASIS VALVE II can accept is 7 Fr. (2.33 mm). Use 7 Fr. (2.33 mm) or smaller sized catheter.
- Handle the valve aseptically.
- Do not inject drugs including oil components such as lipid emulsion, castor oil, interfacial active agent or solubilization agent such as alcohol, through the side tube. It may cause cracks on the stopcock.
- Do not excessively pull, push or bend the side tube and the stopcock.
- Be careful not to cut the side tube when holding it with forceps, or not to cut with scissors and knives.
- Do not turn the three way stop cock more than 180 degrees. The cock could misalign or come off and could cause blood leakage or it could shut down the path for the drug.

Rx ONLY CAUTION: Federal law (U.S.A.) restricts this device to sale by or on the order of a physician.

INSTRUCTIONS

1. Remove the HAEMOSTASIS VALVE II from its packaging carefully.
2. Flush and fill the HAEMOSTASIS VALVE II with heparinized physiological saline solution through the 3-way stopcock.
3. Securely lock the HAEMOSTASIS VALVE II onto the catheter hub.

CAUTION

Be sure the securely lock the HAEMOSTASIS VALVE II onto the catheter hub. Otherwise the HAEMOSTASIS VALVE II can become detached from the catheter hub resulting in blood leakage from the catheter and the possibility of air entering the catheter.

4. Aspirate through the side tubing to remove any air bubbles.
5. Flush the HAEMOSTASIS VALVE II through the side tubing with heparinized physiological saline solution and keep the side tubing completely filled with heparinized physiological saline solution.
6. Insert a guide wire or a catheter through the HAEMOSTASIS VALVE II and advance to the desired position.

CAUTION

- Do not inject any fluid through the side tubing by a power injector.
- Do not flush liquid such as contrast media or heparinized saline solution from side tube while a device such as a dilator or a catheter is inside of the valve.
- Before removing or inserting the catheter through the valve, aspirate blood from the 3-way stopcock to remove any fibrin deposition which may have accumulated in or on the tip of the connected catheters.
- Be careful for rapid aspiration by a syringe. It may cause air to be drawn through the valve.

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