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MEDISAFE FIT C Blood Glucose Meter

Blood glucose monitoring system for self-testing

Instructions for use

Be sure to observe

Caution Carefully follow these instructions.

- Avoid direct sunlight when storing the Test TIPS and the meter. Store the Test TIPS at a suitable temperature (1 to 30°C / 34 to 86°F).
- Before using the MEDISAFE FIT C Blood Glucose Meter, carefully read through this "Instructions for use" and carry out the test in accordance with the instructions from your physician.
- Keep this "Instructions for use" at a place that you can access at any time.
- Carefully read the instructions for use of the medical devices to be used prior to testing.
- Make sure to use MEDISAFE FIT Blood Glucose Test TIP with light blue film-seal. Test TIPS with other color film-seal cannot be used.
- Discuss any questions you might have with your physician prior to your first test.
- Treatment and management of diabetes mellitus should be done under supervision of your physician. Your physician will determine your individual target values for your blood glucose. If you have any doubts about the test results, make sure to consult your physician.
- Do not alter dosage and administration frequency of medication without consulting your physician.
- TERUMO will not be liable for any death, injury, or property damage caused by any use not following the instructions provided.
- MEDISAFE FIT C Blood Glucose Monitoring System is an in-vitro diagnostic system for self-monitoring. Do not use this meter for other purposes, and please obtain latest safety information from your physician or our sales representative.
- Never take decisions on your treatment that are based on error messages.
- There are two different units to indicate blood glucose levels (mg/dL and mmol/L). For each unit, there is a dedicated MEDISAFE FIT C Blood Glucose Meter. The unit your meter displays, cannot be changed.
- Carefully check if your meter is set correctly. Otherwise, you can misinterpret results. In case your meter is set at the wrong unit, please contact our sales representative.

Before use

Diabetes mellitus and self-monitoring of blood-glucose (SMBG)

Diabetes mellitus is a disease resulting in a high glucose level in the blood. Untreated diabetes mellitus has been known to lead to the impairment of the liver, kidney, eyes and nervous system. There are many treatment procedures for diabetes mellitus today, and the first important step is not to panic but consult a physician. It is essential that diabetics measure their own blood glucose level and keep the physician informed about any changes in their blood glucose levels. This is the best way to ensure a more effective course of treatment.


Intended Use

The MEDISAFE FIT C Blood Glucose Monitoring System consisting of a portable instrument and reagents is an in-vitro diagnostic system for the quantitative measurement of glucose concentrations in fresh whole blood. This system is intended for use in self-monitoring of blood-glucose by diabetic persons for management of diabetes mellitus.

How MEDISAFE FIT C Blood Glucose Monitoring System works?

After blood is applied to the MEDISAFE FIT Blood Glucose Test TIP, glucose in the blood sample reacts with the enzyme glucose oxidase on the Test TIP and colours the reagent pad. The intensity of the colour is proportional to the amount of glucose in the blood. The MEDISAFE FIT C Blood Glucose Meter measures this colour concentration and shows the blood glucose level on the display. If you have been using another type of meter - one that provides whole-blood results - you may notice that your test results with the MEDISAFE FIT C Blood Glucose Monitoring System are about 11% higher.


Composition and content of this system

Pictures of the Lancing device and Lancet show TERUMO product examples.

- Make sure to use MEDISAFE FIT Blood Glucose Test TIP with light blue film-seal.
- Check if the contents are complete. If something is missing, please contact our sales representative.
- Blood glucose meter, Test TIP, Lancing device, Lancet and Instructions for use, you can only obtain from our sales representative.
- Lancing device, Lancet and Test TIP are sold separately.
- Do not use if damaged.
- For blood sampling, follow the instructions for use provided with the Lancing device and Lancet.


How to install or replace batteries

- Make sure the orientation of the batteries is correct.
- It is recommended to have back-up batteries available for replacement: two lithium batteries CR2032.
- A new battery lasts around 500 tests or approximately 6 months.
- The clock inside the blood glucose meter stops when batteries are removed for replacement. After replacing batteries, reset date and time.
- Dispose of batteries according to your national regulations.

How to set Time and Date

Make sure to set time and date before first use and after replacing batteries.


How to Use

Step 1 Preparation

- 1 Prepare all necessary materials for blood glucose test.

- Lancet
- Blood glucose meter
- Notebook, pen, pencil, etc.
- Lancing device
- Alcohol swab, tissue, etc. for disinfection
- Test TIP

- 2 Before starting test, wash and dry your hands.


Caution

- Before puncture and drawing blood from the fingertip, wash hands thoroughly with running water before puncture. If blood is drawn from fingertip after touching foods containing sugar such as fruits, blood may be contaminated with sugar causing an incorrect test result. [Disinfection with alcohol swab alone does not sufficiently remove sugar.]

Step 2 Attaching Test TIP

- 1 Press Power button.


Once powered on with a beep, current time and date will be displayed. If anything else appears on the display, see "Error messages and actions to be taken".

Note

- If the time setting screen is displayed, see "How to set Time and Date".


- 3 Peel off the film-seal from the Test TIP case completely.

Caution

- Make sure to use MEDISAFE FIT Blood Glucose Test TIP with light blue film-seal.
- Use Test TIP immediately after removing the film-seal. If the Test TIP is left open after removing the film-seal, it can absorb moisture from the air. This could affect the test results.

- 4 Attach Test TIP case to the edge of the meter firmly.


Push the Test TIP straight until you feel a "click". When the Test TIP is attached, the meter beeps and the screen displays "OK". If anything else appears on the display, see "Error messages and actions to be taken".

Note

- If you do not attach it firmly and correctly, the test results would be affected or an error message is shown.
- Do not slide the eject lever during this procedure.

- 5 Pull the Test TIP case straight off.

Note


- Do not discard the case yet. You will need it for disposing the used Test TIP.

Caution

- Attach and use the Test TIP immediately after removing the film-seal.

- 6 Ready to measure.

Confirm that "OK" is displayed and proceed to Step 3 "Blood sampling". Perform measurement as quickly as possible after attaching the Test TIP.


Note

- If the display is blank, press the Power button.
- If you find any other indications, see "Error messages and actions to be taken".
- For information on the "After Meal" symbol, see "How to use the "After Meal" symbol".

Step 3 Blood sampling

- 1 Choose the puncture site.

Choose the side of fingertip as puncture site.

Note

- Prick a different finger and a different location on your finger each time. Always choosing the same spot may harden the skin and cause calluses.

- 3 Disinfect the puncture site with an alcohol swab, and dry the site thoroughly.


Note

- If the puncture site is wet, blood may not form a sphere, and you may not be able to draw blood.

Caution

- Do not reuse the used Lancet. For single use only. Reusing Lancet may cause blood infection.
- For blood sampling, follow the instructions for use provided with the Lancing device and Lancet.

- 5 Gently squeeze your fingertip to obtain a drop of blood.

Caution

- Blood will clot when it is left in the air. Immediately after pricking, proceed with blood testing.
- To prevent blood from dripping, make sure that the puncture site is facing upward.


Step 4 Testing

- 1 Make sure that "OK" is displayed.

If nothing is displayed, press the Power button and wait until "OK" is displayed.

- 2 Gently apply the tip of the Test TIP to the blood.

Caution


Do not push the Test TIP strongly against your fingertip, the Test TIP may not draw blood sufficiently.

- 4 Gently place the meter aside and wait.

Wipe the puncture site with a tissue and keep it clean.

- 3 When you hear the beep, remove the Test TIP from the blood.


Caution

- Do not touch the Test TIP until the test result is shown in the meter display.
- Stop bleeding (for example by using sterile adhesive tape), if necessary.
- Make sure that blood or devices contaminated with blood, including tissue paper and other items you have used, do not come into contact with any other person.

- 5 Check the test result.

After the beep, the measurement time and test result will be displayed.

*This display is in mg/dL version.

To add the "After Meal" symbol, press the M button. To delete the "After Meal" symbol, press and hold the M button.

If an error message is displayed, see "Error messages and actions to be taken".

- 6 Record the test result in your Notebook.

For details about how to read the test result, see the instructions for use of the MEDISAFE FIT Blood Glucose Test TIP.


Caution

Treatment and management of diabetes mellitus should be done under supervision of your physician. Your physician will determine your individual target values for your blood glucose. If you have any doubts about the test results, make sure to consult your physician. Do not alter dosage and administration frequency of medication without consulting your physician.


Step 5 Cleaning up

- 1 Remove the Test TIP.

Place the empty Test TIP case over the Test TIP.

- 2 Slide the eject lever forward.

- 2 Turning off the power.

Keep pressing Power button for at least 1 second. If the display is blank, this step is not necessary.

Auto power off function

If no operation is performed after the meter is turned on, it will automatically shut down after 2 minutes. (After 5 minutes when "OK" is displayed.)

- 4 Remove Lancet from the Lancing device.

For instructions on how to remove the Lancet, refer to the instructions for use of the Lancing device.

Caution

- Do not reuse the used Lancet. Doing so could cause infection.
- When storing the lancing device, remove the Lancet, otherwise it may become impossible to correctly adjust the puncture depth.

Note

- Be sure to cover the used Test TIP with the TIP case before ejecting. Otherwise blood may be spilled when the Test TIP is ejected.

Caution

- If a fragment of a broken Test TIP remains in the meter, lightly shake the meter or give it a gentle flick with your finger, with the tip attachment facing down, to remove it. Do not use a tool with a sharp tip to remove Test TIP. Otherwise the meter may be damaged and the eject lever may lose mobility.
- Do not remove the Test TIP with your hand directly. Make sure to use the eject lever to remove the Test TIP.
- If you cannot remove the part, contact our sales representative.

- 3 Place the protective cap on the meter.

Note

- If you do not use the protective cap, the inside of the meter can get dirty and test results would be affected.

- 5 Dispose of used Test TIP and Lancet.

Caution

- Dispose of the used Test TIP and used Lancet and any contaminated material as instructed by your physician or in accordance with local regulations for disposal of contaminated materials.

Error messages and actions to be taken

If the following messages are displayed, follow the instructions described below.

Display	Cause	How to deal with it
E-1	• The Test TIP is inserted with an angle. • The Test TIP is not fully inserted.	If you have not yet commenced drawing up the blood, firmly press the Test TIP into place until it clicks. If there is blood on the Test TIP, replace it with a new one and repeat test.
E-2	A used Test TIP is attached to the meter. There may be dirt or dust on the measuring window of the meter.	Replace the Test TIP with a new one and repeat test. Wipe the measuring window with a cotton swab or similar, attach a new Test TIP and repeat test.
E-3	There may be dirt or dust on the measuring window of the meter. The amount of blood drawn up by the Test TIP is insufficient. The meter was shaken and then measurement started before a sufficient amount of blood could be drawn up by the Test TIP. Long time has elapsed after blood sampling.	Remove the Test TIP and inspect the measuring window. If the measuring window is not dirty, the display will disappear. If the display does not disappear, wipe the measuring window with a cotton swab or something similar. Replace the Test TIP with a new one and draw up an appropriate amount of blood (a droplet of approximately 2.5 mm) into the Test TIP to repeat test.
E-4	Eject lever is accidentally pressed during testing, and the Test TIP is loose. You placed the tip into the blood sample twice.	Replace the Test TIP with a new one and draw up an appropriate amount of blood (a droplet of approximately 2.5 mm) into the Test TIP to repeat test.
E-5	Water or disinfecting alcohol are mixed in the blood sample. Other substances than blood were measured (e.g. water).	Adequately dry the finger from which the measurement will be taken, replace the Test TIP with a new one, and repeat test. Make sure that the measurement is only with blood (without any other substances mixed with the blood).
E-6	Test TIP was not used immediately after it was opened. Expired Test TIP was used.	Open a new Test TIP and use it to do test immediately after opening. Replace the Test TIP with one that has not exceeded its expiration date and repeat test.
E-7	The hematocrit value may be outside the measurable range.	If the hematocrit of the blood is above 60% or below 20% the measured value may not be displayed. After checking for other causes of E-3, replace the Test TIP with a new one and repeat test. If an error message is displayed again, contact your physician or our sales representative.
E-8	A measurement error message appears twice in a row.	After rechecking the handling method for E-3, contact your physician or our sales representative.
E-9	The ambient light was too bright and measurement was not performed.	After replacing Test TIP with a new one, point Test TIP toward dark area or move to a location not in direct sunlight or strong lighting and repeat test.

Display	Cause	How to deal with it
"E-6" is displayed.	The Test TIP became detached during measurement.	Replace the Test TIP with a new one and repeat test.
The sun mark is displayed.	The ambient light is too bright to perform measurement. The sun mark is displayed.	If the Test TIP is not attached to the meter, attach it. If the Test TIP is attached, turning the Test TIP toward the shadows or moving to an area out of direct sunlight will cause the display to disappear. If there is no blood on the Test TIP, confirm that the display has disappeared and continue with the measurement. If there is blood on the Test TIP, replace it with a new one and repeat test.
"E-7" is displayed.	The Test TIP became detached during measurement.	Replace the Test TIP with a new one and repeat test.
The temperature mark is displayed.	You tried to perform measurement in a place not within the appropriate temperature (5 to 40°C/41 to 104°F).	Keep the meter out of the carrying case at a location within the appropriate temperature for about 20 minutes. Wait until the error message disappears and repeat test.
"Hi" is displayed.	The blood glucose level is higher than 600 mg/dL (33.3 mmol/L).	Replace the Test TIP with a new one and repeat test. If the same result is displayed again, consult your physician.
"Lo" is displayed.	The blood glucose level is below 20 mg/dL (1.1 mmol/L).	Follow your physician's instructions. Then replace test TIP with a new one and repeat test.
"E14" is displayed.	The meter is out of order.	Please contact your physician or our sales representative.
The battery mark blinks.	Batteries are almost empty.	Replace with new batteries (lithium battery CR2032 X2) as soon as possible.
The battery mark is displayed.	Batteries are empty.	The meter cannot take measurements. Replace with new batteries (lithium battery CR2032 X2) immediately.
The meter cannot be turned off. (Screen example)	The \odot button has not been pressed long enough.	Press the \odot button for more than 1 second. The meter will automatically shut down when not performing any operation. (Auto power off function)
Nothing is displayed.	The power button was not pressed. The batteries have not been installed correctly. The batteries are empty.	Press the \odot button. Remove the batteries and install them correctly. Replace both batteries (lithium battery CR2032 X2) immediately. If the screen remains blank, contact your physician or our sales representative.
The display became blank.	2 minutes have passed without the meter being operated (5 minutes when \odot is displayed).	This is not a malfunction. Press the \odot button to turn the meter on.

Caution
Never take decisions on your treatment that are based on error messages.

How to check previous test results


The meter can store up to 500 test results in its memory.

- When the power is not on, press the \odot button.

Note

- When \odot is displayed, test results cannot be displayed. Operate the meter with the Test TIP removed.

- Press the M button. The most recent test result is displayed.

Note

- To display test results from the oldest to the most recent, press the \odot button.
- Press and hold the M button to move the test results quickly.



- After checking the test results, press and hold the \odot button to turn off the meter.

Note

- The test results of "Lo" (lower than 20 mg/dL or 1.1 mmol/L), and "Hi" (higher than 600 mg/dL or 33.3 mmol/L) are recorded as such.
- When more than 500 tests have been carried out, the oldest test results will be deleted to free up space for the new test results.
- Even when the batteries are removed, the test results are kept in memory.

How to check the average of test results



- When the power is not on, press the \odot button.
- Press the M button. The most recent test result is displayed.
- Press the \odot button. The 7-day average for "After Meal" testing is displayed.

- After checking the test results, press and hold the \odot button to turn off the meter.

How to use the "After Meal" symbol

When testing after a meal, you can record the test result with the "After Meal" symbol. This is available in order to sort results by testing time.


- To delete the "After Meal" symbol, press and hold the M button.

- The "After Meal" symbol cannot be added or deleted from the past test results in the memory.

How to use the alarm

Set the alarm to remind you of the time for testing at the same time each day.

- With the power off, press and hold the M button. "F1" screen appears.
- While "F1" screen is displayed, press the \odot button.
- Press the M button and select the alarm number (1 to 4) displayed in the upper right corner. Press the \odot button to set the chosen alarm number.
- Press the M button and select "ON" to activate the alarm. Pressing the M button switches between "ON" and "OFF".

- Press the \odot button to set the selection. Set the hour, then the minutes.

- Pressing the M button will increase the number.
- Press the \odot button to set the selection.

Note


- Once the minute setting has been made, the meter emits a short beep and the alarm is set.

- To end the alarm setting, hold down the \odot button and turn off the meter.

How to mute the meter

The sounds emitted during operation of the meter and the alarm can be muted.

- With the power off, press and hold the M button. The "F1" screen appears.
- Press the M button. When "F2" screen is displayed, press the \odot button.
- Press the M button and select buzzer "OFF". Pressing the M button will switch between "ON" and "OFF".


- Press the \odot button to set the selection. When the setting is complete, the meter will emit a beep and then turn off.

How to clear test results and reset the meter

All test results are deleted, including the ones done with control solution.

- The buzzer setting and alarm setting will return to the factory setting.
- This operation must be carried out by your physician.

- With the power off, press and hold the M button. The "F1" screen appears.
- Press the M button 4 times. When "F5" is displayed, press the \odot button.

When all test results are deleted, the meter will turn off automatically.

Note


- Once the test results have been deleted, they cannot be restored.
- The time and date do not return to the factory setting.

Maintenance of the meter

Cleaning of the measuring window


If there is dirt or dust on the measuring window of the meter, "E-1" (E-1) or "E-2" (E-2) is displayed and the meter cannot perform measurements.

Wet a new cotton swab with a small amount of water and wipe off dirt or dust along the shape of the measuring window. Then, wipe again with a dry cotton swab.

Cleaning of the meter

If there is dirt on the meter or blood on the tip attachment, wipe it off using a cloth or tissue infused with a small amount of water or disinfecting alcohol.

Caution

- Do not use thinner or benzene on the meter.
- The meter itself is not water-proof. Wring out the fabric or swab to reduce the moisture when cleaning.
- Do not immerse the meter in any liquid.

Meter Operation Check

When performing inspection and maintenance, check that the exterior is not dirty or damaged and that no part of the display is missing.

Item	Frequency	Details
Appearance	Every time	There is no dirt, damage, etc.
Display	Every time	There is no missing icon when the meter is turned on and the backlight is lit.
Self-check	Every time	If \odot is displayed after attaching the Test TIP, the meter is operating normally.

Self-check function during regular tests

The meter performs a self-check automatically every time a test is performed. If \odot is displayed after attaching a Test TIP, the meter is operating as normal. If anything else appears on the display, see "Error messages and actions to be taken".

Control Solution for MEDISAFE FIT C Blood Glucose Monitoring System

MEDISAFE FIT Control Solution verifies that the MEDISAFE FIT C Blood Glucose Meter and MEDISAFE FIT Blood Glucose Test TIP are functioning properly. You do not need to do this control test before testing glucose every time. When you feel Test TIP and meter are not functioning properly, you are able to check the system with the control solution.

- With the power off, press and hold the M button. The "F1" screen appears.
- Press the M button 3 times. When "F4" is displayed, press the \odot button. The "M" mark appears when setting the meter to the control solution testing mode.
- Follow "How to Use" in the instructions for use of MEDISAFE FIT Control Solution to perform control solution testing.
- When the testing is completed, hold the \odot button for more than one second to turn the power off.
- Tolerance is indicated on the MEDISAFE FIT Blood Glucose Test TIP box label. For details about comparing the results of measurements with the tolerance range of Control Solution, see "Test Results" in the instructions for use of MEDISAFE FIT Control Solution.

Waste and Recycle of the meter and batteries

Caution

- The meter may come into contact with blood during use. Please dispose of the meter according to your local regulations. Contact your local distributor or municipality for information about correct disposal.
- Before disposal of the meter, batteries should be removed.

About the battery

- Batteries must not be disposed of within general household waste, but instead, should be collected separately. In this way, the environmental impact associated with disposal is reduced, and the opportunity for recycling. Contact your local distributor or municipality to know the available collection schemes.
- Do not throw batteries into fire. It can cause explosion.

Safety Cautions

Caution

About Blood Glucose Testing

- Do not reuse the Lancet. For single use only. Reusing Lancet may cause blood infection.
- When blood is collected from the fingertip, the hand needs to be washed well with running water before puncture. Dry your hands thoroughly to avoid inaccurate results. Keep the puncture site clean. Otherwise, it may cause blood infection.
- After using the blood sample, wipe the puncture site with a dry and clean tissue paper. Stop bleeding (for example by using sterile adhesive tape), if necessary.
- Make sure that the meter and the Lancing device contaminated with blood, including tissue paper and other items you have used, do not come into contact with any other person.
- In case of suspected hypoglycemia, collect the blood from the fingertip (or the palm) or earlobe. The test results using samples taken from the arm may differ from fingertip or earlobe.
- Always keep the devices needed for testing away from children. There is a risk of children accidentally swallowing the Lancet, batteries, Test TIP and TIP case.
- Do not disassemble or rebuild the meter and Lancing device. It can cause product breakdown or injury.
- Only use meter for measuring blood glucose. Otherwise it can cause product failure or injury.
- Do not use expired Lancets. It may cause failure or injury. The expiration date is printed on each box.
- Do not store the Lancing device and the Lancet at extreme temperatures and humidity.
- If blood is collected from the fingertip immediately after contact with food or the like containing sugar (e.g., fruit), the sugar adhering to the fingertip can mix with blood, possibly leading a falsely high blood glucose value. There is a report that disinfection with alcohol swab is not sufficient in removing the sugar.
- If blood is collected from the fingertip of the patients with reduced peripheral blood flow listed below, the blood glucose value may be lower than the actual value. In such case, blood collected from another site (e.g., venous blood) should be used for testing.
 - Patients with dehydration
 - Patients in shock
 - Patients with disturbed peripheral circulation

About Blood Glucose Testing (Cont.)

- Before testing, Test TIP and the meter should be kept at the testing location minimum twenty (20) minutes before testing. Temperature difference could cause inaccurate test results.
- Testing should be performed indoors in an appropriate setting where the temperature is between 5 - 40°C/41 - 104°F and relative humidity 30 - 85%rh (non-condensing). If the temperature is too low or too high, test results may not be accurate.
- If the blood sample contains more than 5.26 mg/dL (0.30 mmol/L) ascorbic acid, test results may be lower by more than 10%.
- If the blood sample contains more than 85 mg/dL (2.76 mmol/L) glutathione, test results may be lower by more than 10%.
- Test results may be affected if the blood sample contains significant amount of fluoride as a glycolysis inhibitor.

About the meter

- Do not use the meter where it is exposed to strong light such as direct sunlight. It can prevent the meter from testing.
- Use the meter as far away as possible from any equipment generating electromagnetic waves such as cellular phones. The meter may not function correctly in close proximity to equipment of this type.
- Do not drop the meter or expose it to any impact. Do not use the meter if damaged by dropping or exposing to strong impact. In addition, do not place the meter in locations exposed to strong vibration (such as the dashboard of a vehicle).
- Always keep the edge (and measuring window) of the meter clean. When not in use, attach the protective cap. Store the meter within temperature (-10 - 50°C/14 - 122°F) and relative humidity 30 - 95%rh (non-condensing).
- This meter may not be able to be used properly in dry environment or in case there are synthetic fabric nearby due to static electricity.
- When replacing batteries, make sure not to let substances, such as oil or dust, come inside the meter. It may cause product breakdown.
- Prevent foreign matters and fluid from entering the meter. Do not use the meter if damaged by such foreign matters/fluid.
- The clock incorporated in the meter will stop if the batteries are removed. Set date and time again.

About Test TIP

- Do not use a Test TIP if its case or film-seal is open or damaged, or if any of the items is spoiled.
- Do not use a Test TIP after its indicated expiration date has passed. The expiration date is printed on each box of thirty (30) Test TIPs, as well as on the film-seal of the TIP case.
- Store the Test TIP at the suitable temperature (1 to 30°C/34 to 86°F).
- The MEDISAFE FIT Blood Glucose Test TIP is for single use only. Do not reuse.
- When Blood samples are obtained in test tubes (or some other containers), care should be taken that only heparin can be used as anticoagulant. Do not use any other anticoagulants or antiglycolysis reagents.
- Before testing, shake the test tube (or any other container) to mix the blood thoroughly.
- If disinfectant or the like has contaminated the sample, test results may not be accurate.
- There is no significant influence on test results with blood samples having a hematocrit between 20% and 60%. If hematocrit is above 60% (for example, Neonates) or below 20%, test results may be affected.
- Only use MEDISAFE FIT Blood Glucose Test TIP with light blue film-seal.

Specifications of the meter

Name	MEDISAFE FIT C Blood Glucose Meter
Range of Measurement	Blood glucose level of 20 to 600 mg/dL (1.1 to 33.3 mmol/L)
Operating Environment	Avoid direct sunlight. Ambient temperature (5 to 40°C/41 to 104°F), Relative humidity (30 to 85%rh) (non-condensing).
Transport and Storage Environment	Avoid direct sunlight. Ambient temperature (-10 to 50°C/14 to 122°F), Relative humidity (30 to 95%rh) (non-condensing).
Power Source	Lithium batteries CR2032 (X2)
Rated Voltage	3 V
Power Consumption	When meter is turned on: approx. 180 mW When meter is turned off: approx. 100 μ W
EMC Compliance	This product complies with IEC 61326-1: 2012 and IEC 61326-2-6: 2012. CISPR Group 1; Class B
Battery life	Approximately 500 uses or 6 months of use
Memory Capacity	Max. 500 tests (automatically stored)
Clock Accuracy	Within \pm 2 minutes/month
Dimensions	Width: approx. 105 mm, Depth: approx. 33 mm, Height: approx. 21 mm (including protective cap)
Weight	Approx. 50 g (including two batteries)

Specifications and appearance of the product is subject to change without prior notice.

Performance characteristics

① Precision

Repeatability: When standard blood samples were tested with 300 replicates for each concentration across 3 Test TIP lots with 10 MEDISAFE FIT Blood Glucose Meters, repeatability was as follows.

Blood glucose concentration [mg/dL (mmol/L)]	30 - 50 (1.7 - 2.8)	51 - 110 (2.9 - 6.1)	111 - 150 (6.2 - 8.3)	151 - 250 (8.4 - 13.9)	251 - 400 (14.0 - 22.2)
Mean glucose value [mg/dL (mmol/L)]	45 (2.5)	106 (5.9)	143 (7.9)	175 (9.7)	385 (21.4)
Pooled SD [mg/dL (mmol/L)]	2.3 (0.13)	2.3 (0.13)	2.8 (0.16)	3.8 (0.21)	8.3 (0.46)
Pooled CV (%)	5.2	2.2	2.0	2.0	2.2

Intermediate precision: When MEDISAFE FIT Control Solutions were tested on 3 Test TIP lots and 10 MEDISAFE FIT Blood Glucose Meters on each of 10 days, intermediate precision was as follows.

Glucose concentration	Low	Middle	High
Mean glucose value [mg/dL (mmol/L)]	38 (2.1)	113 (6.3)	331 (18.4)
Pooled SD [mg/dL (mmol/L)]	1.3 (0.07)	2.5 (0.14)	7.2 (0.40)
Pooled CV (%)	3.4	2.2	2.2

② System accuracy

System accuracy results for glucose concentration < 100 mg/dL (5.55 mmol/L)

Within \pm 5 mg/dL (Within \pm 0.28 mmol/L)	Within \pm 10 mg/dL (Within \pm 0.56 mmol/L)	Within \pm 15 mg/dL (Within \pm 0.83 mmol/L)
125/156 (80.1%)	156/156 (100%)	156/156 (100%)

System accuracy results for glucose concentration \geq 100 mg/dL (5.55 mmol/L)


Within \pm 5 %	Within \pm 10%	Within \pm 15%
315/444 (70.9%)	426/444 (95.9%)	443/444 (99.8%)

System accuracy results for glucose concentrations between 27 mg/dL (1.48 mmol/L) and 508 mg/dL (28.2 mmol/L)

Within \pm 15 mg/dL or \pm 15% (Within \pm 0.83 mmol/L or \pm 15%)
599/600 (99.8%)

- Reference method: GLUC2 Glucose HK (Hexokinase Method; Roche)
- Regression formula: $< \text{Lot 1} > Y = 1.00x + 2.81$ ($Y = 1.00x + 0.16$)
- $< \text{Lot 2} > Y = 1.03x + 0.25$ ($Y = 1.03x + 0.01$)
- $< \text{Lot 3} > Y = 0.99x + 3.75$ ($Y = 0.99x + 0.21$)
- In brackets is the regression formula of mmol/L unit.

Correlation graph (MEDISAFE FIT Blood Glucose Test TIP vs GLUC2 Glucose HK)

Consensus error grids showed that all glucose concentration results fell within zone A (no effect on clinical action).

③ User Accuracy

A study evaluating glucose values from fingertip capillary blood samples obtained by 100 lay persons showed the following result.

Within \pm 15 mg/dL (\pm 0.83 mmol/L) for glucose concentration below 100 mg/dL (5.55 mmol/L)	Within \pm 15% for glucose concentration at or above 100 mg/dL (5.55 mmol/L)
101/103 (98.1%)	


- Measurement range: 20-600 mg/dL (1.1-33.3 mmol/L)
- Measurement value: Plasma glucose value (whole blood sample; results displayed as plasma glucose value)
- Traceability

The reference method is the glucose analyzer (the glucose oxidase method), and the quality for proofreading of calibrator is glucose standard serum.

List of products

Item	Product Name	Quantity	Product Code
Meter	MEDISAFE FIT C Blood Glucose Meter	1	MS*FR301BM
Test TIP	MEDISAFE FIT Blood Glucose Test TIP	30	MS*FR301BM1
Lancing device	MEDISAFE FINETOUGH	1	MS*GN02
Lancet	MEDISAFE Lancet for FINETOUGH	30	MS*GN4530
Control Solution	MEDISAFE FIT Control Solution	L	MS*CSL
		M	MS*CSM
		H	MS*CSH

About symbols

TERUMO CORPORATION
44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN
TERUMO EUROPE N.V. INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM

TERUMO is a trademark of TERUMO CORPORATION.
®: Registered trademark

©TERUMO CORPORATION OCT. 2015 15H13
MS_FR301_E_500_001