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


Do not use the product, its sterile barrier system or its packaging if damaged or shows any sign of deterioration. Do not use the product, its sterile barrier system or its packaging if damaged or shows any sign of deterioration. Do not use the product, its sterile barrier system or its packaging if damaged or shows any sign of deterioration. Do not use the product, its sterile barrier system or its packaging if damaged or shows any sign of deterioration.

Do not use the product, its sterile barrier system or its packaging if damaged or shows any sign of deterioration. Do not use the product, its sterile barrier system or its packaging if damaged or shows any sign of deterioration. Do not use the product, its sterile barrier system or its packaging if damaged or shows any sign of deterioration. Do not use the product, its sterile barrier system or its packaging if damaged or shows any sign of deterioration.

TERUMO CORPORATION 44-1, 2-CHOME, HATAGAYA, SHBUYA-KU, TOKYO 151-0072, JAPAN MADE IN JAPAN
 TERUMO EUROPE N.V. INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM
 © TERUMO CORPORATION October 2009 08/28

ENGLISH

DESCRIPTION / DESIGNATION

This device is intended for pooling of up to six buffy-coats, for in process filtration and preservation of leukocyte removed platelet concentrates. The integrated platelet storage bag is suitable for storage of platelets at 22±2°C up to 5 days on a platelet agitator.

PRECAUTIONS

This device contains DEHP. Based on animal studies, significant exposure to DEHP may interfere with the normal development of the male reproductive tract. For children, pregnant and nursing woman alternative devices may be appropriate.

INSTRUCTION FOR USE

1. Pooling Procedure

1. Open the packaging.
 2. Close all clamps on the set.
 3. Use a sterile connecting device to connect the platelet medium bag to the platelet medium arm and the buffy-coat bags to the pooling arms.

2. Rinsing Procedure

6. After closing clamp "3", transfer the platelet medium into the buffy-coat bags by using the clamps "1" and "2". Rinse properly the buffy-coat bags and open clamp "3" to pool again into the pooling bag.

3. Centrifugation Procedure

7. At the end of the pooling and rinsing, seal off the tube between the pooling bag and the Y and dispose the emptied bags.

4. Separation and filtration of the platelet concentrate

10. Carefully remove the centrifuged BP-KIT from the centrifugation cup and place the pooling bag into an extractor to apply pressure (Fig. 3).

5. Air removal Procedure from the Platelet Storage Bag and Sample Collection for Quality Control

15. Agitate the platelet storage bag to ensure adequate mixing.
 16. Slightly incline the platelet storage bag to draw the air in the corner closest to the sampling bag.

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