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
Rx Only

Indications For Use

T-Cuff is indicated for use in donor populations during apheresis procedures as an alternative method to maintain pressure on the arm and obtain venipuncture. The T-Cuff operates similar to a pneumatic tourniquet intended to partially restrict blood flow on the upper arm to result in optimal venous access during apheresis collection procedures.

T-Cuff Overview

The T-Cuff is a pneumatic tourniquet. It consists of an inflatable cuff, a pressure release valve, and a bulb for inflation.




Figure 1: T-Cuff

- 1 Cuff:** Consists of a cloth cover and an internal bladder that applies pressure to a donor's arm when inflated
- 2 Bulb:** Allows you to inflate the cuff
- 3 Valve cap:** Allows you to switch mode of operation
- 4 Quick release button:** Allows you to deflate the cuff
- 5 Valve:** Controls cuff pressure

The T-Cuff has two modes of operation; venipuncture mode and collection mode. The mode of operation is changed by adjusting the valve cap.

Venipuncture mode: For use prior to and during venipuncture, this mode allows the T-Cuff to inflate to assist with vein selection. In this mode, the pressure release valve prevents the cuff pressure from exceeding 120 mmHg, with an average maximum of 100 mmHg. The T-Cuff is in venipuncture mode when the valve cap is lowered and the indicator arrows are not aligned.

Collection mode: For use during collection, this mode enables the system to maintain consistent pressure, when assisted by donor squeezing. In this mode, the pressure release valve prevents the cuff pressure from exceeding 40 mmHg. The T-Cuff is in collection mode when the valve cap is raised and the indicator arrows are aligned.

Procedures for Use

- Place the T-Cuff on the donor's upper arm with a snug fit.
 - Position the cuff a minimum of 2 inches above the antecubital area, with the T-Cuff label facing out.
 - Ensure that the cuff is oriented on the donor's arm so that the tubing extends up towards the donor's shoulder to prevent the tubing from interfering with the venipuncture site.

Caution: Do not place the T-Cuff on skin that is compromised or broken.

2. Place the valve in venipuncture mode by rotating the valve cap. Rotate it so that the valve cap is lowered, the indicator arrows are not aligned with each other, and the hexagonal edges do not align as shown in the following figure. When the cap has been rotated far enough to switch modes, you hear an audible click and feel it snap into place.

Figure 2: Venipuncture mode valve position

3. To use the T-Cuff in venipuncture mode:
 - a. Inflate the cuff by squeezing the bulb.
 - b. Identify the venipuncture site.
 - c. Deflate the cuff by pressing and holding the quick release button.

Figure 3: Quick release button

- d. Prepare the venipuncture site according to your center's SOP.
 - e. Inflate the cuff.
 - f. Perform venipuncture.
4. After venipuncture, switch to collection mode by rotating the valve cap. Rotate it so that the valve cap is raised, the indicator arrows align with each other, and the hexagonal edges are aligned as shown in the following figure. When the cap has been rotated far enough to switch modes, you hear an audible click and feel it snap into place.

Figure 4: Collection mode valve position

5. Place the bulb in the donor's hand and instruct the donor to squeeze approximately every 5 seconds throughout the collection procedure.

Cleaning Instructions

Clean the surface of the cuff cover, tubing, bulb and valve as necessary using a damp cloth and mild detergent or 70% isopropyl alcohol wipe. Do not submerge.

To clean the cuff cover separately do the following:

1. Remove the bladder from the cuff cover.
2. Hand wash with a mild detergent and warm water.


When finished cleaning, insert the bladder into the cover opening and close the hook and loop fastener.

Specifications

Table 1: Environmental specifications

Characteristics	Performance
Ambient operating air temperature	15.5 °C to 27.7 °C (60 °F to 82 °F)
Ambient operating humidity	8% to 80% relative humidity (RH), non-condensing

Table 2: Symbols and certifications

Symbol	Definition	Symbol	Definition
	Indicates the medical device manufacturer.		Indicates the date when the medical device was manufactured (or sterilization date, if the product is sterile).
	Indicates the manufacturer's catalog number so that the medical device can be identified.		Indicates that the user should consult the instructions for use.
	Indicates that the product was manufactured in accordance with Annex II of the European Council Directive, 93/42/EEC, as amended.		Indicates the European Authorized Representative of the product when accompanied by the name of the European Authorized Representative.

The valve has been verified to a 3 year life span. The bulb, tubing, and bladder may be replaced as needed. Contact your local sales representative for options regarding replacement parts.

Return of Used Product

If for any reason this product must be returned to Terumo BCT, Inc., a returned goods authorization (an RGA number) is required from Terumo BCT prior to shipping.

IT IS THE RESPONSIBILITY OF THE HEALTH CARE INSTITUTION TO ADEQUATELY PREPARE AND IDENTIFY THE PRODUCT FOR RETURN SHIPMENT.

Please contact your local representative for information regarding returned goods and product complaints.

C E 0086


Terumo BCT, Inc.
10811 W. Collins Avenue
Lakewood, Colorado 80215
USA

USA Phone: +1.877.339.4228

Phone: +1.303.231.4357

USA Fax: +1.866.715.6768

Fax: +1.303.542.5215

Terumo BCT Europe N.V.
Ikarooslaan 41
1930 Zaventem
Belgium

Phone: +32.2.715.05.90

Fax: +32.2.721.07.70

TERUMOBCT.COM

©2017 Terumo BCT, Inc.

2017-03
Part No. 777019-939