TB69A000C2

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TERUMO

TERUFUSION

CE 0123

MD

STERILE EO

VOL

8ml(40°C)
11ml(40°C/2 bar)

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

TERUMO
TERUFUSION

ENGLISH

BLOOD ADMINISTRATION SET for infusion pump

<PRECAUTIONS> • Do not use if protective caps are loose or displaced. • Do not use if unit package is damaged. • Use immediately after opening the unit packaging. • Safely dispose of after use as medical waste following local hospital/health care institution policies to avoid risk of infection. • Ensure that all connections made with the Terufusion blood administration set are secure. • Do not turn upside down or shake the drip chamber during administration. • Disinfect the injection port before use. • We recommend to switch the tubing section used in the pump after maximum 24h to another position, as the flow rate may deviate from the set value due to deformation of the tubing section in the pump after long period of use. • For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device. • This product has been sterilized by ethylene oxide gas.

DO NOT STORE AT EXTREME TEMPERATURES.

• Do not detach the Anti Free Flow clip (AFF clip) from the administration set. • Do not slide up or down the AFF clip while closed. • Make sure to use only TERUMO TERUFUSION MIDPRESS pumps (with AFF function if required).

- ★ a. If the air filter is available, please remove the cap of the filter as showed in Step. 5;
- ★ b. If the AFF clip is available, 1. Make sure that the AFF clip is not closed. 2. Slide the AFF clip down. 3. Make sure that the tube is positioned in the AFF clip as indicated in Fig1-A, and insert the clip into the dedicated aperture of the infusion pump. Then insert the tube properly into the tube guides. (Caution) Read the Instruction manual of the infusion pump before use.

<REPORT OF INCIDENT> If during the use of this device or as a result of its use, a serious incident has occurred, please report it to the manufacturer and/or its authorized representative and to your national authority.

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© : TERUMO CORPORATION 2023/05/Rev.1
All brand names are trademarks or registered trademarks owned by TERUMO CORPORATION, its affiliates, or unrelated third parties.
**TERUMO MEDICAL PRODUCTS (HANGZHOU) CO., LTD.**
AU Sponsor **Terumo Australia Pty Ltd** Macquarie Park NSW 2113 Australia