Ryujin™ Plus OTW

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CE 0086

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1. Bentivoglio LG, et al : Percutaneous Transluminal Coronary Angioplasty (PTCA) in Patients with Relative Contraindications : Results of the National Heart, Lung, and Blood Institute PTCA Registry. Am J Cardiol 1984;53:82C-88C
2. Corbelli J, et al : Percutaneous Transluminal Coronary Angioplasty After Previous Coronary Artery Bypass Surgery. Am J Cardiol 1985;56:398-403
3. Cowley MJ, et al : Efficacy of percutaneous transluminal coronary angioplasty : Technique, patient selection, salutary results, limitations and complications. AMHEARTJ 101:272, 1981
4. Dorros G, et al : Percutaneous Transluminal Coronary Angioplasty : Report of Complications from the National Heart, Lung, and Blood Institute PTCA Registry. Circulation 67, No.4, 1983
5. Gruentzig AR : Percutaneous transluminal coronary angioplasty : six years' experience. AMHEARTJ 107:818, 1984
6. Kelsey SF, et al : Effect of Investigator Experience on Percutaneous Transluminal Coronary Angioplasty. Am J Cardiol 1984;53:56C-64C
7. Kent KM, et al : Long-Term Efficacy of Percutaneous Transluminal Coronary Angioplasty (PTCA) : Report from the National Heart, Lung, and Blood Institute PTCA Registry. Am J Cardiol 1984;53:27C-31C
8. Meier B, et al : Higher balloon dilatation pressure in coronaryangioplasty. AMHEARTJ 107:619, 1984
9. Vandormael MG, et al : Immediate and Short-Term Benefit of Multilesion Coronary Angioplasty : Influence of Degree of Revascularization. J Am Coll Cardiol 1985;6 :983-91
10. Ryan TJ, et al : Guidelines for Percutaneous Transluminal Coronary Angioplasty. JACC Vol.22, No.7 December 1993:2033-54

TERUMO CORPORATION 44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN
MADE IN JAPAN

EC|REP **TERUMO EUROPE N.V.** INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM www.terumo-europe.com

Terumo Australia Pty Ltd Level 4 Building B 11 Talavera Rd Macquarie Park NSW 2113 Australia
T: +61 2 9878 5122 E: cs_australia@terumo.co.jp

®: Registered Trademark
®: Marque Enregistrée
®: Eingetragenes Warenzeichen
®: Marca Registrada
®: Marchio Registrato
®: Geregistreerd Handelsmerk
®: Registrerat Varumärke
®: зарегистрированный товарный знак

Ryujin is a trademark of TERUMO CORPORATION.
Ryujin est une marque de TERUMO CORPORATION.
Ryujin ist eine Marke der TERUMO CORPORATION.
Ryujin es una marca de TERUMO CORPORATION.
Ryujin è un marchio della TERUMO CORPORATION.
Ryujin is een handelsmerk van TERUMO CORPORATION.
Ryujin är ett varumärke tillhörande TERUMO CORPORATION.
Ryujin является товарным знаком компании TERUMO CORPORATION.

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
DESCRIPTION

The Ryujin Plus OTW is a PTCA balloon dilatation catheter with a coaxial system. The tip of the catheter is equipped with a balloon inflatable to a specific diameter and length at recommended pressures. One or two radiopaque markers inside the balloon allow for confirmation of the balloon's position in the patient's vessel under high resolution fluoroscopy. On the shaft, there are two depth markers, one at about 90 cm (for brachial) and the other at 100 cm (for femoral) from the distal end of the catheter. At the proximal end, the catheter has a guide wire port, where a guide wire is to be inserted, and a balloon inflation port used for inflation/deflation of the balloon. The inner lumen is used as a passage of a guide wire, and the outer lumen is used for inflation/deflation of the balloon. The inner lumen cannot be used for angiography or pressure monitoring because of the small diameter. The maximum size guide wire that may be used is 0.014 inches (0.36mm) in diameter. The surface of catheter is partially coated with hydrophilic polymer coating which generates lubricity when wet.

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
ENGLISH

Please read all instructions prior to use.

INDICATIONS

The Ryuji Plus OTW ("dilatation catheter") is intended to be used for percutaneous transluminal coronary angioplasty (PTCA) for the purpose of improving myocardial blood flow in the localized stenotic lesion of the coronary arteries.

CAUTIONS FOR USE

1. Contraindications (patients/conditions in which PTCA must be avoided)

- Lesions in the left main trunk for which no compensation of blood flow by bypass or collateral circulation is available. Failure to observe this warning could result in acute coronary occlusion.
- Patients who had previous coronary artery spasm due to the possibility of acute coronary occlusion.
- Pregnancy or suspected pregnancy. X-ray exposure could damage fetus.

2. Relative contraindications (patients/conditions in which PTCA may carry a higher than usual risk, and should only be attempted if the procedure's benefit outweighs the risk)

- Patients in whom coronary bypass surgery is not applicable. In case of acute-phase ischemic complications; emergency CABG is required.

3. Important safety instructions

WARNINGS

- **Advance the dilatation catheter carefully within the artery and, if any resistance is felt, stop manipulating the dilatation catheter and determine the cause under high resolution fluoroscopy. Continuing to advance the dilatation catheter may result in damage to the vessel and/or separation or laceration of the dilatation catheter. This may necessitate recovery of fragments of the dilatation catheter.**
- **Within the stent strut, advancement or removal or inflation of the dilatation catheter should be done carefully under high resolution fluoroscopy. Failure to take care could result in vascular injury or damage/breakage of the catheter due to stent abrasion, and the balloon bursting below the rated burst pressure.**

CAUTIONS

- Use an inflation/deflation device equipped with an accurate manometer. The balloon may rupture if over-inflation occurs due to inaccurate determination of balloon pressure.
- Administer appropriate anticoagulant and coronary vasodilator to the patient during the PTCA procedure. Carry out appropriate anticoagulant therapy under direction of the physician in charge after completion of the PTCA procedure.
- Do not use agents containing organic solvents or oleaginous contrast media. Contact with these agents may lead to damage of the dilatation catheter and/or rupture of the balloon.
- Always handle catheters with care, and avoid kinks. Do not use if kinked. A kink could damage or break the catheter.
- Operate the catheter with utmost care while performing the kissing balloon technique or parallel wire technique, to avoid entanglement with the accompanying device. If resistance is encountered, remove the catheter and the accompanying device together.
- Choose the appropriate balloon size from the diagnostic site and anatomic aspect.
- After withdrawing the catheter, soak it in physiological saline solution to remove the blood on the surface of the catheter. If the blood is difficult to remove, wipe the catheter once with gauze soaked with physiological saline solution. Before reinserting the same balloon catheter, inspect the entire catheter that neither the lubricity of the surface has decreased, nor any foreign substances are on the catheter. To remove the blood from the guide wire lumen, flush the catheter with heparinized physiological saline solution.
- Refer the SPECIFICATIONS regarding the relation between the diameter of the balloon and the rated burst pressure.

4. Complications

Possible complications of PTCA include, but are not limited to, the following:

- Coronary artery dissection
- Acute myocardial infarction
- Ventricular fibrillation
- Ischaemia due to long-term dilatation
- Arterial perforation
- Arteriovenous fistula
- Palpitation
- Hypotension
- Coronary artery injury
- Unstable angina pectoris
- Intravessel thrombosis
- Arterial rupture
- Infection and pain at the insertion site
- Bradycardia
- Severe arrhythmia
- Cerebral infarction
- Total occlusion of coronary artery or bypass graft.
- Restenosis of coronary artery
- Ischaemia due to spasm
- Distal embolization
- Haematoma
- Nausea and vomiting
- Haemorrhage
- Death

If a complication requiring emergency coronary artery bypass surgery has occurred as the result of PTCA, the mortality of patients who had previous bypass surgery will be higher than that of patients who did not undergo bypass surgery. The long term complications of PTCA remain to be defined.

5. Precautions for application

- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Sterile and non-pyrogenic in an unopened and undamaged unit package. Do not use if the unit package or the product have been damaged or soiled.
- The product should be used immediately after opening the package and be disposed of safely and properly after use.
- The dilatation catheter should be used by a physician who is familiar with, and well trained in, PTCA techniques. Failure to observe this warning could result in failure due to improper operation.
- PTCA procedures should be carried out only at institutions where emergency coronary artery surgery can be performed preparing for severe complications. It is desirable that a cardiovascular surgery team is standing ready during PTCA procedures. Otherwise failure or complications cannot be promptly responded to.
- Do not soak the catheter in sterilizing alcohol or drug solutions containing organic solvents, or wipe the catheter with drugs. Failure to observe this warning could damage or break the catheter or cause loss of lubricity.
- Any advancement after introduction of the dilatation catheter into the vessel should be done under high resolution fluoroscopy.
- The entire operation should be carried out aseptically.

DIRECTIONS FOR USE

CAUTIONS

- Before using the catheter, read the instructions for the accompanying pharmaceutical products and medical devices.
- Before use, confirm that all apparatus including the dilatation catheter is functioning properly. Confirm whether the dilatation catheter is damaged or not, and that the design of the balloon meets the criteria of the procedure and the technique to be used.

1. Preparation of the dilatation catheter

1-1 Carefully remove the catheter from the holder tube.

CAUTIONS

- Remove the catheter in a straight line from the end of the holder tube without bending the catheter.
- Do not remove the catheter by force if resistance is encountered. Forcible removal could cause the balloon not to inflate or deflate.

1-2 After wetting the balloon protection sheath with physiological saline solution, carefully remove the balloon protection sheath and stylet not to damage the balloon part.

CAUTION

- Do not remove the sheath if resistance is felt. Using force may result in damage to the balloon.

1-3 Draw 3 mL of an appropriate contrast media (for example: a 1:1 mixture of contrast media and physiological saline solution) into a 20 mL syringe.

WARNING

Do not use air, gases or liquids other than contrast media to inflate the balloon. In case of leakage from the balloon, such fluids could have serious adverse effects on the patient's health.

1-4 Connect the syringe containing contrast media to the balloon inflation port of the dilatation catheter.

1-5 Holding the syringe with its tip down, aspirate air for 15 seconds.

1-6 While holding the syringe with its tip pointing downwards, inject the contrast media slowly.

1-7 Repeat steps 1-5 and 1-6 several times until contrast media fills the balloon completely.

2. Flushing and Insertion of guide wire

- 2-1 Introduce heparinized physiological saline solution into the catheter through the guide wire port to remove air bubbles.
- 2-2 Visually confirm that the balloon is fully deflated.
- 2-3 Insert a guide wire (not more than 0.014" (0.36 mm) in diameter) through the guide wire port, and advance the guide wire until it reaches the distal end of the dilatation catheter. This manipulation should be carried out with the dilatation catheter held as straight as possible.

WARNING

Before inserting the guide wire through the guide wire port, soak it in physiological saline solution to remove the blood or contrast media on the surface of the guide wire. If the blood or contrast media is difficult to remove, wipe the guide wire once with gauze soaked with physiological saline solution. Inspect the entire guide wire that neither the lubricity of the surface has decreased, nor any foreign substances are on the wire.

Moving the catheter over such residues which are adherent to the guide wire or over a half-wetted wire, may result in the separation or laceration of the dilatation catheter. This may necessitate the recovery of the catheter fragments.

CAUTIONS

- Take care not to kink the distal tip of the dilatation catheter and advance slowly/carefully when inserting the proximal end of the guide wire into the distal tip of the dilatation catheter.
- Carefully insert the guide wire, and use caution to prevent the distal end from damaging the guide wire lumen.
- If the catheter is placed into a bowl of physiological saline solution, carefully coil the shaft to avoid accidental contamination.
- When proximal shaft of the catheter is bent or kinked, do not use the catheter. It may cause the separation of shaft.

3. Connection of an inflation/deflation device equipped with a manometer to the dilatation catheter

- 3-1 Fill an inflation/deflation device equipped with a manometer with the contrast media and expel air from the device.
- 3-2 Attach the inflation/deflation device firmly to the balloon inflation port on the dilatation catheter. To ensure that no air enters the system, the inflation device must be filled adequately with contrast media.

4. Insertion of the dilatation catheter

- 4-1 Insert an introducer sheath into the blood vessel as described in the manufacturer's instruction manual.
- 4-2 Select a guiding catheter that conforms to the label indication, and suitable to the position of the lesion and the patient's anatomy. Flush the guiding catheter with heparinized physiological saline solution before use.

CAUTION

Administer appropriate anticoagulation therapy to the patient before insertion of the guiding catheter.

- 4-3 Position the guiding catheter at the ostium of the desired coronary artery using accepted protocol. Confirm the position of the guiding catheter under high resolution fluoroscopy. After the catheter position is confirmed, administer an appropriate dose of a vasodilator.
- 4-4 Loosen the hemostatic valve of the Y connector attached to the guiding catheter and insert the dilatation catheter with care to avoid kinking.

CAUTION

Make sure the hemostatic valve of the Y connector has been loosened. If tight, the valve will not allow smooth passage of the balloon.

- 4-5 Under high resolution fluoroscopy, advance the dilatation catheter until it reaches a point 2-3 cm proximal to the distal end of the guiding catheter. The depth marker on the shaft will help to confirm how far the catheter has been advanced.
- 4-6 Advance the guide wire into the desired coronary artery under high resolution fluoroscopy. Carry out angiography through the guiding catheter to confirm that the guide wire has crossed the stenotic lesion.

CAUTION

Confirm that the guide wire is correctly inserted into a target vessel by performing contrast radiography from various angles.

- 4-7 Advance the dilatation catheter over the guide wire until the balloon reaches the stenotic lesion.

WARNING

If any resistance is felt, do not advance the guide wire or the dilatation catheter by force. Before proceeding, determine the cause under high resolution fluoroscopy. Advancement by force may result in damage to the vessel and/or laceration or separation of the guide wire or the dilatation catheter. This may necessitate recovery of fragments.

- 4-8 Advance the dilatation catheter to position the balloon at the site of the stenotic lesion with the help of the radiopaque marker, and inflate it at a low pressure of 1-2 atm (101-203 kPa) after tightening the hemostatic valve of the Y connector. Confirm that the balloon is positioned in the centre of the stenotic lesion by checking the resultant unevenness (dumbbell effect).

CAUTION

Do not tighten the hemostatic valve of the Y connector excessively as this may affect the inflation/deflation time and/or kink the catheter shaft.

5. Balloon inflation

- 5-1 Inflate the balloon with appropriate pressure for an appropriate time with the inflation/ deflation device equipped with a manometer; then deflate the balloon.

WARNINGS

- Carefully inflate the balloon under the guidance of high-resolution fluoroscopy, and ensure that compression inflates the balloon. If the balloon does not inflate, do not apply excessive pressure, as this could prevent the balloon from deflating.
- The inflation pressure of the balloon should not exceed the RBP. Pressurization above the RBP may result in rupture of the balloon. The RBP is based on results of in vitro testing. At least 99.9 % of the balloons (with 95 % confidence) will not burst at or below their RBP.
- If a balloon rupture should occur due to pressurization above the RBP, the balloon or its fragments might be released into the vessel, retrieval of which may be needed.
- The short or long term effect of pressurization above the nominal pressure on the coronary arteries is still under investigation.
- Do not inflate the balloon beyond the diameter of the coronary artery proximal or distal to the stenotic lesion.
- Balloon inflation to expand a stent, inside a stent or calcified lesions is combined with a possibility of the balloon rupturing before the RBP is exceeded. Inflate the balloon with due caution.

CAUTION

The balloon may slip out of the lesion when inflated because of the hydrophilic coating. Inflate the balloon carefully under the guidance of high-resolution fluoroscopy so that the balloon does not change position in the lesion.

- 5-2 Pull back the dilatation catheter to withdraw the completely deflated balloon into the guiding catheter after inflation of the balloon, and carry out coronary angiography through the guiding catheter to evaluate the improvement of the stenosis.

CAUTIONS

- Do not move or remove the dilatation catheter before the balloon is deflated completely. Removal of the dilatation catheter should be done after loosening the hemostatic valve of the Y connector.
- While the guide wire is in the vessel, remove the catheter in a straight line along the guide wire.
- Do not remove the catheter if it is bent at the Y connector port. If removal of a bent catheter is attempted, excessive pressure is applied to the area near the guide wire port, possibly damaging or breaking the catheter.

- 5-3 If the improvement of the stenosis is not sufficient, increase the inflation pressure of the balloon gradually to the RBP, or pressurize it repeatedly until no further improvement can be attained. Usually, repeated inflation of the balloon brings about sufficient improvement of the stenosis, which can be confirmed by coronary angiography.

6. Exchange of the dilatation catheter

- 6-1 Loosen the hemostatic valve of the Y connector.

- 6-2 Keep the guide wire in position and remove the dilatation catheter.

- 6-3 Remove the dilatation catheter while maintaining the position of the guide wire in the lesion. Wipe the guide wire surface to avoid problems when inserted in the next catheter.

WARNING

When inserting or exchanging the dilatation catheter, wipe the guide wire once with gauze soaked with physiological saline solution. Inspect the entire guide wire that neither the lubricity of the surface has decreased, nor any foreign substances are on the wire. Moving the catheter over such residues adhered to the guide wire or over a half-wetted wire, may result in the separation or laceration of dilatation catheter. This may necessitate the recovery of the catheter fragments.

CAUTION

Monitor the guide wire position under high resolution fluoroscopy during the exchange.

- 6-4 Insert the next catheter over the proximal end of the guide wire as previously described while maintaining the guide wire position.

CAUTION

Read the manufacturer's instructions when catheters other than the "RyujinPlus OTW" are used.

- 6-5 Follow the directions for use labelled "Insertion of the dilatation catheter" after 4-7 and inflate/exchange dilatation catheters.

7. Removal of the dilatation catheter

After completion of dilatation, deflate the balloon completely and remove the dilatation catheter and guide wire after loosening the hemostatic valve. It is recommended to keep the guide wire in position for a while after the procedure, preparing for any possible unexpected incidents.

8. Temporary storage of the balloon catheter during an operation

Directions for use of balloon protection sheath

CAUTION

Do not reuse the first balloon protection sheath which is on the catheter, after removing it. Failure to observe this warning could cause the balloon to not inflate due to balloon deformation and damage of the shaft.

8-1 Remove the second balloon protection sheath from the compliance sheet.

8-2 Insert the stylet inside the second balloon protection sheath.

8-3 Insert the stylet and the second balloon protection sheath from the tip of the catheter, and carefully cover the balloon with the stylet and the second balloon protection sheath.

CAUTION

Do not forcibly insert the second balloon protection sheath. Lightly twist the balloon and carefully reinsert it.

8-4 When using the catheter, after wetting the second balloon protection sheath at the tip of the catheter in a physiological saline solution, carefully remove the second balloon protection sheath and stylet without damaging the balloon.

CAUTION

Do not remove the second balloon protection sheath by force if resistance is encountered. Forcible removal could cause the balloon to not inflate or deflate.

PRECAUTIONS FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperature, or high humidity during storage.

REFERENCES

Recommended particulars of the techniques, including the positioning of the guiding catheter and the dilatation catheter as well as the inflation cycle, may differ largely with patients, depending on their anatomy and the condition of the lesion. Refer to the following literature for such particulars:

1. Bentivoglio LG, et al : Percutaneous Transluminal Coronary Angioplasty (PTCA) in Patients with Relative Contraindications : Results of the National Heart, Lung, and Blood Institute PTCA Registry. *Am J Cardiol* 1984;53: 82C-88C
2. Corbelli J, et al : Percutaneous Transluminal Coronary Angioplasty After Previous Coronary Artery Bypass Surgery. *Am J Cardiol* 1985;56:398-403
3. Cowley MJ, et al : Efficacy of percutaneous transluminal coronary angioplasty : Technique, patient selection, salutary results, limitations and complications. *AM HEART J* 101:272, 1981
4. Dorros G, et al : Percutaneous Transluminal Coronary Angioplasty : Report of Complications from the National Heart, Lung, and Blood Institute PTCA Registry. *Circulation* 67, No.4, 1983
5. Gruentzig AR : Percutaneous transluminal coronary angioplasty : six years' experience. *AM HEART J* 107:818, 1984
6. Kelsey SF, et al : Effect of Investigator Experience on Percutaneous Transluminal Coronary Angioplasty. *Am J Cardiol* 1984;53:56C-64C
7. Kent KM, et al : Long-Term Efficacy of Percutaneous Transluminal Coronary Angioplasty (PTCA) : Report from the National Heart, Lung, and Blood Institute PTCA Registry. *Am J Cardiol* 1984;53:27C-31C
8. Meier B, et al : Higher balloon dilatation pressure in coronary angioplasty. *AM HEART J* 107:619, 1984
9. Vandormael MG, et al : Immediate and Short-Term Benefit of Multilesion Coronary Angioplasty : Influence of Degree of Revascularization. *J Am Coll Cardiol* 1985;6:983-91
10. Ryan TJ, et al : Guidelines for Percutaneous Transluminal Coronary Angioplasty. *JACC* Vol.22, No.7 December 1993:2033-54

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- Bentivoglio LG, et al : Percutaneous Transluminal Coronary Angioplasty (PTCA) in Patients with Relative Contraindications : Results of the National Heart, Lung, and Blood Institute PTCA Registry. Am J Cardiol 1984;53: 82C-88C
- Corbelli J, et al : Percutaneous Transluminal Coronary Angioplasty After Previous Coronary Artery Bypass Surgery. Am J Cardiol 1985;56:398-403
- Cowley MJ, et al : Efficacy of percutaneous transluminal coronary angioplasty : Technique, patient selection, salutary results, limitations and complications. AM HEART J 101:272, 1981
- Dorros G, et al : Percutaneous Transluminal Coronary Angioplasty : Report of Complications from the National Heart, Lung, and Blood Institute PTCA Registry. Circulation 67, No.4, 1983
- Gruentzig AR : Percutaneous transluminal coronary angioplasty : six years' experience. AM HEART J 107:818, 1984
- Kelsey SF, et al : Effect of Investigator Experience on Percutaneous Transluminal Coronary Angioplasty. Am J Cardiol 1984;53:56C-64C
- Kent KM, et al : Long-Term Efficacy of Percutaneous Transluminal Coronary Angioplasty (PTCA) : Report from the National Heart, Lung, and Blood Institute PTCA Registry. Am J Cardiol 1984;53:27C-31C
- Meier B, et al : Higher balloon dilatation pressure in coronary angioplasty. AM HEART J 107:619, 1984
- Vandormael MG, et al : Immediate and Short-Term Benefit of Multilesion Coronary Angioplasty : Influence of Degree of Revascularization. J Am Coll Cardiol 1985;6:983-91
- Ryan TJ, et al : Guidelines for Percutaneous Transluminal Coronary Angioplasty. JACC Vol.22, No.7 December 1993:2033-54

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The Ruyjin Plus ОТВ ("dilatation catheter") is intended to be used for percutaneous transluminal coronary angioplasty (PTCA) for the purpose of improving myocardial blood flow in the localized stenotic lesion of the coronary arteries.

1. Bentivoglio LG, et al : Percutaneous Transluminal Coronary Angioplasty (PTCA) in Patients with Relative Contraindications : Results of the National Heart, Lung, and Blood Institute PTCA Registry. Am J Cardiol 1984;53: 82C-88C
2. Corbelli J, et al : Percutaneous Transluminal Coronary Angioplasty After Previous Coronary Artery Bypass Surgery. Am J Cardiol 1985;56:398-403
3. Cowley MJ, et al : Efficacy of percutaneous transluminal coronary angioplasty : Technique, patient selection, salutary results, limitations and complications. AM HEART J 101:272, 1981
4. Dorros G, et al : Percutaneous Transluminal Coronary Angioplasty : Report of Complications from the National Heart, Lung, and Blood Institute PTCA Registry. Circulation 67, No.4, 1983
5. Gruentzig AR : Percutaneous transluminal coronary angioplasty : six years' experience. AM HEART J 107:818, 1984
6. Kelsey SF, et al : Effect of Investigator Experience on Percutaneous Transluminal Coronary Angioplasty. Am J Cardiol 1984;53:56C-64C
7. Kent KM, et al : Long-Term Efficacy of Percutaneous Transluminal Coronary Angioplasty (PTCA) : Report from the National Heart, Lung, and Blood Institute PTCA Registry. Am J Cardiol 1984;53:27C-31C
8. Meier B, et al : Higher balloon dilatation pressure in coronary angioplasty. AM HEART J 107:619, 1984
9. Vandormael MG, et al : Immediate and Short-Term Benefit of Multilesion Coronary Angioplasty : Influence of Degree of Revascularization. J Am Coll Cardiol 1985;6:983-91
10. Ryan TJ, et al : Guidelines for Percutaneous Transluminal Coronary Angioplasty. JACC Vol.22, No.7 December 1993:2033-54

TERUMO CORPORATION 44-1, 2-CHOME, HATAGAYA, SHIBUYA-KU, TOKYO 151-0072, JAPAN

MADE IN JAPAN

TERUMO EUROPE N.V. INTERLEUVENLAAN 40, 3001 LEUVEN, BELGIUM www.terumo-europe.com

Terumo Australia Pty Ltd Level 4 Building B 11 Talavera Rd Macquarie Park NSW 2113 Australia

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Ryujin is a trademark of TERUMO CORPORATION.

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© TERUMO CORPORATION 2014-04