2014-04

CE 0086

MÔ TẢ SẢN PHẨM

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ENGLISH	GREEK
FRENCH	RUSSIAN
GERMAN	CZECH
SPANISH	SLOVAK
PORTUGUESE	TURKISH
ITALIAN	LATVIAN
DUTCH	SLOVENIAN
SWEDISH	SERBIAN
DANISH	ROMANIAN
NORWEGIAN	BULGARIAN
FINNISH	UKRAINIAN

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
DESCRIPTION

This dilation catheter is a rapid exchange PTCA balloon dilatation catheter. The tip of the catheter is equipped with a balloon inflatable to a specific diameter and length at recommended pressures. One or two radiopaque markers inside the balloon allow for confirmation of the balloon's position in the patient's vessel under high resolution fluoroscopy. On the shaft, there are two depth markers, one at about 90 cm (for brachial) and the other at 100 cm (for femoral) from the distal end of the catheter. The distal part of the catheter has a coaxial construction. The inner lumen is used to facilitate passage of a guide wire which must not exceed 0.014" (0.36 mm) in diameter, while the other lumen is used for inflation/deflation of the balloon which is attached to the hub fixed to the proximal end of the catheter. The surface of catheter is partially coated with hydrophilic polymer coating which generates lubricity when wet. The flushing needle with a luer port is enclosed with the dilatation catheter for priming the inner lumen at the distal end of the catheter. The CATHETERCLIP and the balloon protection sheath are also provided with the dilatation catheter to aid in handling of the catheter.

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
ENGLISH

INDICATIONS

The Tazuna ("dilatation catheter") is intended to be used for percutaneous transluminal coronary angioplasty (PTCA) for the purpose of improving myocardial blood flow in the localized stenotic lesion of the coronary arteries.

CAUTIONS FOR USE

- Contraindications (patients/conditions in which PTCA must be avoided)**
 - Lesions in the left main trunk for which no compensation of blood flow by bypass or collateral circulation is available. Failure to observe this warning could result in acute coronary occlusion.
 - Patients who had previous coronary artery spasm.
 - Pregnancy or suspected pregnancy. X-ray exposure could damage fetus.
- Relative contraindications (patients/conditions in which PTCA may carry a higher than usual risk, and should only be attempted if the procedure's benefit outweighs the risk)**
 - Patients in whom coronary bypass surgery is not applicable. Emergency CABG is required for acute-phase ischemic complications.
- Important safety instructions**

WARNINGS

- Advance the dilatation catheter carefully within the artery and, if any resistance is felt, stop manipulating the dilatation catheter and determine the cause under high resolution fluoroscopy. Continuing to advance the dilatation catheter may result in damage to the vessel and/or separation or laceration of the dilatation catheter. This may necessitate recovery of fragments of the dilatation catheter.
- Within the stent strut, advancement or removal or inflation of the dilatation catheter should be done carefully under high resolution fluoroscopy. Failure to take care could result in vascular injury or damage/breakage of the catheter due to stent abrasion, and the balloon bursting below the rated burst pressure.

CAUTIONS

- Use an inflation/deflation device equipped with an accurate manometer. The balloon may rupture if over-inflation occurs due to inaccurate determination of balloon pressure.
- Administer appropriate anticoagulant and coronary vasodilator to the patient during the PTCA procedure. Carry out appropriate anticoagulant therapy under direction of the physician in charge after completion of the PTCA procedure.
- Do not use agents containing organic solvents or oleaginous contrast media. Contact with these agents may lead to damage of the dilatation catheter and/or rupture of the balloon.
- Always handle catheters with care, and avoid kinks. Do not use if kinked. A kink could damage or break the catheter.
- Operate the catheter with utmost care while performing the kissing balloon technique or parallel wire technique, to avoid entanglement with the accompanying device. If resistance is encountered, remove the catheter and the accompanying device together.
- Choose the appropriate balloon size from the diagnostic site and anatomic aspect.
- After withdrawing the catheter, it should be placed into a bowl of physiological saline solution. Prior to reinsert the catheter, flush the catheter with heparinized physiological saline solution to remove the blood from the guide wire lumen. Do not reinsert the same catheter if the blood is on the surface.

- Refer to the SPECIFICATIONS regarding the relation between the diameter of the balloon and the rated burst pressure.
- After removing, do not reuse, the balloon protection sheath mounted on the balloon catheter. Failure to observe this warning could cause the balloon not to inflate or deflate due to balloon deformation and damage of the shaft.

4. Complications

Possible complications of PTCA include, but are not limited to, the following:

- Coronary artery dissection
- Acute myocardial infarction
- Ventricular fibrillation
- Ischaemia due to long-term dilatation
- Arterial perforation
- Arteriovenous fistula
- Palpitation
- Hypotension
- Coronary artery injury
- Unstable angina pectoris
- Intravessel thrombosis
- Arterial rupture
- Infection and pain at the insertion site
- Bradycardia
- Severe arrhythmia
- Cerebral infarction
- Total occlusion of coronary artery or bypass graft
- Restenosis of coronary artery
- Ischaemia due to spasm
- Distal embolization
- Haematoma
- Nausea and vomiting
- Haemorrhage
- Death

If a complication requiring emergency coronary artery bypass surgery has occurred as the result of PTCA, the mortality of patients who had previous bypass surgery will be higher than that of patients who did not undergo bypass surgery. The long term complications of PTCA remain to be defined.

5. Precautions for application

- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Sterile and non-pyrogenic in an unopened and undamaged unit package. Do not use if the unit package or the product have been damaged or soiled.
- The product should be used immediately after opening the package and be disposed of safely and properly after use.
- The dilatation catheter should be used by a physician who is familiar with, and well trained in, PTCA techniques.
- PTCA procedures should be carried out only at institutions where emergency coronary artery surgery can be performed preparing for severe complications. It is desirable that a cardiovascular surgery team is standing ready during PTCA procedures.
- Do not soak the catheter in sterilizing alcohol or drug solutions containing organic solvents, or wipe the catheter with drugs. Failure to observe this warning could damage or break the catheter or cause loss of lubricity.
- Any advancement after introduction of the dilatation catheter into the vessel should be done under high resolution fluoroscopy.
- The entire operation should be carried out aseptically.

DIRECTIONS FOR USE

CAUTIONS

- Before using the catheter, read the instructions for the accompanying pharmaceutical products and medical devices.
- Before use, confirm that all apparatus including the dilatation catheter is functioning properly. Confirm whether the dilatation catheter is damaged or not, and that the design of the balloon meets the criteria of the procedure and the technique to be used.

1. Preparation of the dilatation catheter

1-1 Carefully remove the catheter from the holder tube.

CAUTIONS

- Remove the catheter in a straight line from the end of the holder tube without bending the catheter. If the catheter is not removed in a straight line, excessive pressure is applied to the area near the guide-wire port, possibly damaging or breaking the catheter.

- Do not remove the catheter by force if resistance is encountered. Forced removal could cause the balloon not to inflate or deflate.
- 1-2 After wetting the balloon protection sheath with physiological saline solution, carefully remove the balloon protection sheath and stylet not to damage the balloon part. If necessary, dilatation catheter may be coiled and secured using the CATHETERCLIP in accordance with "Directions for use of CATHETERCLIP and balloon protection sheath".

CAUTION Do not remove the sheath if resistance is felt. Using force may result in damage to the balloon.

- 1-3 Draw 3 mL of an appropriate contrast medium (for example: a 1:1 mixture of contrast media and physiological saline solution) into a 20 mL syringe.

WARNING Do not use air, gases or liquids other than contrast media to inflate the balloon. In case of leakage from the balloon, such fluids could have serious adverse effects on the patient's health.

- 1-4 Connect the syringe containing contrast medium to the catheter hub (hereafter referred to as "balloon inflation port").

- 1-5 Holding the syringe with its tip down, aspirate air for 20 ~ 30 seconds.

- 1-6 While holding the syringe with its tip pointing downwards, inject the contrast medium slowly.

- 1-7 Repeat steps 1-5 and 1-6 several times until contrast medium fills the balloon completely.

2. Flushing and Insertion of guide wire

- 2-1 Insert the enclosed flushing needle into the distal tip of the dilatation catheter. Flush with heparinized physiological saline solution in order to remove air bubbles.

CAUTION Take care not to damage the dilatation catheter when inserting the enclosed flushing needle.

- 2-2 Visually confirm that the balloon is fully deflated.

- 2-3 Insert the proximal end of a guide wire (not more than 0.014" (0.36 mm) in diameter) into the distal tip of the catheter. Advance the wire through the guide wire lumen until it protrudes from the guide wire port. Grasp the proximal end of the guide wire and pull it back through the guide wire lumen until guide wire and catheter tip are appropriately positioned.

WARNING Before inserting the dilatation catheter, thoroughly wipe the guide wire with gauze soaked with physiological saline solution to remove any residue of blood or contrast media. Moving the catheter over such residues which are adherent to the guide wire or over a half-wetted wire, may result in the separation or laceration of the dilatation catheter. This may necessitate the recovery of the catheter fragments.

CAUTIONS

- Take care not to kink the distal tip of the dilatation catheter and advance slowly/carefully when inserting the proximal end of the guide wire into the distal tip of the dilatation catheter.

- Carefully insert the guide wire to prevent the proximal end from damaging the guide-wire lumen.
- If the catheter is placed into a bowl of physiological saline solution, carefully coil the shaft to avoid accidental contamination.
- When the proximal shaft of the catheter is bent or kinked, do not use the catheter. It may cause separation of the catheter shaft.

3. Connection of an inflation/deflation device equipped with a manometer to the dilatation catheter

- 3-1 Fill an inflation/deflation device equipped with a manometer with the contrast medium and expel air from the device.

- 3-2 Attach the inflation/deflation device firmly to the balloon inflation port on the dilatation catheter. To ensure that no air enters the system, the inflation device must be filled adequately with contrast medium.

4. Insertion of the dilatation catheter

- 4-1 Insert an introducer sheath into the blood vessel as described in the manufacturer's instruction manual.

- 4-2 Select a guiding catheter that conforms to the label indication, and suitable to the position of the lesion and the patient's anatomy. Flush the guiding catheter with heparinized physiological saline solution before use.

CAUTION Administer appropriate anticoagulation therapy to the patient before insertion of the guiding catheter.

- 4-3 Position the guiding catheter at the ostium of the desired coronary artery using accepted protocol. Confirm the position of the guiding catheter under high resolution fluoroscopy. After the catheter position is confirmed, administer an appropriate dose of a vasodilator.

- 4-4 Insert the dilatation catheter through the hemostatic valve of the Y connector attached to the guiding catheter.

CAUTION Make sure the hemostatic valve of the Y connector has been loosened. If tight, the valve will not allow smooth passage of the balloon.

- 4-5 Under high resolution fluoroscopy, advance the dilatation catheter until it reaches a point 2-3 cm proximal to the distal end of the guiding catheter. The depth marker on the shaft will help confirm how far the catheter has been advanced.

- 4-6 Advance the guide wire into the desired coronary artery under high resolution fluoroscopy. Carry out angiography through the guiding catheter to confirm that the guide wire has crossed the stenotic lesion.

CAUTION Confirm that the guide wire is correctly inserted into a target vessel by performing contrast radiography from various angles.

- 4-7 Advance the dilatation catheter over the guide wire until the balloon reaches the stenotic lesion.

WARNING If any resistance is felt, do not advance the guide wire or the dilatation catheter by force. Before proceeding, determine the cause under high resolution fluoroscopy. Advancement by force may result in damage to the vessel and/or laceration or separation of the guide wire or the dilatation catheter. This may necessitate recovery of fragments.

- 4-8 Advance the dilatation catheter to position the balloon at the site of the stenotic lesion with the help of the radiopaque marker, and inflate it at a low pressure of 1-2 atm (101-203 kPa) after tightening the hemostatic valve of the Y connector. Confirm that the balloon is positioned in the centre of the stenotic lesion by checking the resultant unevenness (dumbbell effect).

CAUTION Do not tighten the hemostatic valve of the Y connector excessively as this may affect the inflation/deflation time and/or kink the catheter shaft.

5. Balloon inflation

- 5-1 Inflate the balloon with appropriate pressure for an appropriate time with the inflation/ deflation device equipped with a manometer; then deflate the balloon.

WARNINGS

- Carefully inflate the balloon under the guidance of high-resolution fluoroscopy. If the balloon does not inflate, do not apply excessive pressure, as this could prevent the balloon from deflating.

- The inflation pressure of the balloon should not exceed the RBP. Pressurization above the RBP may result in rupture of the balloon. The RBP is based on results of in vitro testing. At least 99.9 % of the balloons (with 95 % confidence) will not burst at or below their RBP.
- If a balloon rupture should occur due to pressurization above the RBP, the balloon or its fragments might be released into the vessel, retrieval of which may be needed.
- The short or long term effect of pressurization above the nominal pressure on the coronary arteries is still under investigation.
- Do not inflate the balloon beyond the diameter of the coronary artery proximal or distal to the stenotic lesion.
- Balloon inflation to expand a stent, inside a stent or calcified lesions is combined with a possibility of the balloon rupturing before the RBP is exceeded. Inflate the balloon with due caution.

CAUTION The balloon may slip out of the lesion when inflated because of the hydrophilic coating. Inflate the balloon carefully under the guidance of high-resolution fluoroscopy so that the balloon does not change position in the lesion.

- 5-2 Pull back the dilatation catheter to withdraw the completely deflated balloon into the guiding catheter after inflation of the balloon, and carry out coronary angiography through the guiding catheter to evaluate the improvement of the stenosis.

- CAUTIONS**
- Do not move or remove the dilatation catheter before the balloon is deflated completely. Removal of the dilatation catheter should be done after loosening the hemostatic valve of the Y connector.
 - While the guide wire is in the vessel, remove the catheter in a straight line along the guide wire.
 - Do not remove the catheter if it is bent at the Y connector port. If removal of a bent catheter is attempted, excessive pressure is applied to the area near the guide-wire port, possibly damaging or breaking the catheter.

- 5-3 If the improvement of the stenosis is not sufficient, increase the inflation pressure of the balloon gradually to the RBP, or pressurize it repeatedly until no further improvement can be attained. Usually, repeated inflation of the balloon brings about sufficient improvement of the stenosis, which can be confirmed by coronary angiography.

6. Exchange of the dilatation catheter

- 6-1 Loosen the hemostatic valve of the Y connector.
- 6-2 Grasp the guide wire and the hemostatic valve in one hand and the catheter in the other.
- 6-3 Remove the dilatation catheter while maintaining the position of the guide wire in the lesion. Wipe the guide wire surface to avoid problems when inserted in the next catheter.

WARNING When inserting or exchanging the dilatation catheter, wipe the guide wire once with gauze soaked with physiological saline solution. Inspect the entire guide wire that neither the lubricity of the surface has decreased, nor any foreign substances are on the wire. Moving the catheter over such residues adhered to the guide wire or over an half-wetted wire, may result in the separation or laceration of dilatation catheter. This may necessitate the recovery of the catheter fragments.

CAUTION Monitor the guide wire position under high resolution fluoroscopy during the exchange.

- 6-4 Insert the next catheter over the proximal end of the guide wire as previously described while maintaining the guide wire position.

CAUTION Read the manufacturer's instructions when catheters other than the Tazuna are used.

- 6-5 Follow the directions for use labelled "Insertion of the dilatation catheter" after 4-7 and inflate/exchange dilatation catheters.

7. Removal of the dilatation catheter

After completion of dilatation, deflate the balloon completely and remove the dilatation catheter and guide wire after loosening the hemostatic valve. It is recommended to keep the guide wire in position for a while after the procedure, preparing for any possible unexpected incidents. In order to dispose of the removed dilatation catheter safely and properly, use CATHETERCLIP in accordance with "Directions for use of CATHETERCLIP and balloon protection sheath".

8. Directions for use of CATHETERCLIP and balloon protection sheath


8-1 Directions for use of CATHETERCLIP

1. Remove the CATHETERCLIP from its holder.
2. Shape PTCA dilatation catheter into a single or double loop (Fig. 1).
CAUTION Proceed with care to prevent PTCA dilatation catheter kinking and collapse when forming the loops.
3. Secure the coiled PTCA dilatation catheter with the CATHETERCLIP following the steps below:
 - Hook the tip of the CATHETERCLIP onto the PTCA dilatation catheter (Fig. 2).
 - Fix the PTCA dilatation catheter onto the other end of the CATHETERCLIP (Fig. 3).

WARNING Secure the PTCA dilatation catheter with the CATHETERCLIP at the stiffer, proximal end. Do not use the CATHETERCLIP on the flexible, distal shaft or the PTCA guidewire port of the rapid exchange type PTCA dilatation catheters, it may damage the PTCA dilatation catheter.

4. When removing the CATHETERCLIP from the coiled PTCA dilatation catheter, repeat above steps in reverse (Fig. 3 to Fig. 2).

CAUTION Proceed with care to prevent shaft kinking and collapse while removing the CATHETERCLIP.

8-2 Directions for use of balloon protection sheath

CAUTION It is recommended to use the second balloon protection sheath when the catheter is placed into a bowl of physiological saline solution. Do not reuse, after removing, the balloon protection sheath mounted on the balloon catheter. Failure to observe this warning could cause the balloon to not inflate due to balloon deformation and damage of the shaft.

1. Remove the second balloon protection sheath from the compliance sheet.
2. Insert the stylet inside the balloon protection sheath.
3. Insert the stylet and balloon protection sheath from the tip of the catheter, and carefully cover the balloon with the stylet and balloon protection sheath.

CAUTION Do not forcibly insert the balloon into the balloon protection sheath. Lightly twist the balloon and carefully insert it.

4. When using the catheter, after wetting the balloon protection sheath at the tip of the catheter in a physiological saline solution, carefully remove the balloon protection sheath and stylet without damaging the balloon.

CAUTION Do not remove the balloon protection sheath by force if resistance is encountered. Forcible removal could cause the balloon to not inflate or deflate.

PRECAUTIONS FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperature, or high humidity during storage.

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© TERUMO CORPORATION 2014-04
- Tazuna is a trademark of TERUMO CORPORATION.
- Tazuna ir Terumo Corporation tirguzime