



Runthrough® NS

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CE 0086

<TÊN CỦA TỪNG BỘ PHẦN>

TW34M001-01

TERUMO®

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Terumo Australia Pty Ltd Level 4 Building B 11 Talavera Rd Macquarie Park NSW 2113 Australia

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

INSTRUCTIONS FOR USE

Read the following warnings, precautions, and directions for use carefully.

INDICATION

The RUNTHROUGH NS is indicated to be used to guide dilatation catheters to a lesion for percutaneous transluminal coronary angioplasty (PTCA) for the purpose of improving myocardial blood flow in a stenotic lesion in coronary vessels. Do not use for any other purpose.

DESCRIPTION

The RUNTHROUGH NS is a guide wire for use in PTCA with good conformability and trackability in the vessels. The distal part of the wire has a spring coil with hydrophilic polymer coating which generates high lubricity when wet. The distal portion of the spring coil has a radiopaque spring coil which provides good visibility and shapeability.

Some types of RUNTHROUGH NS have a modified proximal end that enables to attach the RUNTHROUGH NS Extension Wire. Please refer to the product label which explains availability to use Guide Wire Extension system.

CONTRAINDICATIONS

- Patients who had previous coronary artery spasm.
- Pregnant patients or patients who may be pregnant.
- Patients who, because of their condition, are contraindicated for surgical operation.
- Lesions that follow the bypass constructed less than one month ago or false aneurysm.
- Patients who are contraindicated for anti-platelet or anti-coagulation treatment.
- Patients who have a critical allergy to contrast media.

Carefully read the list of CONTRAINDICATIONS in the instructions accompanying the other interventional devices to be used with the RUNTHROUGH NS.

COMPLICATIONS

Possible complications associated with the use of the RUNTHROUGH NS include, but are not limited to, the following:

- Acute myocardial infarction
- Hypotension
- Coronary artery dissection, Arterial perforation, Arterial rupture, Coronary artery injury
- Infectious disease and puncture site complication
- A hemorrhagic complication
- Myocardium ischemia
- Coronary artery spasm
- Bradycardia, Palpitation
- Arteriovenous fistula
- Unstable angina pectoris
- Arrhythmia including the ventricular fibrillation
- Formation of femoral false aneurysm, false aneurysm
- Allergy for medicines
- Total occlusion of coronary artery
- Distal embolization
- Cerebral infarction

Carefully read the COMPLICATIONS in the instructions accompanying the other interventional devices to be used with the RUNTHROUGH NS.

WARNINGS AND CAUTIONS

WARNINGS

Prior to the procedure, carefully examine all equipment to verify its proper function and integrity. Failure to abide by the following warnings might result in damage to the vessel, abrasion of the hydrophilic coating, release of plastic fragments from the RUNTHROUGH NS and/or damage to or breakage/ separation of the wire, that may necessitate retrieval.

- Do not manipulate and/or withdraw the RUNTHROUGH NS through a metal entry needle, a metal dilator, or a metal guide wire inserter. Manipulation and/or withdrawal through a metal entry needle, a metal dilator, or a metal guide wire inserter may result in destruction and/or separation of the RUNTHROUGH NS.
- Do not use the RUNTHROUGH NS with devices which contain metal parts such as Atherectomy catheters. Do not manipulate the RUNTHROUGH NS through the stent struts. Such manipulations may cause the RUNTHROUGH NS hydrophilic polymer coating to peel off, and damage and/or sever the wire.
- The distal 1 cm of the RUNTHROUGH NS is shapeable. Follow the DIRECTIONS FOR USE 1-5) when shaping the RUNTHROUGH NS.
- Manipulate the RUNTHROUGH NS slowly and carefully in the vessel while confirming the behaviour and location of the wire's tip under high resolution fluoroscopy.
- If any resistance is felt or if the tip's behaviour and/or location seems improper, stop manipulating the RUNTHROUGH NS and/or the dilatation catheter and determine the cause by high resolution fluoroscopy. Failure to exercise proper caution may result in bending, kinking, separation of the guide wire's tip, damage to the dilatation catheter, or damage to the vessel.
- A retrieving device, such as a gripper or basket forceps, can only be used after the RUNTHROUGH NS has been removed from the patient's vessel. Using a retrieving device while the RUNTHROUGH NS is in the vessel may cause the RUNTHROUGH NS to break.
- All devices used with RUNTHROUGH NS should be sufficiently wetted prior to and during use with heparinized physiological saline solution. Improper priming of the system may result in damage to or the release of hydrophilic polymer fragments, and/or may prevent smooth movement of the wire.
- Do not use the RUNTHROUGH NS for intervention of or through internal mammary artery. It may result in breakage or separation of the RUNTHROUGH NS.
- Do not use the RUNTHROUGH NS as outer guide wire in case of stenting in the parallel wire technique.
- Do not bend the RUNTHROUGH NS repeatedly at the same point. It may result in deformation, breakage, or separation of the RUNTHROUGH NS.
- If catheter is used with a stopcock, do not operate this stopcock while the RUNTHROUGH NS is in the catheter. It may result in breakage or separation of the RUNTHROUGH NS.
- Manipulate the RUNTHROUGH NS carefully when crossing it through a deployed stent. Stent migration, stent deformation, or breakage or separation of the RUNTHROUGH NS may be caused if the RUNTHROUGH NS is caught by the deployed stent.
- Metal is plated onto the inner core of the RUNTHROUGH NS. Do not use if any disturbance from metal may be expected to MRI or other medical devices.
- If the RUNTHROUGH NS should stick out of the side hole of the catheter, do not remove the RUNTHROUGH NS by force, but remove the RUNTHROUGH NS as one unit with the other ancillary device.
- After replacement of catheter using the RUNTHROUGH NS Extension Wire, make sure that the extension wire is removed before next intervention. Separation of the joint may result from manipulation such as rotation while the RUNTHROUGH NS and the RUNTHROUGH NS Extension Wire are connected.
- The RUNTHROUGH NS Extension Wire is exclusively permitted for attachment to the RUNTHROUGH NS. If another extension wire than the RUNTHROUGH NS Extension Wire is attached (connected), this may result in damage to the RUNTHROUGH NS for compatibility reasons.

CAUTIONS

- The RUNTHROUGH NS should be used by, or under the order of a physician, who is well trained in manipulation and observation of guide wires under high resolution fluoroscopy.
- This product has been sterilized by ethylene oxide gas. For single use only. Do not reuse. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.
- Sterile in an unopened and undamaged unit package. Do not use if the unit package or the RUNTHROUGH NS is damaged or soiled.
- The RUNTHROUGH NS should be used immediately after opening the packaging and be disposed of safely and properly after use, following local regulations for medical waste management.
- The entire procedure should be carried out aseptically.
- When using a drug or a device concurrently with the RUNTHROUGH NS, the operator should have a full understanding of the properties/characteristics of the drug or device so as to avoid damage to the RUNTHROUGH NS.
- Do not reposition the torque device while grasping the RUNTHROUGH NS. It may result in breakage of the RUNTHROUGH NS.

DIRECTIONS FOR USE

1. Preparation of the RUNTHROUGH NS

- 1-1 Remove the RUNTHROUGH NS and the holder from the package.
- 1-2 Before removing the RUNTHROUGH NS, from its holder, inject heparinized physiological saline solution through the hub of the holder with a syringe.
- 1-3 After removing the wire stopper, remove the RUNTHROUGH NS from the holder and inspect the RUNTHROUGH NS prior to use to verify if it is lubricated. If the RUNTHROUGH NS cannot be easily removed from the holder, inject more heparinized physiological saline solution into the holder and try again.
- 1-4 Soak the RUNTHROUGH NS sufficiently in heparinized physiological saline solution.

WARNING

The proximal end of the RUNTHROUGH NS exposes from the holder. Make sure of the direction of the guide wire when inserting the guide wire into the dilatation catheter.

CAUTIONS

- The distal part of the RUNTHROUGH NS is not lubricious unless it is wet. Damage to or the release of hydrophilic polymer fragments may result if it is handled with its surface dry or not sufficiently wet.
- Prior to use, wet the surface of the RUNTHROUGH NS sufficiently. The distal part is not flexible unless it is handled with its surface sufficiently wet.
- The distal part of the RUNTHROUGH NS may be curved due to the hydrophilic polymer when taken out from the holder tube unless it is wet. It can be straightened to the original shape when its surface is wetted sufficiently.

- 1-5 If desired, gently shape the guide wire flexible tip according to standard practices (Fig. 1).

WARNING

Damage to or separation of the guide wire may result if it is handled with its surface dry, heated, held with fingernails, forceps, or some other tool, bent sharply, or bent repeatedly at the same point.

CAUTION

The shapeability of the RUNTHROUGH NS may be insufficient unless it is wetted sufficiently.

- 1-6 If not intended to insert the RUNTHROUGH NS into a dilatation catheter immediately, keep the RUNTHROUGH NS in the holder with heparinized physiological saline solution till the insertion.

CAUTION

When reinserting the RUNTHROUGH NS into the holder, take care not to damage the surface of the guide wire with the edge of the holder.

2. Insertion into a balloon dilatation catheter

CAUTION

Use an appropriate balloon dilatation catheter corresponding to the outer diameter (O.D.) of the RUNTHROUGH NS. Due to the variation of certain catheter tip inner diameters, abrasion of the hydrophilic coating may occur during manipulation. If any resistance is felt during introduction of the guide wire into the catheter, it is advisable to stop using such catheters.

- 2-1 Inspect and prepare the dilatation catheter according to the instructions for use of the balloon dilatation catheter.
- 2-2 If necessary, connect a Y connector to the guide wire port of the balloon dilatation catheter, and prime the balloon dilatation catheter with heparinized physiological saline solution.

WARNING

Prime the balloon dilatation catheter enough so that no air bubbles remain inside. Improper priming may prevent smooth movement of the guide wire.

- 2-3 Make sure that the surface of the RUNTHROUGH NS is wet enough. Insert the RUNTHROUGH NS into the dilatation catheter and advance the guide wire slowly until guide wire and catheter tip are appropriately positioned.

When a rapid exchange type balloon dilatation catheter is used, insert the RUNTHROUGH NS's proximal end into the distal tip of the catheter, advance the guide wire through the guide wire lumen until guide wire and catheter tip are appropriately positioned.

CAUTION

Insertion of the RUNTHROUGH NS into the balloon dilatation catheter while bent may result in damage to the catheter by the proximal end of the guide wire.

- 2-4 Tighten the haemostasis valve of the Y connector while verifying that of the RUNTHROUGH NS's proximal end protrudes and moves smoothly.

CAUTION

Do not manipulate the RUNTHROUGH NS when strongly locked by the haemostasis valve as breakage of the RUNTHROUGH NS may result.

3. System operation

- 3-1 Position the guiding catheter at the desired lesion using accepted protocol.

Insert the dilatation catheter through the haemostasis valve of the Y connector attached to the guiding catheter placed to the desired lesion, and advance the dilatation catheter up to 1-2 cm from the distal end of the guiding catheter under high resolution fluoroscopy.

WARNING

Carefully manipulate the guiding catheter and observe it for undesirable behaviour while the RUNTHROUGH NS and/or the dilatation catheter is in the vessel. Damage to or separation of the guide wire and/or damage to the vessel may result if the guiding catheter slips out of the position and exerts too large load on the guide wire.

- 3-2 When advancing the RUNTHROUGH NS across the lesion, manipulate it slowly and carefully while monitoring the behaviour and the position of the guide wire's tip under high resolution fluoroscopy to prevent the tip from being trapped in the lesion.
- 3-3 Confirm that the position of the RUNTHROUGH NS is firmly inserted into the desired vessel under high resolution fluoroscopy in multi-directions.

WARNINGS

- **When selecting the vessel in which the RUNTHROUGH NS is to be inserted or when advancing the guide wire through the stenosis, do not turn the proximal end of the guide wire three or more turns in succession in the same direction if the guide wire's tip is trapped. Separation of the guide wire may result.**
- **If any resistance to the RUNTHROUGH NS or the dilatation catheter is felt during manipulation or if the shape, behaviour or position of the guide wire's tip seems improper, e.g., in a case where the tip is trapped due to vessel spasm or some other cause or the tip is folded as shown in Fig. 2, stop manipulating the guide wire (and the catheter) and determine the cause carefully by fluoroscopy and take suitable remedial actions. Then remove the guide wire slowly, without turning, to exchange it for a new one. Continuing guide wire or catheter manipulation in the said situations may result in damage to the vessel, damage to or separation of the guide wire's tip, and/or damage to the dilatation catheter.**

- 3-4 Advance the dilatation catheter to the desired lesion over the RUNTHROUGH NS.

4. When exchanging the RUNTHROUGH NS

- 4-1 Loosen the haemostasis valve of the Y connector, and remove the RUNTHROUGH NS slowly while confirming the tip position under high resolution fluoroscopy. Carefully exchange the guide wire so that no air enters the dilatation catheter.
- 4-2 Gently insert new RUNTHROUGH NS into the dilatation catheter in accordance with the above-mentioned items 2-3 and 2-4. In case of using a Rapid Exchange type dilatation catheter, remove the dilatation catheter and exchange the guide wire.

WARNINGS

- **If resistance is felt during the RUNTHROUGH NS insertion, stop manipulation and exchange it for a new one. Manipulation against resistance may result in damage to the vessel due to breakage of the RUNTHROUGH NS or the dilatation catheter.**

- Before advancing the guide wire's tip throughout the dilatation catheter, make sure that the direction of the dilatation catheter's tip is parallel to the vessel. Damage to the vessel may result if the RUNTHROUGH NS is advanced with RUNTHROUGH NS's tip pointing toward the vessel wall or not parallel to the vessel.

5. For temporary storage during the procedure

- 5-1 Remove the blood residues on the RUNTHROUGH NS by rinsing the guide wire in a bowl of heparinized physiological saline solution. Wipe once with sterile gauze moistened with heparinized physiological saline solution if any blood residues still adhere to the guide wire.
- 5-2 Reinsert the RUNTHROUGH NS into the holder and fill with heparinized physiological saline solution.

CAUTIONS

- If the RUNTHROUGH NS moves less smoothly than before, exchange it for a new one.
- Do not wipe the RUNTHROUGH NS with alcohol, or gluconic acid chlorhexidine solution because they may adversely affect the surface of the RUNTHROUGH NS.
- Do not use dry gauze as this may damage the surface coating.

6. When exchanging the balloon dilatation catheter

CAUTION

Carefully read the following directions and the instructions manual accompanying the balloon dilatation catheters.

180 cm guide wire

< Using with an "Over-The-Wire" dilatation catheter >

- Extend the RUNTHROUGH NS by connecting the extension wire, before withdrawing/inserting the catheter.

CAUTION

When withdrawing the Over-The-Wire dilatation catheter in combination with the RUNTHROUGH NS, manipulate the RUNTHROUGH NS slowly and carefully confirming the behaviour of the distal end of the RUNTHROUGH NS by injecting heparinized physiological saline solution or contrast media through the guide wire lumen connected to a syringe or an injector. The RUNTHROUGH NS may advance and may damage the vessel.

< Using with a "Rapid Exchange" dilatation catheter >

- Maintain the position of the RUNTHROUGH NS and carefully remove/insert the dilatation catheter while confirming the position of the guide wire's tip under high resolution fluoroscopy.

300 cm guide wire

< Using with an "Over-The-Wire" dilatation catheter >

- Prior to exchanging the dilatation catheter, wipe the RUNTHROUGH NS with sterile gauze moistened with heparinized physiological saline solution.
- Proceed with exchanging the dilatation catheter, while maintaining the guide wire's position in the vessel under high resolution fluoroscopy.

PRECAUTIONS FOR STORAGE

Avoid exposure to water, direct sunlight, extreme temperatures and high humidity during storage. Store under controlled room temperature.

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© TERUMO CORPORATION 2014-06