Archimede

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24050 Grassobbio (BG) - Italy

ĐĂNG KÝ

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Tiêu chuẩn tham khảo IEC 60601-1 3rd Ed. Medical device safety

IEC 60601-1-2 4th Ed. Electromagnetic compatibility (EMC)

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www.delcon-usa.com

Archimede

Automatic blood components extractor

Instructions for use

REVISION HISTORY

REVISION	DATE	DESCRIPTION
Rev.1.0	September 2021	First release

PREFACE

These Instructions for Use are intended to give information about the automatic blood components separator ARCHIMEDE, its use, and maintenance, including its accessories. For further details on technical assistance, please refer to our Service Manual. Please read these Instructions for Use thoroughly before using the device.

CONTACTS

Manufacturer	Delcon Srl
Italy	Via Zanica 19/F/I 24050 Grassobbio (BG) - Italy
Phone	(+39) 035 4522336
E-mail	info@delcon.it
Web	www.delcon.it
USA	1704 East Boulevard, Suite 101 Charlotte, NC 28203 – USA
Phone	(+1) 704 243 9942
E-mail	info@delcon-usa.com
Web	www.delcon-usa.com

COMPLIANCE

Quality System compliant with the latest ISO 13485 edition.

The product carries the CE marking, in compliance with MDR 2017/745, Regulation EU on medical devices

MANUFACTURER'S RESPONSIBILITIES

The instructions may contain a statement of responsibility regarding the safety and correct use of this, only if device:

- Assembly, setting up, modifying and repair operations must be carried out by authorized personnel.
- The electrical system is proceeded in accordance with the instructions.
- The device is used in accordance with these instructions for use.

Delcon S.r.l. reserves the right to make corrections, modifications, improvements, improvements and other changes to its products, instructions for use and services at any time without notice.

All images in these instructions are for demonstration purposes only and may differ from the final product.

COPYRIGHT

All rights reserved. No part of this document may be reproduced, stored or transmitted in any form without the prior written permission of Delcon S.r.l. Despite all the care taken during the preparation of this document, Delcon S.r.l. assumes no responsibility for the improper use or negligence. Nor is it responsible for any damage caused by the use of the information contained in this document. These same rights are reserved for any person or company involved in the creation and production of these instructions for use. All trademarks mentioned in these instructions are the property of their respective owners.

TABLE OF CONTENTS

INTRODUCTION	1
Symbols and conventions	1
Rating plate data for the machine	3
Special warnings	4
Warnings for personnel safety	5
Notes on electromagnetic emissions	6
Environmental risks and end of life management	10
Warranty	10
OVERVIEW	11
Intended purpose	11
<i>Main features of the device</i>	12
DESCRIPTION OF THE DEVICE	13
Upper section	13
Front section	14
Side panel	14
Back panel	15
Keyboard	16
Backlit screen	17
<i>Main menu</i>	17
<i>List of active procedures</i>	17
<i>Service menu</i>	17
<i>Date and time</i>	18
INSTALLATION	19
Unpacking	19
Where to install the device	20
Installation	20
<i>Mounting the Balance Media, SAG-M Press, Filter.</i>	20
Switching on the device	22
Settings	23
<i>Date and time</i>	23
<i>Procedure settings check</i>	23
Switching off the device	24
HOW TO USE ARCHIMEDE	25
Turning the power on	25
System self-check	26
Main menu	30
Using ARCHIMEDE as a scale	31
Using ARCHIMEDE as a sealer	31

Using ARCHIMEDE as a manual separator	32
SEPARATION PROTOCOLS	33
Questions	34
Meaning of protocol parameters	34
DESCRIPTION OF PROTOCOLS	36
<i>PROTOCOL 1: T&T Triple PPP or PRP + RBC</i>	38
<i>PROTOCOL 2 : Quadruple T&T PPP or PRP BC + RBC</i>	40
<i>PROTOCOL 3: Triple or Quadruple T&B with PPP + BC + RBC Leucocyte Filter</i>	42
<i>PROTOCOL 4: Quadruple or Quintuple T&B PPP + BC + RBC</i>	44
<i>PROTOCOL 5: Separation of PRP to PPP + PRP</i>	46
<i>PROTOCOL 7: PLT + BC with Pool Filter</i>	48
<i>PROTOCOL 8: Washed RBC Separation</i>	50
<i>PROTOCOL 10 [a]: Umbilical Cord Blood Separation</i>	52
<i>PROTOCOL 11 [b]: Component Aliquoting</i>	54
<i>PROTOCOL 12 [c]: PLT + BC of Single Buffy Coat</i>	56
<i>PROTOCOL 13 [d]: T&T RBC with PPP + PPP</i>	58
SCALES CALIBRATION	60
DATA TRANSFER	62
LAN	62
WLAN	62
Updating the firmware and procedure protocols	63
MAINTENANCE	64
Weekly Maintenance	65
<i>Cleaning of the body</i>	65
<i>Cleaning of the electrodes</i>	65
Decontamination	65
<i>Clamps decontamination</i>	66
Biannual maintenance	66
WARNINGS AND ERRORS	67
ARCHIMEDE errors without codes	67
Coding errors and possible solutions	68
COMMERCIAL AND REGULATORY DATA	76
Trade name	76
Manufacturer	76
Regulatory	76
Accessories	77

INTRODUCTION

SYMBOLS AND CONVENTIONS

This table defines the meaning of symbols used in these Instructions for Use, on the documentation, and on the device.

Symbol	Description
	CE Marking
	General warning sign. Read the related documentation to prevent injuries to persons and to the device.
	Refer to instructions for use. "Follow instruction for use"
	Read carefully
	Manufacturer's address
	YYYY-MM-DD Format Manufacturing Date
	Product code
	XXNNNYY serial number format
	MD – Medical Device
	UDI Code (<i>sample</i>): (01) Manufacturer and device identifier (11) Production date (21) Serial number
	This electrical and electronic device cannot be disposed of with domestic waste. It must be disposed of separate from urban waste using the specific waste collection system designated by the local government or authorities.
	Wi-Fi module that complies with the Federal Communications Commission (USA) regulations

Symbol	Description
	Medical device that includes RF transmitters
	Fuses
	Warning: Risk of hand crushing
Operator	The person who uses the device.
Device or ARCHIMEDE	The device described in these instructions for use
IFU	Instructions for use

RATING PLATE DATA FOR THE MACHINE

<p>1</p>	 <div data-bbox="284 454 687 584">  <p>(01) 08033928530429 (11) YYMMDD (21) SSA-YMNNNN</p> </div>		<p>1. UDI Code</p> <p>2. FCC-ID Code Manufacturer</p> <p>3. REF, serial number, and production date</p> <p>4. Electrical Features</p> <p>5. Symbols</p>
<p>2</p>	 <p>via Zanica 19 F/I Grassobbio (BG) 24050 - Italy</p> <div data-bbox="488 618 683 651"> <p>REF MOE1504TER</p> </div> <div data-bbox="488 663 683 696"> <p>SN SSA-YMNNNN</p> </div> <div data-bbox="488 707 683 741"> <p>YYYY-MM-DD</p> </div>	<p>3</p>	
<p>4</p>	<p>100-240 V~ · 50-60Hz · 300VA</p> <div data-bbox="292 831 336 875">  </div> <p>Fuses: 2 x T3,15A 250V IP: 41 FCC ID: R68XPICOW Label version: R1.0</p>		
<p>5</p>	<div data-bbox="284 969 687 1021">  </div> <div data-bbox="371 1088 603 1167"> <p>DELCON <i>Our technology. Your health</i></p> </div>		

SPECIAL WARNINGS

READ CAREFULLY

These instructions for use are intended to be an integral part of the device itself and must always be kept with the product. Before using the device or any accessories, read this publication carefully.

The device must be used only for the purposes for which it is designed, as described in these instructions for use.

ARCHIMEDE uses advanced technology and is therefore fragile. Be careful when handling it including but not limited to dropping or banging against other equipment. Do not sit or step on the device.

Although the device was designed to protect against electromagnetic field interference, the unit must be installed at a certain distance from high emission sources such as UPS, high voltage housings, radio transmitters, and mobile phones (for mobile phones distance should be at least 1.5m or 5ft). Please refer to chapter *Notes on electromagnetic emissions*.

Do not place the device on top of other devices. Additionally, the device must be kept at a distance from life support, first aid, and IVD apparatus. Likewise, the device may not work properly if placed near diagnostic systems or medical-therapeutic treatment units with non-compatible emissions, as per the table in the chapter *Notes on electromagnetic emissions*.

Use only original accessories and replacement parts. The use of non Delcon cables or accessories can cause a decline in electromagnetic compatibility and damage to both ARCHIMEDE and nearby devices.

The device is not meant to be used in the presence of flammable mixtures, such as anesthetics, oxygen, nitrous oxide, or others.

Daily and weekly cleaning must be carried out as described in the chapter maintenance. Always follow the procedures. At least once a year, a Technical Assistance Service technician (or a qualified and authorized personnel) must carry out maintenance on the device.

Only fully trained, authorized Delcon S.r.l. personnel can do service on the device and/or remove the device's external safety protections. Always disconnect the device from its battery or power cord before service and maintenance.

Delcon S.r.l. decline every responsibility for the incorrect operation, action, or modification made to the device by unqualified personnel or by personnel who are not expressly authorized by the manufacturer, and specifically in the following cases:

- For any use of the device outside the intended purpose described in this Instructions for use.
- For damages to persons or objects caused by the incorrect installation or by maintenance carried out by personnel who are not authorized by Delcon S.r.l.
- For electrical or mechanical modifications made during or after the installation, for any reason, without Delcon S.r.l. authorization

NOTE: upon request, Delcon S.r.l. can supply diagrams, component lists, and other documentation to those who are authorized technical personnel.

Adverse events reporting

Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user is established. The report should include at least the name of the device and its code, the serial number, the UDI code, the description of the incident. Please contact:

- DELCON: info@delcon.it or service@delcon.it

WARNINGS FOR PERSONNEL SAFETY

The use of the device is safe unless equipment is not used properly. To avoid injuries, we recommend carefully reading this section and following its instruction. The device must be used by qualified operators who have a basic knowledge of handling and managing blood and its components.

- Never try to pull or move the tubing while the sealing procedure is in progress. If you apply force on the tubing during the sealing procedure, the tube might break and cause the leakage of potentially harmful fluids
- Never touch the electrodes while sealing. While the surge of energy used by the sealer lasts only a few milliseconds, the temperature between the electrodes is high and could burn your fingers
- Never place any metallic or other objects between the electrodes. Doing so could result in physical burns for the user or damage to the unit
- Before using any bag, whether empty or filled with blood or blood components, be sure that the bag is undamaged. Damaged bags can be dangerous for you and the device
- Caution! Watch your fingers. Many parts on the device (such as presses, clamps, sealing heads, and breakers) can exert enough force to injure hands and fingers
- To avoid the risk of electric shock, the equipment must be only connected to a supply mains with protective earth;
- The power cord plug shall match the power receptacle for the country in which the power source is being used. If it does not, contact an authorized Delcon service center
- ARCHIMEDE requires a supply voltage of 100-240V 50-60Hz 230VA. Check if your available power voltage is inside of above declared range”
- Any serious incident that has occurred in relation to the device should be reported to the manufacturer and the Competent Authority of the Member State in which the user is established. The report should include at least the name of the device and its code, the serial number, the UDI code, the description of the incident.

NOTES ON ELECTROMAGNETIC EMISSIONS

Manufacturer's guide and declarations regarding electromagnetic emissions.

CAUTION!

Mobile and portable RF communication devices may interfere with the device during its functioning.

Be careful.

Manufacturer's guide and declarations – Electromagnetic emissions

ARCHIMEDE is intended for use in the electromagnetic environment specified below. It is the customer's or the user's responsibility to ensure that ARCHIMEDE is used in such an environment.

Emission test	Compliance	Electromagnetic environment – Guidance
RF emissions <i>CISPR11</i>	Group 1	ARCHIMEDE uses RF energy only for its internal function. Therefore, its RF emissions are very low and unlikely to affect nearby electronic equipment.
RF emissions <i>CISPR11</i>	Class B	ARCHIMEDE is suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions <i>IEC/EN 61000-3-2</i>	Class A	
Voltage fluctuations / flicker emissions <i>IEC/EN 61000-3-3</i>	Complies	

Manufacturer's guide and declarations – Electromagnetic immunity

ARCHIMEDE is intended for use in the electromagnetic environment specified below. It is the customer's or the user's responsibility to ensure that ARCHIMEDE is used in such an environment.

Immunity test	IEC 60601 Test Level	Compliance	Electromagnetic environment – Guidance
Electrostatic discharge <i>(ESD) IEC 61000-4-2</i>	Contact: $\pm 8\text{kV}$ Air: $\pm 2, 4, 8, 15\text{ kV}$	Complies	Floor covering should be made of wood, concrete or ceramic tile. If the floor is covered with a synthetic material, the relative humidity should be at least 30%.
Fast transient <i>IEC 61000-4-4</i>	Input AC Power Ports: $\pm 2\text{kV}$ Signal Input/Output Ports: $\pm 1\text{kV}$	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Surge Immunity Test <i>IEC 61000-4-5</i>	Input Power Ports: $\pm 0,5, 1\text{ kV}$ ($\pm 2\text{kV}$ Line to Earth) Signal Input/Output: $\pm 0,5, 1, 2\text{ kV}$	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips and interruptions <i>IEC 61000-4-11</i>	Voltage dips: > 95% U_T for 0.5 cycles > 95% U_T for 1 cycle 30% U_T for 25 cycles (50Hz) 30% U_T for 30 cycles (60Hz) > 95% U_T for 250 cycles (50Hz) > 95% U_T for 300 cycles (60Hz)	Complies	Mains power quality should be that of a typical commercial or hospital environment.
Power frequency (50-60 Hz) magnetic field <i>IEC 61000-4-8</i>	3 A/m	Complies	Power frequency magnetic fields should be at levels characteristic of a typical location on a typical commercial or hospital environment.

NOTE: U_T is the AC mains voltage prior to application of the test level.

Immunity Test	IEC 60601 Test Level	Compliance Level
Conducted RF <i>IEC 61000-4-6</i>	Input AC Power Ports: 3V rms outside the ISM band 6V rms in the ISM band 150 kHz to 80 MHz	Complies
Radiated RF <i>IEC 61000-4-3</i>	3V/m 80 MHz to 2,7GHz	Complies

Electromagnetic environment - Guidance

Portable and mobile RF communications equipment should be used no closer to any part of ARCHIMEDE, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter.

Recommended separation distance

$$d = \left[\frac{3.5}{3} \sqrt{P} \right] = 0.28m \quad 0.150 \text{ MHz to } 80 \text{ MHz}$$

$$d = \left[\frac{3.5}{3} \sqrt{P} \right] = 0.28m \quad 80 \text{ MHz to } 800 \text{ MHz}$$

$$d = \left[\frac{7}{3} \sqrt{P} \right] = 0.56m \quad 800 \text{ MHz to } 2,5 \text{ GHz}$$

where *P* is the maximum output power rating of the transmitter in Watts (W) according to the transmitter manufacturer and *d* is the recommended separation distance in meters (m).

Field strength from fixed RF transmitters, as determined by an electromagnetic site survey ^a, should be less than the compliance level in each frequency range. ^b

Interference may occur near equipment marked with the following symbol:

NOTE: At 80 MHz and 800 MHz, the higher frequency range applies

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

^a Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which ARCHIMEDE is used exceeds the applicable RF compliance level above, ARCHIMEDE should be observed, additional measures may be necessary, such as re-orienting or relocating ARCHIMEDE.

^b Over the frequency range 150 kHz to 80 MHz, field strengths should be less than [V1] V/m.

Recommended separation distances between portable and mobile RF communications equipment and ARCHIMEDE

ARCHIMEDE is intended for use in an electromagnetic environment in which radiated RF disturbances are controlled. The customer or the user of ARCHIMEDE can help prevent electromagnetic interference by maintaining a minimum distance between portable and mobile RF communications equipment (transmitters) and ARCHIMEDE recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power or transmitter <i>P_n</i>	Separation distance according to frequency of transmitter <i>m</i>		
	150 kHz to 80 MHz $d \left[\frac{3,5}{3} \sqrt{P} \right]$	80 MHz to 800 MHz $d \left[\frac{3,5}{3} \sqrt{P} \right]$	800 MHz to 2,5 GHz $d \left[\frac{7}{3} \sqrt{P} \right]$
0,01	0,12	0,12	0,23
0,1	0,37	0,37	0,74
1	1,17	1,17	2,33
10	3,69	3,69	7,38
100	11,67	11,67	23,33

For transmitters rated at a maximum output power not listed above, the recommended separation distance *d* in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where *P* is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer.

NOTE: At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE: These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

ENVIRONMENTAL RISKS AND END OF LIFE MANAGEMENT

The device contains materials that, at the end of their life cycle, must be treated as waste in compliance with local laws and regulations.

Specifically, the device contains the following materials and/or components:

- Iron / Steel
- Copper
- Aluminum
- Non-biodegradable plastics
- Fibers glass for printed circuits
- Silicon based electronic components
- ABS
- External surface paint

When this symbol is applied to a product or on its packaging, it means that the product is subject to the Directive 2011/65/EU (ROHS II) on the restriction of the use of certain hazardous substances in electrical and electronic equipment and the Directive 2012/19/EU on waste electrical and electronic equipment (WEEE) and to the directive 2006/66/EC regarding batteries and accumulators and their disposal.

This electric and electronic device cannot be disposed of with normal household waste. It must be handled separately from normal urban waste and sent exclusively to collection centres that have been designated by the relevant local authorities.

To prevent possible health risks for the personnel in charge of disposal and recycling, the user is obligated to remove from the device any trace of potentially dangerous biological, chemical, or radioactive substances.

Illegal disposal by the holder may lead to the application of administrative and penal sanctions as provided for by the national laws in force.

WARRANTY

For information regarding the warranty conditions, contact Delcon S.r.l.

SECTION1.

OVERVIEW

INTENDED PURPOSE

ARCHIMEDE is intended for the extraction of blood components from whole blood collected in single-use bags.

Before blood components extraction, blood bags must be subjected to centrifugation and this process doesn't take place in the device but in specific centrifuges for blood bags located in blood bank laboratories. During centrifugation, the four essential blood components (plasma, platelets, buffy coat and erythrocytes) are separated into different layers according to their specific density.

In no cases the device and its accessories enter in contact with blood.

ARCHIMEDE is intended to be used in blood banks laboratory.

Users (operators) have to be qualified and well-skilled persons with a basic knowledge of the preparation and management of blood components.

Operators are in contact with the machine during:

- bags loading
- magnetic holder positioning (if at least one magnetic holder is required)
- positioning of the profile plate (if a different profile plate has to be mounted)
- activation and setting processes
- cleaning purposes
- service

The operator's involvement is not expected during the ordinary machine process. However, his support is required in case of failure or interruption of machine programs or in case of alerts.

No contact with the patient is expected in the intended purpose.

To function optimally ARCHIMEDE must be placed on a solid, stable surface and it is preferable not to use it near flammable sources and thermal induction.

ARCHIMEDE must be positioned in such a way as to allow the operator to turn it off easily and disconnect it from the mains.

IMPORTANT: To avoid the risk of electric shock, this device must be exclusively connected with a plug with ground protection!

CAUTION!

Position the device on a flat and stable surface or a stable table without wheels, capable to hold at least the weight of the device/m² and comfortable for the operator

Never position the device too close to centrifuges for blood bags. The vibrations generated by the large centrifuges can create operational problems and if they are old type, they can also create an electromagnetic disturbance

The device must be used and placed in blood bank laboratory away from medical devices for life support, first aid medical devices and IVD devices used for diagnosis

Be careful.

READ CAREFULLY!

ARCHIMEDE should be used by trained operators, with basic knowledge of the preparation and handling of blood components. Use is not recommended for operators with mental or physical disabilities.

Operator training does not substitute reading and understanding these instructions for use.

Main features of the device

- The device can be used with certified bags containing up to 600 ml and with standard tubes from 3.5 to 6.0 mm external diameter.
- Never use **frozen** or refrigerated bags: the least flexibility to perform an extraction is not guaranteed.
- The device can store up to 100 programs and handle up to 800 barcodes or RFID tags, with the possibility of format and consistency control, and bidirectional functions. Always remember to control the residual storage capacity and export data to free up space for new records.
- Work can be done with single devices, in stand-alone mode, or connected to a network and in communication with data management software. Bidirectional communication supports external data collection, remote adjustment of the structure and program operating parameters, automatic program selection according to work lists.

SECTION 2.

DESCRIPTION OF THE DEVICE

ARCHIMEDE consists of the following parts:

- An electronic system with an integrated state-of-the-art microcontroller that stores all the procedural data and manages the screen assembly, motor drives, scales, various sensors, clamps, sealing unit, and RFID reader
- A mechanical system with an axial system for handling the movement of the plate and the buffy coat separation system


Figure 1: ARCHIMEDE device

UPPER SECTION

The upper section accommodates the tubes and collection bags. It consists of the following parts:

1. Upper scale, with a range from 0 to 2,000 g and a tolerance of ± 1 g.
2. Photometer for detecting the presence of red blood cells (RBC), consisting of a 540nm optical system and a pair of infrareds (IR) sensors that detect the presence of air bubbles. Two reflective optical systems detect the status of the lid (closed/open) and the presence of a tube.
3. Three clamps, with "Top," "Plasma," and "Buffy coat" sealing heads.
4. One "SAG-M" clamp without a sealing head.
5. One clamp for plasma flow control.
6. "SAG-M" bag holder.
7. Actuator for automatic breaking of the "SAG-M" cannula.

FRONT SECTION

The front section accommodates the bag from which the blood components are being extracted. It consists of the following parts:

8. Mobile extraction plate operated by a stepper motor and two sensors: the first detects the distance between the machine body and the contoured plate, and the second detects the force exerted.
9. Fixed plate equipped with hooks to hang the bag on, depending on the procedure.
10. Primary bag support with a scale from 0 to 2,000 g with a tolerance of ± 1 g.
11. Actuator for automatic breaking of the primary bag cannula.
12. Magnetic Hall effect sensors for automatic detection of the contoured plate.
13. Optical unit to detect the buffy coat level. It consists of 18 IR sensors arranged linearly.
14. "Bottom" clamp with sealing head.

SIDE PANEL



Figure 2: Side panel

The side panel accommodates the side bag and any filter holder. It consists of the following:

15. Lateral scale, with a range from 0 to 2 grams and a tolerance of ± 1 g.
16. Cannula-breaking system, which automatically breaks the cannula of the bag located on the side holder.

BACK PANEL

The back panel consists of the power socket and two connectors that allow data exchange with the barcode reader and the LAN.

Power supply unit that consists of a power socket, on/off switch, and fuse holder.

Antenna for WLAN connection to BloodDonorsLINK.

PS2 connector to connect with a barcode reader.

LAN connector with one yellow and one green LED light.

The yellow LED light does the following:

- Flashes: in case of errors.
- Off: normal operation.

The green LED light is lit up, when there is activity on the LAN.

Note:

Consult the attached documentation.

KEYBOARD

The keyboard has various buttons, some of which are multifunctional.


Figure 3: Keyboard

Stops some functions and allows you to return to the top menu level.

Allow you to move within the menu and to change the value of the function selected.

Allow you to confirm the change or selection made.

Allows you to start the extraction procedure.

Allows you to terminate the current procedure or application.

Allows you to start the function indicated on the screen.

BACKLIT SCREEN

The large 240 × 128 LCD screen is backlit to enhance reading in low-light environments.

Main menu

Scan the bag type or select the procedure

By pressing the  button, you can access the list of installed procedures if the icon is selected (when ARCHIMEDE is turned on, this icon is already selected).

List of active procedures


Figure 4: List of procedures

Use the up and down arrow buttons  to move the selection cursor. Press the  button to start the procedure selected with ✓.

Service menu

To access this menu, press the  button above the  icon and use the arrow buttons  to select the desired item, and confirm by pressing the  button.


Figure 5: Service menu

CALIBRATION: Allows you to calibrate the three ARCHIMEDE scales by selecting the desired weight.

LIST OF PARAMETERS:	Allows you to browse all the numeric values of all installed procedural parameters.
LIST OF OPERATORS:	Allows you to browse all the codes list associated with the operators authorised to use ARCHIMEDE.
CONTRAST	Allows you to set the contrast.
MAINTENANCE:	Allows you to browse the list of active maintenance activities and to perform them (only if connected to BloodDonorsLINK).

Date and time

Even if the screen only sometimes shows the hour and minute, ARCHIMEDE is equipped with an integrated calendar used to send procedural data, including the date and time, to BloodDonorsLINK.

To learn how to set up date and time, consult section 2.5.1.

Note: If the centre uses BloodDonorsLINK, the date and time are automatically updated at the first start up or during the first connection by referring to the date and time on the PC.

SECTION 3.

INSTALLATION

UNPACKING

Carefully remove the device from its packing and check that all the components inside the box are undamaged.

Delcon recommends preserving the original packaging for eventual future usage (service or other features).

The device is supplied with:

- Barcode reader
- Filter support
- Bag mounts
- Power cord
- Air removal system
- Side bag supports and breaker
- Profile plate shapes
- Pressure SAG-M

Check that the electrical power cable is suitable for the electrical socket in the country where the device is to be used. If it is not, contact Delcon S.r.l. or an authorized service center.

Delcon S.r.l. recommends keeping the original box and packaging for possible future use or storage.

CAUTION!

The device weighs approximately 48 Kg without Sag-M press or 54 Kg with Sag-M press: do not move it without help, you could hurt your back, or drop the device.

2 people can easily move the device using the 2 retractable handles found under its base. Each person can take hold of a handle with one hand, and with the other one guide the upper part of the device, to counterbalance the weight.

Carefully take the device out of its packaging together with all the accessory parts.

Check that all device parts are whole.

No modifications to the device are allowed

WHERE TO INSTALL THE DEVICE

CAUTION!

To avoid electrocution, this device must be connected exclusively to electric sockets with earthing.

To achieve best performance in safe conditions, it is important that the device be positioned on work tables meant to support at least 50Kg/m², so as to prevent accidental falls which would be destructive for the device and dangerous for operators.

The work tables must be stable, have no wheels, and be free from mechanical vibrations which can disrupt the scales. If multiple devices are installed on the same table, sufficient space must be left around each single device. Check that table height makes use of the device easy for operators of any stature.

Do not position the device near large floor-mounted centrifuges for bags or in structures which are connected mechanically to large refrigerators. The vibrations transmitted by the motors can create operating problems and, in case of poor-quality or antiquated devices, there can be interfering electromagnetic disturbances.

The work environment needs to be sufficiently well-lit and climatized, while not being exposed to direct sunlight or intense heat sources, such as radiators or stoves. Adequate ventilation must prevent excesses in humidity and concentration of aggressive air-borne agents, be they acidic or alkaline, such as those used for cleaning or disinfecting hospital areas, as well as from air-borne organic solvents, which can make some shiny surfaces opaque.

The device must be positioned at a reasonable distance from hospital equipment which is critical for patients' lives, in so far as it can condition their correct functioning.

CAUTION!

The device as to be position far away from heat sources in order to avoid possible explosion and/or fire.

ARCHIMEDE was not designed to be used in presence of flammable mixture such as anesthetics, oxygen, nitrous oxide or others.

INSTALLATION

Remove ARCHIMEDE from his packaging taking care not to violently bump the scales.

Check that all parts removed from the packaging are not damaged.

ARCHIMEDE kit consists of:

- Side balance panel.
- Filter support.
- 1D barcode reader.
- Power cord.
- Top balance plastic support

Mounting the Balance Media, SAG-M Press, Filter.

Pay attention to scales because even if they are mechanically protected, it is necessary not to hit them violently.

- Top balance (2): Remove the two screws from the top balance stand, place the plate above the top balance stand, and secure it with the two screws removed earlier.
- Side balance (3): Fix the support of the metal side plate with three screws on the support bracket placed on the right side of ARCHIMEDE (the three screws are already present on the bracket).
- SAG-M press: After removing the 5 press fixing screws from the upper right corner of the module that also includes the cannula breaking system and fix it with the five screws removed previously. Connect the cable that handles the signals of the cannula breakage system and bag presence to the CB3 connector, locate on the back of ARCHIMEDE in the upper left.
- Filter carrier: The filter support is magnetic and can be placed in the most suitable position to support the filter for the procedures that provide for it.
- Installation of forks of cannula breakers

DANGER!

Never touch cannula breaker forks, moving parts and cannulas during breaking cycles.

The device uses enough force to harm your fingers.

SWITCHING ON THE DEVICE

First of all, check the effectiveness of the electrical connection's earthing. If the earthing is damaged, both effectiveness and safety can be reduced. In this case, do not turn on the device at all until the problem has been resolved.

CAUTION!

The device requires power of 100-240 VAC at 50-60 Hz, with maximum absorbency of 300 VA

Check that the electrical line has compatible voltage and separation power

Be sure that the main switch on the back of the device is switched off before proceeding to connect to the electric power supply.

Without any bag loaded on the device, with the door closed completely, and without touching the scales, turn on the device with the main switch.

The device will start the initialization process and calibration of all its components. The screen will show the brand and will report what process is under way.

Figure 6: ARCHIMEDE screen

Note: If ARCHIMEDE is not in standalone mode, but is communicating with BlooDonorsLINK on the network, before proceeding with other checks, you must set the LAN/WLAN network parameters. This configuration is described in the Service Manual in the LAN/WLAN configuration section.

The screen should turn on and indicate the firmware model and software version.


Figure 7: "Self-check in progress" screen

After verifying and initialising the various data in the memory, ARCHIMEDE performs all self-checks in all measurement systems. In case of an anomaly, the screen displays the type of problem detected and a possible solution.

If the self-check is passed, the device will move on to the operating menu, otherwise the screen will show appropriate diagnostic messages.


CAUTION!

The start procedure can last some minutes. During start up, do not touch any part of the device. It is dangerous for your fingers, and you could easily falsify the calibration procedure. To improve the device’s stability regarding variations in temperature, and thus the accuracy of scale readings, it is advisable to wait at least 30 minutes before starting the first extraction.

SETTINGS

All device parameters are present at the factory in permanent memory and can be edited using BloodDonorsLINK; but the date and time are managed by a circuit that can operate for more than 60 days, even without a primary power source. However, if ARCHIMEDE is not operating in standalone mode, the date and time are automatically updated when it is turned on, synchronising with BloodDonorsLINK. The date and time must be checked at the first use only if the system is not integrated with BloodDonorsLINK.

Date and time

Using the arrow buttons,  highlight the screen area that indicates the current date and time and press  to enter the editing mode or enter directly by pressing the  button over the clock.

Using the arrow buttons , set the time and date format and confirm the change with the  button. If you do not want to confirm the change, press the  button.

FORMAT	READING
12	0-24
24	0-12 am/pm
dd/mm/yyyy	01/06/2017
mm/dd/yyyy	06/01/2017
yyyy/mm/dd	2017/01/06

Figure 8: List of date and time formats

Procedure settings check

This check allows you to verify the various procedure parameters, and it must be carried out by the application specialist. The procedure parameters and their definitions are described and set in BloodDonorsLINK. They are displayed on ARCHIMEDE screen to allow users, especially those, who work in standalone mode, to verify their accuracy. To access this menu, press the  button

above the  icon and use the arrow buttons  to select the PARAMETERS LIST. Confirm by pressing the  button.

Use the arrow buttons  to select the desired procedure.

DANGER!

Never insert anything other than a single-use bag for processing blood components in the press.

Before loading a bag with blood or blood components in the press, be sure that the bag is whole. Never use clearly damaged bags: this can be harmful both for you and the device.


CAUTION Watch your fingers! Both presses (the main and lateral ones) exert enough force to injure fingers. During processing never insert hands or fingers between the press plates.

After verifying and initializing the various data in memory, ARCHIMEDE will perform self-control of all measurement systems. In case there is any anomaly, the display will display the type of problem detected and the possible remedy.

WARNING!

Be careful not to get your hands on the plate work area

SWITCHING OFF THE DEVICE

Switching off the device is possible at any time using the main switch positioned on the back of the device.

CAUTION!

In case of emergency, power supply for ARCHIMEDE can be interrupted by removing the plug from the electrical socket or by removing the connector for the AC entrance tray, positioned on the back near the main switch.

SECTION 4.

HOW TO USE ARCHIMEDE

The device works through a series of *extraction programs*, a sequence of instructions with associated parameters that describes what the device will do, from start to end. Some programs are pre-configured on the device, but to achieve the best possible within the customer requirement they need to be customized through the device interface or, more effectively, through the BlooDonorsLINK software.


The perfect functioning of ARCHIMEDE depends on the environmental conditions, therefore it is necessary to avoid the following conditions:

- Places with high temperatures, high humidity (see Chapter 13 of these instructions for use)
- Places near open flames
- Places in direct sunlight
- Places where dust accumulates
- Places in the vicinity of flammable substances
- On non-flat and unstable surfaces
- Place ARCHIMEDE so that the minimum distance between the antenna and the operator is 20 cm

Sensors inside the device measure the following conditions: placement of the tubes, buffy coat level, presence of red blood cells (RBC), weights, main bag volume, extraction speed, and extraction force. ARCHIMEDE uses these measurements to properly execute the procedure. During the procedure, if ARCHIMEDE is connected to BlooDonorsLINK, it sends out all sensor data and the data on the extraction procedure in progress, including all data acquired through the barcode reader. If ARCHIMEDE is in standalone mode, the procedural data are stored in first-in, first-out (FIFO) memory and are then automatically sent to BlooDonorsLINK once the connection is established. ARCHIMEDE can save up to 300 procedures before it starts overwriting the oldest one.

TURNING THE POWER ON

When the power is turned on, the entire system performs a self-check. Before powering on ARCHIMEDE, you must remove all bags from their holders and make sure that the plate and bag holders can move freely and are properly positioned.

After connecting the power cable, turn on the main switch, located on the back of the device.

The screen must be backlit, and it must show the word "ARCHIMEDE". If this is not the case, check that the power cable is properly plugged in and that there is voltage in the socket used.

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SYSTEM SELF-CHECK

Every time ARCHIMEDE is powered on, a self-check takes place to check that the system is operating properly.



Figure 9: "Self-check in progress" screen

Actions taken after a negative check are related to the connection mode.


Standalone mode

Connected to BloodDonorsLINK

E²prom memory

Checks the integrity of the calibration data and the procedural parameters.


If the procedural parameters are incorrect, after a warning, the system initialises, resetting them to the factory default values. Calibration of the data is required.

The procedural data are automatically restored. For calibration data, calibration is required.

Clock

Checks the time and date.


If the time and date values are entered in a different format, after confirmation, the system initialises their settings to HH:MM:SS MM/DD/YY.

Synchronisation with the date and time of BloodDonorsLINK takes place at each start up.

Tube sensors

Checks for the presence of tubes and self-calibration of the sensors.


If the calibration values exceed the predefined limits, an error is displayed. If the problem persists, consult the "Troubleshooting" section.

Same as in standalone.

CSU (Central Sealing Unit)

Checks the sealing station and tube presence sensor.


If the CSU is not in the "ready" status and it detects the presence of moisture, an error is displayed. If the problem persists, consult the "Troubleshooting" section.

Same as in standalone mode.

CBC (Cannula-Breaking System)

Checks the cannula-breaking system.


If optical limits are not detected within the predefined time intervals, an error is displayed. If the problem persists, consult the "Troubleshooting" section.

Same as in standalone mode.

Plate motor

Checks the motor movement and limit.


If the limit is not activated within the predefined travel limits, an error is displayed. If the problem persists, consult the "Troubleshooting" section.

Same as in standalone mode.

Force sensor

Checks the force sensor.


If the sensor detects a force higher than the predefined limits, an error is displayed. If the problem persists, consult the "Troubleshooting" section.

Same as in standalone mode.

Buffy coat separator motor

Checks the motor movement and limit.


If the limit is not activated within the predefined travel limits, an error is displayed. If the problem persists, consult the “Troubleshooting” section.

Same as in standalone mode.

IR sensors

Checks and self-calibrates the front infrared (IR) sensors.


If the sensor detects a reading outside of the predefined limits, an error is displayed. If the problem persists, consult the “Troubleshooting” section.

Same as in standalone mode.

HB sensor

Checks and self-calibrates the hemoglobin (HB) sensor.


If the HB sensor has a high dark current or a low sensitivity, an error is displayed. If the problem persists, consult the “Troubleshooting” section.

Same as in standalone mode.

Flow regulator

Checks the actuator movement.


If a problem with the movement of the regulator motor is detected, an error is displayed. If the problem persists, consult the “Troubleshooting” section.

Same as in standalone mode.

Scales

Checks the tare values of all the scales.

If the tare values are not within the predefined limits, an error is displayed. If the problem persists, consult the “Troubleshooting” section.






Same as in standalone mode.

MAIN MENU



There are two different ways to access the various operations:

- With the arrow buttons , you can select a specific icon.
- With the operational buttons , you select the operation above or below the button.

The operations associated with the various icons are as follows:

	Allows you to access the service submenu: calibration, list of parameters, list of operators, maintenance.
	Allows you to log in and log out using the barcode reader.
	Area that displays machine status.
	Area that displays the current time, which also allows access to the time and date editing menu.
	Tare menu and display of the front, upper, and side scales.
	Area that displays the date on which the scales were calibrated.
	Allows you to select manual squeezing.
	Allows you to open and close the selected clamp.
	The arrow appears when sealing is active.
	Allows you to access the procedure selection menu. You can select the procedure manually or by scanning the bag type barcode.

USING ARCHIMEDE AS A SCALE

You can use ARCHIMEDE's three installed load cells as a scale. The maximum measurable weight is 2 kg, with a tolerance of ± 1 g. Before using the system as a scale, you must tare it. Move the cursor using the arrow buttons  until a zero is displayed in the scale area and press the  button to reset the scales. Now move the cursor to another option to avoid pressing the  button to re-run the tare of the scales. Hang the product you want to weigh on the front or side scale, or place it on the upper scale.

The weights are displayed in the following order:


Figure 10: Scale reading area

USING ARCHIMEDE AS A SEALER

You can use ARCHIMEDE's four integrated sealing heads to seal tubes.

Press the  button under the sealing head that you want to use or use the arrow buttons  to highlight the desired clamp.

The sealing heads are located as indicated below:


Figure 11: Sealing heads

To perform sealing with the desired sealing head, press the  button. To open and close the selected clamp without performing a seal, press the  button.

The sealing operation is displayed on the left side of the clamp with the  symbol.

If any errors occur, these messages are displayed in the warning window:

- **Clean the sealing head:** The CSU has detected the presence of moisture. Clean the sealing head and the tube and check for leaks. The same error also indicates that a seal was attempted without inserting the tube.
- **Insert tube:** The tube sensor has detected that the tube is not inserted correctly or that it is not inserted into the clamp.
- **CSU busy:** The CSU is busy with a sealing cycle; try again after a few seconds.

Using ARCHIMEDE as a continuous sealing unit requires **non-intensive use**; otherwise, the sealing unit may overheat and go into thermal protection mode.

USING ARCHIMEDE AS A MANUAL SEPARATOR

ARCHIMEDE can also be used as a manually operated separation system.

In this mode, the following controls are enabled:

- Position sensor that determines the position of the plate.
- Force sensor that stops the plate if the force exceeds the set point; the plate moves again when the force drops below the threshold. The force setting is shown in the related screen.

NOTE: THE FORCE SENSOR IS ACTIVE ONLY WHEN THE PLATE MOVES TOWARD THE MACHINE BODY.

NOTE: DO NOT PUT YOUR HANDS IN THE WORKING AREA OF THE PLATE WHILE IT IS MOVING.

To enter this mode, use the arrow buttons  to highlight the  icon and press the  button. Press the up or down arrow button to move the plate in the desired direction, and use the right or left arrow button to vary the force applied at fixed increments of 5 units.

To stop the plate while maintaining its position, press . To exit and bring the plate to the home position, press the  button.

In this mode, you can also activate the clamps by pressing the  button, and you can perform a seal by pressing the  button.

Clamp 5 is not equipped with a sealing head and therefore cannot be used as a sealer.

In addition, the scales are active and can be used to detect the weights of the hanging bags.

CAUTION!


Caps of the clamps must be removed only for maintenance, while the ordinary use of the device foresees the clamps with their caps.

SECTION 5.

SEPARATION PROTOCOLS

ARCHIMEDE can store up to 18 protocols. The operator can select protocols manually or automatically if the system is connected to BlooDonorsLINK.

- **In manual mode**, to facilitate the selection of protocols in use at the centre, you can configure the system to display only the desired protocols. In this mode, only the protocols that are used frequently will be displayed.

Press the  button to access the list of installed protocols if the icon is highlighted.

Using the arrow buttons , highlight the protocol you want to use and press the  button to start the procedure.


Figure 12: List of active protocols

- **In automatic mode**, you can associate the bag type with the protocol type in BlooDonorsLINK. In this mode, when you scan the barcode of the bag type in the main menu, BlooDonorsLINK sends the type of protocol to ARCHIMEDE based on the available data. If the bag type information is not enough to select the type of protocol, ARCHIMEDE will ask more questions, or you may select the protocols that are associated with that bag type from a smaller list, if there is more than one.

The separation protocols that are currently available are as follows:

PRIMARY

TYPE	NAME
1	PROTOCOL 1: T&T Triple PPP or PRP + RBC
2	PROTOCOL 2 : Quadruple T&T PPP or PRP BC + RBC
3	PROTOCOL 3: Triple or Quadruple T&B with PPP + BC + RBC Leucocyte Filter
4	PROTOCOL 4: Quadruple or Quintuple T&B PPP + BC + RBC
10	PROTOCOL 10 [a]: Umbilical Cord Blood Separation
11	PROTOCOL 11 [b]: Component Aliquoting
13	PROTOCOL 13 [d]: T&T RBC with PPP + PPP

SECONDARY

TYPE	NAME
5	PROTOCOL 5: Separation of PRP to PPP + PRP
7	PROTOCOL 7: PLT + BC with Pool Filter
8	PROTOCOL 8: Washed RBC Separation
12	PROTOCOL 12 [c]: PLT + BC of Single Buffy Coat

QUESTIONS

This feature, which is not available in standalone mode, allows for complete traceability of each procedure. You can set up a series of questions in BloodDonorsLINK to be asked right after selecting the protocol to be used. The answers are checked based on various criteria: length, type of barcode, font type, authorised operator IDs, donor IDs and bag lot validation, LIS/HIS validation, suffix, and prefix.

There are two possible modes that can be set in BloodDonorsLINK:

- **Manual answer-question association:** An arrow is displayed that indicates the code to read. After you scan it with the barcode reader, the barcode data is displayed and the arrow moves to the next question. The operator must read all the barcodes set in the questions sequentially.
- **Automatic answer-question association:** Only the questions and any automatic answers are displayed, and the operator can read any of the desired barcodes. The system then uses the barcode features to recognise which question to associate with the barcode that is read. The operator can read the barcodes in any sequence.

Use the arrow buttons  to move the arrow among questions in manual association. In automatic association, re-read the code to automatically replace it in the correct position.

If the bag has already been processed (the donation ID is already in the database) with a particular protocol, BloodDonorsLINK does not allow separation. An authorised operation can override this.

To cancel the question sequence and return to the main menu, press the  button.

Completing all the questions automatically opens the initial screen of the selected protocol.

MEANING OF PROTOCOL PARAMETERS

ARCHIMEDE uses various parameters during the extraction procedure. These can be fixed (predefined in the protocol) or variable (set by the application specialist via BloodDonorsLINK).

All the variable parameters are pre-set at the factory and are stored in non-volatile memory with a high retention capacity.

The fixed parameters are those that control events with absolute priority, such as time needed to enable various actuators and SAG-M dispensing speeds.

For better management, the variable parameters used in the protocols have been divided into groups based on their meaning. Consult the service manual for a more detailed explanation.

DESCRIPTION OF PROTOCOLS

Some operations are common to all protocols.

The following information is displayed in the initial screen:

Status of the protocol (performing separation or paused)

Clamps: the clamps displayed are those used for separation or for sealing only

HB sensor with its related status (open or closed, with or without tube)

Force detected by the related sensor

Last four digits of the operator code (if previously read)

Number that indicates the type of protocol (left) and the sequential number of the procedure list (right)

Percentage of estimated time for the protocol

Separation time from the breaking of the cannulas to the end of sealing

'Read comments': You can read any comments available on the outside of the fixed front plate, if connected to BlooDonorsLINK

BCL (Buffy coat level)

SET Level of buffy coat set

ACT Actual buffy coat level

STOP

Permanent stop:

Press the  button, then press either the  button or the  button to exit the procedure. ARCHIMEDE displays the sealing page and sends the data to BlooDonorsLINK.

Pause:

Press the  button to pause the procedure, then press the  button or the  button to resume it.

PLEASE NOTE THAT THE FORCE SENSOR IS ACTIVE ONLY WHEN THE PLATE IS MOVING TOWARDS THE MACHINE BODY.

DO NOT PUT YOUR HANDS IN THE WORKING AREA OF THE PLATE WHILE IT IS MOVING.

Sealing

The sealing cycle is linked to the value set in parameters #11 and #12.

Automatic sealing: The sealing cycle of the selected clamp starts automatically.

Manual sealing: Press the  button to seal the enabled lines; otherwise, press the  button and then the  button to continue without sealing.

Note: To better manage the various kits that have multiple tube lengths, ARCHIMEDE provides the option to change the valves used in separation (#40).

Sending data

ARCHIMEDE automatically sends all the events and data related to the procedure to BlooDonorsLINK.

When the procedure is complete, the screen displays the detected weights. If the tubes are not properly arranged, scale readings may be distorted. If this occurs, rearrange the tubes, then weigh the blood components again and send the new weights to BlooDonorsLINK.

Correct weights

After you remove all the tubes and bags, ARCHIMEDE automatically starts a procedure of the same type that we just performed.

You can also start the same type of procedure by pressing the  button. To return to the list of protocols, press the  button.

Incorrect weights

Press one of the function buttons  below the screen to weigh the blood components again. To accept the new values and send them to BlooDonorsLINK, press the  button.

Normal limits

With BlooDonorsLINK, you can specify the acceptable weight of blood components. If control is enabled and the detected weights are not within the set limits, the screen displays the weight of the detected blood components and the limits allowed. Press any button to continue.

Remove the tubes and bags and press the  button to continue with another procedure of the same type; otherwise, press the  button to exit.

PROTOCOL 1: T&T Triple PPP or PRP + RBC

This protocol uses standard triple "top and top" (T&T) bags, with or without SAG-M, and allows the preparation of red blood cells (RBC) and platelet-poor plasma (PPP) or platelet-rich plasma (PRP). It weighs the plasma and RBC.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 13: T&T Triple PPP or PRP + RBC protocol

2. Hang the primary bag on the bag holder on the front panel. The label showing the bag data should face outward. To help with the placement of the SAG-M bag, after detecting the presence of the primary bag, the plate partially closes, and the screen displays the bag icon.
3. Hang the SAG-M bag on the bag supports.
4. Insert the SAG-M tube firmly into clamp 3.
5. Open the HB sensor cover and insert the tube from the primary bag into the sensing area.
6. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
7. Firmly insert the tube coming from the HB sensor into clamp 1.
8. Place the plasma bag on the upper scale (scale 2).
9. Firmly insert the tube of the plasma bag into clamp 2.

If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.

10. Press the  button to start the procedure if it does not start automatically.

If the check for the expected weight of the primary bag has been enabled through BlooDonorsLINK, and the weight detected by ARCHIMEDE is not within the predefined limits, a warning is displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, confirm the weight with the  button or press the  button to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube. The procedure will not start until all tubes are inserted into the active clamps. The related error is also displayed for the undetected bags.

11. When the plate approaches the bag, wait until the force equals the predefined value in parameter #3 for this protocol.

- If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
- If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannulas and the SAG-M-bag, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.
 - If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.
- If cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the collected blood components. If enabled by parameter #6, ARCHIMEDE checks the level of RBC in the primary bag. When the RBC level reaches the selected IR sensor, the screen displays the RBC in the plate area, and the force applied by the primary bag is reduced to 10 kg.

The photometer that detects the presence of RBC now checks, in real time, the flow of the blood components. When the concentration of RBC exceeds the predefined threshold in parameter #8, the RBC is displayed on the sensor icon. If dispensing of excess plasma (#19) is not enabled, the extraction process is blocked, and the procedure moves on to the SAG-M dispensing phase.

12. Air recovery.

- Automatic air recovery (#42):
 - The system sends the air contained in the plasma bag to the selected bag.
- Manual air recovery (#42):
 - If automatic air recovery is not enabled, perform it manually: Place the plasma bag in a vertical position, press the  button to open clamp 3, and manually remove the air by pressing on the bag. When you have finished, press the  button to close clamp 3, and reposition the plasma bag for weight detection. The procedure continues automatically.

13. Remove the tubes and bags.

- To continue with another separation of the same type, press the  button.
- To exit, press the  button.

PROTOCOL 2 : Quadruple T&T PPP or PRP BC + RBC

This procedure uses standard quadruple "top and top" (T&T) bags, with or without SAG-M, and allows the preparation of red blood cells (RBC) and platelet-poor plasma (PPP) or platelet-rich plasma (PRP) and buffy coat (BC) with platelets. It weighs the plasma, the BC, and the RBC.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 14: Quadruple T&T PPP or PRP BC + RBC protocol

2. Hang the primary bag on the bag holders on the front panel. The label showing the bag data should face outward. To help with the placement of the SAG-M bag, after detecting the presence of the primary bag, the plate partially closes, and the screen displays the bag icon.
3. Hang the SAG-M bag on the bag supports.
4. Open the HB sensor cover and insert the tube from the primary bag into the sensing area.
5. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
6. Firmly insert the tube coming from the HB sensor into clamp 1.
7. Place the buffy coat bag on the side scale (scale 3). INVERTED ORDER
8. Insert the buffy coat tube into clamp 2.
9. Place the plasma bag on the upper scale (scale 2).
10. Insert the plasma tube into clamp 3.

If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.

11. Press the  button to start the procedure if it does not start automatically.

If the check for the expected weight of the primary bag has been enabled through BloodDonorsLINK, and the weight detected by ARCHIMEDE is not within the predefined limits, a warning is displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, confirm the weight with the  button or press the  button to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube. The procedure will not start until all tubes are inserted into the active clamps. The related error is also displayed for the undetected bags.

12. When the plate approaches the bag, wait until the force equals the predefined value in parameter #3 for this protocol.

- If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
- If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannulas and the SAG-M bag, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.
 - If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.
- If cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weights of the blood components. ARCHIMEDE dispenses a quantity of plasma into the buffy coat bag that equals the predefined value in parameter #34. If enabled by parameter #6, the system then checks the level of RBC in the primary bag. When the RBC level reaches the selected IR sensor, the screen displays the RBC in the plate area, and the force applied by the primary bag is reduced to 10 kg.

The photometer that detects the presence of the RBC now checks, in real time, the flow of the blood components. When the concentration of RBC exceeds the predefined threshold in parameter #8, the RBC is displayed on the sensor icon. If dispensing of excess plasma (#19) is not enabled, the plasma line is blocked, and the line for the buffy coat opens.

When the amount of buffy coat reaches the predefined value from parameter #22, the system activates the upper separation mechanical system and adjusts the amount of buffy coat that was set by start flowing (#37). The system retracts the plate and dispenses an additional amount of plasma into the buffy coat bag, if this was defined in parameter #36, to clean the buffy coat tube.

If it was defined in parameter #35, before dispensing the SAG-M into the RBC, the buffy coat clamp is kept open to clean the line. The remaining SAG-M is then dispensed.

13. Air recovery.

- Automatic air recovery (#42):
The system sends the air contained in the plasma bag to the selected bag.
- Manual air recovery (#42):
If automatic air recovery is not enabled, perform it manually: Place the plasma bag in a vertical position, press the  button to open clamp 3, and manually remove the air by pressing on the bag. When you have finished, press the  button to close clamp 3, and reposition the plasma bag for weight detection. The procedure continues automatically.

14. Remove the tubes and bags.

- To continue with another separation of the same type, press the  button.

To exit, press the  button.

PROTOCOL 3: Triple or Quadruple T&B with PPP + BC + RBC Leucocyte Filter

This procedure uses triple or quadruple top and bottom (T&B) bags with a leucocyte filter to obtain buffy coat (BC), red blood cells (RBC) and plasma (PPP). It weighs the plasma, BC and RBC.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 15: Triple or Quadruple T&B with PPP + BC + RBC Leucocyte Filter protocol

2. Hang the primary bag on the bag holders on the front panel. The label showing the bag data should face outward. After detecting the presence of the primary bag, the plate partially closes, and the screen displays the bag icon.
3. Open the HB sensor cover and insert the tube from the primary bag into the sensing area.
4. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
5. Place the plasma bag on the upper scale (scale 2).
6. Place the RBC bag on the side scale (scale 3).
7. Firmly insert the tube coming from the HB sensor into clamp 6.
8. Firmly insert the tube coming from clamp 6 into clamp 3.
9. Insert the bottom tube into clamp 4.
10. If the kit is equipped with a leucocyte filter:
 - Hang the SAG-M bag on the bag supports.
 - Insert the filter into the filter holder, if present, and open the manual clamp.

If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps.

The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.

11. Press the  button to start the procedure if it does not start automatically.

If the check for the expected weight of the primary bag has been enabled through BlooDonorsLINK, and the weight detected by ARCHIMEDE is not within the predefined limits, a warning is displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, press the  button to confirm the weight or press the  button to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube, and the procedure will not start until all tubes are inserted into the active clamps. The related error is also shown for the undetected bags.

12. When the plate approaches the bag, wait until the force equals the predefined value in parameter #3 for this protocol.

- If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
- If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannulas and the SAG-M bag, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.
 - If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.
- If automatic cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the blood components. ARCHIMEDE checks the level of buffy coat with the IR optical sensors in clamps 6, 3, and 4, while the distance sensor blocks the extraction process when the volume of buffy coat reaches the predefined value in parameter #18.

When the procedure is complete, the system retracts the plate. If enabled, it dispenses an amount of plasma into the buffy coat bag equal to the predefined value in parameter #36. In addition, if parameter #20 is enabled, the system dispenses the SAG-M.

13. Air recovery.

Automatic air recovery (#42):

The system sends the air contained in the plasma bag to the selected bag.

- Manual air recovery (#42):
 - If automatic air recovery is not enabled, perform it manually: Place the plasma bag in a vertical position, press the  button to open clamp 3, and manually remove the air by pressing on the bag. When you have finished, press the  button to close clamp 3, and reposition the plasma bag for weight detection. The procedure continues automatically.

14. Remove the tubes and bags.

- To continue with another separation of the same type, press the  button.
- To exit, press the  button.

PROTOCOL 4: *Quadruple or Quintuple T&B PPP + BC + RBC*

This procedure uses quadruple or quintuple top and bottom (T&B) bags with an in-line filter to obtain buffy coat (BC), red blood cells (RBC), and two platelet-poor plasma (PPP) products. It weighs the plasma, buffy coat and RBC.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 16: Quadruple or Quintuple T&B PPP + BC + RBC protocol

2. Hang the primary bag on the bag holders on the front panel. The label showing the bag data should face outward. After detecting the presence of the primary bag, the plate partially closes, and the screen displays the bag icon.
3. Open the HB sensor cover and insert the tube from the primary bag into the sensing area.
4. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
5. Place the plasma bag on the upper scale (scale 2).
6. Place the RBC bag on the side scale (scale 3).
7. Place the secondary plasma bag on the upper scale (scale 2).
8. Firmly insert the tube coming from the HB sensor into clamp 6.
9. Firmly insert the tube coming from clamp 6 into clamp 3.
10. Firmly insert the tube of the secondary plasma bag into clamp 2.
11. Insert the bottom tube into clamp 4.
12. If the kit is equipped with a leucocyte filter:
 - Hang the SAG-M bag on the bag supports.
 - Insert the filter into the filter holder, if present, and open the manual clamp.

If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon

and the indicator light turns green. When all the bags and tubes are in place, the text at the top right will change to START and the procedure starts automatically.

13. Press the  button to start the procedure if it does not start automatically.

If the check of the expected weight of the primary bag has been enabled through BlooDonorsLINK, and the weight detected by ARCHIMEDE is not within the predefined limits, a warning will be displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, press the  button to confirm the weight or press the  button to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen shows the clamp number that did not properly detect the tube. The procedure will not start until all tubes are inserted into the active clamps. The related error is also displayed for the undetected bags.

14. When the plate approaches the bag, wait until the force equals the predefined value in parameter #3 for this protocol.

- If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
- If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannulas and the SAG-M bag, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.

- If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.

- If cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the blood components. ARCHIMEDE checks the level of buffy coat with the IR optical sensors in clamps 6, 2, 3 and 4, while the distance sensor blocks the extraction process when the volume of buffy coat reaches the predefined value in parameter #18.

During this phase, if it was enabled with parameter #34, a quantity of plasma is dispensed into the secondary plasma bag.

At the end of the procedure, the system retracts the plate. If enabled, it dispenses an amount of plasma into the buffy coat bag equal to the predefined value in parameter #36. In addition, if parameter #20 is enabled, the system dispenses the SAG-M.

15. Air recovery.

- Automatic air recovery (#42):

The system sends the air contained in the plasma bag to the selected bag.

- Manual air recovery (#42):

If automatic air recovery is not enabled, perform it manually: Place the plasma bag in a vertical position, press the  button to open clamp 3, and manually remove the air by pressing on the bag. When you have finished, press the  button to close clamp 3, and reposition the plasma bag for weight detection. The procedure continues automatically.

16. Remove the tubes and bags.

- To continue with another separation of the same type, press the  button.

- To exit, press the  button.

PROTOCOL 5: Separation of PRP to PPP + PRP

This procedure separates platelet-rich plasma (PRP) from platelet-poor plasma (PPP). It weighs both the platelets and the residue.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 17: Separation of PRP to PPP + PRP protocol

2. Hang the primary bag on the bag holders on the front panel.
 3. Firmly insert the plasma tube into clamp 3.
 4. Place the PPP plasma bag on the upper scale (scale 2).
- If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.
5. Press the  button to start the procedure if it does not start automatically.
- If the check of the expected weight of the primary bag has been enabled through BlooDonorsLINK, and the weight detected by ARCHIMEDE is not within the predefined limits, a warning is displayed that highlights both the limits and the actual weight.
 - After checking the cause of the problem, confirm the weight with  or press  to exit the procedure.

6. When the plate approaches the bag, wait until the force equals the predefined value in parameter #3 for this protocol.
- If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
 - If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannula, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.
 - If clogged tubes or fluid obstructions prevents the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.
 - If cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen displays the current operating phase number and the weight of the collected plasma. ARCHIMEDE checks the volume of the primary bag. When the volume equals the predefined value in parameter #30 or it reaches the distance set in #18, the system blocks the extraction process and retracts to the home position.

7. Remove the tubes and bags.
- To continue with another separation of the same type, press the  button.
 - To exit, press the  button.

PROTOCOL 7: PLT + BC with Pool Filter

This procedure uses of double bags for collecting buffy coat pool platelets with a filter to obtain platelets (PLT) and residual buffy coat (Res). It weighs the platelets and residual buffy coat.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 18: PLT + BC With Pool Filter protocol

2. Hang the primary bag on the bag holders on the front panel. The label that contains the bag data should face outward. After detecting the presence of the primary bag, the plate moves back, and the screen displays the bag icon.
 3. Insert the filter into the appropriate filter holder.
 4. Open the HB sensor cover and insert the tube from the filter into the sensing area.
 5. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
 6. Firmly insert the tube coming from the HB sensor into clamp 1.
 7. Place the PLT bag on the upper scale (2).
- If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right will change to START and the procedure will start automatically.
8. Press the  button to start the procedure if it does not start automatically.

If the check of the expected weight of the primary bag has been enabled through BlooDonorsLINK, and the weight detected by ARCHIMEDE is not within the limits set, a warning is displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, confirm the weight with  or press  to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube and the procedure will not start until all tubes are inserted into the active clamps. The related error is also displayed for the undetected bags.

9. When the plate is approaches the bag, wait until the force equals the predefined value in parameter #3 for this protocol.
 - If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and start the separation cycle.
 - If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannula, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.

- If clogged tubes or fluid obstructions prevents the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.

- If cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the collected PLT. ARCHIMEDE controls the flow by stopping the plate in case of filter obstruction or RBC presence in the HB sensor.

- If the filter is obstructed, press the  button to the next warning.
- If the HB sensor detects RBC and parameter #19 is set to "manual," the system asks whether to proceed with manual recovery or to end the procedure, if a quantity of plasma is set to be dispensed automatically.

10. Air recovery.

- Automatic air recovery (#42):

The system sends the air contained in the plasma bag to the selected bag.

- Manual air recovery (#42):

If automatic air recovery is not enabled, perform it manually: Place the plasma bag in a vertical position, press the  button to open clamp 3, and manually remove the air by pressing on the bag. When you have finished, press the  button to close clamp 3, and reposition the plasma bag for weight detection. The procedure continues automatically.

11. Remove the tubes and bags.

- To continue with another separation of the same type, press the  button.

- To exit, press the  button.

PROTOCOL 8: Washed RBC Separation

This procedure enables separation of the washing fluid or any supernatant from 1 or 2 units of red blood cells (RBC). It weighs the RBC and supernatant.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 19: Washed RBC Separation protocol

2. Hang the primary bag on the bag holders on the front panel. The label that contains the bag data should face outward. After detecting the presence of the primary bag, the plate moves back, and the screen displays the bag icon.
3. Open the HB sensor cover and insert the tube from the primary bag into the sensing area.
4. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
5. Firmly insert the tube coming from the HB sensor into clamp 1.
6. Place the bag for the washing fluid on the upper scale (scale 2).

If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.

7. Press the  button to start the procedure if it does not start automatically.

If the check for the expected weight of the primary bag has been enabled through BloodDonorsLINK, and the weight detected by ARCHIMEDE is not within the limits set, a warning is displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, confirm the weight with  or press  to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube and the procedure will not start until all tubes are inserted into the active clamps. The related error is also displayed for the undetected bags.

8. When the plate is approaches the bag, wait until the force equals the predefined value in parameter #3 for this protocol.
 - If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
 - If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannula, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.
 - If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.
 - If the cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the collected washing fluid.

ARCHIMEDE checks the level of RBC in the primary bag, and when it reaches the IR sensor predefined in parameter #6, the clamp closes by applying a delay equal to parameter #26. If, on the other hand, the HB sensor is active in parameter #6, the clamp closes when the RBC reach the HB sensor, and ARCHIMEDE displays RBC on the sensor cover. The system then opens the plate by dispensing an amount of supernatant equal to the predefined value in parameter #26.

9. Remove the tubes and bags.
 - To continue with another separation of the same type, press the  button.
 - To exit, press the  button.

PROTOCOL 10 [a]: Umbilical Cord Blood Separation

This procedure uses triple top and bottom (T&B) bags to obtain red blood cells (RBC), plasma (PPP), and buffy coat (BC) diluted with plasma added by pressure. It weighs the plasma, BC, and RBC.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 20: Umbilical Cord Blood Separation protocol

2. Hang the primary bag on the bag holders on the front panel. The label that contains the bag data should face outward. After detecting the presence of the primary bag, the plate moves back, and the screen displays the bag icon.
 3. Open the HB sensor cover, and insert the tube from the primary bag into the sensing area.
 4. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
 5. Firmly insert the tube coming from the HB sensor into clamp 6.
 6. Firmly insert the tube coming from clamp 6 into clamp 3.
 7. Place the plasma bag on the upper scale (scale 2).
 8. Place the RBC bag on the side scale (scale 3).
 9. Insert the bottom tube into clamp 4.
- If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.
10. Press the  button to start the procedure if it does not start automatically.

If the check for the expected weight of the primary bag has been enabled through BlooDonorsLINK, and the weight detected by ARCHIMEDE is not within the limits set, a warning is displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, confirm the weight with  or press  to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube. The procedure will not start until all tubes are inserted into the active clamps. The related error is also shown for the undetected bags.

11. When the plate is approaches the bag, wait until the force equals the predefined value in parameter #3 for this protocol.

- If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
- If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannula, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.
 - If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.
- If the cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the blood components. ARCHIMEDE checks the level of buffy coat with the IR optical sensors in clamps 6, 3 and 4, while the extraction process is blocked when the volume of buffy coat reaches the predefined value in parameter #18.

If parameter #37 is greater than parameter #18, the plasma is dispensed into the buffy coat to obtain the desired weight.

12. Remove the tubes and bags.

- To continue with another separation of the same type, press the  button.
- To exit, press the  button.

PROTOCOL 11 [b]: Component Aliquoting

This procedure enables the separation of a blood component into the four available bags. It weighs each bag separately.

Operational sequences:

Parameters #14, #29, and #39 are loaded.

1. If the volumes shown are the desired volumes, press **START** to proceed with the placement of the bags.
Otherwise, if you want to change the volume setting, press the up and down arrow buttons .
2. To change the bag selection, press the left and right arrow buttons .
3. and to disable the bag, reset the volume or press **ESC**.
4. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 21: Component Aliquoting protocol

5. Hang the primary bag on the bag holders on the front panel.
6. Insert the tubes from the bags placed on the upper scale into clamps 1 and 3.
7. Insert the tubes from the bags placed on the side scale into clamp 2.
8. If selected, place bag 1, if it is active, on the upper scale (scale 2).
9. If selected, place bag 2, if it is active, on the upper scale (scale 2).
10. If selected, place bag 3, if it is active, on the upper scale (3).

If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon

and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.

11. Press the  button to start the procedure if it does not start automatically.

If the check for the expected weight of the primary bag has been enabled through BlooDonorsLINK, and the weight detected by ARCHIMEDE is not within the limits set, a warning is displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, confirm the weight with  or press  to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube. The procedure will not start until all tubes are inserted into the active clamps. The related error is also displayed for the undetected bags.

12. When the plate approaches the bag, wait until the force equals the predefined value in parameter #3 of this protocol.
 - If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
 - If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannula, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.

- If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.

- If the cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the various aliquots. ARCHIMEDE checks the weight of the bags. When the predefined volume is reached, the clamp for that aliquot closes. After the last aliquot, the system retracts the plate to the resting position.

13. Remove the tubes and bags.

- To continue with another separation of the same type, press the  button.
- To exit, press the  button.

PROTOCOL 12 [c]: PLT + BC of Single Buffy Coat

This procedure uses double bags for collecting single buffy coat platelets or buffy coat pool platelets with a filter to obtain platelets (PLT) and residual buffy coat (Res). It weighs the platelets and residual buffy coat.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 22: PLT + BC of Single Buffy Coat protocol

2. Hang the primary bag on the bag holders on the front panel. The label that contains the bag data should face outward. After detecting the presence of the primary bag, the plate moves back, and the screen displays the bag icon.
 3. Open the HB sensor cover and insert the tube from the filter into the sensing area.
 4. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
 5. Firmly insert the tube coming from the HB sensor into clamp 1.
 6. Place the PLT bag on the upper scale (scale 2).
If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.
 7. Press the  button to start the procedure if it does not start automatically. If the check for the expected weight of the primary bag has been enabled through BloodDonorsLINK, and the weight detected by ARCHIMEDE is not within the limits set, a warning is displayed that highlights both the limits and the actual weight.
- After checking the cause of the problem, confirm the weight with  or press  to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube. The procedure will not start until all tubes are inserted into the active clamps. The related error is also displayed for the undetected bags.

8. When the plate approaches the bag, wait until the force equals the predefined value for parameter #3 for this protocol.
 - If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
 - If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannula, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.
 - If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.
 - If the cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the collected PLT. ARCHIMEDE controls the flow by stopping the plate in case of RBC presence in the HB sensor. If parameter #19 is set to "manual", you are asked whether to proceed with manual recovery or to end the procedure, if a quantity of plasma is set to be dispensed automatically.

9. Air recovery.
 - Automatic air recovery (#42):
The system sends the air contained in the plasma bag to the selected bag.
 - Manual air recovery (#42):
If automatic air recovery is not enabled, perform it manually: Place the plasma bag in a vertical position, press the  button to open clamp 3, and manually remove the air by pressing on the bag. When you have finished, press the  button to close clamp 3, and reposition the plasma bag for weight detection. The procedure continues automatically.
10. Remove the tubes and bags.
 - To continue with another separation of the same type, press the  button.
 - To exit, press the  button.

PROTOCOL 13 [d]: T&T RBC with PPP + PPP

This procedure uses standard top and top (T&T) double bags to obtain RBC diluted with plasma in a predefined ratio without SAG-M. It weighs the diluted plasma and RBC.

Operational sequences:

1. Place the bags and tubes, avoiding tension and kinks that could create fluid fluctuations and weighing errors.


Figure 23: T&T RBC with PPP + PPP protocol

2. Hang the primary bag on the bag holders on the front panel. The label showing the bag data should face outward.
3. Open the HB sensor cover and insert the tube from the primary bag into the sensing area.
4. Close the lid. The tube must not be stretched so that automatic cannula breaking is possible.
5. Place the plasma bag on the upper scale (scale 2).
6. Firmly insert the tube coming from the HB sensor into clamp 1.
If the tube control option is enabled, ARCHIMEDE uses its optical sensors to check for proper insertion of the tubes into the relevant clamps. The tube is properly inserted when the numeric reference of the clamp is replaced by the  icon and the indicator light turns green. When all the bags and tubes are in place, the text at the top right changes to START and the procedure starts automatically.
7. Press the  button to start the procedure if it does not start automatically.

If the check for the expected weight of the primary bag has been enabled through BloodDonorsLINK, and the weight detected by ARCHIMEDE is not within the limits set, a warning is displayed that highlights both the limits and the actual weight.

- After checking the cause of the problem, confirm the weight with  or press  to exit the procedure.

In addition, if the check for the presence of tubes is enabled, the screen displays the clamp number that did not properly detect the tube. The procedure will not start until all tubes are inserted into the active clamps. The related error is also displayed for the undetected bags.

8. When the plate approaches the bag, wait until the force equals the predefined value for parameter #3 for this protocol.
 - If automatic cannula breaking (#45) is enabled, ARCHIMEDE automatically breaks the cannulas and starts the separation cycle.
 - If automatic cannula breaking (#45) is disabled and the presence of cannulas (#2) is enabled, the message "Break the cannulas" is displayed. Apply force sideways on the upper part of the primary cannula, using the pegs for assistance, and wait for a few seconds for the positive pressure to be released. After detecting the decrease in pressure, ARCHIMEDE automatically starts the separation cycle.
 - If clogged tubes or fluid obstructions prevent the circuit pressure from dropping below the predefined threshold, press the  button to force the start of the cycle.
 - If the cannula breaking (#2) is disabled, ARCHIMEDE automatically starts the separation cycle.

The screen now displays the current operating phase number and the weight of the collected plasma. ARCHIMEDE checks the level of the RBC and closes the procedure when the level reaches the desired value. Then the plate will retract to the home position.

9. Remove the tubes and bags.
 - To continue with another separation of the same type, press the  button.
 - To exit, press the  button.

SECTION 6.

SCALES CALIBRATION

Prior to calibrating the scales, ARCHIMEDE must be turned on for a few minutes for thermal stabilization of the measuring sensors. You must also have certified sample weights.

1. To access calibration, press the  button above the  icon and use the arrow buttons  to select CALIBRATION.
 - Press the  button to confirm.


Figure 24: Service menu

2. Press the arrow buttons  to select the value of the sample weight to use for the calibration.
 - Press the  button to confirm.


Figure 25: Calibration, selection of sample weight

3. To perform the calibration, remove all weight from the front scale.
 - Press the  button to confirm.


Figure 26: Calibration, first phase

4. Place the selected sample weight.
 - Verify that the displayed value is stable.
 - Press Enter to confirm.


Figure 27: Calibration, second phase

5. Repeat the previous operations for the upper scale (scale 2) and the lateral scale (scale 3).


Figure 28: Calibration, display of correction factors

6. Check the calibration factors calculated by the front scale (scale 1), the upper scale (scale 2) and the lateral scale (scale 3).
 - To confirm the factors and store them, press Enter; otherwise, press ESC to restore the previous calibration factors.

SECTION 7.

DATA TRANSFER

ARCHIMEDE can communicate bidirectionally with BlooDonorsLINK in both WLAN and LAN mode. This allows you to transfer data from multiple ARCHIMEDE systems, even those that are a few meters away from the PC where the BlooDonorsLINK server is installed. In the WLAN mode, the maximum possible distance depends on the environment and any disturbances that may be present.

LAN


Figure 29: LAN connection

Every ARCHIMEDE system connected to an IP network must have a unique IP address. This is assigned automatically if the DHCP mode is enabled. To communicate uniformly across a network with multiple units, ARCHIMEDE always sends the serial number to BlooDonorsLINK.

WLAN


Figure 30: WLAN connection

The PC used for the communication must have a wireless 802.11 b/g LAN installed, so that BloodDonorsLINK can identify ARCHIMEDE's unique serial number to communicate with multiple installed systems that are powered on.

Note:

The user must read the accompanying documents to ensure safe operation.

UPDATING THE FIRMWARE AND PROCEDURE PROTOCOLS

The following updates are possible using BloodDonorsLINK:

- Update the internal ARCHIMEDE firmware. It is preferable for the firmware update to be carried out by qualified technical staff who are authorised by the manufacturer.
- Insert new procedures up to a maximum of 18. It is preferable that new procedures be created by application specialists who are familiar with the functions of the various procedural parameters in detail.
- Modify the procedural parameters. Modifications must be performed by the trained operators at the centre or by application specialists to avoid entering values that are not suitable for the smooth operation of the separator.
- Verify the proper operation of all ARCHIMEDE peripherals by calling on qualified technical personnel who are authorised by the manufacturer.

Prior to calibrating the scales, ARCHIMEDE must be turned on for a few minutes for thermal stabilization of the measuring sensors.

SECTION 8.

MAINTENANCE

Regular maintenance involves keeping the device clean (both the body and the electrodes).

WARNING!

Before any cleaning or maintenance, always disconnect the power cord from the power supply.
Don't place the device in water or other liquids, don't spray the device with water or other liquids.

WARNING FOR CLEANING AND DISINFECTION!

Be careful when cleaning or disinfecting ARCHIMEDE. **Some products might contain dangerous ingredients that could damage the equipment.** Damage due to the use of improper cleaners will result in the loss of warranty on the product.

In case of prolonged unused, the device must be properly stored and protected from contaminants.

Products containing the following substances must not be used

Alcohols

- Pure Amyl Alcohol
- Pure Benzyl Alcohol
- Pure Butanol
- Pure Anhydrous Alcohol
- Methyl Alcohol
- Propyl Alcohol

Disinfectants and cleaning solutions

- Peroxide 90 volumes
- Ammonia concentrated
- Bleach > 25% chlorine
- Iodine Tincture

Acceptable solutions for the cleaning and disinfection of the device are:

Alcohols

- Ethanol in low concentrations <10%
- Ethanol < 20% (it can be used, but it must be immediately removed with a clean wet cloth)
- Methyl Alcohol < 20% (but request fast rinsing after application with wet clean cloth)

Disinfectants and cleaning solutions

- Peroxide 40 volumes
- Ammonia, dilute solution <1%
- Benzalkonium chloride <2% even blended with chlorhexidine 0.5–1%
- Formal <40%
- Bleach, with chlorine content <10%
- Merbromin

WEEKLY MAINTENANCE

Cleaning of the body

Wipe the body and the sealing device with a clean cloth dampened with an approved detergent or disinfectant. Avoid dropping or spilling any liquid inside the device. Never use pure alcohol or any other strong solvent to clean the painted body of the sealer. At the end, dry the surface with a clean gauze or cloth.

Cleaning of the electrodes

CAUTION!

Caps of the clamps must be removed only for maintenance, while the ordinary use of the device foresees the clamps with their caps on.

Use a clean gauze wet with alcohol to clean the electrodes frequently; old dirt left in place for a long time becomes very hard due to the effect of the heat making it hard to remove later without damaging the electrodes. At the end, dry the surface with a clean gauze or cloth and wait a few minutes before using the system again.

CAUTION!

Never use any abrasive material like glass-paper or metallic brushes to clean the electrodes and never apply too much force as you may damage the electrodes or compromise their alignment.

After cleaning the equipment, **rinse all parts with water (using a clean, damp cloth)** in order to remove traces of cleaning solutions that could damage the equipment over time.

DECONTAMINATION

Decontamination is an emergency procedure to be carried out only if strictly necessary, for example, if the device is doused in blood or any other possible fluid of human or animal origin which can cause infections on the exposed subject.

CAUTION!

Before starting the decontamination procedure, **always** wear proper personal protective equipment (gloves, glasses, protective masks, lab coats, etc.).

Always refer to your internal decontamination procedures, if any, prior to these instructions for use.

Once fully protected, you can proceed with cleaning with approved detergent (see table above) and subsequent sanitization using a microfiber cloth soaked with commercial hypochlorite solution. After the treatment, always rinse the parts with a damp cloth to remove any residual product.

Never use acetone, trichloroethylene, or other aggressive solvents.

Decontamination **must be** performed every time the device is returned for service or repair, and it is the responsibility of the hospital or the blood center to ensure that no traces of organic material remain on the device

Clamps decontamination

If a tube breaks inside a clamp, the liquid spilt should drain outside the clamp through the special slits at its base. If that's not the case, remove the head by unscrewing the cap and carefully clean the electrodes and the adjacent surfaces; **avoid pouring any liquid in the clamp**: use dampened gauze or cotton wool.

Once done, dry and verify that no residue is left in the clamp, blow strongly to eliminate all residue of gauze or cotton.

Before turning the device on again, wait at least 30 minutes to allow the parts to dry completely. Check clamps' functionality as described in the paragraph "Sealing clamp test", below.

BIANNUAL MAINTENANCE

Every 6 months be sure to:

- check the status of the scales with a 500g calibration weight. Access the *Balance calibration* menu through the *Technical Service* menu.
- check the status of the clamps. Access the *Clamp test* through the *Diagnostic* menu
- check the status of the proportional clamp. Access the *Calibration* menu through the *Technical Service* menu

CAUTION!

In case of prolonged unused, the device must be properly stored and protected from contaminants.

Verify that the symbols affixed to the device does not become illegible after a repeated number of cleaning operations.

WARNINGS AND ERRORS

ARCHIMEDE always monitors the proper operation of the various peripherals. It displays an error code if it finds an issue that cannot be resolved through specially designated algorithms.

ARCHIMEDE ERRORS WITHOUT CODES

Error 1A	THE DEVICE DOES NOT TURN ON
	Cause 1: Mains power is lacking.
	Solution: Check that the mains power is present and that the power cable is properly inserted.
	Cause 2: Fuses are blown.
	Solution: Replace the blown fuses with fuses that have the same characteristics as indicated in the nameplate information.
Error 2A	LOSS OF DATE AND TIME
	Cause 1: The real-time clock (RTC) circuit has been reset due to large network disturbances.
	Solution: Reprogram the date/time and leave ARCHIMEDE on for a few hours. If it is connected to BloodDonorsLINK, the date and time synchronisation occurs automatically when it is powered on.
	Cause 2: The device has been turned off for a long time.
	Solution: Leave the device turned on for at least one hour and try again.
	Cause 3: The RTC backup capacitor failed.
	Solution: Call the authorised technical assistance.
Error 3A	THE DISPLAY SHOWS INCOMPREHENSIBLE CHARACTERS
	Cause 1: Electronic interference.
	Solution: Press the ESC button. Turn ARCHIMEDE off, then on again. Update the firmware via ALNK. Call the authorised technical assistance.
Error 4A	ARCHIMEDE DOES NOT CONNECT TO BLOODONORSLINK
	Cause 1: BloodDonorsLINK is not active or ARCHIMEDE network is not connected.
	Solution: Activate the BloodDonorsLINK server. Connect to ARCHIMEDE WLAN network. Check that the sub-networks of ARCHIMEDE and the BloodDonorsLINK server are correct. Check the firewall settings.
	Cause 2: There is a communication error.
	Solution: Check the signal strength and make sure ARCHIMEDE antenna is screwed in properly and not blocked by metal objects. Move the access point or server closer to ARCHIMEDE. Change the WLAN module transmission channel. Call the authorised technical assistance.

CODING ERRORS AND POSSIBLE SOLUTIONS

Error 1	E²PROM IS BUSY
	Cause 1: E ² prom memory is blocked.
	Solution: Turn ARCHIMEDE off, then on again.
	Cause 2: E ² prom memory is defective.
	Solution: Call the authorised technical assistance.
Error 2	CHECK E²PROM
	Cause 1: E ² PROM has failed.
	Solution: Turn ARCHIMEDE off, then on again.
	Cause 2: E ² prom memory is defective.
	Solution: Call the authorised technical assistance.
Error 3	SPI TIMEOUT
	Cause 1: A peripheral has unexpectedly interrupted communication.
	Solution: Turn ARCHIMEDE off, then on again. Call the authorised technical assistance.
	Cause 2: The CPU or SO card is defective
	Solution: Call the authorised technical assistance.
Error 4	RFID TIMEOUT
	Cause 1: The RFID label is defective.
	Solution: Replace the RFID label.
	Cause 2: The RFID module does not respond.
	Solution: Turn the RFID module off, then on again. Call the authorised technical assistance.
Error 5	RFID COMMUNICATION ERROR
	Cause 1: RFID data is corrupt.
	Solution: Read another RFID label. Turn it off, then on again. Call the authorised technical assistance.
Error 6	BOOT LOADER IS CORRUPT
	Cause 1: The checksum of the boot loader is not correct.
	Solution: Turn ARCHIMEDE off, then on again. Call the authorised technical assistance.
Error 7	FIRMWARE IS CORRUPT
	Cause 1: The checksum of the firmware is not correct.
	Solution: Update the firmware via BlooDonorsLINK.
Error 8	ERROR IN CALIBRATION DATA
	Cause 1: Data is lost or corrupt.
	Solution: Calibrate the scales.
Error 9	VALUES OF THE PARAMETERS OUTSIDE OF THE RANGE

	Cause 1:	An error was shown during the check for the parameters.
	Solution:	Carry out the update procedure through BlooDonorsLINK.
Error 10 (1-4)	A/D CONVERTER IS IN OVERFLOW	
	Cause 1:	The analog tension is not correct.
	Solution:	Check that the scale holders are properly mounted and free of obstructions. Turn ARCHIMEDE off, then on again.
Error 11 (1-4)	A/D CONVERTER TIMEOUT	
	Cause 1:	The A/D convertor is blocked.
	Solution:	Turn ARCHIMEDE off, then on again. Call the authorised technical assistance.
Error 12 (1-4)	UNDERFLOW OF A/D CONVERTER DURING THE TARE OF THE SCALE (X)	
	Cause 1:	The scale holder is improperly positioned.
	Solution:	Check that the holders are properly mounted and free of obstructions. Turn ARCHIMEDE off, then on again.
Error 13 (1-4)	SCALE TARE (X) IS OUTSIDE OF THE LIMITS	
	Cause 1:	Holders are improperly positioned.
	Solution:	Check that the holders are properly mounted and free of obstructions. Turn ARCHIMEDE off, then on again. Call the authorised technical assistance.
Error 14 (1-4)	IN SEPARATION, THE TARE (X) IS OUTSIDE OF THE RANGE	
	Cause 1:	Weight detection circuit $\geq 1,000$ mV.
	Solution:	Remove all unnecessary weight from the scale holder. Call the authorised technical assistance.
Error 15 (1-4)	VALUES OF K SCALE (X) ARE OUTSIDE OF THE RANGE	
	Cause 1:	During the calibration, the system detected that the calibration factor is outside of the limit. The value of the weight used during the calibration does not match the weight that was set.
	Solution:	Repeat the calibration. Call the authorised technical assistance.
Error 16 (1-10)	BRIGHTNESS OF THE FRONT SENSORS	
	Cause 1:	The ambient brightness level is very high.
	Solution:	Decrease the ambient brightness level.
	Cause 2:	The IR sensors are defective.
	Solution:	Call the authorised technical assistance.
Error 17 (1-10)	SENSITIVITY OF THE IR SENSOR IS TOO LOW	
	Cause 1:	The IR sensor detects a difference that is too low between the on and off status. $V_{tr} < 4,500$ mV with the LEDs > 300 mA
	Solution:	Clean the yellow strip and the IR sensors.
	Cause 2:	The IR sensor coupler is faulty.
	Solution:	Call the authorised technical assistance.

Error 18 (1-10)

SENSITIVITY OF THE IR SENSORS IS TOO HIGH

Cause 1: The IR sensor detects a difference that is too high between the on and off status.

Solution: Decrease the ambient brightness.

Cause 2: The IR sensors are faulty.

Solution: Call the authorised technical assistance.

Error 19

SENSITIVITY OF THE HB SENSOR IS TOO LOW

Cause 1: The HB sensor does not detect a sufficient difference between the status with the LED on and off.

Solution: Clean the sensors.

Cause 2: The green sensor is faulty.

Solution: Call the authorised technical assistance.

Error 20

SENSITIVITY OF THE HB SENSOR IS TOO HIGH

Cause 1: The HB sensor detects a difference that is too high between the status with the LED on and off.

Solution: Call the authorised technical assistance.

Error 21

TIME OF PLATE OPENING > 30 SECONDS

Cause 1: A mechanical obstruction is preventing the plate from reaching the open or closed position. Time > 30 seconds.

Solution: Remove all obstructions. Clean the scroll bars.

Cause 2: The open position sensor is faulty.

Solution: Call the authorised technical assistance.

Error 22

TIME OF PLATE CLOSURE > 30 SECONDS

Cause 1: Mechanical obstructions are preventing the plate from reaching the closed position. Time > 30 seconds.

Solution: Remove all obstructions. Clean the scroll bars.

Cause 2: The force sensor is defective.

Solution: Call the authorised technical assistance.

Error 23

BLADE POSITIONING TIME > 10 SECONDS

Cause 1: Mechanical obstructions are preventing the blade from reaching the open or closed position. Time > 10 seconds.

Solution: Remove all obstructions. Clean the blade assembly.

Cause 2: BC position sensor is faulty.

Solution: Call the authorised technical assistance.

Error 24

MOVEMENT OF PLATE TO THE HOME POSITION > 30 SECONDS

Cause 1: Mechanical obstructions are preventing the plate from reaching the home position. Time > 30 seconds.

Solution: Remove all obstructions. Clean the scroll bars.

Cause 2: The home position sensor is faulty.

Solution: Call the authorised technical assistance.

Error 25

MOVEMENT OF PLATE TO THE IR CAL POSITION > 10 SECONDS

Cause 1: Mechanical obstructions are preventing the plate from reaching the home position. Time > 10 seconds.

Solution: Remove all obstructions. Clean the scroll bars.

Cause 2: The IR calibration position sensor is faulty.

Solution: Call the authorised technical assistance.

Error 26

IN ZERO POSITION SENSOR IS OVER RANGE

Cause 1: Mechanical obstructions are preventing the plate from reaching the home position.

Solution: Remove all obstructions.

Cause 2: The distance sensor is faulty.

Solution: Call the authorised technical assistance.

Error 27

WITH THE PLATE FREE, THE FORCE SENSOR HAS DETECTED A VALUE > 1

Cause 1: Mechanical obstructions are preventing the correlation of the distance sensor values with those of the force.

Solution: Remove all obstructions.

Cause 2: The force sensor is defective.

Solution: Call the authorised technical assistance.

Error 28

CALIBRATION POSITION DISTANCE WITH TENSION > 4600 mV

Cause 1: Mechanical obstructions are preventing the plate from reaching the calibration position.

Solution: Remove all obstructions.

Cause 2: The linear sensor is faulty.

Solution: Call the authorised technical assistance.

Error 29 (1-4)

SEALER HAS DIRTY ELECTRODES

Cause 1: The sealer has dirty or wet heads.

Solution: Clean and dry the sealing head that has an error. Press the  button to open the clamp that has an error, reinsert the tube, and press the  button to repeat the sealing procedure. If it is in the middle of the procedure, press the  button to continue.

Cause 2: The CSU or a sealing head is defective.

Solution: Call the authorised technical assistance.

Error 30

SEALING UNIT IS OCCUPIED

Cause 1: The CSU cable is disconnected.

Solution: Check and connect the CPU-CSU cable.

Cause 2: The CSU is defective.
Solution: Call the authorised technical assistance.

Error 31

TUBE IS IN THE CLAMP

Cause 1: The self-check detected the presence of a tube.

Solution: Remove all tubes and repeat the self-check.

Cause 2: The tube sensor is faulty.

Solution: Call the authorised technical assistance.

Error 32 (1-4)

OFFSET LOAD CELL HAS EXCEEDED THE LIMITS

Cause 1: During the self-check, a weight or a bag was placed on the load cell.

Solution: Remove all weights or bags and repeat the self-check.

Cause 2: The load cell is faulty.

Solution: Call the authorised technical assistance.

Error 33

CLAMP 6 MOTOR TIMEOUT

Cause 1: The clamp actuator movement is blocked.

Solution: Remove tubes or any obstructions from the actuator and try again.

Cause 2: The position sensor is defective.

Solution: Call the authorised technical assistance.

Error 34

THE NUMBER OF STEPS TO CARRY OUT A COMPLETE CYCLE OF CLAMP 6 IS NOT CORRECT

Cause 1: Clamp actuator movement is blocked.

Solution: Remove any tubes or obstructions and try again.

Cause 2: The motor is defective.

Solution: Call the authorised technical assistance.

Error 35 (1-3)

K CELL IS OVER RANGE

Cause 1: During the initial self-test, an error was detected in the calibration factors.

Solution: Recalibrate the scales.

Error 36

K SENSOR FORCE IS OVER RANGE

Cause 1: During the initial self-test, an error was detected in the calibration factor.

Solution: Recalibrate the force sensor.

Error 37

WEIGHT OF THE PLASMA-PLT IS INSUFFICIENT

Cause 1: The separated plasma volume is less than what is required.

Solution: Check the value that was set.

Cause 2: Scale is not calibrated.

Solution: Calibrate the scale.

Error 38

WEIGHT OF THE PLASMA IS INSUFFICIENT

Cause 1: Weight of the separated plasma is insufficient.

Solution: Check that the cannulas are broken and that there are no obstructions, and then try again.

Error 39

WEIGHT OF THE BUFFY COAT IS INSUFFICIENT

Cause 1: The volume of buffy coat collected is less than what was set.

Solution: Check the actual weight.

Cause 2: The scale is not calibrated.

Solution: Calibrate the scale (see scale calibrations).

Error 40

WEIGHT OF THE SAG-M IS INSUFFICIENT

Cause 1: The volume of SAG-M is insufficient to complete the procedure.

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Cause 2: This scale is not calibrated.

Solution: Calibrate the scale (see the 'Calibration of Scales' section).

Error 41

CHECK THE SAG-M FLOW

Cause 1: The SAG-M tube is kinked or obstructed.

Solution: Place the tube and remove any obstructions.

Cause 2: Clamp 5 is closed.

Solution: Unlock clamp 5.

Cause 3: The scale is not calibrated.

Solution: Calibrate the scale (see the "Calibration of Scales" section).

Call the authorised technical assistance.

Error 42

END OF RUN ACHIEVED WITH SAG-M IS INSUFFICIENT

Cause 1: The SAG-M bag is missing. With the plate completely open, the system did not detect any SAG-M dispensing.

Solution: Place the SAG-M bag.

Cause 2: The fully open sensor is faulty.

Solution: Call the authorised technical assistance.

Error 43

RBC IS NOT DETECTED

Cause 1: The detection threshold value is incorrect.

Solution: Check parameter 8.

Cause 2: The RBC activation distance value is incorrect.

Solution: Check parameter 18.

Cause 3: The RBC volume is insufficient.

Solution: Repeat the procedure.

Cause 4: The HB sensor is faulty.

Solution: Call the authorised technical assistance.

Error 44

CHECK THE FLOW

Cause 1: The filter is blocked.
Solution: Try again with a new tube set.
Cause 2: The tubes are improperly positioned.
Solution: Reposition the tubes.
Cause 3: The scale is not calibrated.
Solution: Calibrate the scale (see the "Calibration of Scales" section).

Error 45

RBC IS DETECTED

Cause 1: RBC is detected.
Solution: Check the actual quantity of supernatant and dispense it manually if it is insufficient. Press the  button to continue or the  button to exit.

Error 46

UNEXPECTED WEIGHT

Cause 1: The scale detected a change due to an external factor.
Solution: Remove any foreign objects from the scale.
Cause 2: The scale is unstable.
Solution: Check that the scale holder is not touching the plate.
Call the authorised technical assistance.

Error 47

NOT USED

Error 48

PROCEDURE WAS INTERRUPTED BY THE OPERATOR

Cause 1: The operator pressed the  button and confirmed the interruption by pressing the  button.
Solution: Repeat the procedure.

Error 50

DISTANCE-WEIGHT ERROR

Cause 1: The plate is in a position where the weight detected should be higher.
Solution: Check for any obstructions that could result in a false reading of the weight. Check the parameter settings for the bag used.

Error 51

WEIGHT-DISTANCE ERROR

Cause 1: The plate is in a position where the weight detected should be lower.
Solution: Check for any obstructions that could result in a false reading of the weight. Check the parameter settings for the bag used.

Error 52

PLATE TIMEOUT

Cause 1: The distance detected has not changed in the past 18 seconds.
Solution: Check for any obstructions in the tubes and for the proper breaking of the cannulas.

Error 53 (1-6)

TUB SENSOR BRIGHTNESS IS HIGH

Cause 1: The ambient brightness detected by the tube sensor is too high.

Solution: Decrease the ambient brightness level.
Replace the tube sensor (x).

Error 54

NO BAGS ARE PRESENT

Cause 1: One or more bags were not detected.

Solution: Check whether the bags are properly positioned.

Error 55

TUBE (X) NOT PRESENT

Cause 1: The tube (x) is not in the proper position or the bags were not detected.

Solution: Check the placement of the tube (x).

Error 56 (1-4)

TIMEOUT OF THE CANNULA-BREAKING MOTOR.

Cause 1: A mechanical obstruction is preventing the movement of the system.

Solution: Remove the obstruction.

Cause 2: The sensor of the cannula-breaking system is defective.

Solution: Call the authorised technical assistance.

Error 57

TIMEOUT OF ACTIVATION OF THE CANNULA BREAKING (X).

Cause 1: The cannula-breaking system (x) did not achieve the end run between the proper time intervals.

Solution: Turn ARCHIMEDE off, then on again. Call the authorised technical assistance.

Error 58

SAG-M SENSOR ERROR.

Cause 1: The IR detection system of the SAG-M bag is defective.

Solution: Check whether the ambient light is too high, then clean the red optical protection screen.

Call the authorised technical assistance.

COMMERCIAL AND REGULATORY DATA

TRADE NAME

Name and Product Code: ARCHIMEDE = MOE1504TER

MANUFACTURER

DELCON S.r.l.

Via Zanica 19/F/1

24050 Grassobbio (BG) - Italy

REGULATORY

Protection class 1 (according to IEC 60601-1)

European regulation MDR 2017/745

Reference standard IEC 60601-1 3rd Ed. Medical device safety

IEC 60601-1-2 4th Ed. Electromagnetic compatibility (EMC)

ACCESSORIES

ARCHIMEDE may include the following optional accessories:

- SAG-M mechanical press assembly with breaker
- RFID kit
- Omnidirectional barcode reader
- Customer certificated weight 1000g
- Profile plates N°2, 3, 5
- Splash guard
- Filter holder

And the following service accessories:

- Service certificated weight kit 500g calibration mass ($\pm 0.01\%$)
- Buffy coat press calibrator
- Service tool to keep the cover open
- Service tool to mount force plate
- Force tool

APPENDIX 3.

TECHNICAL DATA

Dimensions	435 x 500 x 425H mm
Weight	~ 48 kg
Power supply	100-240 VAC
Line frequency	50÷60 Hz
Power Consumption	90 VA max
Fuses	2 x T3,15A
Number of protocols:	18
Number of weighing scales	3
Weighing range	0~2000 ml or g
Weighing accuracy	±2 g
Weighing maximum load	5000 g
Number of clamps	6
Number of sealing heads	4
Optoelectronic sensors	2 in each clamp 1 in Haemoglobin sensor 18 in frontal press
Max collection time:	20 min
Internal memory:	Up to 999 records
Export memory:	MMC / SD cartridges and USB memory stick
Display and user interface:	Backlighted graphic display 240 x 128mm: membrane keyboard
Environmental conditions for use	Temperature: +5°÷ +45°C Humidity: 0 ÷ 80% non-condensing Pressure: 500 ÷ 1060 hPa
Environmental conditions for storage and transportation:	Temperature: -20°÷ +70°C Humidity: 20 ÷ 95% non-condensing Pressure: 500 ÷ 1060 hPa

Delcon S.r.l.

via Zanica 19/F/1
24050 Grassobbio (BG) - Italy
(+39) 035 452 23 36

info@delcon.it
service@delcon.it

www.delcon.it

Delcon USA

1704 East Boulevard, Suite 101
Charlotte, NC 28203
(+1) 704 243 9942

info@delcon-usa.com

www.delcon-usa.com