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
IMUGARD III-RC

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ENGLISH

DESCRIPTION / DESIGNATION

- This device is indicated for the removal of leukocytes from a single unit of packed red blood cells or whole blood.
- Sterilized by ethylene oxide. Sterile in an unopened and undamaged unit package.
- Non-pyrogenic fluid path.

PRECAUTIONS

- This device contains DEHP. Based on animal studies, significant exposure to DEHP may interfere with the normal development of the male reproductive tract. For children, pregnant and nursing women, alternative devices may be appropriate.
- Do not use if the unit package or the product has been damaged or soiled.
- Dispose of safely after single use to avoid risk of infection.
- Do not store at extreme temperatures and humidity. Avoid direct sunlight.
- Read the following instructions carefully before using.
- For single use only. Do not reuse. Do not resterilize. Do not reprocess. Reprocessing may compromise the sterility, biocompatibility and functional integrity of the device.

INSTRUCTIONS FOR USE

CAUTION • Strict aseptic technique should be observed at all times.

Attach blood bag

- Close the RED "C" clamp.
- Move the roller clamp to just below the drip chamber and close it.
- Confirm that both air vents are closed.
- Open the peel tab of the outlet port on the unit to be filtered.
- Remove the protector from the spike located at the inlet end of the filter tubing and insert the spike into the outlet port on the unit to be filtered.

CAUTION Do not hold the tubing part of the spike while inserting it into the outlet port. Instead, hold the hub part of the spike firmly. This will prevent loosening of tubing part from the hub resulting in leakage. (Fig. 1)

- Hang the unit of blood to be filtered on an I.V. stand.
- Attach the venipuncture needle to the luer lock connector located at the end of the tubing (leave the needle cap on the needle).

Prime filter

- To start filtration, open the RED "C" clamp which is located just below the spike going into the blood bag, and open the GREEN air vent which is located on the filter housing. (Fig. 2)
- Let the blood fill the filter housing, and when it reaches the GREEN air vent located at the top of the housing, close the GREEN air vent firmly.
- Invert the drip chamber and open the roller clamp. (Fig. 2A)
- Once the drip chamber has filled half way, return it to the upright position. (Fig. 2B)
- Continue the priming process until the air in the tubing below the drip chamber and the air in the needle has been removed.
- Close the roller clamp.

Transfuse blood

- Perform the venipuncture.
- Open the roller clamp. Use the roller clamp to adjust the flow rate (the filter is designed to deliver approximately 15 drops per ml).
- After the blood bag containing the unfiltered blood has emptied completely, open the air vent which is located just above the filter on the inlet tubing. This will allow the blood remaining in the inlet tubing to drain into the filter.
- Once the inlet tubing has emptied of blood, open the GREEN air vent which is located on the filter housing. This will allow the blood remaining in the filter, the outlet section of filter tubing, and in the drip chamber to be transfused. (Fig. 3)
- After the drip chamber has emptied, close the roller clamp.
- Close the RED "C" clamp and both air vents.
- Remove the needle from the patient.

Fig.1
Euk. 1

Hub, Partie Centrale, Hub, Nεφρά διέξυ (αυ σπικην), Partie Central, Tubing part, Partie tubulaire, Schlauchteil, Tubular, Tubo, Slangegedeelte, Slangdelen, Τμήμα του αούλίου, Filter housing, Chambre du filtre, Filteergehuse, Alojamiento del filtro, Corpo del filtro, Filter behuizing, Φίλτρον, Filtrationskammer, Cámaras de goeteo, Drip chamber, Chambre de goutte à goutte, Tropfkammer, Cámara de goteo, Camera di gocciolamento, Druppelkamer, Στραγγιστικός θαλάμος, Roller clamp, Clamp a rullo, Rolliem, Rolklammer, Pinza de rodillo, Καλιφόρας ρυθμιστός

Fig.2
Euk. 2

Prime filter, Amorrage, Fullen, Cebado, Priming, Vullen, Fyllning, Fyllning, Fyllning, Fyllning, Fyllning

Fig.3
Euk. 3

Blood recovery, Restitucion du sang, Blutrückführung, Recuperacion de la sangre, Recuperar sangue, Opvangen van het bloed, Bloedruyete, Ανάκτηση αίματος

Fig.2A
Euk. 2A

Fig.2B
Euk. 2B

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© TERUMO CORPORATION May 2013 13868