

Insulin


Our full line of insulin syringes offers patients with diabetes comfort and accuracy for their daily self-injections.

- Unsurpassed needle sharpness and smooth plunger movement help make injections more comfortable.
- Bold, clear scale markings help ensure accurate dosing.
- A wide range of sizes.

Insulin Syringes

Terumo Standard Insulin Syringes

Smooth movement of plunger and minimized dead space alleviates patient discomfort. Semi-transparent material and clear prints on the syringe allow users to inject exact amount of insulin needed. Extra thin wall technology for maximum flow rate. Reduced injection time provides virtually pain free experience. The needle surface is siliconized to minimized penetration resistance. Barrel, plunger and cap are made with PP, Latex free. Produced using SUS 304 stainless steel tube. Comply with ISO 8537 standard.



PERMANENTLY ATTACHED NEEDLE – virtually eliminates needle separation and insulin leakage.

SYNTHETIC GASKET – provides smoother plunger action so that less injection pressure is necessary. Gasket will not react with insulin when prefilling your syringe.

CLEARLY CALIBRATED SCALE MARKINGS – for precise dosage measurements.

CLEAR BARREL – enhances visibility of insulin level for accurate dosage control.

NEEDLE SHIELD – color-matched to U-100 insulin concentration, snaps on for a secure fit...very important when you are prefilling your syringe.

31G NEEDLE – the thinnest needle available with Terumo sharpness for improved injection comfort

TAPERED HUB, FLUSH-FIT PLNGER GASKET – eliminates dead space and facilitates removal of air bubbles. Measure dose on the front gasket ring

Syringe specification

Syringe type	Barren Length (mm)	Outer Diameter (mm)	Inner Diameter (mm)	Cap Length (mm)
1.0mL	81.4	6.73	4.5	26
0.5mL	75.1	5.50	3.36	26

Needle Specification

Gauge (G)	Outer Diameter (mm)	Inner Diameter (Min: mm)	Length (mm)
31	0.254 - 0.267	0.114	8mm (5/16")
30	0.298 - 0.320	0.113	

Product code

SS*10M3108KS
 SS*10M3008KS
 SS*05M3108KS
 SS*05M3008KS

Descriptions

31G x 8mm (5/16") Graduations: 0.02cc, 2 units of insulin
 30G x 8mm (5/16") Graduations: 0.02cc, 2 units of insulin
 31G x 8mm (5/16") Graduations: 0.01cc, 1 unit of insulin
 30G x 8mm (5/16") Graduations: 0.01cc, 1 unit of insulin

Packed: 100/box
Sterilized by
Ethylene oxide

Made in: Korea



Terumo Vietnam Medical Equipment Co., Ltd.
 Head office
 14th & 16th Floor, Geleximco Building, 36 Hoang
 Cau Street, O Cho Dua Ward, Dong Da District
 Hanoi, Vietnam
 Tel: +84 24 3936 1643
 Fax: +84 24 3936 1641

Terumo Vietnam Medical Equipment Co., Ltd.
 HCMC Branch office
 Room 505, Hoang Anh Safomec Building, 7/1
 Thanh Thai, Ward 14, District 10
 HCM City, Vietnam
 Tel: +84 28 3866 9263

TVME-2310-003