


SAFEED[®] Extension Tube



A wide product line for infusion line extension applications

Luer lock connectors and Luer slip connectors are available in 10 different lengths and diameters. SAFEED extension tubes can be connected to an infusion set or 3-way-stop-cock; they can also be used for medicine infusion from a syringe pump.

Available sizes

Item code	Tip	Capacity	Tubing length	Inner diameter	Outer diameter
SF•ET2022	Luer	2.0mL	25cm	3.1mm	5.0mm
SF•ET2022L	Luer lock	2.0mL	25cm	3.1mm	5.0mm
SF•ET0222L22	Luer	0.25mL	25cm	1.1mm	2.6mm
SF•ET1725	Luer	1.7mL	50cm	2.1mm	3.8mm
SF•ET3825	Luer	3.8mL	50cm	3.1mm	5.0mm
SF•ET2527	Luer	2.5mL	75cm	1.1mm	2.6mm
SF•ET2527L	Luer lock	2.5mL	75cm	1.1mm	2.6mm
SF•ET5527L	Luer lock	5.5mL	75cm	3.1mm	5.0mm
SF•ET0727L22	Luer lock	0.75mL	75cm	1.1mm	2.6mm
SF•ET152EL22	Luer lock	1.5mL	150cm	1.1mm	2.6mm

Tubing material: PVC Sterilization: EOG Unit (pcs.): Unit box: 50 pcs. Unit (pcs.): Shipping carton: 400 pcs.



Terumo Corporation
2-44-1 Hatagaya Shibuya-ku Tokyo 151-0072 Japan
<https://www.terumo.co.jp/>

Terumo Singapore Pte. Ltd. (Yangon Branch)
Room no. 1220, 12th Floor, Myanmar Centre Business Suites, Myanmar Centre Tower 2, 192 Kaba Aye Pagoda Road, Bahan Township, Yangon, Myanmar
TEL: +95-1-934-5065 (Ext.1228)

Terumo Vietnam Medical Equipment Co., Ltd. Head Office
14th and 16th Floor, Geleximco Building, 36 Hoang Cau Street, O Cho Dua Ward, Dong Da District, Hanoi Vietnam.
TEL: +84-24-3936-1643 / +84-24-3936-1644
FAX: +84-24-3936-1641
<https://terumo.com.vn/en/>

PT. Terumo Indonesia
Wisma KEIAI 5th Floor JL. Jend. Sudirman KAV. 3 Jakarta 10220 Indonesia M

Terumo (Thailand) Co., Ltd.
No.88 The PARQ Building, Unit 8W9-16, 8th Floor, Ratchadaphisek Road, Klongtoey Subdistrict, Klongtoey District, Bangkok 10110 Thailand MAP
TEL: +66 2111 2133

Distributor

Terumo Vietnam Medical Equipment Co., Ltd.

Head office
14th & 16th Floor, Geleximco Building, 36 Hoang Cau Street, O Cho Dua Ward, Dong Da District Hanoi, Vietnam
Tel: +84 24 3936 1643
Fax: +84 24 3936 1641

Terumo Vietnam Medical Equipment Co., Ltd.

HCMC Branch office
Room 505, Hoang Anh Safomec Building, 7/1 Thanh Thai, Ward 14, District 10, HCM City, Vietnam
Tel: +84 28 3866 9263

Terumo Marketing Philippines Inc. 3203 West Tower, Pse Centre Exchange Road, Ortigas Center Pasig City Philippines

Terumo Australia Pty Limited
Level 4, Building A 5 Talavera Road, Macquarie Park NSW 2113 Australia

Essen Technology (Beijing) Co., Ltd.
Building 1, NO.A12 Tianzhu West Road, Tianzhu Airport Industrial Zone A, Shunyi District, Beijing 101318, PR China

Terumo Korea Corporation
23rd Fl., GT Tower, 411 Seocho-daero, Seocho-gu Seoul, 06615 Korea

Terumo China (Hong Kong) Ltd. 5/F 80 Gloucester Road Wanchai Hong Kong

Terumo Medical (Shanghai) Co., Ltd.
Rm.901, Zhongshan Expo Plaza 666 Huaihai Road West, Changning Shanghai 200052 China

Terumo Medical Products (Hangzhou) Co., Ltd.
M4-9-5, Economic & Technological Development Zone Hangzhou 310018 China