

# Kim lều tĩnh mạch **SURFLO**<sup>®</sup> I.V. CATHETER

*Kim lều tĩnh mạch ngoại biên **SURFLO** của Terumo đượ thiết kế và sản xuất theo tiêu chuẩn chất lượng toàn cầu nghiêm ngặt, bao gồm tiệt trùng bằng chùm tia điện tử và tiêu chuẩn quản lý chất lượng ISO/hệ thống mã màu đặc biệt, giúp sử dụng dễ dàng, chắc chắn và an toàn hơn.*




SURFLO đượ biết đến nhờ **ĐỘ SẮC BÉN**

**Mới**

### **Đầu catheter dạng SL<sup>(\*)</sup>**

Kim lều tĩnh mạch SURFLO đượ làm từ chất liệu ethylene tetrafluorethylene (**ETFE**), đem lại sức cản thấp khi thâm nhập. Đầu ca-tê-te thiết kế ôm chặt thân kim giúp thâm nhập qua da dễ dàng, giảm sức cản thấp khi thâm nhập, giảm đau hơn.

(\*)SL=dạng thuôn, mượt và sức cản thấp



**Mới**

### **Tiệt trùng bằng chùm tia điện tử**

Tiệt trùng bằng chùm tia điện tử đã đượ chấp nhận rộng rãi như là tiêu chuẩn của ngành công nghiệp bởi vì đạt đượ tính hiệu quả và thân thiện môi trường.

### **Sản xuất theo tiêu chuẩn toàn cầu**

Kim lều tĩnh mạch SURFLO đượ sản xuất theo tiêu chuẩn ISO 9002 đạt chứng chỉ hệ thống quản lý chất lượng và bảng mã màu theo tiêu chuẩn ISO giúp xác định nhanh chóng và chính xác.



## Dòng sản phẩm

Cỡ Catheter	Mã màu	Mã sản phẩm	Độ dài catheter	Đường kính ngoài	Đường kính trong	Cỡ kim	Lưu lượng mL/phút
18G	Xanh lá cây	SR+OT1832C	1 1/4" (32mm)	1.30mm	0.95mm	20G	100
		SR+OT1851C	2" (51mm)				90
20G	Hồng	SR+OT2032C	1 1/4" (32mm)	1.10mm	0.80mm	22G	60
22G	Xanh da trời	SR+OT2225C	1" (25mm)	0.85mm	0.60mm	24G	35
24G	Vàng	SR+OT2419C	3/4" (19mm)	0.67mm	0.47mm	27G	15
14G	Cam	SR+OX1451C	2" (51mm)	2.17mm	1.73mm	16G	305
		SR+OX1464C	2 1/2" (64mm)				300
16G	Xám	SR+OX1632C	1 1/4" (32mm)	1.70mm	1.30mm	18G	200
		SR+OX1651C	2" (51mm)				190
		SR+OX1664C	2 1/2" (64mm)				180
18G	Xanh lá cây	SR+OX1832C	1 1/4" (32mm)	1.30mm	0.95mm	20G	100
		SR+OX1851C	2" (51mm)				90
		SR+OX1864C	2 1/2" (64mm)				85
20G	Hồng	SR+OX2025C	1" (25mm)	1.10mm	0.80mm	22G	65
		SR+OX2032C	1 1/4" (32mm)				60
		SR+OX2051C	2" (51mm)				55
22G	Xanh da trời	SR+OX2225C	1" (25mm)	0.85mm	0.60mm	24G	35
24G	Vàng	SR+OX2419C	3/4" (19mm)	0.67mm	0.47mm	27G	15

Số công bố: 20000318/PCBB-BYT  
 Chủ sở hữu số công bố: Công ty TNHH Thiết bị Y tế Terumo Việt Nam  
 Tầng 14, 16, tòa nhà GELEXIMCO, 36 Hoàng Cầu, P. Ô Chợ Dừa, Q. Đống Đa, Hà Nội, Việt Nam  
 Hãng sản xuất: Terumo (Philippines) Corporation  
 Sản xuất tại: Philippines

Số văn bản công khai hình thức và nội dung quảng cáo: 009/23-CKQC-Terumo



### Công ty TNHH Thiết bị y tế Terumo Việt Nam

Distributor

Trụ sở  
 Tầng 14 và 16, tòa nhà Geleximco, Số 36 Hoàng Cầu, P. Ô Chợ Dừa, Q. Đống Đa, TP. Hà Nội, Việt Nam Tel: +84 24 3936 1643  
 Fax: +84 28 3936 1641

### Chi nhánh Công ty TNHH Thiết bị y tế Terumo Việt Nam tại TP. Hồ Chí Minh

Phòng 505, tòa nhà Hoàng Anh Safomec, 7/1 Thành Thái, Phường 14, Quận 10, TP. Hồ Chí Minh, Việt Nam Tel: +84 28 3866 9263



### TERUMO CORPORATION

44-1, 2-chome, Hatagaya, Shibuya-ku, Tokyo 151-0072, Japan  
 Phone: 81-3-3374-8111 Fax: 81-3-3374-8198

**TERUMO (PHILIPPINES) CORPORATION**  
 124 East Man Avenue, Laguna Technopark Bldg., Laguna, Philippines  
 Phone: 63-49-541-2111 Fax: 63-49-541-2121

**TERUMO MARKETING PHILIPPINES, INC.**  
 Unit W-1508 The Philippines Stock Exchange Centre West Tower, Exchange Road Ortigas Center, Pasig City, Philippines  
 Phone: 63-2-632-1674 Fax: 63-2-632-7966

Other offices: •New Jersey •Miami •Mexico •São Paulo •Brussels •Rome •Madrid •London •Göteborg •Vienna •Paris •Frankfurt

**TERUMO CORPORATION Singapore Branch**  
 300 Beach Road, #32-00 The Concourse, Singapore 199566  
 Phone: 65-6-291-3603 Fax: 65-6-291-2696

**TERUMO (THAILAND) CO., LTD.**  
 #1206, 12th Floor 54 B.B. Building Asok Road (Sukhumvit 21), Klongtoey, Bangkok 10110, Thailand  
 Phone: 66-2-260-7020 Fax: 66-2-260-7019

**TERUMO (THAILAND) CO., LTD. Hanoi Representative Office**  
 International Centre 17 Ngo Duyen str., Unit 05, 8th Floor Hanoi, Vietnam  
 Phone: 84-4-936-1643 Fax: 84-4-936-1841

**P. T. TERUMO INDONESIA**  
 Wisma Kyoo Prince 11th Floor J.L. Jend. Sudirman Kav. 3, Jakarta 10220, Indonesia  
 Phone: 62-21-672-0271 Fax: 62-21-672-4072

**TERUMO CORPORATION Kuala Lumpur Branch**  
 Kelana Business Centre, Unit 301 Block A, No. 57 Jalan SS7/2, Kelana Jaya, 47301 Petaling Jaya, Selangor Darul Ehsan, Malaysia  
 Phone: 60-3-7492-9688 Fax: 60-3-7492-9689

**TERUMO CORPORATION Australian Branch**  
 Level 4, Bld B 11 Talavera Rd, Macquarie Park, New South Wales 2113, Australia  
 Phone: 61-2-9678-6122 Fax: 61-2-9678-5085

**TERUMO CORPORATION Taipei Branch**  
 7A-1, No. 170 Tun-Hua North Road Taipei, Taiwan, R.O.C.  
 Phone: 886-2-2545-1250 Fax: 886-2-2545-1251

**TERUMO KOREA CORPORATION**  
 1st Fl. Kyungnam Bldg., 833-48, Yeoksam-dong, Kangnam-gu, Seoul 135-836, Korea  
 Phone: 82-2-585-8226 Fax: 82-2-585-9224

**TERUMO CORPORATION Hong Kong Branch**  
 13/F, King's Commercial Centre, 25 King's Road, North Point, Hong Kong  
 Phone: 852-2866-0811 Fax: 852-2829-0451

**TERUMO CORPORATION Shanghai Office**  
 Room 81, Wan Zhong Commercial Building, 1300 Yan An Road (West), Shanghai 200050, China  
 Phone: 86-21-62406593 Fax: 86-21-62117257

**TERUMO CORPORATION Beijing Office**  
 Room 411, Building B, Jin Fu Commercial Plaza No. 13, Eastern Beijingfang Fu Wai Street, Xicheng District, Beijing 100037, China  
 Phone: 86-10-68360571 Fax: 86-10-68362150

**TERUMO CORPORATION Guangzhou Office**  
 A511 Floor, Chin Lin Plaza, No. 13-35 Pan Fu Road, Guangzhou 510180, China  
 Phone: 86-20-81366681 Fax: 86-20-81366145

**TERUMO CORPORATION Chennai Branch**  
 Alexander Square, 4th Floor, 34 & 35 Sarda Patel Road, Gurgaon, Chennai 600 052, India  
 Phone: 91-44-22300634 Fax: 91-44-22300622

**TERUMO CORPORATION Dubai Branch**  
 Arab Tower, Off. No 1308, P.O. Box 20291 Dubai, U.A.E.  
 Phone: 971-4-2212220 Fax: 971-4-2213300