



SURFLO®

I.V. CATHETER

Terumo SURFLO® I.V. CATHETER is designed and manufactured under stringent global quality standards, including electron beam sterilization and ISO quality management/color code systems, for easier, surer and safer use.




SURFLO is known for its **SHARPNESS**.

NEW

SL-type* catheter tip

The SURFLO I.V. CATHETER is made of ethylene tetrafluorethylene (ETFE), which always provides the appropriate penetration resistance. The tip design offers smooth insertion, reducing penetration resistance for lessened pain.

*SL = Smooth tapering and Low resistance



NEW

Electron Beam Sterilization

Electron Beam Sterilization has gained widespread acceptance as the industry standard because it is both effective and environmentally friendly.

Manufacturing under global standards

The SURFLO I.V. CATHETER is manufactured under the ISO 9002 certified quality management system and color-coded in accordance with ISO standard for quick and accurate identification.



Product range

| Catheter gauge | Color code | Product code | Catheter length | Catheter external dia. | Catheter internal dia. | Needle gauge | Flow rate MI/min |
|----------------|-------------|--------------|-----------------|------------------------|------------------------|--------------|------------------|
| 18G | Deep Green | SR+OT1832C | 1 1/4" (32mm) | 1.30mm | 0.95mm | 20G | 100 |
| | | SR+OT1851C | 2" (51mm) | | | | 90 |
| 20G | Pink | SR+OT2032C | 1 1/4" (32mm) | 1.10mm | 0.80mm | 22G | 60 |
| 22G | Deep Blue | SR+OT2225C | 1" (25mm) | 0.85mm | 0.60mm | 24G | 35 |
| 24G | Yellow | SR+OT2419C | 3/4" (19mm) | 0.67mm | 0.47mm | 27G | 15 |
| 14G | Orange | SR+OX1451C | 2" (51mm) | 2.17mm | 1.73mm | 16G | 305 |
| | | SR+OX1464C | 2 1/2" (64mm) | | | | 300 |
| 16G | Medium Grey | SR+OX1632C | 1 1/4" (32mm) | 1.70mm | 1.30mm | 18G | 200 |
| | | SR+OX1651C | 2" (51mm) | | | | 190 |
| | | SR+OX1664C | 2 1/2" (64mm) | | | | 180 |
| 18G | Deep Green | SR+OX1832C | 1 1/4" (32mm) | 1.30mm | 0.95mm | 20G | 100 |
| | | SR+OX1851C | 2" (51mm) | | | | 90 |
| | | SR+OX1864C | 2 1/2" (64mm) | | | | 85 |
| 20G | Pink | SR+OX2025C | 1" (25mm) | 1.10mm | 0.80mm | 22G | 65 |
| | | SR+OX2032C | 1 1/4" (32mm) | | | | 60 |
| | | SR+OX2051C | 2" (51mm) | | | | 55 |
| 22G | Deep Blue | SR+OX2225C | 1" (25mm) | 0.85mm | 0.60mm | 24G | 35 |
| 24G | Yellow | SR+OX2419C | 3/4" (19mm) | 0.67mm | 0.47mm | 27G | 15 |



TERUMO CORPORATION

44-1, 2-chome, Hatagaya, Shibuya-ku, Tokyo 151-0072, Japan
Phone: 81-3-3374-8111 Fax: 81-3-3374-8196

TERUMO (PHILIPPINES) CORPORATION

124 East Main Avenue, Laguna Technopark Bldg., Laguna, Philippines
Phone: 63-49-541-2111 Fax: 63-49-541-2121

TERUMO MARKETING PHILIPPINES, INC.

Unit W-1503 The Philippines Stock Exchange Centre West Tower, Exchange Road Ortigas Center, Pasig City, Philippines
Phone: 63-2-632-1674 Fax: 63-2-632-7950

Other offices: •New Jersey •Miami •Mexico •Sao Paulo •Brussels •Rome •Madrid •London •Göteborg •Vienna •Paris •Frankfurt

® Registered Trademark
Manufactured by TERUMO (PHILIPPINES) CORPORATION.
© TERUMO CORPORATION, June 2003. Printed in Japan.

TERUMO CORPORATION Singapore Branch

300 Beach Road, #32-06 The Concourse, Singapore 199555
Phone: 65-6-291-3903 Fax: 65-6-291-2686

TERUMO (THAILAND) CO., LTD.

#1206, 12th Floor 54 B.S. Building Asok Road (Bukhumi 21), Kwaeng Klongtoeyus, Khet Wattana, Bangkok 10110, Thailand
Phone: 66-2-260-7020 Fax: 66-2-260-7019

TERUMO (THAILAND) CO., LTD. Hanoi Representative Office

International Centre 17 Ngo Duyen str., Unit 05, 6th Floor Hanoi, Vietnam
Phone: 84-4-936-1643 Fax: 84-4-936-1641

P. T. TERUMO INDONESIA

Wisma Kyosai Prince 11th Floor J.L. Jend. Sudirman Kav. 3, Jakarta 10220, Indonesia
Phone: 62-21-572-6071 Fax: 62-21-572-4072

TERUMO CORPORATION Kuala Lumpur Branch

Kolomo Business Centre, Unit 301 Block A, No. 97, Jalan SS17/2, Kelana Jaya, 47301 Petaling Jaya, Selangor Darul Ehsan, Malaysia
Phone: 60-3-7492-9888 Fax: 60-3-7492-9669

Distributor

Terumo Vietnam Medical Equipment Co., Ltd.

Head office

14th & 16th Floor, Geleximco Building, 36 Hoang Cau Street, O Cho Dua Ward, Dong Da District Hanoi, Vietnam

Tel: +84 24 3936 1643

Fax: +84 24 3936 1641

Terumo Vietnam Medical Equipment Co., Ltd. HCMC Branch office

Room 505, Hoang Anh Safomec Building, 7/1 Thanh Thai, Ward 14, District 10, HCM City, Vietnam

Tel: +84 28 3866 9263

TERUMO CORPORATION Australian Branch

Level 4, Bld B 11 Talavera Rd, Macquarie Park, New South Wales 2113, Australia
Phone: 61-2-9878-5122 Fax: 61-2-9878-5085

TERUMO CORPORATION Taipei Branch

7A-1, No. 170 Tun-Hua North Road Taipei, Taiwan, R.O.C.
Phone: 886-2-2545-1250 Fax: 886-2-2545-1251

TERUMO KOREA CORPORATION

151 P. Kyungnam Bldg., 830-45, Yeoksam-dong, Kangnam-gu, Seoul 135-836, Korea
Phone: 82-2-965-9225 Fax: 82-2-965-9224

TERUMO CORPORATION Hong Kong Branch

13/F, King's Commercial Centre, 25 King's Road, North Point, Hong Kong
Phone: 852-2666-0811 Fax: 852-2629-0451

TERUMO CORPORATION Shanghai Office

Room 8J, Wan Zhong Commercial Building, 1303 Yan An Road (West), Shanghai 200060, China
Phone: 86-21-62406590 Fax: 86-21-62117257

TERUMO CORPORATION Beijing Office

Room 411, Building B, Jin Fu Commercial Plaza No. 13, Eastern Beiyangfang Fu Wai Street, Xi Cheng District, Beijing 100037, China
Phone: 86-10-68360574 Fax: 86-10-68362150

TERUMO CORPORATION Guangzhou Office

A511 Room, QinLin Plaza, No. 13-35 Pan Fu Road, Guangzhou 510180 China
Phone: 86-20-81399681 Fax: 86-20-81386145

TERUMO CORPORATION Chennai Branch

Alexander Square, 4th Floor, 34 & 35 Sardar Patel Road, Gundy, Chennai 600 052, India
Phone: 91-44-22300634 Fax: 91-44-22300622

TERUMO CORPORATION Dubai Branch

Arbif Tower, Off. No.1306, P.O. Box 20291 Dubai, U.A.E.
Phone: 971-4-2212220 Fax: 971-4-2213330